

PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

| STANDARD OPERATING PROCEDURE | | | |
|------------------------------------|------------------------|--|--|
| Department: Quality Control | SOP No.: | | |
| Title: Validation Of Calculator | Effective Date: | | |
| Supersedes: Nil | Review Date: | | |
| Issue Date: | Page No.: | | |

1.0 OBJECTIVE:

To lay down a procedure for validation of Digital Calculator.

2.0 SCOPE:

This SOP is applicable to all digital calculator use in Quality Control Department at Macleods Pharmaceuticals Limited, Nalagarh.

3.0 RESPONSIBILITY:

Officer, Executive – Quality Control Department

Head – Quality Control Department

4.0 DEFINITION (S):

NA

5.0 PROCEDURE:

5.1 **Identification to Calculator**

5.1.1 Calculators existing in department shall be identified by specific alpha numeric pattern as given

EQ/QCD/SC/01...02...03

5.1.2 Ensure that all the digit of a calculator is visible and printed matter on buttons are clearly legible.



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

| STANDARD OPERATING PROCEDURE | | | |
|------------------------------------|------------------------|--|--|
| Department: Quality Control | SOP No.: | | |
| Title: Validation Of Calculator | Effective Date: | | |
| Supersedes: Nil | Review Date: | | |
| Issue Date: | Page No.: | | |

5.2 **Procedure for calculator Validation**

5.2.1 Check the accuracy of the calculator by carrying out the following set of calculation.

| Set No. | Data | | | | Mean | Sample Standard Deviation | Sample Variance | |
|------------|-------|-------|-------|-------|-------|---------------------------------|----------------------|-------------------------|
| 1 | 8 | 9 | 10 | 11 | 12 | 10 | 1.581138 | 2.50 |
| 2 | 0.8 | 0.9 | 1.0 | 1.1 | 1.2 | 1.0 | S x 10 ⁻¹ | 2.50 x 10 ⁻² |
| 3 | 0.08 | 0.09 | 0.10 | 0.11 | 0.12 | 0.10 | S x 10 ⁻² | 2.50 x 10 ⁻⁴ |
| 4 | 0.008 | 0.009 | 0.010 | 0.011 | 0.012 | 0.010 | S x 10 ⁻³ | 2.50 x 10 ⁻⁶ |
| 5 | 78 | 79 | 80 | 81 | 82 | 80 | S | 2.50 |
| 6 | 708 | 709 | 710 | 711 | 712 | 710 | S | 2.50 |
| 7 | 7008 | 7009 | 7010 | 7011 | 7012 | 7010 | S | 2.50 |

- 5.2.2 Record the observation as per Annexure I and match with the above figures.
- 5.2.3 Carry out the similar calculation by another chemist and record the observation.
- 5.2.4 Record the observation made by both the chemist. The acceptance criteria for accepting the calculator for laboratory calculation are that "No difference up to 6th digit after decimal should be observed".

6.0 **ABBREVIATION** (S):

QCD – Quality Control Department

SOP - Standard Operating Procedure

CAL - Calculator

S – Standard deviation

7.0 REFERENCE (S):

Geenfield and Diday, The statistician, 29: 33-35 (1980)



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

| STANDARD OPERATING PROCEDURE | | | |
|------------------------------------|------------------------|--|--|
| Department: Quality Control | SOP No.: | | |
| Title: Validation Of Calculator | Effective Date: | | |
| Supersedes: Nil | Review Date: | | |
| Issue Date: | Page No.: | | |

8.0 ANNEXURE (S):

Annexure I - Validation Datasheet for Calculator

9.0 REVISION CARD:

| S.No. | Revision no. | Revision date | Details of revision | Reason (s) for revision |
|-------|--------------|---------------|---------------------|-------------------------|
| | | | | |