



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

**RISK ASSESSMENT
REPORT BY FMEA**

| | |
|-----------------------------------|---|
| Product/System/Equipment | AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE |
| Risk Assessment Report No. | |
| Report Date | |



PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

TABLE OF CONTENTS

| S.No. | Description | Page No. |
|--------------|---|-----------------|
| 1.0 | Introduction | 4 |
| 2.0 | Objective | 4 |
| 3.0 | Scope | 4 |
| 4.0 | Risk Assessment Approach | 4 |
| 5.0 | Responsibility | 4 |
| 6.0 | Reference Documents | 5 |
| 7.0 | Risk Ranking Parameters | 5-7 |
| 8.0 | Acceptance Criteria for risk assessment by FMEA | 7 |
| 9.0 | Risk assessment as per FMEA | 8-11 |
| 9.1 | Review of Risk assessment as per FMEA after action taken. | 12 |
| 10.0 | Risk Control Measures | 13 |
| 11.0 | Summary and Conclusion Report for Risk Assessment | 14 |
| 12.0 | Final Report Approval | 15 |



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

DOCUMENT APPROVAL:

This risk analysis study for the preapproval of report by following:

| Responsibility | Department | Name | Signature | Date |
|----------------|-------------------|------|-----------|------|
| Prepared by | Quality assurance | | | |
| Reviewed by | Production | | | |
| | Quality control | | | |
| | Engineering | | | |
| | Store | | | |
| | Quality assurance | | | |
| Approved by | Head-QA | | | |



PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

1.0 Introduction

The “automatic air jet bottle air and vacuum cleaning machine” is inbuilt with turn table for smooth transfer of container to the cleaning section. The automatic air jet bottle air and vacuum cleaning machine works the cleaning of bottle started through air and transferred of cleaned bottle through conveyor for filling material with assurance of product safety.

2.0 Objective

Objective of this report is to assess the risk associated with the equipment Automatic air jet bottle air and vacuum cleaning machine in pre assessment in the manufacturing facility of Cepha Block at, in line with the guidance of the Risk Management manual of and ICH Q9.

3.0 Scope

The scope of this document is limited to the design, installation, operation, performance and safety of equipment “Automatic air jet bottle air and vacuum cleaning machine” and define its failure mode at pre assessment in the manufacturing facility of Cepha Block at

4.0 Risk assessment approach

Risk assessment is carried out as per FMEA (Failure mode effects analysis) method.

5.0 Responsibility

Quality Assurance
Engineering
Production
Quality Control
Store

6.0 Reference Documents

1. ICH Q9-Quality Risk Management
2. guidance on Risk assessment.



PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

Background

..... is intended to start manufacturing of powder solid oral facility at Risk assessment is a part of corporate quality assurance. Pre Quality Risk assessment of “Automatic air jet bottle air and vacuum cleaning machine” is done to check the system is capable of providing quality product throughout the life cycle of the drug product.

7.0 RISK RANKING PARAMETERS

7.1 Rating parameters for Severity

| Effect | Scale | Description |
|-------------|-------|--|
| No effect | 1 | No effect on output |
| Very slight | 2 | Customer not annoyed |
| Slight | 3 | Slight |
| Minor | 4 | Minor effect on performance |
| Moderate | 5 | Moderate effect on performance |
| Significant | 6 | Partial failure but operable |
| Major | 7 | Product performance severely affected, but some operability and safe |
| Extreme | 8 | Very dissatisfied, product inoperable but safe |
| Serious | 9 | Potentially hazardous effect, time-dependent failure |
| Hazardous | 10 | Hazardous effect, safety related sudden failure |

7.2 Rating parameters for Occurrence

| Occurrence | Scale | Description |
|-----------------|-------|---|
| Almost never | 1 | Failure unlikely; history shows no failures |
| Remote | 2 | Rare number of historical failure |
| Very Slight | 3 | Very few failures likely |
| Slight | 4 | Few failures likely |
| Low | 5 | Occasional number of failures likely |
| Medium | 6 | Medium number of failures likely |
| Moderately High | 7 | Moderately high number of failures likely |



PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

| Occurrence | Scale | Description |
|----------------|-------|-------------------------------------|
| High | 8 | High number of failures likely |
| Very High | 9 | Very high number of failures likely |
| Almost certain | 10 | Failure almost certain |

7.3 Rating parameters for Detection control

| Detection | Scale | Description |
|-----------------|-------|--|
| Almost certain | 1 | Proven detection methods with high reliability |
| Very High | 2 | Proven detection methods available |
| High | 3 | Detection tools have high chance of detecting methods |
| Moderately High | 4 | Almost certain not to detect failure |
| Medium | 5 | Detection tools have moderate chance of detecting defect |
| Low | 6 | Detection tools have a low chance of detecting failure |
| Slight | 7 | Detection tools may not detect failure |
| Very Slight | 8 | Detection tools will probably not detect failure |
| Remote | 9 | Detection tools most likely will not detect failure |
| Impossible | 10 | Failure not detected |

Note: Individual contributory factor for each potential failure mode shall be rated. Other scale parameters may also be selected based on the process.

8.0 ACCEPTANCE CRITERIA FOR RISK ASSESSMENT BY FMEA

Acceptance criteria for FMEA are as follows:

| S.No. | RPN Rating | RPN Category | Action Status |
|-------|------------|--------------|----------------|
| 1. | ≥ 76 | Critical | CAPA Required |
| 2. | 51 to 75 | Major | CAPA Required |
| 3. | 26 to 50 | Moderate | CAPA Required |
| 4. | Up to 25 | Minor | Not applicable |



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

9.0 PRE-RISK ASSESSMENT AS PER FMEA:

Name of facility/Utility/Equipment/Process/Operation: Automatic air jet bottle air and vacuum cleaning machine

| S.No. | Potential Failure Mode | Potential effect (s) of failure | Severity (S) | Potential cause/ Mechanism of failure | Occurrence (O) | Current Control | Detection (D) | RPN (S x O x D) | Recommended action | Responsibility and TCD | Action Results | | | | |
|-------|---|--|--------------|--|----------------|--|---------------|-----------------|--|------------------------------|----------------|----------|------------|-----------|---------|
| | | | | | | | | | | | Action taken | Severity | Occurrence | Detection | New RPN |
| 1 | Required Area (floor, Temperature, RH, Differential pressure) not proper for the Automatic air jet bottle air and vacuum cleaning machine | Area will not be suitable for proper functioning of Equipment. | 6 | No or less clarity of the product requirement and machine functionality. | 3 | Approved layout is in place with dimensions & required environmental condition | 3 | 54 | Care has to be taken during Area Qualification | Engineering, QA, Production, | | | | | |
| 2 | Required parameter not defined in URS. / URS not proper for system | Systems not receive suitable for proper output of quality with all parameter as per specification. Affect the product quality. | 4 | No or less clarity of the product requirement and machine functionality. | 3 | Preparation of URS before procurement of equipment is in place with all pre-specified parameter. | 2 | 24 | Current control measures are adequate | NA | NA | NA | NA | NA | NA |



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

| S.No. | Potential Failure Mode | Potential effect (s) of failure | Severity (S) | Potential cause/ Mechanism of failure | Occurrence (O) | Current Control | Detection (D) | RPN (S x O x D) | Recommended action | Responsibility and TCD | Action Results | | | | |
|-------|---|--|--------------|--|----------------|--|---------------|-----------------|---------------------------------------|------------------------|----------------|----------|------------|-----------|---------|
| | | | | | | | | | | | Action taken | Severity | Occurrence | Detection | New RPN |
| 3 | Required utilities (compressed air, electricity) are not available | Machine will not function as expected. | 7 | No or less clarity of the product requirement and machine functionality with respect to utility requirement. | 2 | URS is in place for system with all predefined requirement of utility like electricity, compressed air. | 1 | 14 | Current control measures are adequate | NA | NA | NA | NA | NA | NA |
| 4 | Wrong machine selection in terms of Dimension, capacity and output. | Installation will be affected if dimension is not considered. Output will also get affected if capacity is not considered. | 6 | No or less clarity of the machine. | 2 | URS is in place for dimension , capacity of the Automatic air jet bottle air and vacuum cleaning machine | 1 | 12 | Current control measures are adequate | NA | NA | NA | NA | NA | NA |



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

| S.No. | Potential Failure Mode | Potential effect (s) of failure | Severity (S) | Potential cause/ Mechanism of failure | Occurrence (O) | Current Control | Detection (D) | RPN (S x O x D) | Recommended action | Responsibility and TCD | Action Results | | | | |
|-------|---|--|--------------|---|----------------|---|---------------|-----------------|---------------------------------------|------------------------|----------------|----------|------------|-----------|---------|
| | | | | | | | | | | | Action taken | Severity | Occurrence | Detection | New RPN |
| 5 | MOC and machine contact parts , not meeting GMP requirement | Not meeting GMP requirements and product get affected. | 7 | No or less clarity of the machine contact part and MOC. | 3 | URS is in place for MOC (contact part should be of SS316 or 316L and non contact parts will be of SS304 and machine contact parts to fulfill GMP requirements. | 1 | 21 | Current control measures are adequate | NA | NA | NA | NA | NA | NA |
| 6 | Equipment not received with the process safety measures. | Accident may happen. | 9 | No or less clarity about equipment safety measures. | 2 | Requirement of Safety measures like Earthing, Emergency stop button, safetyvalve defined in URS. | 1 | 18 | Current control measures are adequate | NA | NA | NA | NA | NA | NA |



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

| S.No. | Potential Failure Mode | Potential effect (s) of failure | Severity (S) | Potential cause/ Mechanism of failure | Occurrence (O) | Current Control | Detection (D) | RPN (S x O x D) | Recommended action | Responsibility and TCD | Action Results | | | | |
|-------|-----------------------------|---------------------------------|--------------|--|----------------|---|---------------|-----------------|---|------------------------|----------------|----------|------------|-----------|---------|
| | | | | | | | | | | | Action taken | Severity | Occurrence | Detection | New RPN |
| 7 | No dedicated air filter bag | Product contamination occurs | 8 | No or less clarity about equipment and product Quality | 2 | Decision is taken to use dedicated filter bag | 2 | 32 | Product dedicated filter bag to be procured | Production | | | | | |



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

Corrective Action: *(an extra sheet can be used if space is insufficient)*

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

(Sign/Date)



PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

11.0 SUMMARY AND CONCLUSION REPORT FOR RISK ASSESSMENT

Summary:

.....

.....

.....

.....

.....

.....

.....

.....

.....

Conclusion:

.....

.....

.....

.....

.....

.....

.....

.....

.....



PRE RISK ANALYSIS FOR AUTOMATIC AIR JET BOTTLE AIR AND VACUUM CLEANING MACHINE

12.0 FINAL REPORT APPROVAL:

The final report shall be signed after identifying all the risks and critical control parameters. All the reports or documents have been attached to the respective report (if applicable).

Signature in the block below indicates that all the control measures taken are documented and have been reviewed and found to be acceptable.

| Department | Name | Designation | Signature | Date |
|-------------------|------|-------------|-----------|------|
| Quality assurance | | | | |
| Production | | | | |
| Quality control | | | | |
| Engineering | | | | |
| Store | | | | |
| Head-QA | | | | |