



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

**PERFORMANCE QUALIFICATION
REPORT
FOR
SUPER HEATED WATER SPRAY
STERILIZER**

| | |
|------------------------------|----------------------|
| EQUIPMENT ID No. | |
| LOCATION | LOADING AREAK |
| DATE OF QUALIFICATION | |
| SUPERSEDES REPORT No. | NIL |



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

REPORT CONTENTS

| S.No. | TITLE | PAGE No. |
|--------------|---|-----------------|
| 1.0 | REPORT PRE APPROVAL | 3 |
| 2.0 | OBJECTIVE | 4 |
| 3.0 | SCOPE | 4 |
| 4.0 | RESPONSIBILITY | 4 |
| 5.0 | EQUIPMENT DETAILS | 5 |
| 6.0 | PQUALIFICATION REQUIREMENTS: | 6-7 |
| 7.0 | TESTS & CHECKS | 8-36 |
| 8.0 | CHECKLIST OF ALL TESTS & CHECKS | 37 |
| 9.0 | DOCUMENTS ATTACHED | 38 |
| 10.0 | NON-COMPLIANCE | 38 |
| 11.0 | DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY | 38 |
| 12.0 | CHANGE CONTROL, IF ANY | 38 |
| 13.0 | REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY) | 39 |
| 14.0 | CONCLUSION | 39 |
| 15.0 | RECOMMENDATION | 39 |
| 16.0 | ABBREVIATION | 40 |
| 17.0 | REPORT POST APPROVAL | 41 |



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

1.0 REPORT PRE APPROVAL:

PREPARED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|--|------|-----------|------|
| OFFICER / EXECUTIVE (QUALITY ASSURANCE) | | | |

REVIEWED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|--|------|-----------|------|
| OPERATING MANAGER (QUALITY ASSURANCE) | | | |
| HEAD (ENGINEERING) | | | |
| HEAD (PRODUCTION) | | | |

APPROVED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|-----------------------------|------|-----------|------|
| HEAD (QUALITY ASSURANCE) | | | |



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

2.0 OBJECTIVE:

- To compile the Validation report carried out as per Protocol for Super Heated Water Spray Sterilizer installed in the LVP Line Loading Area.

3.0 SCOPE:

- The Report covers all aspects of Performance Qualification for the Super Heated Water Spray Sterilizer installed in the LVP Line Loading Area.

4.0 RESPONSIBILITY:

- The Validation Group, comprising of a representative from each of the following Departments, shall be responsible for the overall compliance of this Report:

| DEPARTMENTS | RESPONSIBILITIES |
|--------------------------|---|
| Quality Control | <ul style="list-style-type: none">• Preparation of Reports and submission to Quality Assurance Department.• To conduct Validation activity as per the Approved Protocol.• To provide analytical support for validation activity. |
| Quality Assurance | <ul style="list-style-type: none">• Review of Performance Qualification Report.• To compile and approval of report.• To monitor all Validation Activities and ensure the Validation is carried out as per the Protocol.• To review Report for completeness and Technical Accuracy. |
| Production | <ul style="list-style-type: none">• Review of Performance Qualification Report.• To co-ordinate and support Performance qualification Activity. |
| Engineering | <ul style="list-style-type: none">• Review of Performance Qualification Report.• To co-ordinate and support Validation Activity.• Responsible for Trouble shooting during execution (If Occurs). |



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

5.0 EQUIPMENT DETAILS:

| | |
|------------------------------|--|
| Equipment Name | Super Heated Water Spray Sterilizer |
| Equipment ID. No. | |
| Size | 1750 DIA X 4500 LG mm |
| Chamber volume | 10800 Liters |
| Manufacturer's Name | M/s Machine Fabrik |
| Supplier's Name | M/s Machine Fabrik |
| Place of Installation | LVP Line Loading Area |



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.0 PRE – QUALIFICATION REQUIREMENTS:

6.1 Training Record of Validation Team:

- All the persons involved in the execution of qualification activity including the persons of outside agencies must be trained in all aspects of the qualification activity including the test methodology, acceptance criteria and safety precautions to be followed during working at service floor.

Verify the training records and attached

6.2 Calibration of Test Instruments: Calibration of all the instruments used for qualification should be mentioned along with Calibration Certificates.

| S. No. | Name of Test Instrument | Date of Last Calibration | Next Due on | Status | Availability of Calibration Certificate | Verified By (QA) Sign/Date |
|--------|-------------------------|--------------------------|-------------|--------|---|----------------------------|
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |

Inference:

.....
.....
.....

Reviewed By:
(Manager QA)
(Sign & Date)



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.3 Biological Indicator Detail :

| BI Detail | Observation | Verified By (QA) Sign/date |
|-------------------------------------|--------------------|---|
| Name of Biological Indicator | | |
| Code Number | | |
| Lot Number | | |
| Spore Population | | |
| Z Value | | |
| D Value | | |
| Manufacturing Date | | |
| Expiry Date | | |

Inference:

.....
.....
.....

Reviewed By:
(Manager QA)
(Sign & Date)



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.0 TESTS & CHECKS:

6.1 EMPTY CHAMBER HEAT DISTRIBUTION STUDY

| | |
|---------------|--|
| Name of Cycle | Heat Distribution Study (For Empty Chamber) |
|---------------|--|

6.1.1 OBSERVATION OF CYCLE SET PARAMETER

| Parameters | Set Value | Cycle Observed Value | | |
|---------------------------|-----------|----------------------|----------|----------|
| | | Cycle -01 | Cycle-02 | Cycle-03 |
| Cycle Start Date | | | | |
| Cycle Start Time | | | | |
| Add water in | 30 Sec. | | | |
| Initial H/E Exhaust | 03 min | | | |
| Set Point 1 | 95.0°C | | | |
| Set Point 2 | 100.0°C | | | |
| Set Point 3 | 105.0°C | | | |
| Rate 1 | 5.0°C | | | |
| Rate 2 | 4.0°C | | | |
| Rate 3 | 2.0°C | | | |
| Sterilization Temperature | 108.0 °C | | | |
| Sterilization Time | 60 min. | | | |
| Control Band | 0.2 °C | | | |



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

| Parameters | Set Value | Cycle Observed Value | | |
|---------------------------------|-----------|----------------------|----------|----------|
| | | Cycle -01 | Cycle-02 | Cycle-03 |
| Overshoot Temperature | 110.0 °C | | | |
| Sterilization Stop Temperature | 107.5 °C | | | |
| Sterilization Reset Temperature | 107.0 °C | | | |
| H/E Exhaust Delay Time | 3 min. | | | |
| H/E Cooling Exhaust | 3 min. | | | |
| Slow Cooling & Temperature | 85 °C | | | |
| Cooling End Temperature | 50 °C | | | |
| H/E Drain Time | 5 min. | | | |
| Process End Pressure | 0.030 Bar | | | |
| Cycle End Date & Time | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date.....

Inference:

.....
.....

Reviewed By
(Manager QA)
Sign/Date:



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.1.2 SUMMARY DETAIL OF STERILIZATION PROCESS

| S. No. | Critical variables | CYCLE | | | | | |
|--------|-------------------------------|----------|----------|----------|----------|----------|----------|
| | | Internal | | | External | | |
| | | Cycle-01 | Cycle-02 | Cycle-03 | Cycle-01 | Cycle-02 | Cycle-03 |
| 01 | Date | | | | | | |
| 02 | Set sterilization temperature | | | | | | |
| 03 | Time process start | | | | | | |
| 04 | Sterilization start Time | | | | | | |
| 05 | Sterilization End Time | | | | | | |
| 06 | Cycle end time | | | | | | |
| 07 | Location of sensor | | | | | | |
| 08 | Sensor no at Cold point | | | | | | |
| 09 | Equilibrium Time | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....
.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.4 HEAT PENETRATION STUDY MINIMUM LAOD

| | |
|---------------|------------------------|
| Name of Cycle | Heat Penetration Study |
|---------------|------------------------|

6.4.1 OBSERVATION OF CYCLE SET PARAMETER

| Parameters | Set Value | Cycle Observed Value | | |
|--------------------------------|----------------------|----------------------|----------|----------|
| | | Cycle-01 | Cycle-02 | Cycle-03 |
| Cycle Start Date | | | | |
| Cycle Start Time | | | | |
| Add water in | 30 Sec. | | | |
| Initial H/E Exhaust | 03 min | | | |
| Set Point 1 | 95.0 ⁰ C | | | |
| Set Point 2 | 100.0 ⁰ C | | | |
| Set Point 3 | 105.0 ⁰ C | | | |
| Rate 1 | 5.0 ⁰ C | | | |
| Rate 2 | 4.0 ⁰ C | | | |
| Rate 3 | 2.0 ⁰ C | | | |
| Sterilization Temperature | 108.0 ⁰ C | | | |
| Sterilization Time | 60 min. | | | |
| Control Band | 0.2 ⁰ C | | | |
| Overshoot Temperature | 110.0 ⁰ C | | | |
| Sterilization Stop Temperature | 107.5 ⁰ C | | | |



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

| Parameters | Set Value | Cycle Observed Value | | |
|---------------------------------|-----------|----------------------|----------|----------|
| | | Cycle-01 | Cycle-02 | Cycle-03 |
| Sterilization Reset Temperature | 107.0 °C | | | |
| H/E Exhaust Delay Time | 3 min. | | | |
| H/E Cooling Exhaust | 3 min. | | | |
| Slow Cooling & Temperature | 85 °C | | | |
| Cooling End Temperature | 50 °C | | | |
| H/E Drain Time | 5 min. | | | |
| Process End Pressure | 0.030 Bar | | | |
| Cycle End Date & Time | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....
.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.4.2 SUMMARY DETAIL OF STERILIZATION PROCESS

| Sr. NO. | Critical variables | CYCLE | | | | | |
|---------|-------------------------------|----------|----------|----------|----------|----------|----------|
| | | Internal | | | External | | |
| | | Cycle-01 | Cycle-02 | Cycle-03 | Cycle-01 | Cycle-02 | Cycle-03 |
| 01 | Date | | | | | | |
| 02 | Set sterilization temperature | | | | | | |
| 03 | Time process start | | | | | | |
| 04 | Sterilization start Time | | | | | | |
| 05 | Sterilization End Time | | | | | | |
| 06 | Cycle end time | | | | | | |
| 07 | Location of sensor | | | | | | |
| 08 | Sensor no at Cold point | | | | | | |
| 09 | Equilibrium Time | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....

.....

Reviewed By
(Manager QA)
Sign/Date:



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.4.3 OBSERVATION REPORT OF BIOLOGICAL INDICATOR & CHEMICAL INDICATOR.

| Status of Biological Indicator | | | | | |
|--------------------------------|-------------|----------|-------------|-----------|-------------|
| Cycle-I | | Cycle-II | | Cycle-III | |
| Sr. No. | Observation | Sr. No. | Observation | Sr. No. | Observation |
| 01 | | 01 | | 01 | |
| 02 | | 02 | | 02 | |
| 03 | | 03 | | 03 | |
| 04 | | 04 | | 04 | |
| 05 | | 05 | | 05 | |

| Status of Chemical Indicator | | | | | |
|------------------------------|-------------|----------|-------------|-----------|-------------|
| Cycle-I | | Cycle-II | | Cycle-III | |
| Sr. No. | Observation | Sr. No. | Observation | Sr. No. | Observation |
| 01 | | 01 | | 01 | |

Verified By
(Quality Assurance)
Sign/Date.....



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.4.4 F₀ CALCULATION

(a) Numerical F₀ Value:

The actual observations obtained during the heat penetration study at different temperature sensing locations are compiled in the table and the observed temperature shall be subjected for calculation of F₀ values at that particular location. The lethality factor calculations are done by using the following formula and the computed (during the sterilization period) are given in the following table.

$$F_0 = dt \sum 10^{(T-121.1)/Z} \quad \text{_____ (a)}$$

$$F_0 = dt \sum (\text{Sum of lethality factors})$$

- Where,**
- dt:** Time interval between successive temperature measurements (in min)
 - T:** Observed temperature at that particular time
(as per the actual temperatures recorded)
 - Z:** change in the heat resistance of *Bacillus Subtilis* spores as temperature is changed (10⁰C or as mentioned in COA)

(b) F₀ Value for Biological Indicators:

The biological F₀ value for biological indicator strip exposed during the sterilization can be calculated as follows.

$$F_0 = D_{121} (\log A - \log B) \quad \text{..... (b)}$$

$$F_0 = \text{.....}(\log \text{.....} - \log \text{.....})$$

$$F_0 =$$

$$F_0 =$$

- Where,**
- D₁₂₁:** D Value of the Biological Indicator at 121⁰C
 - A:** Experimental Biological indicator concentration or spore population
 - B:** Desired level of sterility (SAL- 10⁰)

(c) Desired Spore log reduction:

Calculate the desired reduction in spore population by using the formula-

$$SLR_{\text{desired}} = \log A - \log SAL_{\text{desired}} \quad \text{----- (c)}$$

$$SLR_{\text{desired}} =$$

$$SLR_{\text{desired}} =$$

- Where,**
- A:** Experimental population of Biological Indicator
 - SAL:** Desired level of Sterility (10⁰)

(d) Actual Spore log reduction:

Calculate actual reduction in spore population by using the formula-

$$SLR_{\text{Actual}} = F_0 / D_{121} \quad \text{----- (d)}$$

Where,

$$F_0 \quad : \quad \text{Minimum Calculated } F_0 \text{ value}$$

$$D_{121} \quad : \quad \text{D value of the Biological Indicator at } 121^0\text{C}$$



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

| | | | | | | | | |
|--|--|--|--|--|--|--|--|--|
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY STERILIZER

| | | | | | | | | |
|--|--|--|--|--|--|--|--|--|
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY STERILIZER

| | | | | | | | | |
|--|--|--|--|--|--|--|--|--|
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:
.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.5 HEAT PENETRATION STUDY MAXIMUM LAOD

| | |
|---------------|------------------------|
| Name of Cycle | Heat Penetration Study |
|---------------|------------------------|

6.5.1 OBSERVATION OF CYCLE SET PARAMETER

| Parameters | Set Value | Cycle Observed Value | | |
|--------------------------------|----------------------|----------------------|----------|----------|
| | | Cycle-01 | Cycle-02 | Cycle-03 |
| Cycle Start Date | | | | |
| Cycle Start Time | | | | |
| Add water in | 30 Sec. | | | |
| Initial H/E Exhaust | 03 min | | | |
| Set Point 1 | 95.0 ⁰ C | | | |
| Set Point 2 | 100.0 ⁰ C | | | |
| Set Point 3 | 105.0 ⁰ C | | | |
| Rate 1 | 5.0 ⁰ C | | | |
| Rate 2 | 4.0 ⁰ C | | | |
| Rate 3 | 2.0 ⁰ C | | | |
| Sterilization Temperature | 108.0 ⁰ C | | | |
| Sterilization Time | 60 min. | | | |
| Control Band | 0.2 ⁰ C | | | |
| Overshoot Temperature | 110.0 ⁰ C | | | |
| Sterilization Stop Temperature | 107.5 ⁰ C | | | |



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

| Parameters | Set Value | Cycle Observed Value | | |
|---------------------------------|-----------|----------------------|----------|----------|
| | | Cycle-01 | Cycle-02 | Cycle-03 |
| Sterilization Reset Temperature | 107.0 °C | | | |
| H/E Exhaust Delay Time | 3 min. | | | |
| H/E Cooling Exhaust | 3 min. | | | |
| Slow Cooling & Temperature | 85 °C | | | |
| Cooling End Temperature | 50 °C | | | |
| H/E Drain Time | 5 min. | | | |
| Process End Pressure | 0.030 Bar | | | |
| Cycle End Date & Time | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:
.....
.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.5.2 SUMMARY DETAIL OF STERILIZATION PROCESS

| Sr. NO. | Critical variables | CYCLE | | | | | |
|---------|-------------------------------|----------|----------|----------|----------|----------|----------|
| | | Internal | | | External | | |
| | | Cycle-01 | Cycle-02 | Cycle-03 | Cycle-01 | Cycle-02 | Cycle-03 |
| 01 | Date | | | | | | |
| 02 | Set sterilization temperature | | | | | | |
| 03 | Time process start | | | | | | |
| 04 | Sterilization start Time | | | | | | |
| 05 | Sterilization End Time | | | | | | |
| 06 | Cycle end time | | | | | | |
| 07 | Location of sensor | | | | | | |
| 08 | Sensor no at Cold point | | | | | | |
| 09 | Equilibrium Time | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....

.....

Reviewed By
(Manager QA)
Sign/Date:



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.5.3 OBSERVATION REPORT OF BIOLOGICAL INDICATOR & CHEMICAL INDICATOR.

| Status of Biological Indicator | | | | | |
|--------------------------------|-------------|----------|-------------|-----------|-------------|
| Cycle-I | | Cycle-II | | Cycle-III | |
| Sr. No. | Observation | Sr. No. | Observation | Sr. No. | Observation |
| 01 | | 01 | | 01 | |
| 02 | | 02 | | 02 | |
| 03 | | 03 | | 03 | |
| 04 | | 04 | | 04 | |
| 05 | | 05 | | 05 | |
| 06 | | 06 | | 06 | |
| 07 | | 07 | | 07 | |
| 08 | | 08 | | 08 | |
| 09 | | 09 | | 09 | |
| 10 | | 10 | | 10 | |
| 11 | | 11 | | 11 | |
| 12 | | 12 | | 12 | |
| 13 | | 13 | | 13 | |
| 14 | | 14 | | 14 | |
| 15 | | 15 | | 15 | |
| 16 | | 16 | | 16 | |
| 17 | | 17 | | 17 | |
| 18 | | 18 | | 18 | |
| 19 | | 19 | | 19 | |
| 20 | | 20 | | 20 | |

| Status of Chemical Indicator | | | | | |
|------------------------------|-------------|----------|-------------|-----------|-------------|
| Cycle-I | | Cycle-II | | Cycle-III | |
| Sr. No. | Observation | Sr. No. | Observation | Sr. No. | Observation |
| 01 | | 01 | | 01 | |
| 02 | | 02 | | 02 | |
| 03 | | 03 | | 03 | |
| 04 | | 04 | | 04 | |

Verified By
(Quality Assurance)
Sign/Date.....



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.5.4 Fo CALCULATION

(a) Numerical Fo Value:

The actual observations obtained during the heat penetration study at different temperature sensing locations are compiled in the table and the observed temperature shall be subjected for calculation of Fo values at that particular location. The lethality factor calculations are done by using the following formula and the computed (during the sterilization period) are given in the following table.

$$F_0 = dt \sum 10^{(T-121.1)/Z} \quad \text{_____ (a)}$$

$$F_0 = dt \sum (\text{Sum of lethality factors})$$

- Where,**
- dt:** Time interval between successive temperature measurements (in min)
 - T:** Observed temperature at that particular time
(as per the actual temperatures recorded)
 - Z:** change in the heat resistance of *Bacillus Subtilis* spores as temperature is changed (10⁰C or as mentioned in COA)

(b) Fo Value for Biological Indicators:

The biological Fo value for biological indicator strip exposed during the sterilization can be calculated as follows.

$$F_0 = D_{121} (\log A - \log B) \quad \text{..... (b)}$$

$$F_0 = \text{.....}(\log \text{.....} - \log \text{.....})$$

$$F_0 =$$

$$F_0 =$$

- Where,**
- D₁₂₁:** D Value of the Biological Indicator at 121⁰C
 - A:** Experimental Biological indicator concentration or spore population
 - B:** Desired level of sterility (SAL- 10⁰)

(c) Desired Spore log reduction:

Calculate the desired reduction in spore population by using the formula-

$$SLR_{\text{desired}} = \log A - \log SAL_{\text{desired}} \quad \text{----- (c)}$$

$$SLR_{\text{desired}} =$$

$$SLR_{\text{desired}} =$$

- Where,**
- A:** Experimental population of Biological Indicator
 - SAL:** Desired level of Sterility (10⁰)

(d) Actual Spore log reduction:

Calculate actual reduction in spore population by using the formula-

$$SLR_{\text{Actual}} = F_0 / D_{121} \quad \text{----- (d)}$$

Where,

$$F_0 \quad : \quad \text{Minimum Calculated } F_0 \text{ value}$$

$$D_{121} \quad : \quad \text{D value of the Biological Indicator at 121}^0\text{C}$$



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

| | | | | | | | | |
|--|--|--|--|--|--|--|--|--|
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:
.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

| | | | | | | | | |
|--|--|--|--|--|--|--|--|--|
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:

Inference:
.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY STERILIZER

| | | | | | | | | |
|--|--|--|--|--|--|--|--|--|
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.6.1 OBSERVATIONS OF BOTTLES LEAK TEST

| Sr. No | Parameter | CYCLE | | | | | | | | | | | |
|--------|---------------------|------------|-------------|--------------|-------------|------------|-------------|--------------|-------------|------------|-------------|--------------|-------------|
| | | Cycle-01 | | | | Cycle-02 | | | | Cycle-03 | | | |
| | | Trolley -I | Trolley -II | Trolley -III | Trolley -IV | Trolley -I | Trolley -II | Trolley -III | Trolley -IV | Trolley -I | Trolley -II | Trolley -III | Trolley -IV |
| 01 | No. of Leak Bottles | | | | | | | | | | | | |
| 02 | Status | | | | | | | | | | | | |
| 03 | Checked By | | | | | | | | | | | | |
| 04 | Verified By | | | | | | | | | | | | |
| 05 | Remark | | | | | | | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date.....

Inference:

.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:

.....



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.7 BOTTLES LEAK TEST CHALLENGE STUDY MAXIMUM LOAD

BOTTLES LEAK TEST DETAILS

| | |
|---------------------------|--|
| Equipment Name | |
| Equipment Make | |
| Equipment Location | |
| Equipment ID No. | |

| Name of Cycle | | Bottles Leak Test Challenge | | |
|-----------------------------|-------------------|-----------------------------|-----------------|------------------|
| Date | | | | |
| Parameters | Set Value | Cycle Observed Value | | |
| | | Cycle-I | Cycle-II | Cycle-III |
| Cycle Started Date | | | | |
| Cycle Started Time | | | | |
| Cycle Completed Date | | | | |
| Cycle Completed Time | | | | |
| Vacuum | -0.700 Bar | | | |
| Vacuum Band | 0.050 Bar | | | |
| Vacuum hold time | 10 Min. | | | |
| Pressure | -0.200 Bar | | | |
| No of pulses | 3 Nos. | | | |
| Process end pressure | -0.030 Bar | | | |
| Actual Leakage | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....
.....

Reviewed By
(Manager QA)
Sign/Date:



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

6.7.1 OBSERVATIONS OF BOTTLES LEAK TEST

| Sr. No. | Critical variables | CYCLE | | | | | | | | | | | |
|---------|---------------------|------------|-------------|--------------|-------------|------------|-------------|--------------|-------------|------------|-------------|--------------|-------------|
| | | Cycle-01 | | | | Cycle-02 | | | | Cycle-03 | | | |
| | | Trolley -I | Trolley -II | Trolley -III | Trolley -IV | Trolley -I | Trolley -II | Trolley -III | Trolley -IV | Trolley -I | Trolley -II | Trolley -III | Trolley -IV |
| 01 | No. of Leak Bottles | | | | | | | | | | | | |
| 02 | Status | | | | | | | | | | | | |
| 03 | Checked By | | | | | | | | | | | | |
| 04 | Verified By | | | | | | | | | | | | |
| 05 | Remark | | | | | | | | | | | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:.....

Inference:

.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:

.....



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

7.0 CHECKLIST OF ALL TESTS AND CHECKS:

| TESTS OR CHECKS | EXECUTED [Y/N] | CHECKED BY (SIGN & DATE) QA | COMMENT |
|---|---------------------------|--|----------------|
| Heat Distribution Study | | | |
| Minimum Loaded Chamber Heat Penetration Studies with Biological Indicator Placement For 100 ml LDPE bottles. (3960 Bottles) | | | |
| Maximum Loaded Chamber Heat Penetration Studies with Biological Indicator Placement For 100 ml LDPE bottles. (15840 Bottles) | | | |
| Bottle Leakage Test (Minimum Load) (3960 Bottles) | | | |
| Bottle Leakage Test (Maximum Load) (15840 Bottles) | | | |

Compiled By
(Quality Assurance)
(Sign & Date).....

Inference:

Reviewed By:
(Manager QA)
(Sign & Date).....



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

8.0 DOCUMENTS ATTACHED:

- Raw data of Microbiological Analysis
- Calibration Certificates for Data Logger & Temperature Sensor.
- COA of Biological Indicator.
- Data Logger Printouts.
- Super Heated Water Spray Sterilizer PLC Printouts.

9.0 NON COMPLIANCE:

.....
.....
.....
.....
.....

10.0 DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY:

.....
.....
.....
.....

11.0 CHANGE CONTROL, IF ANY:

.....
.....
.....
.....



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

12.0 REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY) :

.....
.....
.....
.....

13.0 CONCLUSION:

.....
.....
.....
.....

14.0 RECOMMENDATION:

.....
.....
.....
.....



**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

15.0 ABBREVIATIONS:

| | | |
|------|---|--------------------------------------|
| QA | : | Quality Assurance |
| QC | : | Quality Control |
| No. | : | Number |
| Ltd. | : | Limited |
| FDA | : | Food and Drug Administration |
| CFR | : | Code of Federal Regulations |
| CQA | : | Corporate Quality Assurance |
| GMP | : | Good Manufacturing Practices |
| cGMP | : | Current Good Manufacturing Practices |
| WHO | : | World Health Organization |
| SOP | : | Standard Operating Procedure |



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**PERFORMANCE QUALIFICATION REPORT FOR SUPER HEATED WATER SPRAY
STERILIZER**

16.0 REPORT POST APPROVAL:

PREPARED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|--|------|-----------|------|
| OFFICER / EXECUTIVE (QUALITY ASSURANCE) | | | |

REVIEWED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|--|------|-----------|------|
| OPERATING MANAGER (QUALITY ASSURANCE) | | | |
| HEAD (ENGINEERING) | | | |
| HEAD (PRODUCTION) | | | |

APPROVED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|-----------------------------|------|-----------|------|
| HEAD (QUALITY ASSURANCE) | | | |