



STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation and Calibration of Ultra-Violet Chamber	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

1.0 OBJECTIVE:

To lay down a procedure for Operating & Calibration of Ultra-Violet Chamber.

2.0 SCOPE:

This SOP is applicable for Operation and Calibration of Ultra-Violet Chamber in the Quality Control Department.

3.0 RESPONSIBILITY:

Officer, Executive - Quality Control Department
Head - Quality Control Department

4.0 DEFINITION(S):

NA

5.0 PROCEDURE:

Make: Camag, Model: LAAAG

5.1 Operation:

- 5.1.1 Ensure that instrument is clean and free from dust.
- 5.1.2 Place the developed TLC plate in the chamber.
- 5.1.3 Close the Chamber.
- 5.1.4 Select the desired wavelength by switching on the respective lamp (254nm, 365nm).
- 5.1.5 Observe the developed plate through the eye viewer.
- 5.1.6 Switch 'OFF' the illuminating lamp before removing the plate
- 5.1.7 Record the details in the equipments/instrument log book for results & breakdown.

5.2 Calibration:

- 5.2.1 Prepare 0.2% w/v solution of sodium salicylate in 95% ethanol (for 365 nm) and 0.04% solution of sodium salicylate (for 254 nm) in 95 % ethanol.
- 5.2.2 Apply both conc. (5 µl) on plate coated with silica gel G.
- 5.2.3 Operate the instrument as operating procedure.



STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation and Calibration of Ultra-Violet Chamber	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

- 5.2.4 Examine the spot in position normal to radiation.
- 5.2.5 The lamp should be capable of revealing without doubt a standard spot of sodium salicylate with a diameter of about 5 mm on a chromatographic plate coated with silica gel G.
- 5.2.6 Frequency – Quarterly
- 5.2.7 If instrument is out of calibration, affix “UNDER MAINTENANCE” label on the instrument and call for service engineer.
- 5.2.8 Fill the calibration status on metallic calibration label of the instrument, record the calibration results in annexure-I
- 5.2.9 Note calibration activity in the instrument log book

5.3 Cleaning:

Clean the Instrument properly with cotton cloth.

5.4 Precaution:

- 5.4.1 UV Radiations are harmful to eyes and skin therefore never looks directly at the lightened UV Lamp or tube.
- 5.4.2 Also never insert your hand or escape your body surface inside the UV chamber when the lamps are illuminated.
- 5.4.3 Do not rapidly turn /off the switch as it will shorten the life of the lamp.

6.0 ABBREVIATION(S):

QCD - Quality Control Department
SOP - Standard Operating Procedure
TLC - Thin layer chromatography
μl – Microlitre

7.0 REFERENCE(S):

IP-1996

8.0 ANNEXURE(S):

Annexure – I Calibration Record of Ultra violet Chamber



PHARMA DEVILS
QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation and Calibration of Ultra-Violet Chamber	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

9.0 REVISION CARD:

S.No.	REVISION No.	REVISION DATE	DETAILS OF REVISION	REASON (S) FOR REVISION



PHARMA DEVILS
QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation and Calibration of Ultra-Violet Chamber	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

ANNEXURE I

Calibration Record: Ultra Violet Chamber			
Location		Model No.	
Make	Camag	Identification No.	
Calibration Done On		Calibration Due on	
Reference SOP No.:			Page No. : 4 of 1

S.No.	Test	At nm	Observation
1.	0.2 % w/v solution of sodium salicylate in 95 % ethanol	365	
2.	0.04 % w/v solution of sodium salicylate in 95 % ethanol	254	
LIMIT: The lamp should be capable of revealing clearly a standard spot of sodium salicylate with a diameter of about 5 mm on a chromatographic plate coated with silica gel G.			
Remarks: The Instrument Calibration complies / Does Not Comply.			

Calibrated By :	Checked By :
Date :	Date :