

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Preparation & Standardization of Volumetric solution	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

1.0 OBJECTIVE:

To lay down a procedure for Preparation and Standardization of volumetric solution.

2.0 SCOPE:

This SOP is applicable for Preparation and Standardization of volumetric solution use in the Quality Control Laboratory.

3.0 RESPONSIBILITY:

Executive, Officer – Quality Control. Head – Quality Control

4.0 **PROCEDURE** :

4.1 **Preparation of volumetric solution:**

- 4.1.1 Prepare the volumetric solution as per the General test procedure or respective pharmacopoeia.
- 4.1.2 Always use purified water or Milli Q water whenever intended for use.
- 4.1.3 Use analytical grade reagents.
- 4.1.4 Record the details of volumetric solution preparation in 'Preparation log of volumetric solution' (as per Annexure I).
- 4.1.5 Assign batch number to the volumetric solution in the following manner.

VS-XXXX /YY-00, Where,

- VS : Volumetric solution
- _ : Dash

XXXX : Serial number

- / : Slash
- YY : Last two digits of the current year.
- _ : Dash
- 00 : Initial Standardization number
- 01 : Re-standardization number
- 4.1.6 If the solution is not clear and contains suspended particles then filter the solution and standardise.



STANDARD OPERATING PROCEDURE

STANDARD OF ERATING TROCEDURE			
Department: Quality Control	SOP No.:		
Title: Preparation & Standardization of Volumetric solution	Effective Date:		
Supersedes: Nil	Review Date:		
Issue Date:	Page No.:		

4.1.7 If the solution is still hazy/turbid then Record the discard details in the Remarks column of 'Preparation log of volumetric solution'.

(As per Annexure- I).

4.2 Standardization of Volumetric solution:

- 4.2.1 Standardize the volumetric solution as per the respective procedure given in GTP or the respective pharmacopoeia.
- 4.2.2 Use primary standards of NIST grade or A. R. grade reagent for standardization.
- 4.2.3 Dry the primary standards or A. R. grade reagents before use, as mentioned in the respective GTP, pharmacopoeia or reference document.
- 4.2.4 Keep separate A. R. grade reagent for standardization purpose.
- 4.2.5 The solution shall be standardized in triplicate and mean of the Molarity/ Normality value shall be calculated. Calculate the precision of all three values it should not be more than \pm 0.2 % (as per Annexure -II).
- 4.2.6 Record the details of 'Volumetric solution preparation and standardization" in a respective validated calculation excel sheet or as per format given in annexure-II.
- 4.2.7 The actual molarity of the volumetric solution should not differ by \pm 10% of the prescribed strength and the variation of Molarity/Normality between three readings should not be more than 1.0 %.
- 4.3 An authorized person shall check all the preparation and standardization details and shall ensure that the solution is prepared as per the standard procedure and sign the record.
- 4.4 If the actual Molarity/Normality differ from ± 10 % of the prescribed strength then discard the solution and prepare a new fresh solution. Record the discard details in the Remarks column of 'Preparation log of volumetric solution' (as per Annexure I).
- 4.5 If the Variation between three Molarities/ Normality's exceeds 1.0 %, then re-standardize the solution and determine the Molarity/Normality. If again it exceeds 1.0 % then discard the same and prepare a new fresh solution. Record the discard details in the Remarks column of 'Preparation log of volumetric solution' (as per Annexure I).
- 4.6 Alcoholic and Iodine solutions shall be standardised immediately before use.
- 4.7 Ensure the expired solutions are discarded after validity period is over and record the same in the log.



STANDARD OPERATING PROCEDURE			
Department: Quality Control	SOP No.:		
Title: Preparation & Standardization of Volumetric solution	Effective Date:		
Supersedes: Nil	Review Date:		
Issue Date:	Page No.:		

4.8 A daily monitoring of solutions on the shelf shall be carried out by concerned personnel for verifying the validity of the solutions. In case any solution is expiring on that day or becomes turbid, same shall be removed from the shelf and destroyed and a new fresh solution shall be prepared.

4.9 Storage of volumetric solution:

4.9.1 Store all the volumetric solutions in dry and cleaned glass bottles as per the storage conditions given in the respective GTP or pharmacopoeia. Containers other than glass shall be used if prescribed in the pharmacopoeia (e.g. LDPE container).

Note: Always store light sensitive solutions in amber coloured glass bottles.

- 4.9.2 Use the prepared solution within one month or as indicated on label.
- 4.9.3 During usage if any change is observed in the physical properties of the solution, then discard the same and record in the log.
- 4.9.4 Store all volumetric solutions at $25 \pm 2^{\circ}$ C.

4.10 Labeling:

4.10.1 Paste a label with all details on the bottle as per Annexure–III.

5.0 ANNEXURE (S):

Annexure - I: Preparation log of volumetric solutionAnnexure - II: Volumetric solution preparation and standardization record.Annexure - III: Specimen Label for volumetric solution.

6.0 **REFERENCE** (S):

USP/BP/Ph. Eur./IP

SOP: Preparation, Approval, Distribution control, revision and destruction of Standard operating Procedure (SOP).

7.0 ABBREVIATION (S)/ DEFINITION (S):

- A.R. Analytical Reagent
- GTP General Test Procedure



STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Preparation & Standardization of Volumetric solution	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

NIST - National Institute of Standards and Technology

- USP : United State Pharmacopoeia
- BP : British Pharmacopoeia
- Ph. Eur. : European Pharmacopoeia
- IP : Indian Pharmacopoeia

REVISION CARD

S. No.	REVISION No.	REVISION DATE	DETAILS OF REVISION	REASON (S) FOR REVISION	REFERENCE CHANGE CONTROL No.
1	00			New SOP	



STANDARD OPERATING PROCEDURE

SOP No.:

Page No.:

Effective Date:

Review Date:

Department:

Title: Preparation & Standardization of Volumetric solution

Supersedes: Nil

Issue Date:

ANNEXURE I PREPARATION LOG OF VOLUMETRIC SOLUTION

S. No	Date of Preparation	Name of Volumetric solution	B. No. / Lot No.	Quantity Prepared	Valid upto	Prepared by	Remarks



	STANDARD	OPERATING PROC	EDURE	
Department:		SOP No.:		
Title: Preparation & Sta	andardization of Volumetri	ic solution	Effective Date:	
Supersedes: Nil			Review Date:	
Issue Date:			Page No.:	
VOLUMETRIC	ANN C SOLUTION PREPARA	NEXURE II ATION & STANDAR	DIZATION RECORD	
VOLUMETRI	IC SOLUTION PREPAR	RATION & STANDA	RDIZATION RECORD	
			Reference SOP	
Name of Volumetric s	olution			
Reference GTP No.				
Batch No. of Volumet	ric solution			
1.0 Preparation:				
Balance ID: EQ/QCD	/			
Reagent :				
Batch No.:	Grade :	Mfr :		
Procedure :				
Prepared By/ date:		Checked By / Dat	e:	
2.0 Calculation for S	Standardization:			



STANDARD OPERATING PROC	EDURE
Department:	SOP No.:
Title: Preparation & Standardization of Volumetric solution	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

STANDARDIZATION					
Balance	ID:				
Primary	v Standard:	1			
Batch N	0.:	Grade	:	Mfr :	
Primary	standard dryi	ng procedure :			
1 st	Titration	2 nd Tit	ration		3 rd Titration
Weight (g)	Titer volume (mL)	Weight (g)	Titer volume (mL)	Weight (g)	Titer volume (mL)
Calculation (1):		Calculation (2):		Calculation (3):	
Molarity/Normality: Molarity/Norma			lity:	Molarity/Norm	ality:
Variatio	Variation (NMT 1.0 %) : Precision :(NMT 0.2 %) :				
Average Molarity / Normality :(± 10 % of the Prescribed Strength)					
Due dat	e of Re-standar	dization:			
Standard	lized By/Date:		Checked By/Da	ate:	



STANDARD OPERATING PROCEDURE			
Department:	SOP No.:		
Title: Preparation & Standardization of Volumetric solution	Effective Date:		
Supersedes: Nil	Review Date:		
Issue Date:	Page No.:		

ANNEXURE III SPECIMEN LABLE FOR VOLUMETRIC SOLUTION

Name of Volumetric Solution	
Date of Preparation	
Sign /Date	
Valid up to	
Date of standardization	
Due date of Restandardization	
Initial Standardisation	Restandardisation
B.No:	B.No.
Normality/Molarity:	Normality/Molarity:
Sign/Date:	Sign/Date:
Storage condition	At $25^{\circ}C \pm 2^{\circ}C$
	· · ·