

### PHARMA DEVILS ENGINEERING DEPARTMENT

Title: Preventive Maintenance of Deduster			
SOP No.:	Revision No.:	00	
Effective Date:	Supersedes No.	Nil	
<b>Review Date:</b>	Page No.	1 of 3	

#### 1.0 OBJECTIVE

1.1 To describe a procedure for the preventive maintenance of deduster.

#### 2.0 SCOPE

2.1 This procedure applies to the Engineering department.

#### 3.0 **RESPONSIBILITY**

3.1 Engineering

#### 4.0 ACCOUNTABILITY

4.1 Plant Head

#### 5.0 **REFERENCE** (S)

5.1 In-house.

#### 6.0 **PROCEDURE**

- 6.1 Procedure for mechanical work:
  - 6.1.1 Switch off the power supply to the machine and put "under maintenance" label.
  - 6.1.2 Check locking arrangement for rotor to motor.
  - 6.1.3 Check compressed air pipeline holes for chocking if found chocked, clean the same.
  - 6.1.4 Check the wheels for any wear and tear.
- 6.2 Procedure for electrical work:
  - 6.2.1 Check and tighten all the electrical connections.
- 6.3 Procedure for trial run:
  - 6.3.1 Switch on the power supply.
  - 6.3.2 Run the deduster according to SOP.
  - 6.3.3 Check working of the deduster.
  - 6.3.4 Check & record the observations as per Annexure-1.
- 6.4 Procedure for handover of machines:
  - 6.4.1 Remove under maintenance label & affix preventive maintenance status tag and inform the concerned department.
- 6.5 Frequency for preventive maintenance is monthly.

#### 7.0 HISTORY

7.1 Details are given below.

SOP No.	REASON FOR CHANGE	EFFECTIVE DATE



# PHARMA DEVILS

ENGINEERING DEPARTMENT

Title: Preventive Maintenance of Deduster			
SOP No.:	Revision No.:	00	
Effective Date:	Supersedes No.	Nil	
Review Date:	Page No.	2 of 3	

8.0 ABBREVIATIONS: The abbreviations used in the SOP are as follows:

- 8.1 SOP - Standard Operating Procedure
- 8.2 No. - Number
- 8.3 -Quality Assurance QA



## PHARMA DEVILS

ENGINEERING DEPARTMENT

Title: Preventive Maintenance of Deduster			
SOP No.:		<b>Revision No.:</b>	00
Effective Date:		Supersedes No.	Nil
Review Date:		Page No.	3 of 3

#### ANNEXURE-I

#### PREVENTIVE MAINTENANCE CHECK LIST FOR DEDUSTER

Equipment ID No.\_\_\_\_\_

Date: \_\_\_\_\_

S.No.	Particulars	Done /Not Done	Remarks
1	Check and tight all electrical connections.		
2	Check pressure pulley, worm gear.		
3	Check alignment and pressure setting.		
4	Check dust collector's pipe.		
5	Check gearbox oil level & if required top up with servo 90 oil.		
6	Check all pressure gauges.		
7	Others.		

Done by (Sign & date) Checked by (Sign & date)