

PHARMA DEVILS ENGINEERING DEPARTMENT

Title: Preventive Maintenance of Fluidized Bed Dryer

| SOP No.: | Revision No.: | 00 |
|-----------------|----------------|--------|
| Effective Date: | Supersedes No. | Nil |
| Review Date: | Page No. | 1 of 3 |

1.0 **OBJECTIVE**

1.1 To describe a procedure for the preventive maintenance of fluidized bed dryer.

2.0 SCOPE

2.1 This procedure applies to the Engineering department.

3.0 **RESPONSIBILITY**

3.1 Engineering

4.0 ACCOUNTABILITY

4.1 Plant Head

5.0 **REFERENCE (S)**

5.1 In-house.

6.0 **PROCEDURE**

- 6.1 Procedure for mechanical work:
 - 6.1.1 Switch off the power supply of the machine and put "under maintenance" label.
 - 6.1.2 Check all the dampers movement, if they are not functioning smoothly, and then lubricate them. Tighten all the nuts and bolts.
 - 6.1.3 Check the steam coil fins for clogging and physical damage.
 - 6.1.4 Check the dutch weave mesh for any damage.
 - 6.1.5 Check for the smooth operation of all the pneumatic cylinders, and check chambersealing gasket.
 - 6.1.6 Clean and lubricate the wheels of fluidized bed dryer bowl trolley, the blower motor and check for the bolt of impellers.
 - 6.1.7 Check for any abnormal sound or vibration and grease the bearings.
 - 6.1.8 Check the "V" belt and motor pulley condition. If required, replace.
- 6.2 Procedure for electrical work:
 - 6.2.1 Open the panel covers & remove the dust from the panel by compressed air.
 - 6.2.2 Tight all wire connections and check the motor connections.
 - 6.2.3 Check operation of pneumatically operated solenoid valves.
 - 6.2.4 Check PID control valve of container seals (Up & Down bolts)
 - 6.2.5 Clean all the sensors with compressed air.
- 6.3 Procedure for trial run:
 - 6.3.1 Switch on the power supply.
 - 6.3.2 Run the machine according to SOP.
 - 6.3.3 Set the required air pressure.
 - 6.3.4 Check the operation of electrical safety: -

6.3.4.1 Higher current



PHARMA DEVILS ENGINEERING DEPARTMENT

Title: Preventive Maintenance of Fluidized Bed Dryer

| SOP No.: | Revision No.: | 00 |
|-----------------|----------------------|--------|
| Effective Date: | Supersedes No. | Nil |
| Review Date: | Page No. | 2 of 3 |

- 6.3.4.2 Safety pneumatics
- 6.3.4.3 Torn bag detector
- 6.3.4.4 Steam coil
- 6.4.4 Check interlocking of blower with exhaust damper.
- 6.4.5 Check & record the observations.

6.5 Procedure for handover of the machines:

- 6.6.1 Remove "under maintenance" label & affix the preventive maintenance status tag and inform the concerned department.
- 6.7 Frequency of preventive maintenance is monthly.

7.0 HISTORY

7.1 Details are given below.

| SOP No. | REASON FOR CHANGE | EFFECTIVE DATE |
|---------|-------------------|----------------|
| | | |
| | | |

8.0 ABBREVIATIONS: The abbreviations used in the SOP are as follows:

- 8.1 SOP Standard Operating Procedure
- 8.2 No. Number
- 8.3 QA Quality Assurance



PHARMA DEVILS ENGINEERING DEPARTMENT

| Title: Preventive Maintenance of Fluidized Bed Dryer | | | |
|--|--|----------------------|--------|
| SOP No.: | | Revision No.: | 00 |
| Effective Date: | | Supersedes No. | Nil |
| Review Date: | | Page No. | 3 of 3 |
| | | | |

ANNEXURE-I

PREVENTIVE MAINTENANCE CHECK LIST FOR FLUIDIZED BED DRYER

Equipment ID No._____

| S.No. | Particulars | Done/Not Done | Remarks |
|------------------|---|---------------|---------|
| MECHANICAL WORK | | | |
| 1 | Switch off the power supply of the machine and put "under maintenance" label. | | |
| 2 | Check all the dampers movement | | |
| 3 | Check the steam coil fins. | | |
| 4 | Check all the pneumatic cylinders. | | |
| 5 | Checking of 0.3 micron prefilter. | | |
| 6 | Checking of limit switch. | | |
| 7 | Checking of heat exchanger coils. | | |
| 8 | Checking of broken bag detector. | | |
| 9 | Check the "V" belt and motor pulley condition | | |
| 10 | Others. | | |
| ELECTRICAL WORKS | | | |
| 1 | Tight all wire connections | | |
| 2 | Check the motor connections | | |
| 3 | Clean all the sensors | | |
| 4 | Check operation of pneumatically operated solenoid valves. | | |

Done by (Sign & date) Checked by (Sign & date)

Date: _____