

PHARMA DEVILS ENGINEERING DEPARTMENT

Title: Preventive Maintenance of Tablet Inspection Machine

SOP No.:	Revision No.:	00
Effective Date:	Supersedes No.	Nil
Review Date:	Page No.	1 of 3

1.0 OBJECTIVE

1.1 To describe a procedure for the preventive maintenance of Tablet inspection machine.

2.0 SCOPE

2.1 This procedure applies to the Engineering department. .

3.0 **RESPONSIBILITY**

3.1 Engineering

4.0 ACCOUNTABILITY

4.1 Plant Head

5.0 **REFERENCE** (S)

5.1 In-house.

6.0 **PROCEDURE**

- 6.1 Procedure for mechanical work:
 - 6.1.1 Switch off the power supply of the machine and put "under maintenance" label.
 - 6.1.2 Clean and lubricate the motor bearing.
 - 6.1.3 Check belt tension and if required tighten the same.
 - 6.1.4 Check foot paddle whether, working or not.
- 6.2 Procedure for electrical work:
 - 6.2.1 Clean the panel with blower.
 - 6.2.2 Check the operation of the limit switch.
- 6.3 Procedure for trial run:
 - 6.3.1 Switch on the power supply.
 - 6.3.2 Run the machine as per the SOP.
 - 6.3.3 Check & record the observations as per Annexure-1.
- 6.4 Procedure for handover of machines:
 - 6.4.1 Remove "under maintenance" label & affix preventive maintenance status tag and inform to the concerned department.
- 6.5 Frequency of preventive maintenance is once in two months.

7.0 HISTORY

7.1 Details are given.





Title: Preventive Maintenance of Tablet Inspection Machine			
SOP No.:	Revision No.:	00	
Effective Date:	Supersedes No.	Nil	
Review Date:	Page No.	2 of 3	

SOP No.	REASON FOR CHANGE	EFFECTIVE DATE

8.0 ABBREVIATIONS: The abbreviations used in the SOP are as follows:

- 8.1 SOP Standard Operating Procedure
- 8.2 No. Number
- 8.3 QA Quality Assurance



PHARMA DEVILS ENGINEERING DEPARTMENT

Title: Preventive Maintenance of Tablet Inspection Machine

SOP No.:	Revision No.:	00
Effective Date:	Supersedes No.	Nil
Review Date:	Page No.	3 of 3

ANNEXURE-I

PREVENTIVE MAINTENANCE CHECK LIST FOR TABLET INSPECTION BELTS

Equipment ID No._____

Date: _____

S.No.	Particulars	Done/Not Done	Remarks
1.	Check electrical connections.		
2.	Check moving parts of belts if required do the lubrication of the parts.		
3.	Checking of the conveyer belt.		

Done by (Sign & date) Checked by (Sign & date)