

PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE		
Department: Quality Control	SOP No.:	
Title: Storage, Backup And Retrievibility Of Analytical Data (HPLC)	Effective Date:	
Supersedes: Nil	Review Date:	
Issue Date:	Page No.:	

1.0 **OBJECTIVE**:

To lay down a procedure for storage, Backup & Retrievibility of analytical data (HPLC).

2.0 SCOPE:

This SOP is applicable for storage, Backup & Retrievibility of analytical data (HPLC).

3.0 **RESPONSIBILITY:**

Executive, Officer – Quality Control Department

Head – Quality Control Department

4.0 **DEFINITION(S)**:

NA

5.0 PROCEDURE:

For Chemstation Software:

5.1 Weekly Backup:

- 5.1.1 Take weekly backup of HPLC data from "D drive" to "E drive".
- 5.1.2 For selection of HPLC data from "D drive " following path are given

D drive/Data/select files/folder, which you want for backup.

- 5.1.3 Copy the file/folder and save it into "E drive"
- 5.1.4 Take the snap shot of data transfer from D drive to E drive and record this activity in Annexure I.

5.2 Monthly Backup:

- 5.2.1 Insert the backup CD into the CDROM.
- 5.2.2 Select the files/folder in "E drive" which has to be restored.
- 5.2.3 After restore the data in CD, delete the data from "E drive"
- 5.2.4 Store the CD in CD case and keep the same in with proper labeling as per following pattern.

e.g. CD/MM./YY /XX

Where,



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE			
Department: Quality Control	SOP No.:		
Title: Storage, Backup And Retrievibility Of Analytical Data (HPLC)	Effective Date:		
Supersedes: Nil	Review Date:		
Issue Date:	Page No.:		

MM- Month

YY -Instrument ID No.

XX- Year....07,08,09

- 5.2.5 We can also restore the data in C.D. in between month if our system slow due to excess capacity of data.
- 5.2.6 Record the above activity in Annexure I.

6.0 ABBREVIATION (S):

QCD - Quality Control Department

SOP - Standard Operating Procedure

CD - Compact Disc

7.0 REFERENCE(S):

NA

8.0 ANNEXURE(S):

Annexure- I: Details of weekly/Monthly Backup of analytical Data

9.0 REVISION CARD:

S.No.	REVISION No.	REVISION DATE	DETAILS OF REVISION	REASON (S) FOR REVISION