



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

<b>EQUIPMENT ID. No.</b>	
<b>LOCATION</b>	<b>Entry A/L -2 For Manufacturing Area FFS Line</b>
<b>DATE OF QUALIFICATION</b>	
<b>SUPERSEDE PROTOCOL No.</b>	<b>NIL</b>



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**PROTOCOL CONTENTS**

<b>S.No.</b>	<b>TITLE</b>	<b>PAGE NO.</b>
<b>1.0</b>	<b>PRE-APPROVAL</b>	<b>3</b>
<b>2.0</b>	<b>OBJECTIVE</b>	<b>4</b>
<b>3.0</b>	<b>SCOPE</b>	<b>4</b>
<b>4.0</b>	<b>RESPONSIBILITY</b>	<b>5</b>
<b>5.0</b>	<b>EQUIPMENT DETAILS</b>	<b>6</b>
<b>6.0</b>	<b>SYSTEM DESCRIPTION</b>	<b>6</b>
<b>7.0</b>	<b>PRE-QUALIFICATION REQUIREMENTS</b>	<b>7</b>
<b>8.0</b>	<b>CRITICAL VARIABLES TO BE MET</b>	<b>8</b>
<b>9.0</b>	<b>REFERENCES</b>	<b>16</b>
<b>10.0</b>	<b>DOCUMENTS TO BE ATTACHED</b>	<b>16</b>
<b>11.0</b>	<b>DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY</b>	<b>17</b>
<b>12.0</b>	<b>CHANGE CONTROL, IF ANY</b>	<b>17</b>
<b>13.0</b>	<b>REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY)</b>	<b>17</b>
<b>14.0</b>	<b>CONCLUSION</b>	<b>18</b>
<b>15.0</b>	<b>RECOMMENDATION</b>	<b>18</b>
<b>16.0</b>	<b>ABBREVIATIONS</b>	<b>19</b>
<b>17.0</b>	<b>POST APPROVAL</b>	<b>20</b>



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**1.0 PROTOCOL PRE – APPROVAL:**

**INITIATED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>OFFICER/EXECUTIVE (QUALITY ASSURANCE)</b>			

**REVIEWED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>HEAD (PRODUCTION)</b>			
<b>HEAD (ENGINEERING)</b>			

**APPROVED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>HEAD (QUALITY ASSURANCE)</b>			



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**2.0 OBJECTIVE:**

- To provide documented evidence for the Installation Qualification of Dynamic Garment Cabinet for FFS Line.
- To confirm that the equipment and its components are installed as per the Specifications mentioned in the design qualification document and other requirements given by supplier.

**3.0 SCOPE:**

- The Scope of this installation qualification Protocol cum Report is limited to qualification of Dynamic Garment Cabinet (Make: Chempharm Industries India Ltd.) to be installed in Entry A/L -2 for Mfg Area of FFS Line.
- This document provides all the relevant information related to specification, installation checks and acceptance criteria to be required to perform Installation qualification activity of Dynamic Garment storage Cabinet.



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**4.0 RESPONSIBILITY:**

The Validation Group, comprising of a representative from each of the following departments, shall be responsible for the overall compliance of this Protocol cum Report:

<b>DEPARTMENTS</b>	<b>RESPONSIBILITIES</b>
<b>Quality Assurance</b>	<ul style="list-style-type: none"><li>• Initiation, Approval and Compilation of the Installation Qualification Protocol cum Report.</li><li>• Co-ordination with Production and Engineering to carryout Installation Qualification.</li><li>• Monitoring of Installation Qualification Activity.</li><li>• Post Approval of Installation Qualification Protocol cum Report after Execution.</li></ul>
<b>Production</b>	<ul style="list-style-type: none"><li>• Review &amp; Pre-Approval of Installation Qualification Protocol cum Report.</li><li>• To Co-ordinate and support for Execution of Qualification study as per Protocol.</li><li>• Post Approval of Installation Qualification Protocol cum Report after Execution.</li></ul>
<b>Engineering</b>	<ul style="list-style-type: none"><li>• Review &amp; Pre-Approval of Installation Qualification Protocol cum Report.</li><li>• Co-ordination, Execution and technical support in Dynamic Garment cabinet Installation Qualification Activity.</li><li>• Responsible for Trouble Shooting (if occurs during execution).</li><li>• Post Approval of Installation Qualification Protocol cum Report after Execution.</li></ul>



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**5.0 EQUIPMENT DETAILS:**

<b>Equipment Name</b>	Dynamic Garment Cabinet
<b>Equipment</b>	
<b>Manufacturer's Name</b>	Chempharm Industries India Ltd.
<b>Model</b>	CP-GSC-3' X 1.5' X 7'
<b>Supplier's Name</b>	Chempharm Industries India Ltd.
<b>Location of Installation</b>	Entry A/L -2 for Mfg Area FFS Line

**6.0 SYSTEM DESCRIPTION:**

Dynamic Garment storage cabinet is used to maintain Class 100 through HEPA filter having an efficiency of 99.99% down to 0.3 $\mu$ , with a velocity of 90 $\pm$ 20% FPM, at its face to remove atmosphere contaminants from air and maintain garments in Class 100 environment.

Dynamic Garment storage cabinet consists of HEPA filter with an efficiency of 99.99% down to 0.3 $\mu$  with permitted pressure drop. The system is equipped with a motor blower assembly and Pre-filter & fresh air filter to suck air from atmosphere and to pass it through HEPA filter.



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**7.0 PRE – QUALIFICATION REQUIREMENTS:**

**7.1 Verification of Documents:**

- Executed and approved design qualification document.
- instrumentation diagram
- Technical specification of equipment.
- Calibration certificate of components.
- Certificate of material of construction of components.

**7.1.1 Procedure:**

- Verify the above mentioned documents for availability, completeness and approval status.
- If any deviation is observed the same has to be recorded giving reasons for deviation and deviation should be approved by authorized person.
- Approved Drawings and supporting documents would form a part of the IQ Protocol cum report.

**7.1.2 Acceptance Criteria:**

- All the documents should be available, complete and approved by respective authorities.



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**8.0 CRITICAL VARIABLES TO BE MET:**

**8.1 General Checks and Location Suitability:**

<b>INSTALLATION CHECKS</b>	<b>ACCEPTANCE CRITERIA</b>	<b>OBSERVATION</b>	<b>OBSERVED BY (ENGINEERING) SIGN/DATE</b>
Leveling	Should be properly balanced and leveled		
Edges of parts	Metal parts should be properly grind without any sharp edges		
Welding of Joints	Welding of joints should be without any welding burrs		
Place of Installation	Entry A/L -2 for Mfg. Area, FFS Line		
Room Condition	General working condition		
Working space around the equipment	Should be sufficient for easy operation, cleaning, sanitation and maintenance		

**Checked By  
(Production)  
Sign/Date: .....**

**Verified By  
(Quality Assurance)  
Sign/Date: .....**

**Inference:**

.....  
.....  
.....

**Reviewed By  
(Manager QA)  
Sign/Date: .....**





**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**8.2 Equipment Verification:**

<b>INSTALLATION CHECKS</b>	<b>ACCEPTANCE CRITERIA</b>	<b>OBSERVATION</b>	<b>OBSERVED BY (ENGINEERING) SIGN/DATE</b>
<b>Equipment</b>	Dynamic Garment Cabinet		
<b>Model</b>	CP-GSC-3' X 1.5' X 7'		

**ELECTRICAL INSTALLATION:**

Electricity	Voltage	220-230 V		
	Phases	3 Phase		
	Power consumption	350 W		
	Frequency	50 -60 Hz		
Electrical connections have been provided and secured.	Should be provided & secured			
All components in the panel are properly secured	Should be properly secured			
All terminals are tightened	Should be tightened			
Earthing connection to control panel & equipment	Earthing connection to control panel & equipment should be provided.			

**Checked By  
(Production)**

**Sign/Date:** .....

**Verified By**

**(Quality Assurance)**

**Sign/Date:** .....

**Inference:**

.....  
.....  
.....

**Reviewed By**

**(Manager QA)**

**Sign/Date:** .....



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**8.3 Installation Checks:**

S.No.	SPECIFICATION	OBSERVATION	OBSERVED BY (ENGINEERING) SIGN/DATE
1.	Check the proper mechanical installation of Dynamic Garment cabinet		
2.	Check the proper electrical installation of Dynamic Garment cabinet		
3.	Check the parts are working properly		
4.	Check the equipment is free from any defects		
5.	Check the finishing of product contact parts		

**Checked By  
(Production)**

**Sign/Date:** .....

**Verified By**

**(Quality Assurance)**

**Sign/Date:** .....

**Inference:**

.....  
.....  
.....

**Reviewed By**

**(Manager QA)**

**Sign/Date:** .....



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**8.4 MOC Verification List:**

<b>COMPONENTS</b>	<b>ACCEPTANCE CRITERIA</b>	<b>OBSERVED BY (ENGINEERING) SIGN/DATE</b>
<b>Body</b>	SS304	
<b>Dore Hinge</b>	SS304	
<b>Blower</b>	SS	
<b>HEPA</b>	Micro Fiber Glass	
<b>Fresh Air &amp; Exhaust Filter</b>	Al Expended+3HDPE+Al Expended	
<b>PAO Port</b>	SS	
<b>Switch</b>	SS	
<b>Indicator</b>	STD	
<b>Hanging Hook Set</b>	SS	
<b>Hanging pipe</b>	SS	

**Checked By  
(Production)  
Sign/Date: .....**

**Verified By  
(Quality Assurance)  
Sign/Date: .....**

**Inference:**

.....  
.....  
.....

**Reviewed By  
(Manager QA)  
Sign/Date: .....**



PHARMA DEVILS

INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET

PROTOCOL No.:

8.5 EQUIPMENT VERIFICATION

CRITICAL VARIABLE	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY (ENGINEERING) SIGN/DATE
<b>Body</b>			
Manufacturer	Chempharm Industries		
Type	100 % Exhaust Dynamic		
Overall Dimension	1050 X 450 X 2105 mm		
Capacity	190 CFM		
Static Pressure	30 mm of water		
MOC	SS 304 Sheet of 1.0 mm		
Surface Finish	Hair Line Finish		
Door	Double Wall Sandwich Doors—Double Door		
Door Hinge	SS 304 ,06 Nos.		
Door's Glass	125X 875 X5 mm – 4 Nos.		
<b>Blower Assembly</b>			
Make	Air Scanner		
MOC of Blower	SS		
MOC of impeller	Aluminium		
RPM	1350 RPM		
Motor Capacity	1/3 HP - single phase		
<b>HEPA Filter</b>			
Make	Chempharm		
Type	Minipleat		
Size	313 X 783 X 69 mm		
Quantity	01 No.		
Media	Micro Fiber Glass		
Efficiency	99.99% down to 0.3 $\mu$ ,		
Filter class	H-14		
<b>Pre-filter</b>			
Make	Chempharm Industries		



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
 PROTOCOL CUM REPORT  
 FOR  
 DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

CRITICAL VARIABLE	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY (ENGINEERING) SIGN/DATE
Size	176 X 196 X 20 mm		
Quantity	02 No.		
Type	Box type		
Efficiency	90% Down to 5 $\mu$		
Media	Al Expended+3HDPE+Al Expended		
<b>Fresh air filter</b>			
Make	Chempharm Industries		
Size	270 X 510 X 50 mm		
Quantity	01 No.		
Type	Box type		
Efficiency	90% Down to 5 $\mu$		
Media	Al Expended+3HDPE+Al Expended		
<b>Magnehelic Gauge</b>			
Make	Dwyer		
Range	HEPA filter: 0-50 mm of		
<b>Accessories</b>			
UV Light	Make	Philips	
	Watt	15 W	
	Quantity	1 No.	
Hour meter	Make	Nishant	
	Quantity	1 No.	
Switch	Make	Roma	
	Qty.	03 Nos.	
	voltage	5/15 Amp.	
Tube Light	Make	Philips	
	Watt	14 W	
	Quantity	1 No.	
PAO Port	SS		



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

<b>CRITICAL VARIABLE</b>	<b>ACCEPTANCE CRITERIA</b>	<b>OBSERVATION</b>	<b>OBSERVED BY (ENGINEERING) SIGN/DATE</b>
Door Handle With	STD. – 01 No.		
Levelling Screw	SS – 04 Nos.		
Indicator	Make : Laptron		
Hanging Hook Set	SS --- 03 NOS.		
Hanging pipe	SS (19 mm dia)		
<b>Electrical Supply</b>			
Power Supply	220-230 V AC/ 50-60 Hz		
Power Consumption	350 Watts		

**Checked By  
(Production)**

**Sign/Date:** .....

**Verified By**

**(Quality Assurance)**

**Sign/Date:** .....

**Inference:**

.....  
.....  
.....

**Reviewed By**

**(Manager QA)**

**Sign/Date:** .....



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**8.6 Safety:**

<b>CHECKS</b>	<b>ACCEPTANCE CRITERIA</b>	<b>OBSERVATION</b>	<b>OBSERVED BY ENGINEERING SIGN/DATE</b>
Well embedded equipment	For proper sifting		
Electrical wiring and Earthing	Electrical wiring should be as per approved drawings. Double external earthing to control machine (panel and motors).		
<b>Start ON/OFF switch:</b> To stop the process immediately	Should be provided for safety		

**Checked By  
(Production)  
Sign/Date: .....**

**Verified By  
(Quality Assurance)  
Sign/Date: .....**

**Inference:**

.....  
.....  
.....

**Reviewed By  
(Manager QA)  
Sign/Date: .....**



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**9.0 REFERENCES:**

**The Principle Reference is the following:**

- Validation Master Plan
- Schedule-M – “Good Manufacturing Practices and Requirements of Premises, Plant and Equipment for Pharmaceutical Products.”
- WHO Essential Drugs and Medicines Policy, QA of Pharmaceuticals, Vol-2 – Good Manufacturing Practices and Inspection.

**10.0 DOCUMENTS TO BE ATTACHED:**

- Technical details for Equipment Requirement with Engineering Drawings.
- Certificate of MOC
- Calibration certificates





**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**11.0 DEVIATION FROM PRE-DEFINED SPECIFICATION IF, ANY:**

.....  
.....  
.....  
.....  
.....  
.....

**12.0 CHANGE CONTROL, IF ANY:**

.....  
.....  
.....  
.....  
.....  
.....

**13.0 REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY ):**

.....  
.....  
.....  
.....  
.....  
.....  
.....



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**14.0 CONCLUSION:**

.....

.....

.....

.....

.....

.....

**15.0 RECOMMENDATION:**

.....

.....

.....

.....

.....

.....

.....

.....



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**16.0 ABBREVIATIONS:**

%	:	Percent
μ	:	Micron
AC	:	Alternate current
cGMP	:	Current Good Manufacturing Practice
CQA	:	Corporate Quality Assurance
EU	:	European Union
FFS	:	Form Fill Seal
FPM	:	Feet per minute
GA	:	General Arrangement
HEPA	:	High Efficiency Particulate Air
HP	:	Horse Power
Hz	:	Horse Power
IB	:	Injection block
ID.	:	Identification
IQ	:	Installation qualification
Ltd.	:	Limited
mm	:	Millimeter
MOC	:	Material of Construction
Nos.	:	Number
PAO	:	Poly Alpha olefin
Pvt.	:	Private
QA	:	Quality Assurance
RPM	:	Rotation per minute
SS	:	Stainless Steel
SSG	:	Sterile Storage Cabinet
UV	:	Ultra Violet
V	:	voltage
W	:	Watt



**PHARMA DEVILS**

**INSTALLATION QUALIFICATION  
PROTOCOL CUM REPORT  
FOR  
DYNAMIC GARMENT CABINET**

**PROTOCOL No.:**

**17.0 PROTOCOL POST APPROVAL:**

**INITIATED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>OFFICER/EXECUTIVE (QUALITY ASSURANCE)</b>			

**REVIEWED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>HEAD (PRODUCTION)</b>			
<b>HEAD (ENGINEERING)</b>			

**APPROVED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>HEAD (QUALITY ASSURANCE)</b>			