



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**PERFORMANCE QUALIFICATION  
REPORT  
FOR  
VIAL SEALING MACHINE**

<b>EQUIPMENT ID. No.</b>	
<b>LOCATION</b>	<b>Vial Capping Room</b>
<b>DATE OF QUALIFICATION</b>	
<b>SUPERSEDES PROTOCOL No.</b>	<b>NIL</b>



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**REPORT CONTENTS**

<b>S.No.</b>	<b>TITLE</b>	<b>PAGE No.</b>
<b>1.0</b>	<b>Pre-Approval</b>	<b>3</b>
<b>2.0</b>	<b>Objective</b>	<b>4</b>
<b>3.0</b>	<b>Scope</b>	<b>4</b>
<b>4.0</b>	<b>Responsibility</b>	<b>5</b>
<b>5.0</b>	<b>Equipment Details</b>	<b>6</b>
<b>6.0</b>	<b>Pre-Qualification Requirement</b>	<b>6</b>
<b>7.0</b>	<b>Tests &amp; Checks</b>	<b>7</b>
<b>8.0</b>	<b>Checklist Of All Tests And Checks</b>	<b>10</b>
<b>9.0</b>	<b>Documents To Be Attached</b>	<b>11</b>
<b>10.0</b>	<b>Non Compliance</b>	<b>11</b>
<b>11.0</b>	<b>Deviation From Pre-Defined Specification, If Any</b>	<b>11</b>
<b>12.0</b>	<b>Change Control, If Any</b>	<b>11</b>
<b>13.0</b>	<b>Review Inclusive of Follow Up Action, If Any</b>	<b>12</b>
<b>14.0</b>	<b>Conclusion</b>	<b>12</b>
<b>15.0</b>	<b>Recommendations</b>	<b>12</b>
<b>16.0</b>	<b>Abbreviations</b>	<b>13</b>
<b>17.0</b>	<b>Post Approval</b>	<b>14</b>



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**1.0 PRE – APPROVAL:**

**INITIATED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>OFFICER/EXECUTIVE (QUALITY ASSURANCE)</b>			

**REVIEWED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>HEAD (PRODUCTION)</b>			
<b>HEAD (QUALITY CONTROL)</b>			
<b>HEAD (ENGINEERING)</b>			

**APPROVED BY:**

<b>DESIGNATION</b>	<b>NAME</b>	<b>SIGNATURE</b>	<b>DATE</b>
<b>HEAD (QUALITY ASSURANCE)</b>			



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**2.0 OBJECTIVE:**

- To provide documented evidence that the Equipment is performing consistently, repeatedly and reproducibly within its established operating range and the results of all the test parameters meet the pre-defined acceptance criteria.
- To confirm the suitability of the Standard Operating Procedures for all routine activities associated with the system.

**3.0 SCOPE:**

- The scope of this report is limited for qualification of Vial Sealing Machine installed in Vial Capping Room.
- This report provides all the relevant information of the performance qualification activity, In-process observations and analytical data of testing of collected samples.



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**4.0 RESPONSIBILITY:**

The Validation Group, comprising of a representative from each of the following departments, shall be responsible for the overall compliance of this Report:

<b>DEPARTMENTS</b>	<b>RESPONSIBILITIES</b>
<b>Quality Assurance</b>	<ul style="list-style-type: none"><li>• Preparation, Authorization, Approval and Compilation of the Performance Qualification Report.</li><li>• Co-ordination with Quality Control, Production and Engineering to carryout Performance Qualification Activity.</li><li>• Monitoring of Performance Qualification.</li></ul>
<b>Production</b>	<ul style="list-style-type: none"><li>• Review of Report.</li><li>• To co-ordinate and support Performance Qualification Activity.</li></ul>
<b>Quality Control</b>	<ul style="list-style-type: none"><li>• Review of Report.</li><li>• Analytical Support (Microbiological Testing/Analysis)</li></ul>
<b>Engineering</b>	<ul style="list-style-type: none"><li>• Reviewing of qualification protocol for correctness, completeness and technical excellence.</li><li>• Responsible for trouble shooting (if occurred during execution).</li><li>• Maintenance &amp; preventive maintenance as per schedule.</li></ul>



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**5.0 EQUIPMENT DETAILS:**

<b>Equipment Name</b>	Vial Sealing Machine
<b>Equipment</b>	.....
<b>Manufacturer's Name</b>	Aegis Pharma Tech
<b>Supplier's Name</b>	Aegis Pharma Tech
<b>Location of Installation</b>	Vial Capping Room

**6.0 PRE – QUALIFICATION REQUIREMENTS:**

Verification for availability, completeness and approval status of all the required relevant documents shall be done and observations shall be recorded in the performance qualification report.

- Executed and approved Design Qualification document.
- Executed and approved Installation Qualification document.
- Executed and approved Operational Qualification document.
- SOP for Operation & Cleaning of Vial Sealing Machine.
- SOP for Preventive Maintenance of Vial Sealing Machine.



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**7.0 TESTS AND CHECKS:**

**7.1 Verification of Documents:**

Record the observations for documents in the below mentioned table.

S.No.	Document Name	Document/SOP No.	Completed (Yes/No)	Checked By (Engineering) Sign/Date	Verified By (QA) Sign/Date
1.	Executed and approved Design Qualification document				
2.	Executed and approved Installation Qualification document				
3.	Executed and approved Operational Qualification document				
4.	PQ Protocol approved				
5.	SOP for Operation & Cleaning of Vial Sealing Machine.				
6.	SOP for Preventive Maintenance of Vial Sealing Machine.				

**Checked By  
(Production)**

**Sign/Date:** .....

**Verified By  
(Quality Assurance)**

**Sign/Date:** .....

**Inference:**

.....  
.....  
.....

**Reviewed By  
(Manager QA)**

**Sign/Date:** .....



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

PROTOCOL No.:

**7.2 Capping Quality:**

Date of Test	
Vial Size	

Cycle .....				
Parameter	Low Speed ( )			
Sample after .....(min.)				
Vials No.	1	2	3	4
Proper Sealing (Vacuum Leak Test)				
Proper Crimping				
Breakage of Flip off Seals				
Cut on Aluminium Seals				
Aluminium seal Rotation				
Parameter	Optimum Speed ( )			
Proper Sealing (Vacuum Leak Test)				
Proper Crimping				
Breakage of Flip off Seals				
Cut on Aluminium Seals				
Aluminium seal Rotation				





**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

Parameter	High Speed ( )			
Proper Sealing (Vacuum Leak Test)				
Proper Crimping				
Breakage of Flip off Seals				
Cut on Aluminium Seals				
Aluminium Seal Rotation				

**Checked By  
(Production)  
Sign/Date: .....**

**Verified By  
(Quality Assurance)  
Sign/Date: .....**

**Inference:**

.....

.....

.....

.....

.....

**Reviewed By  
(Manager QA)  
Sign/Date: .....**



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**8.0 CHECKLIST OF ALL TESTS & CHECKS:**

This checklist is provided to ensure that all tests or checks required for this protocol have been executed.

Tests or Checks	Executed (Yes/No)	Remarks
Verification of DQ, IQ & OQ & Other Documents		
Verification of Machine Performance		

**Checked By  
(Production)  
Sign/Date: .....**

**Verified By  
(Quality Assurance)  
Sign/Date.....**

**Inference:**

.....

.....

.....

.....

.....

**Reviewed By  
(Manager QA)  
Sign/Date: .....**



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**9.0 DOCUMENTS TO BE ATTACHED:**

- Operation and Maintenance Manual.
- Copy of SOP's.
- Any Other Relevant Documents.

**10.0 NON COMPLIANCE:**

.....  
.....  
.....  
.....  
.....  
.....

**11.0 DEVIATION FROM PREDEFINED SPECIFICATION IF, ANY:**

.....  
.....  
.....  
.....  
.....  
.....  
.....  
.....  
.....

**12.0 CHANGE CONTROL, IF ANY:**

.....  
.....  
.....  
.....  
.....  
.....



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**13.0 REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY):**

.....  
.....  
.....  
.....  
.....  
.....

**14.0 CONCLUSION:**

.....  
.....  
.....  
.....  
.....  
.....  
.....  
.....

**15.0 RECOMMENDATION:**

.....  
.....  
.....  
.....  
.....  
.....  
.....  
.....



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**16.0 ABBREVIATIONS:**

WHO	:	World Health Organization
FDA	:	Food and Drug Administration
CFR	:	Code of Federal Regulations
cGMP	:	Current Good Manufacturing Practices
EU	:	European Union
QC	:	Quality Control
DQ	:	Design Qualification
IQ	:	Installation Qualification
OQ	:	Operational Qualification
PQ	:	Performance Qualification
SOP	:	Standard Operating Procedure
NLT	:	Not Less Than
KW	:	Kilo watt
SS	:	Stainless Steel
ID.	:	Identification
mm	:	Mili meter
MCB	:	Miniature Circuit Breaker
ID	:	Inner Diameter



**PHARMA DEVILS**

**PERFORMANCE QUALIFICATION REPORT  
FOR  
VIAL SEALING MACHINE**

**PROTOCOL No.:**

**17.0 POST APPROVAL:**

**INITIATED BY:**

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER/EXECUTIVE (QUALITY ASSURANCE)			

**REVIEWED BY:**

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			
HEAD (QUALITY CONTROL)			
HEAD (ENGINEERING)			

**APPROVED BY:**

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			