



STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Daily Temperature Monitoring of Walk-In Type Incubator	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

1.0 OBJECTIVE:

To lay down a Procedure for Daily Temperature Monitoring of Walk-In Type Incubator.

2.0 SCOPE:

This SOP is applicable for Daily Temperature Monitoring of Walk-In Type Incubator ($32.5^{\circ}\text{C} \pm 2.5^{\circ}\text{C}$ and $22.5^{\circ}\text{C} \pm 2.5^{\circ}\text{C}$) in Microbiology Section of Quality Control.

Equipment I.D.

3.0 RESPONSIBILITY:

Operating Person - Microbiology

4.0 ACCOUNTABILITY:

Head – QC

5.0 ABBREVIATIONS:

Ltd	Limited
No.	Number
QA	Quality Assurance
QC	Quality Control
SOP	Standard Operating Procedure

6.0 PROCEDURE:

6.1 Manual Temperature Monitoring:

6.1.1 Observed the temperature after pressing hold for minimum and maximum on digital screen of Digital Temperature Meter. After observed the temperature Digital temperature meter shall be reset immediately.

6.1.2 Record the Observations in **Annexure-I**, Titled “**Incubator Temperature Monitoring Record**”.

6.1.3 Record the Temperature of walk-in type incubator twice in a working day:

6.1.4 If temperature goes out of limit with respect to set temperature, acknowledge the alarm and inform to engineering department for rectification of problem.

6.1.5 If alarm observed continues more than one hour same shall be addressed through lab incident.

6.1.6 Temperature Limit: $22.5^{\circ}\text{C} \pm 2.5^{\circ}$ and $32.5^{\circ}\text{C} \pm 2.5^{\circ}\text{C}$.




STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Daily Temperature Monitoring of Walk-In Type Incubator	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

6.2 Online Temperature Monitoring:

6.2.1 Online temperature monitoring is applicable for walk-in type incubator.

6.2.2 Open the SCADA Online temperature monitoring system and display show as below



6.2.3 Select the login and click on “LOGIN” and display show as below






STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Daily Temperature Monitoring of Walk-In Type Incubator	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

6.2.4 Select the “user” in user name column and enter the password than click on “OK’ and display show as below.



6.2.5 After login click on “MAIN PAGE” and display show as below




6.2.6 Select the incubator which to be print and click on Incubator I.D and display show as below




STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Daily Temperature Monitoring of Walk-In Type Incubator	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:



6.2.7 Click on “REPORT” and display show as below after that select the date and time which to be print and click on “QUERY REPORT”.

6.2.8 Select the data 24 hours for print with print interval of every hour.



6.2.9 Click on “PRINT REPORT” and display show as below



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Daily Temperature Monitoring of Walk-In Type Incubator	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

Print dialog box showing settings for HP LaserJet 1022 printer, ISO A4 paper size, and Landscape orientation. The 'Print Current' button is highlighted.

5/25/2020 11:00...	22.9
5/25/2020 12:00...	22.5
5/25/2020 1:00...	22.5
5/25/2020 2:00...	22.6

HOME

Monday, May 25, 2020 9:00:00 AM

Tuesday, May 26, 2020 9:00:00 AM

QUERY REPORT

PRINT REPORT

6.2.10 Click on “Print Current” for print and display show as below

DELTA SCADA :- DEVLOPMENT SOFTWARE

DATE/TIME	ROOM ID (20°C-25°C)
MINIMUM TEMP	21.2
MAXIMUM TEMP	23
5/25/2020 10:00...	22.7
5/25/2020 11:00...	22.9
5/25/2020 12:00...	22.5
5/25/2020 1:00...	22.5
5/25/2020 2:00...	22.6

HOME

Monday, May 25, 2020 9:00:00 AM

Tuesday, May 26, 2020 9:00:00 AM

QUERY REPORT


PRINT REPORT

6.2.11 After print click on “HOME” and display show as below



STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Daily Temperature Monitoring of Walk-In Type Incubator	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:



6.2.12 Print the temperature report of incubators.

6.2.13 Temperature monitoring printout report shall be verified and reviewed with sign & date.

6.2.14 During temperature monitoring printout report verification, if any temperature variation with respect to temperature limits observed at single point; justification shall write on temperature monitoring printout.

6.2.15 If temperature variation with respect to temperature limits observed at multiple times; same shall be address through lab incident.

6.2.16 Temperature Limit: 22.5°C ± 2.5° for PC/INS/S-39/05 and 32.5°C ± 2.5°C.

6.2.17 In case of malfunction of SCADA online temperature monitoring system; manual temperature Monitoring shall be performed.

7.0 REFERENCES:

Not Applicable

8.0 ANNEXURES:

ANNEXURE No.	TITLE OF ANNEXURE	FORMAT No.
Annexure-I	Incubator Temperature Monitoring Record	

ENCLOSURES: SOP Training Record.



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Daily Temperature Monitoring of Walk-In Type Incubator	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

9.0 DISTRIBUTION:

- Controlled Copy No. 01 Quality Assurance
- Controlled Copy No. 02 Microbiology
- Master Copy Quality Assurance

10.0 REVISION HISTORY:

CHANGE HISTORY LOG

Revision No.	Change Control No.	Details of Changes	Reason for Change	Effective Date	Updated By

