

MULTI MIX MANUFACTURING PLANT

PROTOCOL No.:

INSTALLATION QUALIFICATION PROTOCOL CUM REPORT FOR MULTI MIX PLANT

EQUIPMENT ID NO.	
LOCATION	Manufacturing Line
DATE OF QUALIFICATION	
SUPERSEDES PROTOCOL NO.	NIL



PROTOCOL No.:

MULTI MIX MANUFACTURING PLANT

PROTOCOL CONTENTS

S.No.	TITLE	PAGE No.
1.0	Protocol pre –approval	3
2.0	Objective	4
3.0	Scope	4
4.0	Responsibility	5
5.0	Equipment details	6
6.0	System description	6
7.0	Pre-qualification requirements	7
8.0	Critical variable to be met	8
9.0	References	25
10.0	Documents to be attached	25
11.0	Deviation from pre-defined specification, if any	25
12.0	Change control, if any	25
13.0	Review (inclusive of follow up action, if any)	26
14.0	Conclusion	26
15.0	Recommendation	26
16.0	Abbreviations	27
17.0	Post approval	28



PROTOCOL No.:

MULTI MIX MANUFACTURING PLANT

1.0	PRE -APPROVAL:
	INITIATED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER/EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			



PROTOCOL No.:

MULTI MIX MANUFACTURING PLANT

2.0 OBJECTIVE:

- To carry out the Installation Qualification of Multi-Mix Plant with Model MP 500 procured Propack Technologies Pvt. Ltd. to be used for Manufacturing of Creams/Gel & Multi mix.
- To confirm that the equipment and its components are as per the Specifications and Installed as per the approved Design and complies with cGMP practices.
- To prove that each Operation proceeds as per the Design Specification and the tolerances prescribed there in the document, are the same at utmost transparency.
- To ensure that there is sufficient information available to enable the equipment to be operated and maintained safely, effectively and consistently.

3.0 SCOPE:

- To verify the critical dimensions of the unit and record Serial Numbers/Model number of critical components.
- To verify that the correct hardware has been installed, system initializes correctly.
- To record the as-built drawing numbers of equipment drawing, P & ID and circuit diagram.
- To Calibrate Temperature and Pressure measurements of Control System, Recorder, Gauges and displays.



PROTOCOL No.:

MULTI MIX MANUFACTURING PLANT

4.0 **RESPONSIBILITY:**

The Validation Group, Comprising of a representative from each of the following Departments, shall be responsible for the overall compliance of this Protocol.

DEPARTMENTS	RESPONSIBILITIES	
Quality Assurance	 Preparation, Review, Approval and Compilation of the Installation Qualification Protocol cum Report. Co-ordination with Production and Engineering to carryout Installation Qualification. Monitoring of Installation Qualification Activity. 	
Production	 Review of Protocol cum Report. Execution of Installation Qualification. 	
Engineering	 Review of Protocol cum Report. To co-ordinate and support Installation Qualification Activity. Calibration of Process Instruments. 	

5.0 EQUIPMENT DETAILS:

Equipment Name	Multi-Mix Manufacturing Plant	
Equipment ID		
Manufacturer's Name	Propack Technologies Pvt. Ltd.	
Supplier's Name	Propack Technologies Pvt. Ltd.	
Model	MP 500	
Location of Installation	Manufacturing Line	



FOR

MULTI MIX MANUFACTURING PLANT

PROTOCOL No.:

6.0 SYSTEM DESCRIPTION:

To design and manufacture multi mix plant for processing of Multi mix/Cream/Gels/Lotion as per product safety, cGMP guideline and to provide assurance that the equipment is manufactured as per the URS and it complies with the scope of supply.

- 1. Multi mixer manufacturing vessel
- 2. Wax phase vessel
- 3. Transfer pump
- 4. Electric control panel
- 5. Vacuum pump
- 6. Utility system
- 7. Batch storage vessel working platform
- 8. Homogenizer
- 9. Meter in jump

Multi Mixer manufacturing vessel:

It consists of cylindrical shell and jacketed vessel. It is fitted with the top mounted SS 316 shaft with anchor having baffles and Teflon scrappers moving in a clockwise direction. One more baffles system is mounted in the inner side of the vessel. The vessel is provided with pressure release vent, safety valve rupture disc, gauge and a temperature sensor with digital display. The vessel is provided with bottom homogenizer and unloading of finished product to storage vessel using lobe pump. The vessel is also provided with steam and cooling water to the jacketed tank. The vessel is also provided with light glass, sight glass, charge hole and hand hold on top dished end.

High speed homogenizer is installed at the manufacturing vessel. It is a Silverson type homogenizer and consists of slit sleeve type SS 316 blade and rotates at 2800 RPM.

Wax phase Vessel:

It is fitted with bottom mounted stirrer coupled to SS 316 shaft with agitator, pressure gauge, vent valve, safety valve rupture disc, and a temperature sensor with digital display. It is provided with bottom outlet connected to manufacturing vessel through a conical filter having SS mesh screen of 100# filter of melted waxes. It is also provided with the steam supply to the jacket.



PROTOCOL No.:

MULTI MIX MANUFACTURING PLANT

Utility system:

A utility pendant is provided to bring the utility lines from the service floor to the platform so as to run the utility line below the platform.

There is a manual mode of operation for manufacturing plant-400 kg. For manual mode selector switches are provided on control panel to control the parameter.

Water inlet : 1" dia. TC flanged end.
Water outlet : 1" dia. TC flanged end.
Cooling water inlet : 1" dia. TC flanged end.
Cooling water outlet : 1" dia. TC flanged end.

7.0 PRE-QUALIFICATION REQUIREMENTS:

The results of any tests should meet the limits and acceptance criteria specified in the test documents. Any deviations or issues should be rectified and documented prior to IQ commencing.

7.1 System Pre-requisites:

		COMPLETED	CHECKED BY	VERIFIED BY
S.No.	DESCRIPTION OF PRE-REQUISITE	(YES/NO)	(PRODUCTION)	(QA)
			(SIGN & DATE)	(SIGN & DATE)
	Verify that the DQ of the Multi-Mix Plant			
	has been executed and approved.			
1.	DQ Protocol Document No.:			



8.0

INSTALLATION QUALIFICATION PROTOCOL CUM REPORT FOR

PROTOCOL N	No.:
------------	------

MULTI MIX MANUFACTURING PLANT

8.1 GENERAL CHECKS AND LOCATION SUITABILITY:

CRITICAL VARIABLE TO BE MET:

INSTALLATION CHECKS	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY (ENGINEERING) (SIGN & DATE)
Grouting and Mounting of			
Equipment			
Equipment should be			
Properly Balanced & Leveled			
All The Metal Parts Should	Eng. / Production / QA		
Be Properly Grounded	to Certify		
Without any Sharp			
Edges.			
Welding of Joints Without			
Any Welding Burrs			
Place of Installation	Manufacturing Line		
Room Condition	Temp NMT 25°C		
	RH NMT 55%		
Illumination	NLT 300 Lux.		
Working space around the	Should be sufficient for		
Equipment	easy operation,		
	cleaning, sanitation and		
	maintenance		

Checked By	Verified By
Production	Quality Assurance
Sign & Date	Sign & Date
Informaci	
	•••••••••••••••••••••••••••••••••••••••
•••••	
•••••	•••••••••••••••••••••••••••••••••••••••
	Reviewed By:
	Manager QA
	Sign & Date



8.2

INSTALLATION QUALIFICATION PROTOCOL CUM REPORT FOR

MULTI MIX MANUFACTURING PLANT

PROTOCOL No.:

INSTALLATION CHECKS VERIFICATION:

S.No.	TEST PARAMETER	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
A.	Multi Mixing ma	nufacturing vessel:		
1.	Check any physical damage to multi mix vessel.	There should no damage to the mixing vessel		
2.	Check specification of multi mixing vessel.	Working capacity: 625 lt. Gross capacity: 750 lt. Min. Batch size: 200 lt. Shell thickness: 6 mm cyl. Insulation : 2 mm Thk. Legs : SS 304 3 Nos. Jacket thickness: 5 mm Bottom thickness: 7 mm conical. Inside finish: mirror finish Outside finish: matt finish		
3.	Check the specification of agitator motor	5 HP / internal agitator / 415 V / 1440 RPM FLP Make: Hindustan		
4.	Check the specification of gear Box	GEAR BOX, TYPE- W, GEAR BOX SIZE -110, REDUCTION RATIO -20:1, FRAME TYPE- 112, MOUNTING TYPE-U MAKE :BONIFIGLIOLI		
5.	Check the specification of product sensor	Make: radix Type: PT – 100 MOC: SS 316 L End connection: 1/2" BSP Temperature range (0 – 200)		
6.	Check the specification of product outlet valve	Make : Cipriani Ø38mm dia., SS 316 L outlet Pneumatic operated redial diaphragm valve		



FOR

MULTI MIX MANUFACTURING PLANT

PROTOCOL No.:

S.No.	TEST PARAMETER	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
7.	Check the specification of water Inlet.	make: Inoxpa Size: 25 mm dia. With pneumatic operated butterfly valve		
8.	Check the specification of CIP.	make: Inoxpa Size: 40 mm dia. With pneumatic operated butterfly valve MOC: SS 316 L		
9.	Check the specification of ingredient suction valve.	make: Inoxpa Size : 25 mm dia. TC Type : Flexible hose & pneumatic operated butterfly valve		
10.	Check the specification of vent filter valve	make: Inoxpa Size: 25 mm dia. TC end Type: pneumatic operated butterfly valve		
11.	Check the specification of safety valve.	make : Inoxpa Size : 25 mm dia. TC		
12.	Check the specification of pressure gauge	Make: Baumer MOC: SS 304 Range: 0- 10 kg/cm ² End conn-1/4" BSP Dial Size: 2.5"		
13.	Check the specification of compound gauge	Make: Baumer MOC: SS 316 L Range: 760mmHg to 10 kg/cm ² End conn-3/8" BSP Dial Size: 4"		
14.	Check the specification of recirculation	make: Inoxpa Size: 38 mm dia. TC end Type: pneumatic operated divert valve		
15.	Check the specification of wax inlet	make: Inoxpa Size: 25 mm dia. TC end Type: pneumatic operated butterfly valve		



PROTOCOL No.:

S.No.	TEST PARAMETER	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
16.	Check the specification of Spray ball	Make : Atul Ø50mm dia., End conn-3/4" BSP SS 316 L		
17.	Check the specification of additive hopper	make : Inoxpa Size : 38 mm dia. TC end Type : manual butterfly valve		
18.	Steam inlet / cooling outlet	Size : 25 mm TC with pneumatic		
19.	Steam outlet / cooling inlet	Size : 25 mm TC with pneumatic piston valve Make : Spirex		
В.	Wax Phase vesse	el		
1.	Check any physical damage to Wax phase vessel.	There should no damage to the mixing vessel		
2.	Check specification of wax phase vessel.	Working capacity: 350 lt. Gross capacity: 420 lt. Min. batch size: 105 lt. Shell thickness: 5 mm Vessel top: 2 mm thick SS 316 L with flat top 1/3 openable Insulation: 2 mm thick Legs: made of SS 304 – 04 Nos. Jacket thickness: 4 mm Vessel inside finish: mirror finish Outside finish: matt finish		
3.	Check the specification of motor	Top entry saw cutter type Make : Hindustan 2 HP 960 RPM / 415 V FLP		



PROTOCOL No.:

S.No.	TEST PARAMETER	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE	
4.	Check the specification of product temperature sensor	Make: radix Type: PT – 100 MOC: SS 316 L End connection: 1/2" BSP Temperature range (0 – 200)			
5.	Check the specification of product outlet valve	Make : Cipriani Size : 25 mm Type : Pneumatic operated redial diaphragm valve			
6.	Check the specification of safety valve.	make : Inoxpa Size : 25 mm dia. TC			
7.	Check the specification of air vent on jacket	Size : 25 mm dia. TC end			
8.	Check the specification of hot steam inlet	25 mm TC with pneumatic operated piston			
9.	Check the specification of hot steam outlet	25 mm TC with pneumatic operated piston valve Make: Spirex			
10.	Check the specification of Spray ball	Make : Atul End conn-3/4" BSP SS 316 L			
11.	Check the specification of pressure gauge	Make: Baumer MOC: SS 304 Range: 0- 10 kg/cm ² End conn-1/4" BSP			
C.	C. Water phase heating vessel				
1.	Check any physical damage to water phase heating vessel.	There should no damage to the water phase vessel			



PROTOCOL No.:

S.No.	TEST PARAMETER	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
2.	Check specification of multi mixing vessel.	Working capacity: 350 lt. Gross capacity: 420 lt. Min. batch size: 105 lt. Shell thickness: 5 mm Vessel top: 2 mm thick SS 316 L with flat top 1/3 open able Insulation: 2 mm thick Legs: made of SS 304 – 04 Nos. Jacket thickness: 4 mm Vessel inside finish: mirror finish Outside finish: matt finish		
3.	Check the specification of motor	Top entry marine type impeller Make : Hindustan 2 HP 960 RPM / 415 V FLP		
4.	Check the specification of product sensor	Make: radix Type: PT – 100 MOC: SS 316 L End connection: 1/2" BSP Temperature range (0 – 200)		
5.	Check the specification of product outlet valve	Make : Cipriani Size : 25 mm bottom valve Type : Pneumatic operated redial diaphragm valve		
6.	Check the specification of CIP.	make: Inoxpa Size: 50 / 25 mm With pneumatic operated butterfly valve		
7.	Check the specification of safety valve.	make : Inoxpa Size : 25 mm dia. TC		
8.	Check the specification of air vent on jacket	Size : 25 mm dia. TC end		
9.	Check the	Make : Atul		



PROTOCOL No.:

S.No.	TEST PARAMETER	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
	specification of Spray ball	End conn-3/4" BSP SS 316 L		
10.	Check the specification of hot steam inlet	25 mm TC with pneumatic operated piston valve Make: Spirex		
11.	Check the specification of hot steam outlet	25 mm TC with pneumatic operated piston valve Make: Spirex		
12.	Check the specification of pressure gauge	Make: Baumer MOC: SS 304 Range: 0- 10 kg/cm ² End conn-1/4" BSP		
13.	Air vent	Size: 25 mm TC end MOC: 304		
D.	Working Storage	e vessel		
1.	Check any physical damage to working storage vessel.	There should no damage to the working storage vessel.		
2.	Check specification of multi mixing vessel.	Working capacity: 600 lt. Gross capacity: 690 lt. Min. batch size: 180 lt. Shell thickness: 2 mm flat Legs: made of SS 304 – 04 Nos. Vessel inside finish: mirror finish Outside finish: matt finish		
3.	Check the specification of CIP.	make : Inoxpa Size : 50 / 25 mm With pneumatic operated butterfly valve		
4.	Check the specification of product outlet valve valve Type: Pneumatic operated redial diaphragm valve			
5.	Check the	Make : Cipriani		



PROTOCOL No.:

S.No.	TEST PARAMETER	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
	specification of product inlet valve	Size: 40 mm TC end Type: Pneumatic operated butterfly valve		
6.	Check the specification of Spray ball	Make : Atul End conn-3/4" BSP SS 316 L		
E.	Transfer Pump ((Lobe Pump)		
1.	Check any physical damage to transfer Pump	There should no damage to the transfer pump		
2.	Check specification of Transfer Pump	Type: Twin Lobe MOC: SS 316 L Inlet / outlet-25 mm TC Motor: 2.0 HP / 415V / 760 RPM Make: Hindustan		
F.	Electrical Contro	ol panel		
1.	Check any physical damage to electrical control panel	There should no damage to the electrical control panel		
2.	Check specification of Transfer Pump	Make : PTPL		
G.	Vacuum Pump			
1.	Check any physical damage to Vacuum Pump	There should no damage to the Vacuum Pump		
2.	Check specification of Transfer Pump	Type: Water Ring Size: 3 HP / 415V / 2800 RPM Make: Hindustan 1015M566260		
H.	Working Platfor	m		
1.	Check any	There should no damage to the		



PROTOCOL No.:

S.No.	TEST PARAMETER	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
	physical damage to Working platform	working platform		
2.	Check specification of Transfer Pump	Size: 4000 X 1000 X 1000 mm MOC: SS3 304		
I.	Meter In Pump Check any	There should no damage to the		
1.	physical damage to meter in pump	meter in pump		
2.	Check specification of meter in pump	Type: Reciprocating Plunger Type, volume adjustable Motor: 1.0 HP / 415V / 1440 RPM Make: Hindustan Size: 25 mm TC 1015M56504		
J.	Homogenizer			
1.	Check any physical damage to Homogenizer	There should no damage to the Homogenizer		
2.	Check specification of meter in pump	Motor: 1.0 HP 2800/ RPM / 440V Mounting- bottom of the vessel in conical portion Make: Hindustan 1015M56504		
3.	VFD , POWER – 10hp , 50Hz FLP Voltage : AC 3 phase Make Siemens			



PROTOCOL No.:

Checked By Production Sign & Date	Verified By Quality Assurance Sign & Date
Inference:	
	Reviewed By Manager QA Sign & Date



PROTOCOL No	١.	
-------------	----	--

MULTI MIX MANUFACTURING PLANT

8.3 UTILITIES REQUIRED:

PARAMETERS	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
Electricity	3 Phase Plus Earthing, 5 Wire Line Up To The Panel Board Terminal. Voltage- 440 ± 10% V -18.5 HP Frequency- 50 ± 3% Hz		
Steam	Shall be properly connected and Identified		
Cooling Water	Shall be properly connected and Identified		
Compressed Air	6.0 kg/cm ²		
Vacuum	6-8 bar		

8.4 MATERIAL OF CONSTRUCTION:

S.No.	COMPONENTS	МОС	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
1.	All contact Parts	AISI SS 316		
2.	Shell	AISI SS 316		
3.	Bottom	AISI SS 316		
4.	Тор	SS 316		
5.	Flange	AISI SS 316L		
6.	Gaskets	Food Grade		
7.	Jacket	SS 304		
8.	Legs	SS 304		
9.	Insulation	Fiber Wool (Mineral Wool Asbestos Free		



PROTOCOL No.:

S.No.	COMPONENTS	MOC	OBSERVATION	OBSERVED BY ENGINEERING SIGN & DATE
10.	Transfer Pump (LOBE PUMP)	SS316 (Product Contact		
		Parts)		
11.	Meter-in Pump	SS 316 (Product Contact		
		Part)		
12.	 Working Platform: Square Pipe Frame Work with Top SS Dimple Sheet Ladder Railing is provided on all three sides of the Platform. Legs (Round Pipe Legs are provided) Platform (GMP Standard Paint Free) 	SS 304		
13.	 Batch Storage Vessel: Shell 16 SWG (Cylindrical type) Bottom 16 SWG (Conical type welded with shell) Top 16 (Loose type Lid) Legs with PU Castor wheel with Bracket 	SS 316 (All Contact Parts) SS 304		

Checked By	Verified By
Production	Quality Assurance
Sign & Date	Sign & Date
Inference:	
•••••	••••••
	Daviewed Dv
	Reviewed By
	Manager QA
	Sign & Date



PROTOCOL No.:

MULTI MIX MANUFACTURING PLANT

8.5 INSTALLATION CHECKS:

S.No.	SPECIFICATION	OBSERVATION	OBSERVED BY (ENGINEERING) (SIGN & DATE)
1.	Verify that major components are		
	protected from shock and there is no		
	physical damage		
2.	Check the proper mechanical		
	installation of multi mix plant		
3.	Check the proper electrical		
	installation of multi mix plant		
4.	Verify that all piping and electrical		
	connection have been done		
	according to the drawings		
5.	Verify that the required electrical		
	connections are tight, weather		
	proofed and grounded.		
6.	Equipment identification nameplate		
	is visible		
7.	Unit are installed on foundation and		
	secured in place as per		
	manufacturer recommendations		
8.	All access ports are examined and		
	cleared of any debris.		
9.	Wiring diagram is glued of taped to		
	inside section of control panel		
10.	Verify that there is sufficient room		
	for servicing provided		
11.	Check the equipment is free from		
	any defects		
12.	Check the finishing of product	`	
	contact parts		
13.	Verify that the 'As Built' Drawings		



PROTOCOL No.:

S.No.	SPECIFICATION	OBSERVATION	OBSERVED BY (ENGINEERING) (SIGN & DATE)
	are Complete and represent the		
	design concept		
14.	There should not be any loose		
	fasteners		
15.	There should not be any loose		
	electrical connection		
16.	There should not be any damage		

Checked By	Verified By		
Production	Quality Assurance		
Sign & Date	Sign & Date		
Inference:			
••••••	••••••		
••••••			
	Reviewed By		
	Manager QA		
	Sign & Date		



MULTI MIX MANUFACTURING PLANT

PROTOCOL No.:

8.6 SAFETY:

CHECKS	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING (SIGN & DATE)
Pressure Switches	For product safety and Gasket		
	safety		
Electrical wiring and	Electrical wiring should be as		
Earthing	per approved drawings. Double		
	external earthing to control		
	machine (panel and motors)		
	and operator should be		
	provided.		
Motor overload relay – The	Should be provided For Motor		
switchgear shall trip if	safety		
overloaded			
Emergency off: To stop the	Should be provided For		
process immediately	equipment and operator safety		
Safety clamps: For holding	For equipment and product		
product container and acting	safety		
as a jumper for earth			
continuity of PC and RC			
MCB	MCB is provided so that when		
	there is an overload in current		
	or any short circuit then the		
	MCB trips.		
Mechanical Guard	Mechanical guard for all		
	rotating parts.		
Joints	Welding of joints without any		
	welding burrs		
Metal Parts	All the metal parts should be		



PROTOCOL No.:

CHECKS	ACCEPTANCE CRITERIA	OBSERVATION	OBSERVED BY ENGINEERING (SIGN & DATE)
	properly grounded without any		
	sharp Edges.		
Leveling And Balancing	Equipment should be properly		
	Balanced & Leveled		
Cooling arrangement	Cooling arrangement provided		
	for motor shaft and seal with		
	special flow switch to detect		
	water, so only if cooling water		
	is circulated then only		
	homogenizer will operate, this		
	ensures long life of the		
	equipment.		
Operating Panel	Is provided on platform of the		
	Equipment for operator safety.		
Dimple sheet	Platform to ensure proper grip		
	during walking and railing is		
	provided as safety feature.		
Checked By Production Sign & Date		Verified By Quality Assurand Sign & Date	ce

Production Sign & Date	Quality Assurance Sign & Date		
Inference:			
•••••••••••••••••••••••••••••••••••••••	••••••		
	Reviewed By Manager QA Sign & Date		



PROTOCOL No.:

MULTI MIX MANUFACTURING PLANT

8.7 Spare parts: List of spare parts to be attached (if any)

9.0 REFERENCES:

- Validation Master Plan
- Schedule-M "Good Manufacturing Practices and Requirements of Premises, Plant and Equipment for Pharmaceutical Products."
- WHO Essential Drugs and Medicines Policy, QA of Pharmaceuticals, Vol-2 Good Manufacturing Practices and Inspection.

10.0 DOCUMENTS TO BE ATTACHED:

- Technical details for Equipment Requirement with Engineering Drawings.
- Certificate of MOC
- Calibration certificates

11.0	DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY:
	•••••••••••••••••••••••••••••••••••••••
4.0	
12.0	CHANGE CONTROL, IF ANY:
13.0	REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY):

PHARMA DEVILS

SWG

cGMP

:

HP

Hz

VFD

RPM

FLP

PTPL

CIP

RH

MMP

P & ID

V

INSTALLATION QUALIFICATION PROTOCOL CUM REPORT

FOR

MULTI MIX MANUFACTURING PLANT

PR	വ	Γ	C	ΩI	ן . ז	No	•

PHARI	MA DEVILS		WIGHT WITA WIANGFACTORING I LANT	
14.0	CONCLUS	SION:		
	•••••	•••••	••••••	•••••
	•••••	•••••	••••••	•••••
	•••••	•••••	••••••	•••••
	•••••	•••••	••••••	•••••
15.0	RECOMM	IENDA'	ΓΙΟΝ:	
	•••••	•••••	••••••	•••••
	•••••	•••••		•••••
	•••••	•••••		•••••
	•••••	•••••	••••••	•••••
		^		
16.0	ABBREV	'ATIO	N:	
	WHO	:	World health organization	
	MOC	:	Material of Construction	
	MCB	:	Miniature circuit Breaker	
	SS	:	Stainless Steel	
	PU		Poly urethane	

Slandered Wire gauge

Variable frequency Drive

Revolution per minute

Horse power

Flame proof

Cleaning in place

Relative Humidity

Hertz

Volt

Current Good Manufacturing Practice

Propack Technologies private limited

Multi mix manufacturing plant

Piping and Instrumentation diagram



PRO	TO	\mathbf{COL}	No.:
-----	----	----------------	------

MULTI MIX MANUFACTURING PLANT

17.0 POST APPROVAL INITIATED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			