



STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Endotoxin Challenge Test	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

1.0 OBJECTIVE:

To lay down a procedure for Endotoxin Challenge Test.

2.0 SCOPE:

This SOP is applicable for Endotoxin Challenge Test in Microbiology Section of Quality Control Area.

3.0 RESPONSIBILITY:

Officer / Executive – Microbiology

4.0 ACCOUNTABILITY:

Head – QC

5.0 ABBREVIATIONS:

CSE	Control Standard Endotoxin
No.	Number
Ltd.	Limited
LRW	LAL Reagent Water
LAL	Limulus Amebocyte Lysate
QA	Quality Assurance
QC	Quality Control
SOP	Standard Operating Procedure

6.0 PROCEDURE:

6.1 MATERIAL AND INSTRUMENTS:

- 6.1.1 Limulus Amebocyte Lysate Reagent.
- 6.1.2 Endotoxin Indicator Vial (100000 EU/Vial) or as per received from vendor.
- 6.1.3 LAL Reagent Water.
- 6.1.4 Ampoules or Vial according to requirement.
- 6.1.5 Depyrogenated Dilution Tubes (12 x 75 mm, 16 x 100mm)
- 6.1.6 Depyrogenated LAL Assay Tube (10 x 75 mm)
- 6.1.7 Micropipette with Pyrogen free tip (20-200 µl)
- 6.1.8 Micropipette with Pyrogen free tip (100-1000 µl)
- 6.1.9 Vortex Mixer.
- 6.1.10 Heating Block.



STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Endotoxin Challenge Test	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

6.2 Switch “ON” the Heating Block and set the temperature at $37^{\circ} \pm 1^{\circ}\text{C}$.

6.3 PREPARATION OF CHALLENGE VIALS:

6.3.1 Reconstitute the challenge vial of Endotoxin with 1ml LRW or as per vendor COA to yield 100000 EU/ml and vortex according to manufacturers instructions.

6.3.2 Transfer 0.1ml aliquot into ampoules or vials used for endotoxin challenge test.

6.3.3 Keep the ampoules/vial in Laminar Air Flow hood for overnight. Each ampoule/vial now contains 10,000 EU. Mark the above prepared ampoules/vial as 1 to 10 numbers.

6.3.4 Keep at least 1 ampoules/vial as positive control (do not expose through oven/tunnel).

6.3.5 Mark the remaining ampoules/vial as NPC.

6.3.6 Expose these ampoules/vials to appropriate location in DHS/Tunnel as per depyrogenation cycle.

6.4 DILUTION OF POSITIVE CONTROL AMPOULE / VIAL:

6.4.1 Reconstitute the ampoules / vial with 1ml LRW and vortex vigorously for 5 minutes and each subsequent dilution for 3-4 minutes.

6.4.2 Now the concentration of Endotoxin in the PPC will be 10000 EU/ml.

6.4.3 Prepare 1:100 dilution of the above to obtain 100 EU/ml.

6.4.4 Prepare 1:100 dilution of the above to obtain 1 EU/ml.

6.4.5 From the above 1 EU/ml preparation, prepare a two fold dilution series upto 2λ , λ , $\lambda/2$, $\lambda/4$. where λ =Labelled Lysate sensitivity, if $\lambda=0.125$ EU/ml

DILUTION TABLE

Sr. No.	Endotoxin	LRW	Endotoxin Con.(EU/ML)
1.	10000 EU/Vial	1 ml	10000 EU/ml
2.	0.1 ml of 10000 EU/ml	9.9 ml	100 EU/ml
3.	0.1 ml of 100 EU/ml	9.9 ml	1 EU/ml 8λ
4.	0.5 ml of 1 EU/ml	0.5 ml	0.5 EU/ml 4λ
5.	0.5 ml of 0.5 EU/ml	0.5 ml	0.25 EU/ml 2λ
6.	0.5 ml of 0.25 EU/ml	0.5 ml	0.125 EU/ml λ



STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Endotoxin Challenge Test	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

7.	0.5 ml of 0.125 EU/ml	0.5 ml	0.06 EU/ml $\lambda/2$
8.	0.5 ml of 0.06 EU/ml	0.5 ml	0.03 EU/ml $\lambda/4$

6.5 DILUTION OF NPC AMPOULES / VIAL :

- 6.5.1** Reconstitute each of the ampoules/vial with 1ml LRW and vortex vigorously for 5 minutes.
- 6.5.2** It is assumed that three log reduction is achieved after exposure of ampoules/vial in oven/Tunnel, the Endotoxin concentration in the vial is 10 EU/ml.
- 6.5.3** Prepare 1:10 dilution of each ampoules/vial to obtain 1 EU/ml.
- 6.5.4** Further prepare 1:8 dilution to obtain 0.125 EU/ml.

DILUTION TABLE

Tube No.	Endotoxin Indicator	LAL reagent Water	Endotoxin Concentration
1.	0.1 ml of 10 EU/ml	0.9 ml	1 EU/ ml
2.	0.1 ml of 1 EU/ml	0.7 ml	0.125 EU/ml

6.6 LAL TEST PROCEDURE:

- 6.6.1** Test the two fold dilution series prepared from the positive controls ampoules/vial in duplicate.
- 6.6.2** Test the 0.125 EU/ml dilutions prepared from each of ampoules/vial.
- 6.6.3** Test should be carried out in clean depyrogenated 10x 75 mm assay tubes only.
- 6.6.4 Procedure for Positive Control:**

S. No.	Dilutions	CSE Dilution Used	LRW	Lysate in μ l	No. of Replicates
1.	2 λ	100 μ l of 2 λ	–	100 μ l	2
2.	λ	100 μ l of λ	–	100 μ l	2
3.	$\lambda/2$	100 μ l of $\lambda/2$	–	100 μ l	2
4.	$\lambda/4$	100 μ of $\lambda/4$	–	100 μ l	2
5.	Negative Water Control (NWC)	–	100 μ l	100 μ l	2

- 6.6.5 Procedure for Negative Product Control (Challenged ampoules/vial):**



STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Endotoxin Challenge Test	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

Location No.	Dilutions	Endotoxin Indicator Dilution (0.125 EU/ml)	Lysate in μ l	No. of Replicates
1.	NPC	100 μ l	100 μ l	2
2.	NPC	100 μ l	100 μ l	2
3.	NPC	100 μ l	100 μ l	2
4.	NPC	100 μ l	100 μ l	2
5.	NPC	100 μ l	100 μ l	2
6.	NPC	100 μ l	100 μ l	2
7.	NPC	100 μ l	100 μ l	2
8.	NPC	100 μ l	100 μ l	2
9.	NPC	100 μ l	100 μ l	2
10.	NPC	100 μ l	100 μ l	2

6.7 CALCULATION:

Log Reduction = Log value of recovered Endotoxin from positive control – Log value of recovered sample from heat treated sample

Recovered EU/ml from heat treated sample (X) = Reciprocal of last dilution \times λ

Recovered EU/ml from Positive Control (Y) = Reciprocal of last dilution \times λ

6.8 INTERPRETATION OF RESULTS/ACCEPTANCE CRITERIA:

6.8.1 Test results are valid if recovery of Endotoxin in unexposed vials is within a two fold dilution of the labeled claim.

6.8.2 The depyrogenation cycle is considered as successfully validated if there is more than 3 log reduction is achieved in challenge Endotoxin vials exposed into Oven/ tunnels at specified place.

6.8.3 For a valid Depyrogenation cycle, the PPC must be positive and NPC's must be negative indicating a greater than 3-log reduction of endotoxin.

6.9 Record the "Endotoxin Challenge Test Record in **Annexure – I**.



PHARMA DEVILS
MICROBIOLOGY DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Endotoxin Challenge Test	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

ANNEXURE-I
ENDOTOXIN CHALLENGE TEST REPORT

Report. No.			
Date of Vial Exposure		No. of Exposed Vial	
Date of Testing		Date of Release	
Tunnel/oven Location		Tunnel/oven ID No.	
Performed By		Date of Analysis	
Shift		BET Kit	
Heating Block Temperature	37°C ± 1°C	Manufacturer	
Incubation Started at		Incubation Time	60 ± 2 Minutes
Incubation Completed at			

REAGENTS DETAILS:

Reagent Details	Lysate	CSE	LRW
Lot No.			
Sensitivity/Potency			
Date of Reconstitution /Opening			
Use Before			
Expiry Date			
Manufacturer			

DILUTION PREPARATION FOR HEAT TREATED VIAL (10 EU/ VIAL ASSUMED)

S.No.	Test Dilution	Test	LRW	Endotoxin Concentration
-	-	10 EU/Vial or Ampoule Assumed	1 ML	10 EU/ml
1.	1:10	0.1ml	0.9 ml	1 EU/ml
2.	1:8	0.1ml	0.7 ml	0.125 EU/ml
Reciprocal of last dilution		10 X 8 =80		

DILUTION PREPARATION FOR POSITIVE CONTROL VIAL (10,000 EU/ VIAL)

S.No.	Test Dilution	Test	LRW	Endotoxin Concentration
-	-	10,000 EU/Vial or Ampoule	1ML	10,000 EU/ml
1.	1:10	0.1ml	0.9ml	10,00 EU/ml
2.	1:10	0.1ml	0.9ml	100 EU/ml
3.	1:10	0.1ml	0.9ml	10 EU/ml



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Endotoxin Challenge Test	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

4.	1:10	0.1ml	0.9ml	1 EU/ml
5.	1:2	0.5ml	0.5ml	0.5Eu (4λ) EU/ml
6.	1:2	0.5ml	0.5ml	0.25Eu (2λ) EU/ml
7.	1:2	0.5ml	0.5ml	0.125Eu (λ) EU/ml
Reciprocal of last dilution		10 x 10 x 10 x 10 x 2 x 2 x 2 = 80,000		

CALCULATION OF LOG REDUCTION:

Vial No.	Recovered EU/ml = Reciprocal of last dilution x λ		Log reduction = Recovered EU from positive control – Recovered EU from Heat treated vial
	Heat treated vial	Positive control	
1.			
2.			
3.			
4.			
5.			
6.			
7.			
8.			
9.			
10.			
11.			
12.			

PREPARATION OF TEST SOLUTION:

Solution	Tube No.	Product Dilution (1:80)	LAL Water	LAL Reagent	Total Volume
Negative Product Control	01	100 μl	----	100 μl	200 μl
	02	100 μl	----	100 μl	200 μl
Negative Water Control	01	----	100 μl	100 μl	200 μl
	02	----	100 μl	100 μl	200 μl

OBSERVATIONS:

Results of Negative Product Control & Negative Water Control:

	Tube No.	Observation Vial No.											
		1	2	3	4	5	6	7	8	9	10	11	12
Negative Product Control	1												
	2												
Negative Water Control	1												
	2												

Results of Positive Control:



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Endotoxin Challenge Test	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

Tube No.	Endotoxin Dilution				Negative Water Control
	2λ	λ	$\lambda/2$	$\lambda/4$	
01					
02					

+Ve: Gel Formation

-Ve: No Gel Formation

Remark: The depyrogenation cycle of Tunnel complies (three log reduction)/does not comply for Endotoxin challenge Test

Microbiologist:

Checked By:

Date:

Date: