



**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Production	<b>SOP No.:</b>
<b>Title:</b> Operation and Cleaning of Planetary Mixer	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

**1.0 OBJECTIVE:**

To lay down a Procedure for Operation and Cleaning of Planetary Mixer.

**2.0 SCOPE:**

This SOP is applicable to Operation & Cleaning of Planetary Mixer used in Soft Gelatin Capsule section.

**3.0 RESPONSIBILITY:**

Officer / Executive Production

**4.0 ACCOUNTABILITY:**

Head Production

**5.0 ABBREVIATIONS:**

BMR	Batch Manufacturing Record
IPA	Isopropyl Alcohol
Ltd.	Limited
PLM	Planetary Mixer
Pvt.	Private
QC	Quality Control
SOP	Standard Operating Procedure

**6.0 PROCEDURE:**

**6.1 OPERATION:**

**6.1.1** Ensure that the equipment & area is cleaned.

**6.1.2** Take the line clearance from QA as per "Line Clearance" SOP and enter the details in BMR.

**6.1.3** Affix the status label on the machine.

**6.1.4** Ensure that the temperature and humidity of area is within limit as specified in BMR.

**6.1.5** Switch 'ON' the PLM and lift the stirrer by turning the Knob switch anticlockwise to up position.

**6.1.6** Load the medicament in PLM Bowl.

**6.1.7** Take the PLM Bowl below stirrer and adjust the position.

**6.1.8** Take down the stirrer by turning the knob switch clockwise direction to down position.

**6.1.9** Switch 'ON' the stirrer by pressing the slow/fast push button and set the time on timer as



**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Production	<b>SOP No.:</b>
<b>Title:</b> Operation and Cleaning of Planetary Mixer	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

specified in BMR.

**6.1.10** Switch 'ON' the vacuum pump and open the vacuum valve upto 650 to 715 mm/Hg for de-aeration of medicament or as defined in BMR.

**6.1.11** After completion of mixing switch 'OFF' the stirrer by pressing "Red" push button.

**6.1.12** Open the vacuum valve and release the vacuum.

**6.1.13** After completion of process lift the stirrer by turning the Knob switch anticlockwise to up position.

**6.1.14** Take out the PLM bowl from the PLM.

**6.1.15** Enter the operation details in "Machine Utilization Record"

**6.2 CLEANING:**

**6.2.1 For Product Changeover (Type B Cleaning):**

**6.2.1.1** Ensure that mains switch is "OFF".

**6.2.1.2** Dismantle the stirrer and PLM Bowl of PLM and take to washing area.

**6.2.1.3** Clean it with 2 % v/v of Extran MA-02 solution.

**6.2.1.4** Wash with sufficient quantity of potable water.

**6.2.1.5** Finally rinse with sufficient quantity of purified water.

**6.2.1.6** Wipe the stirrer and PLM Bowl with dry lint free cloth and finally mop with 70% v/v IPA solution.

**6.2.1.7** Put the stirrer in a SS tray.

**6.2.1.8** Officer/Executive QA shall check visually & collect the swab sample with intimation slip cum analysis report and send to QC for analysis.

**6.2.1.9** Use the PLM after receiving swab test intimation slip cum analysis report from QC showing negative identification.

**6.2.1.10** Assemble the stirrer in PLM.

**6.2.1.11** Clean the remaining part of PLM with the help of wet cloth followed by lint free cloth.

**6.2.1.12** Mop the PLM with 70% v/v IPA solution.

**6.2.1.13** Affix a label as "Cleaned" on the PLM.

**6.2.1.14** Enter the cleaning details in "Machine Utilization Record"

**6.2.2 For Batch Changeover (Type A Cleaning) :**



**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Production	<b>SOP No.:</b>
<b>Title:</b> Operation and Cleaning of Planetary Mixer	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

**6.2.2.1** Clean the PLM by mopping with dry lint free cloth.

**6.2.2.2** Mop the PLM with 70% v/v IPA solution.

**6.2.2.3** Affix the status label on the PLM.

**6.2.2.4** Enter the cleaning details in “Machine Utilization Record”.

**6.2.3 Frequency of Cleaning:**

**6.2.3.1 For B Type Cleaning:**

- (a) Product to Product Changeover.
- (b) After 5 batches of the same product.
- (c) If cleaned equipment is kept idle more than 72 hours.
- (d) If Dirty equipment is kept idle more than 24 hours.
- (e) After any Maintenance of Product Contact Parts.
- (f) Changeover of one Batch to Next Batch of the same Product with Descending Potency.

**6.2.3.2 For A Type Cleaning:**

- (a) Changeover from one Batch to Next Batch of the same Product with Same Potency.
- (b) Changeover from one batch to next Batch of the same Product with Higher Potency.

**7.0 ANNEXURES:**

Not Applicable

**ENCLOSURES:** SOP Training Record

**8.0 DISTRIBUTION:**

- Controlled Copy No. 01      Quality Assurance
- Controlled Copy No. 02      Production
- Master Copy                      Quality Assurance

**9.0 REFERENCES:**

Not Applicable.



**PHARMA DEVILS**

PRODUCTION DEPARTMENT

**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Production	<b>SOP No.:</b>
<b>Title:</b> Operation and Cleaning of Planetary Mixer	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

**10.0 REVISION HISTORY:**

**CHANGE HISTORY LOG**

<b>Revision No.</b>	<b>Change Control No.</b>	<b>Details of Changes</b>	<b>Reason for Change</b>	<b>Effective Date</b>	<b>Updated By</b>