



## STANDARD OPERATING PROCEDURE

<b>Department:</b> Microbiology	<b>SOP No.:</b>
<b>Title:</b> Procedure for Operation and Calibration of Micropipette	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

### 1.0 OBJECTIVE:

To lay down a procedure for Operation and Calibration of Micropipette.

### 2.0 SCOPE:

This SOP is applicable for Operation and Calibration of Micropipette in Microbiology Section of Quality Control Laboratory.

### 3.0 RESPONSIBILITY:

Officer / Executive - Microbiology

### 4.0 ACCOUNTABILITY:

Head – QC

### 5.0 ABBREVIATIONS:

µl	Micro liter
ml	Milliliter
No.	Number
QA	Quality Assurance
QC	Quality Control
SOP	Standard Operating Procedure
RSD	Relative Standard Deviation

### 6.0 PROCEDURE:

- 6.1 Calibrate the Micropipettes for a volume of 20 to 200 µl 100 to 1000 µl and 1000 to 10000 µl.
- 6.2 Set the desired volume using the Volume-Setting Knob e.g. for 20 to 200 µl set 50 µl, 100 µl & 200 µl and for 100 µl to 1000 µl set 100 µl ,500 µl, 1000 µl and 1000 µl to 10000 µl set 1000 µl ,5000 µl, 10000 µl.
- 6.3 Carefully fit the tip into tip cone and pick up the Water for Injections.
- 6.4 Weigh the empty Vial (A).
- 6.5 Pick up the set water volume and add into the vial after weigh of empty vial (A) at temperature  $25\pm 2^{\circ}\text{C}$ .
- 6.6 Observe the weight of vial with water (B).
- 6.7 Remove set tip with pipettes.
- 6.8 Calculate the exact weight of water  $M = B - A$ .
- 6.9 Divide the value of M with the D i.e. 0.997 gm/ml (D is specific density of water at  $25\pm 2^{\circ}\text{C}$ ) to obtain actual volume in ml /µL.





**PHARMA DEVILS**  
MICROBIOLOGY DEPARTMENT

**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Microbiology	<b>SOP No.:</b>
<b>Title:</b> Procedure for Operation and Calibration of Micropipette	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

**9.0 REFERENCES:**  
Not Applicable.

**10.0 REVISION HISTORY:**

**CHANGE HISTORY LOG**

<b>Revision No.</b>	<b>Change Control No.</b>	<b>Details of Changes</b>	<b>Reason for Change</b>	<b>Effective Date</b>	<b>Updated By</b>



**PHARMA DEVILS**  
MICROBIOLOGY DEPARTMENT

**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Microbiology	<b>SOP No.:</b>
<b>Title:</b> Procedure for Operation and Calibration of Micropipette	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

**ANNEXURE – I**  
**CALIBRATION OF MICROPIPETTE (20-200 µL)**

<b>Date of Calibration</b>		<b>Next Due Date</b>	
<b>Micropipette Make</b>		<b>Micropipette ID</b>	
<b>Volume Range</b>			

**OBSERVATION TABLE:**

S.No.	Empty Vial wt. (A)	Volume Pipette in µl	Wt. of Vial & Water (B)	Wt. of Water (M=B-A)	Density of Water at 25°C	Actual Volume (V=M/D)
<b>For 50 µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	
<b>For 100 µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	
<b>For 200 µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	

RSD should be not more then 2%.

**LIMITS:**

VOLUME IN µL	TOLERANCE LIMIT IN µL	
	MINIMUM	MAXIMUM
50	49.5	50.5
100	99.0	101.0
200	198.0	202.0

**Remarks:** The Calibration of Micropipettes Complies / Doesn't Comply as per requirement.

**Done by:**  
**Date:**

**Checked By:**  
**Date:**

**Approved By:**  
**Date:**



**PHARMA DEVILS**  
MICROBIOLOGY DEPARTMENT

**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Microbiology	<b>SOP No.:</b>
<b>Title:</b> Procedure for Operation and Calibration of Micropipette	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

**ANNEXURE – II**  
**CALIBRATION OF MICROPIPETTE (100-1000 µL)**

<b>Date of Calibration</b>		<b>Next Due Date</b>	
<b>Micropipette Make</b>		<b>Micropipette ID</b>	
<b>Volume Range</b>			

**OBSERVATION TABLE:**

S.No.	Empty Vial wt. (A)	Volume Pipette in µl	Wt. of Vial & Water (B)	Wt. of water (M=B-A)	Density of Water at 25°C	Actual Volume (V=M/D)
<b>For 100 µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	
<b>For 500µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	
<b>For 1000µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	

**RSD should be not more then 2%.**

**LIMITS:**

VOLUME IN µL	TOLERANCE LIMIT IN µL	
	Minimum	Maximum
100	99.0	101.0
500	495.0	505.0
1000	990.0	1010.0

**Remarks:** The Calibration of Micropipettes Complies / Doesn't Comply as per requirement.

**Done By:**  
**Date:**

**Checked By:**  
**Date:**

**Approved By:**  
**Date:**



**PHARMA DEVILS**  
MICROBIOLOGY DEPARTMENT

**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Microbiology	<b>SOP No.:</b>
<b>Title:</b> Procedure for Operation and Calibration of Micropipette	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

**ANNEXURE – III**  
**CALIBRATION OF MICROPIPETTE (1000-10000 µL)**

<b>Date of Calibration</b>		<b>Next Due Date</b>	
<b>Micropipette Make</b>		<b>Micropipette ID</b>	
<b>Volume Range</b>			

**OBSERVATION TABLE:**

S.No.	Empty Vial wt. (A)	Volume Pipette in µl	Wt. of Vial & Water (B)	Wt. of water (M=B-A)	Density of Water at 25°C	Actual Volume IN µL (V=M/D)
<b>For 1000 µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	
<b>For 5000 µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	
<b>For 10000 µl</b>						
1.						
2.						
3.						
					<b>Average</b>	
					<b>% RSD</b>	

RSD should be not more than 2%.

**LIMITS:**

VOLUME IN µL	TOLERANCE LIMIT IN µL	
	Minimum	Maximum
1000	990.0	1010.0
5000	4950.0	5050.0
10000	9900.0	10100.0

**Remarks:** The Calibration of Micropipettes Complies / Doesn't Comply as per requirement.

**Done By:**  
**Date:**

**Checked By:**  
**Date:**

**Approved By:**  
**Date:**