



**PERFORMANCE QUALIFICATION PROTOCOL
FOR
SINGLE HEAD AUTOMATIC POWDER FILLING
MACHINE**

PROTOCOL No.:

**PERFORMANCE QUALIFICATION
PROTOCOL
FOR
SINGLE HEAD AUTOMATIC POWDER
FILLING MACHINE**

EQUIPMENT ID. No.	
LOCATION	Ointment Section
DATE OF QUALIFICATION	
SUPERSEDE PROTOCOL No.	NIL



PHARMA DEVILS

**PERFORMANCE QUALIFICATION PROTOCOL
FOR
SINGLE HEAD AUTOMATIC POWDER FILLING
MACHINE**

PROTOCOL No.:

PROTOCOL CONTENTS

S.No.	TITLE	PAGE No.
1.0	Protocol approval	3
2.0	Objective	4
3.0	Scope	4
4.0	Responsibility	4
5.0	Equipment details	5
6.0	System description	5
7.0	Reason for qualification	5
8.0	Site of study	5
9.0	Frequency of Qualification	5
10.0	Pre-qualification requirement	6
11.0	Tests & checks	7-8
12.0	Check list for all test & checks	8
13.0	References	8
14.0	Documents to be attached	8
15.0	Non compliance	8
16.0	Deviation from pre-defined specification, if any	9
17.0	Change control, if any	9
18.0	Abbreviations	9



PHARMA DEVILS

**PERFORMANCE QUALIFICATION PROTOCOL
FOR
SINGLE HEAD AUTOMATIC POWDER FILLING
MACHINE**

PROTOCOL No.:

1.0 PROTOCOL APPROVAL:

PREPARED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER/EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OPERATING MANAGER (QUALITY ASSURANCE)			
HEAD (QUALITY CONTROL)			
HEAD (ENGINEERING)			
HEAD (PRODUCTION)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			



**PERFORMANCE QUALIFICATION PROTOCOL
FOR
SINGLE HEAD AUTOMATIC POWDER FILLING
MACHINE**

PROTOCOL No.:

2.0 OBJECTIVE:

- To carry out the Performance Qualification of Manufacturing Tank used for manufacturing of liquid Preparation.
- To Provide Documented Verification that the Equipment as connected with ancillary system is suitable for indented purpose and produced product as per pre defined Acceptance Criteria

3.0 SCOPE:

- The scope of this qualification protocol is limited to qualification of Manufacturing Tank Installed **Ointment Section**.

4.0 RESPONSIBILITY:

The Validation Group, comprising of a representative from each of the following departments shall be responsible for the overall compliance of this Protocol cum Report:

DEPARTMENTS	RESPONSIBILITIES
Quality Assurance	<ul style="list-style-type: none"> • Preparation, Review, and Authorization of Performance Qualification Protocol. • Co-ordination with Production and Engineering to carryout Performance Qualification Activity..
Production	<ul style="list-style-type: none"> • Review & Approval of Performance Qualification Protocol. • To Co-ordinate and support for execution of Operational Qualification study as per Protocol.
Quality Control	<ul style="list-style-type: none"> • Review & Pre Approval of Performance Qualification Protocol. • Analytical Support (Microbiological Testing / Chemical Analysis)
Engineering	<ul style="list-style-type: none"> • Review & Pre Approval of Performance Qualification Protocol. • To co-ordinate and support Performance Qualification Activity.

 PHARMA DEVILS	PERFORMANCE QUALIFICATION PROTOCOL FOR SINGLE HEAD AUTOMATIC POWDER FILLING MACHINE	PROTOCOL No.:
--	--	----------------------

5.0 EQUIPMENT DETAILS:

Equipment Name	Single Head Automatic Powder Filling Machine
Equipment	
Manufacturer's Name	D.M. Engineering Co.
Model	NA
Location of Installation	Ointment Section

6.0 SYSTEM DESCRIPTION

Modal SHPF-1 Single Head Powder Filling Machine is designed to fill powder in pharmaceutical bottles. The basic machine has a fabricated frame with SS cladded table top. The table is fitted with SS covers all around.it houses Electrical Motor, Gear Box, Electrical Panel Box and Control station at convenient place. The Conveyor has adjustable size to match with bottle size.

Empty sterile bottles are conveyor fed via turn table to the star wheel which conveys bottles from conveyor to the filling head for till the desired quantity of powder.

The model SHPF-1 can handle a bottle size range from 22 mm diameter to 85 mm diameter with the appropriate change parts. Bottle height adjustable from 50 mm to 200 mm. filling range starts from 01 gms. To 1000 gms. Can be accommodated with the appropriate change parts.

“DMEC” Automatic Powder Filling is designed for accurate volumetric filling of dry syrup powder, granule substances in quantities from 01 gm to 1000 gm per fills. Any type of container like tin, jar, bottles, bags, pouches, cardboard drums can be utilized for filling the powder. General accuracy of filling ranging between 1.5 % to 2.5 % depending upon the density of powder and quality of powder with control humidity and temperature of the room.

7.0 REASON FOR QUALIFICATION:

- New equipment installed in Ointment Section.
- After any major breakdown or after major modification
- Relocation of Equipment.

8.0 SITE OF STUDY:

- Ointment Section

9.0 FREQUENCY OF QUALIFICATION

- After Every Two years as per Validation Master Plan.

 PHARMA DEVILS	PERFORMANCE QUALIFICATION PROTOCOL FOR SINGLE HEAD AUTOMATIC POWDER FILLING MACHINE	PROTOCOL No.:
--	--	----------------------

10.0 PRE – QUALIFICATION REQUIREMENTS:

Verification for availability, completeness and approval status of all the required relevant documents shall be done and observations shall be recorded in the performance qualification report.

10.1 Verification of Documents:

Record the observations for documents in the below mentioned table.

S.No.	Document Name	Completed Yes/No	Checked By Engineering Sign/Date	Verified By QA Sign/Date
1.	Executed and approved Design Qualification cum report.			
2.	Executed and approved Installation Qualification cum report.			
3.	Executed and approved Operational Qualification cum report.			
4.	PQ Protocol approved.			
5.	SOP for Operation & Cleaning of Single Head Automatic powder Filling Machine.			
6.	SOP for Preventive Maintenance Single Head Automatic powder Filling Machine.			

10.2 Training Record of Validation Team:

- All the persons involved in the execution of qualification activity must be trained in all aspects of the qualification activity including the test methodology, acceptance criteria and safety precautions to be followed during working.
- Verify the training records and record the details in table mentioned in performance qualification report.

11.0 TESTS & CHECKS:

11.1 Equipment Volumetric Capacity Test and Uniformity Of Solution:

11.1.1 Objective:

- The purpose of this test is to demonstrate that Vessel Capacity provided is as agreed with the Equipment supplier and meeting User Requirement and solution prepared is homogeneous as seen visually and active contents are uniform.

 PHARMA DEVILS	PERFORMANCE QUALIFICATION PROTOCOL FOR SINGLE HEAD AUTOMATIC POWDER FILLING MACHINE	PROTOCOL No.:
--	--	----------------------

11.1.2 Equipment / Instrument Used:

- Purified Water
- Calibrated Vessel / QC equipment to measure required quantity for charging Water
- Sodium chloride. (0.9%)

11.1.3 Method Applied

- **Minimum Capacity (400 Ltr) at Minimum & Maximum RPM:** Charge 400 liters to 2000 Liter Tank of purified water using calibrated cylinder/ vessel.
- Add NaCl 3.6 Kg (0.9%) to 400 Ltr. charged vessel.
- Operate the equipment as per SOP on operation of manufacturing vessel.
- Agitate the mixture for 10 min. at Minimum RPM (400).
- After the completion of cycle take 100 ml of sample from top and bottom & send to QC lab for assay.
- Above procedure shall be repeated at Maximum RPM (950).

- **Maximum Capacity (2000 Ltr) at Minimum & Maximum RPM:** Charge 2000 liters to 2000 Liter Tank of purified water using calibrated cylinder/ vessel.
- Add NaCl 18.0 Kg (0.9%) to 2000 Ltr. charged vessel.
- Operate the equipment as per SOP on operation of manufacturing vessel.
- Agitate the mixture for 10 min. at Minimum RPM (400).
- After the completion of cycle take 100 ml of sample from top and bottom & send to QC lab for assay.
- Above procedure shall be repeated at Maximum RPM (950).

Acceptance Criteria:

Uniformity of solution

- **For 400 Ltr. And 2000 Ltr:** Assay of NaCl should be between 0.882% W/V – 0.918% W/V.
- Relative standard deviation shall not be more than 2.0 %.
- Equipment should run trouble free without any problems after charging material up to working volume i.e. 2000 Ltr.

 PHARMA DEVILS	PERFORMANCE QUALIFICATION PROTOCOL FOR SINGLE HEAD AUTOMATIC POWDER FILLING MACHINE	PROTOCOL No.:
--	--	----------------------

12.0 CHECKLIST OF ALL TESTS & CHECKS

S.No.	Name of Test or Check
1.	Equipment Volumetric Minimum Capacity Test and Verification of Uniformity of Mixing at Minimum RPM
2.	Equipment Volumetric Minimum Capacity Test and Verification of Uniformity of Mixing at Minimum RPM
3.	Equipment Volumetric Maximum Capacity Test and Verification of Uniformity of Mixing at Maximum RPM
4.	Equipment Volumetric Minimum Capacity Test and Verification of Uniformity of Mixing at Maximum RPM

13.0 REFERENCES:

- Validation Master Plan.
- Schedule - M – “Good Manufacturing Practices and Requirements of Premises, Plant and Equipment for Pharmaceutical Products.”
- WHO Essential Drugs and Medicines Policy, QA of Pharmaceuticals, Vol-2. Good Manufacturing Practices and Inspection.
- SOP for “Operation & Cleaning of Single Head Automatic powder Filling Machine.”

14.0 DOCUMENTS TO BE ATTACHED:

- Test Report from QC lab
- Any other Relevant Documents.

15.0 NON COMPLIANCE:

- In case of any Non-compliance observed during PQ, same shall be handled through SOP for Handling of Non-Compliance.

16.0 DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY

- In case of any deviation observed during PQ, same shall be handled through SOP for Handling of Deviation.

17.0 CHANGE CONTROL, IF ANY

- If any change is required during PQ, same shall be handled through SOP for Change Management.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION PROTOCOL
FOR
SINGLE HEAD AUTOMATIC POWDER FILLING
MACHINE**

PROTOCOL No.:

18.0 ABBREVIATIONS:

cGMP	:	Current Good Manufacturing Practices
DQ	:	Design Qualification
ID.	:	Identification
IQ	:	Installation Qualification
LTD.	:	Limited
No.	:	Number
OQ	:	Operational Qualification
PPQ	:	Performance Qualification Protocol
PVT	:	Private
S.S	:	Stainless Steel
SOP	:	Standard Operating Procedure