



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING
MACHINE**

EQUIPMENT ID No.	
LOCATION	'Q' BLOCK
DATE OF QUALIFICATION	
SUPERSEDES PROTOCOL No.	



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

REPORT CONTENTS

S.No.	TITLE	PAGE No.
1.0	REPORT PRE APPROVAL	3
2.0	OBJECTIVE	4
3.0	SCOPE	4
4.0	RESPONSIBILITY	5
5.0	EQUIPMENT DETAILS	5
6.0	SYSTEM DESCRIPTION	6
7.0	PRE-QUALIFICATION REQUIREMENT	6
8.0	TESTS & CHECKS	8
9.0	CHECKLIST OF ALL TESTS & CHECKS	14
10.0	DOCUMENTS TO BE ATTACHED	14
11.0	NONCOMPLIANCE	14
12.0	DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY	14
13.0	CHANGE CONTROL, IF ANY	14
14.0	REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY)	15
15.0	CONCLUSION	15
16.0	RECOMMENDATION	15
17.0	ABBREVIATIONS	15
18.0	REPORT POST APPROVAL	16



**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

1.0 REPORT PRE- APPROVAL:

INITIATED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER / EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OPERATING MANAGER (QUALITY ASSURANCE)			
HEAD (PRODUCTION)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

2.0 OBJECTIVE:

- To provide documented evidence that the Equipment is performing as per the parameter defined in Performance Qualification and that it gives result as per the predetermined acceptance criteria.
- To demonstrate that the system will operate reproducibly and consistently within its operating range.
- To confirm the suitability of the Standard Operating Procedures for all routine activities associated with the system.

3.0 SCOPE:

This Report is applicable for performance qualification of Airjet Bottle Cleaning machine installed in Bottle Washing Area.



**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

4.0 RESPONSIBILITY:

The Qualification team, comprising of a representative from each of the following departments, shall be responsible for the overall compliance of this Protocol.

DEPARTMENTS	RESPONSIBILITIES
Quality Assurance	<ul style="list-style-type: none">• Preparation, Review and Approval of the Performance Qualification Protocol.• Co-ordination with Quality Control, Production and Engineering to carryout Performance Qualification Activity.• Provide Training to qualification team.
Production	<ul style="list-style-type: none">• To co-ordinate and support Performance Qualification Activity.
Engineering	<ul style="list-style-type: none">• To provide the required Utility and Engineering support.• Responsible for trouble shooting (if occurred during execution).• Maintenance & preventive maintenance as per schedule.

5.0 EQUIPMENT DETAILS:

Equipment Name	Airjet Bottle Cleaning machine
Equipment ID.	
Manufacturer's Name	Bhavani Engineering
Location of Installation	Bottle washing Machine



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

6.0 SYSTEM DESCRIPTION:

The Automatic Airjet Bottle Air and Vacuum Cleaning Machine is compact unit totally made of SS structure with height adjustment legs, are provided to adjust the machine height and highly efficient machine with elegant look. This multifunctional multi featured machine meets the GMP requirements of washing for glass and plastic Bottles. The machine requires manual loading and automatic unloading of Bottles

PROCESS:

A Machine is inbuilt with Turntable for smooth transfer of round Bottle and suite is provided for noncircular Bottle to the Cleaning Section. This machine works with the principal of back pressure of container. Bottle feed by the Turntable or suit to reach to the star wheel, which transfers the Bottle to the cleaning section one after the. Bottle reaches to the cleaning section, air nozzles starts to flush the air inside towards the Bottle and simultaneously vacuum will suck the particles, disturbed by the air.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

7.0 PRE-QUALIFICATION REQUIREMENTS:

7.1 Training Record of Validation Team:

- All the persons involved in the execution of Qualification Protocol must be trained in all aspects of the qualification activity including the test methodology, acceptance criteria and safety precautions to be followed during working at service floor.

7.2 SYSTEM PRE-REQUISITES:

Verify that the SOP for Operating, Cleaning and Preventive Maintenance of the Airjet Bottle Cleaning machine has been prepared.

S. No.	Document Name	Completed (Yes/No)	Checked By (Engineering) Sign/Date	Verified By (QA) Sign/Date
1.	DQ Protocol approved			
2.	IQ Protocol approved			
3.	OQ Protocol approved			
4.	PQ Protocol approved			
5.	SOP for Operation & Cleaning of Airjet Cleaning Machine			
6.	SOP for Preventive Maintenance Airjet Cleaning Machine			

**Checked By
Production
Sign / Date:** _____

**Verified By
Quality Assurance
Sign / Date:** _____

Inference:.....
.....
.....

**Reviewed By
Manager QA
Sign / Date:** _____



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

8.0 TESTS & CHECKS:

8.1 Machine Speed Optimization:

8.1.1 Trial No.: 01:

Date of Test		Equipment ID	
Total Bottles taken for test		Bottle Size	
Parameter	Low Speed ()	Optimum Speed()	High Speed ()
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Total rejection			

Checked By
Production
Sign/Date: _____

Verified By
Quality Assurance
Sign/Date: _____

Inference:

.....
.....
.....

Reviewed By
Manager QA
Sign / Date: _____



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

Trial No.: 02

Date of Test		Equipment ID	
Total Bottles taken for test		Bottle Size	
Parameter	Low Speed ()	Optimum Speed()	High Speed ()
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Total rejection			

Checked By
Production
Sign/Date: _____

Verified By
Quality Assurance
Sign/Date: _____

Inference:
.....
.....
.....

Reviewed By
Manager QA
Sign/Date: _____



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

8.1.2 Trial No.: 03

Date of Test		Equipment ID	
Total Bottles taken for test		Bottle Size	
Parameter	Low Speed ()	Optimum Speed()	High Speed ()
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Sample after(min)			
Machine jam			
Bottle Breakage			
Rejection			
Total rejection			

Checked By
Production
Sign/Date: _____

Verified By
Quality Assurance
Sign/Date _____

Inference:

.....
.....
.....

Reviewed By
Manager QA
Sign/Date: _____



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

8.2 Air Pressure Test :

8.2.1 Trial no.01:

Date of test		Product name	
Batch No.		Pack Size	

At Air Pressure 1.5 kg / cm² :

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

At Air Pressure 2.0 kg / cm²:

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

At Air Pressure 2.5 kg / cm²:

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

Checked By
Production
Sign / Date: : _____

Verified By
Quality Assurance
Sign / Date: : _____

Inference:.....
.....
.....

Reviewed By
Manager QA
Sign / Date: : _____



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

8.2.2 Trial No.02:

Date of test		Product name	
Batch No.		Pack Size	

At Air Pressure 1.5 kg / cm²:

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

At Air Pressure 2.0 kg / cm²:

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

At Air Pressure 2.5 kg / cm²:

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

Checked By
Production
Sign / Date: _____

Verified By
Quality Assurance
Sign / Date: _____

Inference:.....
.....
.....

Reviewed By
Manager QA
Sign / Date: _____



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

8.2.3 Trial No.03:

Date of test		Product name	
Batch No.		Pack Size	

At Air Pressure 1.5 kg / cm²:

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

At Air Pressure 2.0 kg / cm²:

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

At Air Pressure 2.5 kg / cm²:

Date	Time	Qty. Taken	Observation			Checked By	Verified By
			Foreign Particle	Fabrics material	Paper Piece		

Checked By
Production
Sign / Date: _____

Verified By
Quality Assurance
Sign / Date: _____

Inference:.....
.....
.....

Reviewed By
Manager QA
Sign / Date: _____



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

9.0 CHECKLIST OF ALL TESTS & CHECKS:

The following table lists the number of tests / samples to be carried out & comments on the sample record sheet.

TESTS OR CHECKS	EXECUTED [Y/N]	COMMENT
Machine Speed Optimization		
Air Pressure Test		

10.0 DOCUMENTS TO BE ATTACHED:

- Any Other Relevant Documents

11.0 NON COMPLIANCE:

.....

.....

.....

.....

12.0 DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY:

.....

.....

.....

.....

13.0 CHANGE CONTROL, IF ANY:

.....

.....

.....

.....

.....

14.0 REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY):

.....

.....

.....

.....

.....



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

15.0 CONCLUSION:

.....
.....
.....
.....
.....

16.0 RECOMMENDATION:

.....
.....
.....
.....

17.0 ABBREVIATIONS:

- WHO : World Health Organization
- FDA : Food and Drug Administration
- CFR : Code of Federal Regulations
- GMP : Good Manufacturing Practices
- QA : Quality Assurance
- SOP : Standard Operating Procedure
- mm : Millimetre
- Amp. : Ampere
- DQ : Design Qualification
- IQ : Installation Qualification
- OQ : Operational Qualification
- PQ : Performance Qualification



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
AIRJET BOTTLE CLEANING MACHINE**

PROTOCOL No.:

18.0 REPORT POST APPROVAL:

INITIATED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER / EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OPERATING MANAGER (QUALITY ASSURANCE)			
HEAD (PRODUCTI ON)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			