



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

EQUIPMENT ID No.	
LOCATION	
DATE OF QUALIFICATION	
SUPERSEDES PROTOCOL No.	NIL



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

REPORT CONTENTS

S.No.	TITLE	PAGE No.
1.0	Report Pre Approval	3
2.0	Objective	4
3.0	Scope	4
4.0	Responsibility	5
5.0	Equipment Details	5
6.0	System Description	6
7.0	Reason For Qualification	6
8.0	Site of Study	6
9.0	Frequency of Qualification	6
10.0	Pre-Qualification Requirement	7
11.0	Tests & Checks	8
12.0	Checklist of All Tests & Checks	14
13.0	References	14
14.0	Documents To Be Attached	15
15.0	Noncompliance	15
16.0	Deviation From Pre-Defined Specification, If Any	15
17.0	Change Control, If Any	15
18.0	Review (Inclusive Follow up ,Action)	15
19.0	Conclusion	16
20.0	Recommendation	16
21.0	Abbreviations	17
22.0	Report Post Approval	18



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

1.0 PROTOCOL APPROVAL:

INITIATED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER / EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			
HEAD (QUALITY CONTROL)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

2.0 OBJECTIVE:

To provide documented evidence that the Equipment is performing as per the parameter defined in operational qualification and that it gives result as per the predetermined acceptance criteria.

To demonstrate that the system will operate reproducibly and consistently within its operating range.

To confirm the suitability of the Standard Operating Procedures for all routine activities associated with the system.

The document also provides the observed and obtained values indicating compliance to the PQ Protocol.

3.0 SCOPE:

The Protocol covers all aspects of Performance Qualification for the Premix Filling Machine.

This Protocol will define the methods and documentation used to qualify the Premix Filling Machine for PQ.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

4.0 RESPONSIBILITY:

The Qualification team, comprising of a representative from each of the following departments, shall be responsible for the overall compliance of this Protocol.

DEPARTMENTS	RESPONSIBILITIES
Quality Assurance	<ul style="list-style-type: none">• Preparation, Review, Approval and Compilation of the Performance Qualification Protocol.• Co-ordination with Quality Control, Production and Engineering to carryout Performance Qualification Activity.• Monitoring of Performance Qualification.
Production	<ul style="list-style-type: none">• Review of Protocol.• To co-ordinate and support Performance Qualification Activity.
Engineering	<ul style="list-style-type: none">• Reviewing of qualification Protocol for correctness, completeness and technical excellence.• Responsible for trouble shooting (if occurred during execution).• Maintenance & preventive maintenance as per schedule.
Quality Control	<ul style="list-style-type: none">• Review of Performance Qualification report.• Approval of report post approval.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

5.0 EQUIPMENT DETAILS:

Equipment Name	Premix Filling Machine
Equipment ID.	
Manufacturer's Name	
Supplier's Name	
Location of Installation	Cephalosporin block

6.0 SYSTEM DESCRIPTION:

The Pre-Mix Filling packing machines are nozzle filling machines. They are specially designed for filling particular quantity of the solution, suspension and other liquid phase in the bottle. The machine is ideal for the free flow liquid.

7.0 REASON FOR QUALIFICATION:

After completion of the Operational Qualification of the equipments, it is imperative to perform the Performance Qualification. The study will establish that the parameters are followed, critical variables are under control and the quality of the output is, as desired.

8.0 SITE OF STUDY:

Dry Powder Injection.

9.0 FREQUENCY OF QUALIFICATION:

- Once in two year time period.
- After any major breakdown or after major modification.
- After Change of Location.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

10.0 PRE-QUALIFICATION REQUIREMENTS:

10.1 SYSTEM PRE-REQUISITES:

Verify that the DQ / IQ / OQ of the Premix Filling Machine have been executed and approved.

Verify that the SOP for Operating, Cleaning and Preventive Maintenance of the Premix Filling Machine has been prepared.

S. NO.	DESCRIPTION OF PRE-REQUISITE	COMPLETED (YES / NO)	CHECKED BY ENGINEERING SIGN / DATE	VERIFIED BY QA SIGN / DATE
1.	Verify that the DQ / IQ / OQ of the Premix filling machine has been executed and approved.			
2.	DQ Protocol Document No.:			
3.	IQ Protocol Document No.: OQ Protocol Document No.:			
4.	SOP of "Operation and Cleaning of Premix Filling Machine"			



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

C. PRODUCT NAME:

BATCH NUMBER:

MFG DATE:

DATE:

No. of Bottles	Initial			Middle			Final		
	Full Hopper	Half Hopper	Optimum Speed	Full Hopper	Half Hopper	Optimum Speed	Full Hopper	Half Hopper	Optimum Speed
1.									
2.									
3.									
4.									
5.									
6.									
7.									
8.									
9.									
10.									

**Checked By
Production
Sign / Date: _____**

**Verified By
Quality Assurance
Sign / Date: _____**

Inference:.....
.....
.....

**Reviewed By
Manager QA
Sign / Date: _____**



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

C. PRODUCT NAME:

BATCH NUMBER:

MFG DATE:

DATE:

OBSERVATIONS:

No. of Bottles	Initial			Middle			Final		
	Full Hopper	Half Hopper	Optimum Speed	Full Hopper	Half Hopper	Optimum Speed	Full Hopper	Half Hopper	Optimum Speed
1.									
2.									
3.									
4.									
5.									
6.									
7.									
8.									
9.									
10.									

**Checked By
Production
Sign / Date: _____**

**Verified By
Quality Assurance
Sign / Date: _____**

Inference.....
.....
.....

**Reviewed By
Manager QA
Sign / Date: _____**



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

12.0 CHECKLIST OF ALL TESTS & CHECKS:

The following table lists the number of tests / samples to be carried out & comments on the sample record sheet.

TESTS OR CHECKS	EXECUTED [Y/N]	COMMENT
Average fill Volume and Fill variation of fill bottle		

13.0 REFERENCES:

The Principle Reference is the following:

- Validation Master Plan
- Schedule M: "Good Manufacturing Practices and Requirements of Premises, Plant and Equipment for Pharmaceutical Products."
- WHO Essential Drugs and Medicines Policy, QA of Pharmaceuticals, Vol-2: Good Manufacturing Practices and Inspection.

The following references are used to give addition guidance:

- FDA/ISPE Baseline Pharmaceutical Engineering Guide-Volume 5: Commissioning and Qualification Guide, First Edition / March 2001.
- Code of Federal Regulations (CFR), Title 21, Part 210, Current Good Manufacturing Practice (cGMP) in Manufacturing, Processing, Packing, or Holding of Drugs, Beta. April 1, 1998.
- Code of Federal Regulations (CFR), Title 21, Part 211, Current Good Manufacturing Practice (cGMP) for Finished Pharmaceuticals, April 1, 1998.
- EU Guide to Good Manufacturing Practice, Part 4, 1997.
- European Commission's working party on control of medicines and inspections document, Validation Master Plan, Design Qualification, Installation & Operational Qualification, Non Sterile Process Validation, Cleaning Validation, October 1999.
- GMP Guide, Validation of Automated Systems in Pharmaceutical Manufacture, Version 4.0, December 2001.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

14.0 DOCUMENTS TO BE ATTACHED:

- Any Other Relevant Documents

15.0 NON COMPLIANCE:

.....

.....

.....

16.0 DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY:

.....

.....

.....

17.0 CHANGE CONTROL, IF ANY:

.....

.....

.....

18.0 REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY):

.....

.....

.....

19.0 CONCLUSION:

.....

.....

.....

20.0 RECOMMENDATION:

.....

.....

.....



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

21.0 ABBREVIATIONS:

Sr.	:	Senior
No.	:	Number
gm	:	gram
BSS	:	British Standard Sieve
BMR	:	Batch Manufacturing Record
WHO	:	World Health Organization
FDA	:	Food and Drug Administration
CFR	:	Code of Federal Regulations
GMP	:	Good Manufacturing Practices
SOP	:	Standard Operating Procedure
PRF	:	Premix Filling Machine
°C	:	Degree Centigrade
mm	:	Millimeter
Amp.	:	Amper
DQ	:	Design Qualification
IQ	:	Installation Qualification
OQ	:	Operational Qualification
PQ	:	Performance Qualification



PHARMA DEVILS

**PERFORMANCE QUALIFICATION
REPORT
FOR
PREMIX FILLING MACHINE**

PROTOCOL No.:

22.0 REPORT POST APPROVAL

INITIATED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER / EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			
HEAD (QUALITY CONTROL)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			