



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

**PERFORMANCE QUALIFICATION
REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

EQUIPMENT ID. NO.	
LOCATION	Liquid Filling Area
DATE OF QUALIFICATION	
SUPERSEDES REPORT NO.	NIL



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

REPORT CONTENTS

S.No.	TITLE	PAGE No.
1.0	Report –Pre Approval	3
2.0	Objective	4
3.0	Scope	4
4.0	Responsibility	5
5.0	Equipment Details	6
6.0	Pre-Qualification Requirement	7-8
7.0	Tests & Checks	9-15
8.0	Check list for all test & checks	16
9.0	Documents To Be Attached	17
10.0	Non Compliance	17
11.0	Deviation From Pre–Defined Specification, If Any	17
12.0	Change Control, If Any	17
13.0	Review Inclusive of Follow Up Action, If Any	18
14.0	Conclusion	18
15.0	Recommendations	18
16.0	Abbreviations	19
17.0	Report Post - Approval	20



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

1.0 REPORT PRE – APPROVAL:

PREPARED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER/EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OPERATING MANAGER (QUALITY ASSURANCE)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			



**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

2.0 OBJECTIVE:

- To provide documented evidence that the Equipment is performing consistently, repeatedly and reproducibly within its established operating range and the results of all the test parameters meet the pre-defined acceptance criteria.
- To confirm the suitability of the Standard Operating Procedures for all routine activities associated with the system.

3.0 SCOPE:

- The score of this report is limited for qualification of Vertical Laminar Air Flow installed in the **Liquid Filling Area**.
- This report provides all the relevant information of the performance qualification activity, In-process observations and analytical data of testing of collected samples.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

4.0 RESPONSIBILITY:

The Validation Group, comprising of a representative from each of the following departments, shall be responsible for the overall compliance of this Report:

DEPARTMENTS	RESPONSIBILITIES
Quality Assurance	<ul style="list-style-type: none">• Preparation Review & Authorization and Compilation of the Performance Qualification Report.• Co-ordination with Quality Control, Production and Engineering to carryout Performance Qualification Activity.• Monitoring of Performance Qualification.• Post Approval of Performance qualification report After Execution.
Production	<ul style="list-style-type: none">• Review of Performance Qualification Report.• To co-ordinate and support Performance Qualification Activity.• Post Approval of Performance qualification report After Execution.
Quality Control	<ul style="list-style-type: none">• Analytical Support (Microbiological Testing/Analysis)
Engineering	<ul style="list-style-type: none">• Review of Performance Qualification Report for correctness, completeness and technical excellence.• Responsible for trouble shooting (if occurred during execution).• Maintenance & preventive maintenance as per schedule.• Post Approval of Performance qualification report After Execution.
External Qualification Agency (if Applicable)	<ul style="list-style-type: none">• Performance of qualification activity as per Protocol.



**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

PHARMA DEVILS

5.0 EQUIPMENT DETAILS:

Equipment Name	Vertical Laminar Air Flow
Equipment	
Manufacturer's Name	
Model	
S.No.	
Supplier's Name	
Location of Installation	Liquid Filling Area

6.0 PRE – QUALIFICATION REQUIREMENTS:

Verification for availability, completeness and approval status of all the required relevant documents shall be done and observations shall be recorded in the performance qualification report.

6.1 Verification of Documents:

Record the observations for documents in the below mentioned table.

S.No.	Document Name	Document/SOP No.	Completed (Yes/No)	Checked By (Engineering) Sign/Date	Verified By (QA) Sign/Date
1.	Executed and approved Design Qualification document				
2.	Executed and approved Installation Qualification document				
3.	Executed and approved Operational Qualification document				
4.	SOP for operation & Cleaning of Laminar Air Flow				
5.	SOP for Preventive Maintenance of Laminar Air Flow				

Inference:

.....

**Reviewed By
(Manager QA)
(Sign & Date):**



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

6.2 Training Record of Validation Team:

- All the persons involved in the execution of Qualification Protocol must be trained in all aspects of the qualification activity including the test methodology, acceptance criteria and safety precautions to be followed during working at service floor.

6.3 Calibration of Test Instruments:

- Calibration of all the instruments used for qualification should be mentioned along with Calibration Certificates.

S. No.	Name of Test Instrument	Date of Last Calibration	Next Due on	Status	Availability of Calibration Certificate	Verified By (QA) Sign/Date
1.						
2.						
3.						
4.						
5.						
6.						
7.						
8.						

Inference:

.....
.....
.....

**Reviewed By
(Manager QA)
(Sign & Date):.....**



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

7.0 TESTS AND CHECKS:

7.1 AIR VELOCITY MEASUREMENT:

Instrument Name	
Make	
Model / Type	
Calibration Date	
Calibration Due Date	
Calibration Certificate Attached	

OBSERVATION AND RESULTS:

Date	Area	Filter ID	Air Velocity (In Ft. /Min.)					Acceptance Criteria	Average Air Velocity (Ft/Min)
			LOCATION						
			1	2	3	4	5		
								The Average measured clean air velocity should be 90±20 % ft/min at 6 inches downstream from the filter face	

Checked By:
(Engineering)
Sign & Date:.....

Verified By:
Quality Assurance)
Sign & Date:.....

Inference:

.....
.....
.....

Reviewed By
Manager QA)
Sign & Date) :.....



**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

7.2 HEPA FILTER INTEGRITY TEST (PAO TEST) REPORT

Instrument Name	
Make	
Calibration Date	
Calibration Due Date	
Calibration Certificate Attached	
Up Steam Concentration	

TEST RESULTS:

Date	Area Name	HEPA Id. / S. No.	Acceptance Criteria	Observation (% of Leakage)
			The PAO penetration / leak through HEPA filters should not be greater than 0.01% of the upstream PAO concentration.	

Checked By:
(Engineering)
Sign & Date:.....

Verified By:
(Quality Assurance)
Sign & Date:.....

Inference:

.....
.....
.....

Reviewed By:
(Manager QA)
(Sign & Date):.....



**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

7.3 DIFFERENTIAL PRESSURE RECORD

7.3.1 Across HEPA Filter:

Pressure Differential	Across HEPA - Filter
Magnehelic Gauge ID No.	
Date of Calibration	
Calibration due date	
Acceptance Criteria	

Date	Name of Equipment	Observation											
		Time 00 - 04 Hr.		Time 04 - 08 Hr.		Time 08 - 12 Hr.		Time 12 - 16 Hr.		Time 16 - 20 Hr.		Time 20 - 00 Hr.	
		Time	mm of WC	Time	mm of WC	Time	mm of WC	Time	mm of WC	Time	mm of WC	Time	mm of WC

Checked By:
(Engineering)
Sign & Date:.....

Verified By:
(Quality Assurance)
Sign & Date:.....

Inference:
.....
.....
.....

Reviewed By:
(Manager QA)
(Sign & Date):.....



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

7.3.2 Across Pre Filter:

Pressure Differential	Across Pre- Filter
Magnehelic Gauge ID No.	
Date of Calibration	
Calibration due date	
Acceptance Criteria	

Date	Name of Equipment	Observation											
		Time 00 - 04 Hr.		Time 04 - 08 Hr.		Time 08 - 12 Hr.		Time 12 - 16 Hr.		Time 16 - 20 Hr.		Time 20 - 00 Hr.	
		Time	pa	Time	pa	Time	pa	Time	pa	Time	pa	Time	pa

(Checked By:
(Engineering)
Sign & Date:.....

Verified By:
(Quality Assurance)
Sign & Date:.....

Inference:

.....
.....
.....

Reviewed By:
(Manager QA)
(Sign & Date) :.....



**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

7.4 NON – VIABLE PARTICLE COUNT

Name of equipment :
Particle Counter Id. :
Date of Calibration :
Due on Calibration :
Make :

Date	Area /Location	Observation			
		In operation		At Rest	
		≥0.5μ	≥5.0μ	≥0.5μ	≥5.0μ


Checked By:
(Engineering)
Sign & Date:.....

Verified By:
(Quality Assurance)
Sign & Date:.....

Inference:

.....

Reviewed By:
(Manager QA)
(Sign & Date):.....

 PHARMA DEVILS	PERFORMANCE QUALIFICATION REPORT FOR VERTICAL LAMINAR AIR FLOW	PROTOCOL No.:
--	---	----------------------

7.7 AIR FLOW PATTERN:

Date of Testing		Make/Model	
Instrument Name		Calibration Date	
Instrument ID.		Calibration Due Date	

Area	Air Flow Pattern Should Be Moving In Downward Direction	The Air Flow Pattern Shall Be From Supply Air to Return Filter	Visibility of Smoke Generated (Yes/No)

Checked By:
(Engineering)
Sign & Date:.....

Verified By:
(Quality Assurance)
Sign & Date:.....

Inference:

.....

Reviewed By:
(Manager QA)
(Sign & Date):.....

8.0 CHECK LIST OF ALL TESTS & CHECKS

S.No.	Name of Test or Check	Execution (Yes/No.)	Remark	Verified By (Sign & Date)
1.	Air Velocity Measurement			
2.	HEPA Filter Integrity Test (PAO Test) Report			
3.	Differential Pressure Record			
4.	Non – Viable Particle Count			
6.	Environmental Monitoring - (Settle Plate Method)			
7.	Environmental monitoring (Air Sampling Method)			
8.	Air Flow Pattern Test			

Inference:

.....

.....

.....

.....

Reviewed By:
(Manager QA)
(Sign & Date) _____



**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

9.0 DOCUMENTS TO BE ATTACHED:

- Report of QC (Micro) Analysis
- Calibration Certificate of Test Instrument.
- Any Other Relevant Document

10.0 NON COMPLIANCE:

.....

.....

.....

.....

.....

.....

11.0 DEVIATION FROM PREDEFINED SPECIFICATION IF, ANY:

.....

.....

.....

.....

.....

.....

.....

.....

12.0 CHANGE CONTROL, IF ANY:

.....

.....

.....

.....

.....

.....

.....

.....



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

16.0 ABBREVIATIONS:

CFM	:	Cubic flow Metter
DQ	:	Design Qualification
IQ	:	Installation Qualification
LTD.	:	Limited
mm	:	Millimetre
No.	:	Number
OQ	:	Operational Qualification
PAO	:	Poly alpha olefin
PVT	:	Private
SOP	:	Standard Operating Procedure



**PERFORMANCE QUALIFICATION REPORT
FOR
VERTICAL LAMINAR AIR FLOW**

PROTOCOL No.:

17.0 REPORT POST-APPROVAL:

PREPARED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER/EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OPERATING MANAGER (QUALITY ASSURANCE)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			