



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

**PERFORMANCE QUALIFICATION
REPORT
FOR
VIAL OPTICAL INSPECTION
MACHINE**

EQUIPMENT ID. No.	
LOCATION	Packing Hall
DATE OF QUALIFICATION	
SUPERSEDES PROTOCOL No.	NIL



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

REPORT CONTENTS

S.No.	TITLE	PAGE No.
1.0	Pre-Approval	3
2.0	Objective	4
3.0	Scope	4
4.0	Responsibility	5
5.0	Equipment Details	6
6.0	Pre-Qualification Requirement	6
7.0	Tests & Checks	7-11
8.0	Checklist Of All Tests And Checks	12
9.0	Documents To Be Attached	13
10.0	Non Compliance	13
11.0	Deviation From Pre-Defined Specification, If Any	13
12.0	Change Control, If Any	13
13.0	Review Inclusive of Follow Up Action, If Any	14
14.0	Conclusion	14
15.0	Recommendations	14
16.0	Abbreviations	15
17.0	Post Approval	16



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

1.0 PRE – APPROVAL:

INITIATED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER/EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

2.0 OBJECTIVE:

- To provide documented evidence that the Equipment is performing consistently, repeatedly and reproducibly within its established operating range and the results of all the test parameters meet the pre-defined acceptance criteria.
- To confirm the suitability of the Standard Operating Procedures for all routine activities associated with the system.

3.0 SCOPE:

- The scope of this report is limited for qualification of Vial Optical Inspection Machine installed in Packing Hall.
- This report provides all the relevant information of the performance qualification activity, In-process observations.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

4.0 RESPONSIBILITY:

The Validation Group, comprising of a representative from each of the following departments, shall be responsible for the overall compliance of this Report:

DEPARTMENTS	RESPONSIBILITIES
Quality Assurance	<ul style="list-style-type: none">• Preparation, Authorization, Approval and Compilation of the Performance Qualification Report.• Co-ordination with Quality Control, Production and Engineering to carryout Performance Qualification Activity.• Monitoring of Performance Qualification.
Production	<ul style="list-style-type: none">• Review of Report.• To co-ordinate and support Performance Qualification Activity.
Engineering	<ul style="list-style-type: none">• Reviewing of qualification protocol for correctness, completeness and technical excellence.• Responsible for trouble shooting (if occurred during execution).• Maintenance & preventive maintenance as per schedule.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

5.0 EQUIPMENT DETAILS:

Equipment Name	Vial Optical Inspection Machine
Equipment
Manufacturer's Name	Ambica Pharma Machines Private Limited
Model	AVIN - 240
Supplier's Name	Ambica Pharma Machines Private Limited
Location of Installation	Packing Hall

6.0 PRE – QUALIFICATION REQUIREMENTS:

Verification for availability, completeness and approval status of all the required relevant documents shall be done and observations shall be recorded in the performance qualification report.



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

7.0 TESTS AND CHECKS:

7.1 VERIFICATION OF MACHINE SPEED:

Date of Test	
Vial Size	

Parameter	Trial-I	Trial-II	Trial-III
	No. of Vial in One Minutes	No. of Vial in One Minutes	No. of Vial in One Minutes
At 10 (Speed knob Setting)			
At 20 (Speed knob Setting)			
At 30 (Speed knob Setting)			
At 40 (Speed knob Setting)			
At 50 (Speed knob Setting)			
At 60 (Speed knob Setting)			
At 70 (Speed knob Setting)			
At 80 (Speed knob Setting)			
At 90 (Speed knob Setting)			
At 100 (Speed knob Setting)			

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:

Inference:

.....
.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

7.2 LUX LEVEL :

Date of Test

Parameter	Table No.-I	Table No.-II	Table No.-III	Table No.-IV
Lux Level (NLT 2200)				

**Checked By
(Production)**

Sign/Date:

**Verified By
(Quality Assurance)**

Sign/Date:

Inference:

.....
.....
.....

**Reviewed By
(Manager QA)**

Sign/Date:



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

7.3 CHALLENGE TEST:

Date of Test	
Vial Size	

Parameter		Observation	
		Visual Inspector-I	Visual Inspector-II
At 40 (Speed knob Setting)	Trial-I Total 400 Vials Good -350 Vials Rejected-50 Vials (Sealing-15, Cracked-21, Empty-14)		
	Trial-II Total 400 Vials Good -342 Vials Rejected-58 Vials (Sealing-17, Cracked-22, Empty-15)		
	Trial-III Total 400 Vials Good -346 Vials Rejected-54 Vials (Sealing-11, Cracked-26, Empty-17)		
At 50 (Speed knob Setting)	Trial-I Total 500 Vials Good -443 Vials Rejected-57 Vials (Sealing-21, Cracked-17, Empty-19)		
	Trial-II Total 500 Vials Good -449 Vials Rejected-51 Vials (Sealing-20, Cracked-14, Empty-17)		
	Trial-III Total 500 Vials Good -427 Vials Rejected-73 Vials (Sealing-24, Cracked-29, Empty-20)		
At 60 (Speed knob Setting)	Trial-I Total 600 Vials Good -518 Vials Rejected-82 Vials (Sealing-31, Cracked-26, Empty-25)		
	Trial-II Total 600 Vials		



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

Parameter		Observation	
		Visual Inspector-I	Visual Inspector-II
At 70 (Speed knob Setting)	Good -520 Vials Rejected-80 Vials (Sealing-30, Cracked-22, Empty-28)		
	Trial-III Total 600 Vials Good -517 Vials Rejected-83 Vials (Sealing-33, Cracked-28, Empty-22)		
	Trial-I Total 700 Vials Good -609 Vials Rejected-91 Vials (Sealing-37, Cracked-13, Empty-41)		
	Trial-II Total 700 Vials Good -602 Vials Rejected-98 Vials (Sealing-36, Cracked-14, Empty-48)		
	Trial-III Total 700 Vials Good -611 Vials Rejected-89 Vials (Sealing-30, Cracked-41, Empty-18)		

**Checked By
(Production)**

Sign/Date:

Verified By

(Quality Assurance)

Sign/Date:

Inference:

.....
.....
.....

Reviewed By

(Manager QA)

Sign/Date:



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

8.0 CHECKLIST OF ALL TESTS & CHECKS:

This checklist is provided to ensure that all tests or checks required for this protocol have been executed.

Tests or Checks	Executed (Yes/No)	Remarks
Verification of Machine Speed.		
Lux Level		
Challenge Test		

**Checked By
(Production)
Sign/Date:**

**Verified By
(Quality Assurance)
Sign/Date:**

Inference:
.....
.....
.....

**Reviewed By
(Manager QA)
Sign/Date:**



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

9.0 DOCUMENTS TO BE ATTACHED:

- Any Other Relevant Documents.

10.0 NON COMPLIANCE:

.....
.....
.....
.....
.....
.....

11.0 DEVIATION FROM PREDEFINED SPECIFICATION IF, ANY:

.....
.....
.....
.....
.....
.....
.....
.....

12.0 CHANGE CONTROL, IF ANY:

.....
.....
.....
.....
.....



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

13.0 REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY):

.....
.....
.....
.....
.....
.....

14.0 CONCLUSION:

.....
.....
.....
.....
.....
.....
.....
.....

15.0 RECOMMENDATION:

.....
.....
.....
.....
.....
.....
.....
.....



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

16.0 ABBREVIATIONS:

WHO	:	World Health Organization
FDA	:	Food and Drug Administration
CFR	:	Code of Federal Regulations
cGMP	:	Current Good Manufacturing Practices
EU	:	European Union
CQA	:	Corporate Quality Assurance
QC	:	Quality Control
DQ	:	Design Qualification
IQ	:	Installation Qualification
OQ	:	Operational Qualification
PQ	:	Performance Qualification
SOP	:	Standard Operating Procedure
NLT	:	Not Less Than
KW	:	Kilo watt
SS	:	Stainless Steel
ID.	:	Identification
mm	:	Mili meter
MCB	:	Miniature Circuit Breaker
ID	:	Inner Diameter



PHARMA DEVILS

**PERFORMANCE QUALIFICATION REPORT
FOR
VIAL OPTICAL INSPECTION MACHINE**

PROTOCOL No.:

17.0 POST APPROVAL:

INITIATED BY:

DESIGNATION	NAME	SIGNATURE	DATE
OFFICER/EXECUTIVE (QUALITY ASSURANCE)			

REVIEWED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (PRODUCTION)			
HEAD (ENGINEERING)			

APPROVED BY:

DESIGNATION	NAME	SIGNATURE	DATE
HEAD (QUALITY ASSURANCE)			