



# PHARMA DEVILS

PRODUCTION DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Title:** Cleaning and Sterilization of Filters

<b>SOP No.:</b>		<b>Department:</b>	Production	
		<b>Effective Date:</b>		
<b>Revision No.:</b>	00	<b>Revision Date:</b>		
<b>Supersede Revision No.:</b>	Nil	<b>Page No.:</b>	1 of 5	

### 1.0 OBJECTIVE:

To lay down a procedure for cleaning and sterilization of filters.

### 2.0 SCOPE:

This SOP is applicable for cleaning and sterilization of filters used in Production Department.

### 3.0 RESPONSIBILITY:

Officer / Executive – Production

### 4.0 ACCOUNTABILITY:

Head – Production

### 5.0 ABBREVIATIONS:

%	Percent
IPA	Isopropyl Alcohol
Ltd.	Limited
No.	Number
Pvt.	Private
QA	Quality Assurance
SOP	Standard Operating Procedure
WFI	Water for Injection

### 6.0 PROCEDURE:

#### 6.1 PRECAUTIONS:

- 6.1.1 Handle the Filter carefully to avoid any possible damage / extraneous contamination.
- 6.1.2 Never remove the filters from line when it is in under operation.
- 6.1.3 Filter should be clean with their respective cleaning agent.
- 6.1.4 When observed any physical damage, filter should be replaced by newer one and discard the damaged filter.

#### 6.2 CLEANING AND STERILIZATION OF CARTRIDGE FILTERS:

- 6.2.1 Dismantle the filter from their housing and transfer it in the unit preparation area.
- 6.2.2 Clean the filter housing with WFI.
- 6.2.3 Take out the Filter & Dip into 70% IPA solution for 20 minutes.
- 6.2.4 Clean the filter with WFI by flushing of WFI in case of hydrophilic Cartridge only.



**STANDARD OPERATING PROCEDURE**

**Title:** Cleaning and Sterilization of Filters

<b>SOP No.:</b>		<b>Department:</b>	Production	
		<b>Effective Date:</b>		
<b>Revision No.:</b>	00	<b>Revision Date:</b>		
<b>Supersede Revision No.:</b>	Nil	<b>Page No.:</b>	2 of 5	

**6.2.5** Dry the filter by flushing of Compressed Air.

**6.2.6** Re-assembled the cleaned filter in their respective housing, if Cartridge filter used in Aseptic Area or Direct Contact in Product, filter should be sterilized before installation.

**6.2.7** Perform the Sterilization process of the Cartridge filters as per respective SOP.

**6.3 CLEANING AND STERILIZATION OF AUTOCLAVE AIR VENT FILTER:**

**6.3.1** Remove the filter from autoclave chamber vent port and transfer to Unit Preparation area for cleaning and attached another sterilized filter on the same position.

**6.3.2** For cleaning of vent filter, dip the filter in 70% IPA for 20 minute.

**6.3.3** After that dry it by flushing of compressed Air.

**6.3.4** Check physical condition of the filter visually, if it is physically damaged then discard it.

**6.3.5** Perform the Sterilization process of the vent filters as per respective SOP.

**6.4 CLEANING AND STERILIZATION OF CAPSULE FILTERS:**

**6.4.1** Remove the filter from their line and transfer it in Filter Cleaning Area.

**6.4.2** For cleaning of Capsule filter, dip the filter in 70% IPA for 20 minute.

**6.4.3** In case of hydrophilic Capsule clean the filter by flushing of WFI.

**6.4.4** After that dry it by flushing of compressed Air.

**6.4.5** Check physical condition of the filter visually, if it is physically damaged then discard it.

**6.4.6** Perform the Sterilization process of the vent filters as per respective SOP.

**6.5** Record the filter cleaning, sterilization and sterilization cycle details of other than product filter in **Annexure-I, "Filter Cleaning and Sterilization Record"**.

**Note:** Cartridge filter and Capsule Filter should be replaced by newer one when it is failure in BPT/WIT/ FFT or after completion of recommended Sterilization Cycle.

**6.6 FREQUENCY:**

**6.6.1 Hydrophilic Filter** – After Operation/Wherever required.



# PHARMA DEVILS

PRODUCTION DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Title:** Cleaning and Sterilization of Filters

<b>SOP No.:</b>		<b>Department:</b>	Production	
		<b>Effective Date:</b>		
<b>Revision No.:</b>	00	<b>Revision Date:</b>		
<b>Supersede Revision No.:</b>	Nil	<b>Page No.:</b>	3 of 5	

**6.6.2 Hydrophobic Filter** –15± 4 days/Wherever required.

### 7.0 ANNEXURES:

ANNEXURE No.	TITLE OF ANNEXURE	FORMAT No.
Annexure-I	Filter Cleaning and Sterilization Record	
Annexure-II	Integrity Test Value of the Filters	

**ENCLOSERS:** SOP Training Record.

### 8.0 DISTRIBUTION:

- Controlled Copy No. 01      Quality Assurance
- Controlled Copy No. 02      Production
- Master Copy                      Quality Assurance

### 9.0 REFERENCES:

Not Applicable

### 10.0 REVISION HISTORY:

#### CHANGE HISTORY LOG

Revision No.	Change Control No.	Details of Changes	Reason for Change	Effective Date	Updated By





# PHARMA DEVILS

PRODUCTION DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Title:** Cleaning and Sterilization of Filters

<b>SOP No.:</b>		<b>Department:</b>	Production
		<b>Effective Date:</b>	
<b>Revision No.:</b>	00	<b>Revision Date:</b>	
<b>Supersede Revision No.:</b>	Nil	<b>Page No.:</b>	5 of 5

### ANNEXURE – II INTEGRITY TEST VALUE OF THE FILTERS

Category	Filter size	Pore size	Manufacturer	Catalogue No.	Test to be performed		MOC	Max. sterilization cycle
					Test pressure	Test value		
Hydrophilic	10 inch	0.45 $\mu$ + 0.2 $\mu$	Sartorius		3200 mbar	Bubble point should not found below test value	PES	25
	30 inch	0.45 $\mu$ + 0.2 $\mu$	Sartorius		3200 mbar		PES	25
	10 inch	0.2 $\mu$	Pall		2760 mbar		PVDF	25
	10 inch	0.2 $\mu$	Millipore		3100 mbar		PVDF	25
	30 inch	0.2	Millipore		3100 mbar		PVDF	25
	10 inch	0.2 $\mu$	MDI		3440 mbar		PES	15
Hydrophobic	5 inch	0.2 $\mu$	MDI		1520 mbar	Bubble point should not found below test value	PTFE	80
	10 inch	0.2 $\mu$	Pall		1100 mbar		PTFE	80
	5 inch	0.2 $\mu$	Pall		1380 mbar		PTFE	80
	10 inch	0.2 $\mu$	Pall		1380 mbar		PTFE	80
	5 inch	0.2 $\mu$	Millipore		1100 mbar		PTFE	80
	10 inch	0.2 $\mu$	Millipore		1100 mbar		PTFE	80
Hydrophobic	5 inch	0.2 $\mu$	Sartorius		700 mbar		PTFE	80