



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Filtration of Bulk Solution

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	1 of 3

1.0 OBJECTIVE:

To lay down a Procedure for Filtration of Bulk Solution.

2.0 SCOPE:

This SOP is applicable for Filtration of Bulk Solution in Production Area.

3.0 RESPONSIBILITY:

Officer / Executive Production

4.0 ACCOUNTABILITY:

Head Production

5.0 ABBREVIATIONS:

Ltd.	Limited
SOP	Standard Operating Procedure
No.	Number
QA	Quality Assurance
GM	General Manager
BMR	Batch Manufacturing Record
PVT.	Private

6.0 PROCEDURE:

- 6.1 Clean the Cartridge Filter Housing containing 1.2 μ , 0.6 μ and 0.22 μ Cartridge Filters respectively.
- 6.2 Assemble the Cartridge Filter and housing and check for proper fixing of "O – Ring" to ensure complete tightness of the Filtration assembly.
- 6.3 Perform Bubble Point Test of Cartridge Filter to ensure its integrity before Sterilization.
- 6.4 Perform Sterilization of assembled Cartridge filter housing containing Cartridge Filter of 1.2 μ , 0.6 μ and 0.22 μ .
- 6.5 Release the air from filter prior to start up of the filtration operation.
- 6.6 Start the Batch Filtration after release of the Batch from the QC and as per instructions from Manufacturing Chemist.
- 6.7 After completion of filtration, perform the Bubble Point Test of the Cartridge Filter to confirm the Post Filtration Integrity.
- 6.8 Record the Filtration related details at specified place in respective BMR.



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Filtration of Bulk Solution

SOP No.:		Department:	Production	
		Effective Date:		
Revision No.:	00	Revision Date:		
Supersede Revision No.:	Nil	Page No.:	2 of 3	

7.0 ANNEXURES:

ANNEXURE No.	TITLE OF ANNEXURE	FORMAT No.
Annexure-I	Bulk Solution Filtration Record	

ENCLOSURES: SOP Training Record

8.0 DISTRIBUTION:

- Controlled Copy No.1 Quality Assurance
- Controlled Copy No.2 Production
- Master Copy Quality Assurance

9.0 REFERENCES:

Not Applicable.

10.0 REVISION HISTORY:

CHANGE HISTORY LOG

Revision No.	Change Control No.	Details of Changes	Reason for Change	Effective Date	Updated By

