



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	1 of 12

1.0 OBJECTIVE:

To lay down a procedure for Operation and Cleaning of Hi-Cart Cartoning Machine.

2.0 SCOPE:

This SOP is applicable for Operation and Cleaning of Hi-Cart Cartoning Machine (Make **Pampac**) used.

3.0 RESPONSIBILITY:

Operating Person-Production

4.0 ACCOUNTABILITY:

Head Production

5.0 ABBREVIATIONS:

BPR Batch Packing Record
HMI Human Machine Interface
IPA Iso Propyl Alcohol
No. Number
QA Quality Assurance
SOP Standard Operating Procedure

6.0 PROCEDURE:

6.1 Precaution:

- 6.1.1 Line Clearance shall be taken by Production (Packing) Officer / Executive & given by QA Officer / Executive prior to startup of the Hi – Cart Machine.
- 6.1.2 Do not rush while running the Machine to avoid any accidental injury or Packing of substandard Vial / Cartons / Leaflet.
- 6.1.3 Impression of stereo should be taken on respective format and check all the text matter carefully.
- 6.1.4 For verification of MRP of particular product refer the MRP mentioned in the SAP at location 1100 in case of Loan License Product (LL) OR 9700 in case of Principle to Principle (P2P) as well as the MRP mentioned in the Process Order & Conversion note.
- 6.1.5 A special attention needed to check the letter looking similar like O & 0, I & 1, P & R, O & Q, I & L, C & G, E & F. For reference refer the **Annexure No.-I** titled as “**Identification of Similar Looking Letters**”.

6.1.6 Verification of Functional Unit:

Operate the machine as per SOP and Check for the following functions of the machine.



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	2 of 12

S. No.	Station	Specification	Verification Procedure
1.	Product Feeding System	To Transfer the Bottle Into the Product Pocket at Rated Speed	Check the Bottle Dropped without Damage
2.	Prefolded Leaflet transfer system	The Leaflet as per Requirements Transfer Through Turret to Leaflet Clamps Provided on leaflet Chain.	Check the Leaflet Smooth Transfer of Leaflet Without Damage.
3.	Carton Loading System	To Stack /Load the Carton as per Change Part into Carton Magazine & Carton to be Set With the Help of Adjustable Guides	Check Smooth Movement of Stacked Cartons in Carton Magazine for Easy Pickup
4.	Carton Pick Up Transfer System	To Pick Up the Carton From Carton Magazine into Carrier Finger	Carton are Picked and Placed In The Carton Chain Pockets Properly without Damages
5.	Product Pushing System	To Insert the Product into Carton along with Leaflets	Ensure that Smooth Insertion of Product & Leaflets into Carton Without Damage
6.	Carton Closing System	To Close (Lock) the Cartons with the help of Tuck in Pusher & all Guide to be set with the help of Handle and Digital Ganter	Check the Carton are Getting Closed (locked) properly without damage.
7.	Carton rejection System	To Reject the Carton which are Having No Product.	Check the carton and gets rejected as per specification mentioned along side

6.1.7 Operating Access Level: The operating access level and their access authorization is as below:



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	3 of 12

Level	Level Name	Access Authorization
Level -1	Operator	Home Screen, Function keys, Optional Keys, Alarm/Events, Monitor, Shift array, PLC Input, PLC Output, ECR Monitor.
Level -2	Supervisor	Home Screen ,Function keys , Optional Keys, Machine Settings, cam, Pockets setting, Counter, timers, Blister Drop, Alarm/Event, Monitor ,Shift array. PLC Input Output, ECR Monitor,
Level-3	Administrator	Set all above, System Setting +Project Setting, Print Information

6.1.8 Detail of Alarm and Interlocking in the Machine:

S.No.	Description of test	Message displayed on HMI	Result of test
1.	Emergency switch pressed on HMI	Emergency stop at console	Machine should be stopped immediately.
2.	Push Emergency switch pressed on discharge end	Emergency stop at discharge end	Machine should be stopped immediately.
3.	Emergency switch at pressed at infeed	Emergency stop at machine infeed	Machine should be stopped immediately
4.	Air pressure is less or below the set value	Low air pressure	Machine should be stopped immediately
5.	Carton vacuum key off from HMI	Carton vacuum off	Machine should be stopped immediately
6.	Leaflet vacuum key off from HMI	Leaflet vacuum off	Machine should be stopped immediately

6.2 Operation:

- 6.2.1 Switch "ON" the Control Switch of the Machine.
- 6.2.2 Open the Compressed Air Valve & adjust the Compressed Air Pressure to 6.0 to 6.5 Kg/ cm².
- 6.2.3 Select login on HMI and enter user and password.
- 6.2.4 Turn the knob on control panel to SET mode to set the machine.
- 6.2.5 Feed the required Unit Cartons in the Carton Magazine.



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	4 of 12

- 6.2.6 Feed the required Leaflet in the Leaflet Magazine as per applicability.
- 6.2.7 Fix the Rubber Stereo in Printing Block with the help of both side adhesive tape.
- 6.2.8 Feed screw capped & checked labeled Vials in the M/c hopper.
- 6.2.9 Select the function key option in HMI.
- 6.2.10 Select the Machine Home, Turn Table, Conveyor, Carton Vacuum, Leaflet Vacuum and Inch in the Display Panel.
- 6.2.11 Set the machine as per required pack size and carton coding setting.
- 6.2.12 Turn the knob of the Set Mode to the Run Mode and press the "Machine On".
- 6.2.13 Check the Batch Coding details on Unit Carton & get it approved from Production Officer & QA officer before running the Machine.
- 6.2.14 After cartoning the vials, unit cartons comes out on conveyor belt where Air Jet has been provided to displace empty Carton if it comes.
- 6.2.15 Check the individual packed carton for Batch Coding details and any damage.
- 6.2.16 Reject and segregate any damaged carton and carton having improper batch coding detail. For reference refer the **Annexure No.-II** titled as "**Reference examples of rejected cartons due to improper coding**".
- 6.2.17 Segregate the rejected carton and keep into lock and key rejection bin.
- 6.2.18 After completion of packing destroy the rejected carton and also consider during reconciliation of the same.
- 6.2.19 Keep the good Packed Cartons in Shovel for shrink packing.
- 6.2.20 Remove the vials and insert from rejected carton and pack into the good carton.



STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	5 of 12



Figure- 1: Hi-Cart Machine

6.3 Sensor Detail and Its Functional Role: The sensor detail and its functional role is mentioned as below:

S.No.	Sensor	Functional role of sensor
1.	Continuous product sensor	This sensor senses the availability of product continuously in the infeed conveyor before star wheel. If gap between vials is more or any vial comes horizontally in the in feed conveyor the machine will be stopped.
2.	Star wheel loading sensor	This sensor senses the overloading of the loading Star Wheel and in case of over loading the machine will be stopped.
3.	Vertical product check	This sensor senses the direction of vial vertical or horizontal in the pocket. If any vial comes in vertical direction the machine will be stopped.
4.	Product check	This sensor senses the availability of product in the product pocket. If product is not available in the pocket then leaflet will not picked up for the particular pocket.
5.	Leaflet low level check	This sensor senses the level of leaflet in the leaflet magazine. If the level of leaflet comes upto low level sensor lever then machine will be stopped.
6.	Leaflet check	This sensor senses the availability of leaflet for particular product. If leaflet is not available then carton will not picked up for the particular product.
7.	Carton low level check	This sensor senses the level of carton in the carton magazine. If



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	6 of 12

S.No.	Sensor	Functional role of sensor
		the level of carton comes upto low level sensor lever then machine will be stopped.
8.	Carton check	This sensor senses the availability of carton for particular product. If carton is not available for any particular product then product will get diverted without insertion. Means the in such case the pusher for this particular pocket will not move & product will remain in the pocket and product will be dropped into the rejection box at discharge end.
9.	Half product insertion check	If product is not inserted completely into carton then machine will be stopped.
10.	Hand wheel sensor	If during machine operation hand wheel comes out from its position then machine will be stopped to avoid any damage.
11.	Arm reference sensor	This sensor gives command to the arm for rotation.
12.	Carton flap home sensor	This sensor gives command to close the flap of carton.
13.	Carton clutch sensor	This sensor senses the overloading of the carton chain clutch and in case of overloading the machine will be stopped.
14.	Empty carton check	This sensor senses the availability of product in the folded carton. Any carton without product will be rejected and dropped into rejection box at discharge end.
15.	Pusher loading sensor	This sensor senses the overloading on the product pusher and in case of over loading the machine will be stopped.
16.	Machine home sensor	This sensor denotes the homing status of machine.

6.4 Challenge test:

6.4.1 Perform the challenge test of below sensor during operation of Hi-Cart machine. Action and observation is also mentioned in the below table.

S.No.	Sensor	Action during challenge test	Observation
1.	Vertical product check	Keep one product vertically before sensing	Machine should be stopped immediately
2.	Product check	Manually take out one product from its pocket before sensing	Leaflet should not picked up for this particular pocket
3.	Leaflet low level check	Remove leaflet from magazine up to the low level sensor lever	Machine should be stopped immediately
4.	Leaflet check	Remove one leaflet manually	Carton should not picked up for



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	7 of 12

S.No.	Sensor	Action during challenge test	Observation
		before sensing	the particular pocket
5.	Carton low level check	Remove carton from magazine up to the low level sensor lever	Machine should be stopped immediately
6.	Carton check	Remove carton manually from sucker arm before sensing	Product should get diverted without insertion
7.	Empty carton check	Remove one product before insertion into carton	Carton should get rejected at discharge end

6.4.2 Frequency of Challenge Test: Perform the challenge test initially before start the machine operation and towards the end of operation.

6.4.3 Record the challenge test and its observation as per **Annexure No-III** titled as **“Challenge test record during operation of Hi-Cart machine”**.

6.4.4 Give the justification if observation is found **Not OK**.

6.4.5 Attach the filled annexure in the respective BPR.

6.4.6 Record the operation details in respective **“Equipment’s /Instruments Usage Cleaning Log”** as per **SOP**.”.

6.5 Cleaning:

6.5.1 After completion of Packing Operation remove the excess Packing Material from the Packing line and Machine.

6.5.2 Collect the Rejected Cartons in separate crate or Poly Bag.

6.5.3 Count the Rejected Cartons and record in the BPR.

6.5.4 Put the Status Label **“TO BE CLEANED”** on the Machine.

6.5.5 Transfer the Finished Goods to the Finished Goods Store.

6.5.6 Remove the Stereos from the Machine and hand over the Production Officer/Executive for verification and destruction in presence of QA Officer / Executive.

6.5.7 Run the Machine for cleaning.

6.5.8 Use Compressed Air to clean the inner side of the Machine.

6.5.9 Check the Machine and Area for remnant of the Processed Batch.

6.5.10 Clean the surface of the Machine with the 70% IPA.

6.5.11 Clean the area by Moping.

6.5.12 Put the Status Label **“CLEANED”** on the Machine.



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production	
		Effective Date:		
Revision No.:	00	Revision Date:		
Supersede Revision No.:	Nil	Page No.:	8 of 12	

6.5.13 Record the maintenance/ cleaning details in respective “Equipment’s /Instruments Usage Cleaning Log” as per **SOP No. HPD/078**.

7.0 ANNEXURES:

ANNEXURE No.	TITLE OF ANNEXURE	FORMAT No.
Annexure-I	Identification of Similar Looking Letters	
Annexure-II	Reference sample of rejected cartons due to improper coding	
Annexure-III	Challenge test record during operation of Hi-Cart machine	

ENCLOSURES: SOP Training Record.

8.0 DISTRIBUTION:

- Controlled Copy No. 01 Quality Assurance
- Controlled Copy No. 02 Production (I-Block)
- Controlled Copy No. 03 Production (Q-Block, Three Piece Section)
- Master Copy Quality Assurance

9.0 REFERENCES:

- SOP: Product Packing Operation in Parenteral blocks.

10.0 REVISION HISTORY:

CHANGE HISTORY LOG

Revision No.	Change Control No.	Detail of Changes	Reason for Change	Effective Date	Updated By



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:	00	Revision Date:	
Supersede Revision No.:	Nil	Page No.:	9 of 12

ANNEXURE –I

IDENTIFICATION OF SIMILAR LOOKING LETTERS

LETTER	SIMILAR LOOKING LETTER
C	G
E	F
I (I)	1 (One)
L	I (I)
O (O)	0 (Zero)
O (O)	Q (Q)
P	R



PHARMA DEVILS

PRODUCTION DEPARTMENT


STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production	
		Effective Date:		
Revision No.:	00	Revision Date:		
Supersede Revision No.:	Nil	Page No.:	10 of 12	


ANNEXURE –II

REFERENCE SAMPLE OF REJECTED CARTONS DUE TO IMPROPER CODING




“O” “Zero”


Differentiation between “O” and “0” (zero)




Incomplete printing of batch detail




Incomplete printing of batch detail



Incomplete printing of batch detail



Excess spreading of ink



Rubbing on printed batch detail




PHARMA DEVILS

PRODUCTION DEPARTMENT


STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production	
		Effective Date:		
Revision No.:	00	Revision Date:		
Supersede Revision No.:	Nil	Page No.:	11 of 12	



Faded printing of batch detail



Printing of batch detail in out of printing zone



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE

Title: Operation and Cleaning of Hi-Cart Cartoning Machine

SOP No.:		Department:	Production
		Effective Date:	
Revision No.:		Revision Date:	
Supersede Revision No.:		Page No.:	12 of 12

ANNEXURE –III

CHALLENGE TEST RECORD DURING OPERATION OF HI-CART MACHINE

Date	Product Name	B. No	Time	Challenge Test Observation(OK/Not OK)							Done By (Prod)	Checked By (Prod)	Verified By (IPQA)
				Vertical product check	Product check	Leaflet low level check	Leaflet check	Carton low level check	Carton check	Empty carton check			

- **Frequency of Challenge test:** Initially before start the machine operation and towards the end of operation.
- Give Justification if observation is found **Not OK**.