



Title: Operation of Purified Water Circulation Loop

| | | | |
|------------------------|--|-----------------------|--------|
| SOP No.: | | Revision No.: | 00 |
| Effective Date: | | Supersedes No. | Nil |
| Review Date: | | Page No. | 1 of 7 |

1.0 OBJECTIVE:

To lay down a procedure for Operation of Purified Water Circulation Loop.

2.0 SCOPE:

This SOP is applicable for Operation of Purified Water Circulation Loop.

3.0 RESPONSIBILITY:

Officer / Executive – Engineering

4.0 ACCOUNTABILITY:

Head – Engineering

5.0 DEFINATION:

5.1 Purified water is water that has been mechanically filtered or processed to remove impurities and make it suitable for use.

6.0 PROCEDURE:

6.1 PRE STARTUP CHECKS:

6.1.1 Check the Water Level of PW make up.

6.1.2 Check the Air Pressure should be 5 Kg / cm² for Pneumatic Valve Operation.

6.2 OPERATION:

6.2.1 Switch “ON” the Main Control of PW distribution system.

6.2.2 ‘Log in’ screen shall appear, enter the user name and password.

6.2.3 Press the PW distribution button, the PW distribution screen shall appear.

6.2.4 Press the AUTO button, then the START and STOP button shall appear.

6.2.5 Press the START button, the PW circulation pump shall start.

6.2.6 Switch ON the main panel of WFI distribution system.

6.2.7 ‘Log in’ screen shall appear, enter the user name and password.

6.2.8 Press the AUTO button, then the START and STOP button shall appear.

6.2.9 Press the START button, the WFI circulation pump shall start.



Title: Operation of Purified Water Circulation Loop

| | | | |
|------------------------|--|-----------------------|--------|
| SOP No.: | | Revision No.: | 00 |
| Effective Date: | | Supersedes No. | Nil |
| Review Date: | | Page No. | 2 of 7 |

6.2.10 Observe parameters mentioned in following table.

| S.No. | Parameter | Range |
|-------|-------------------|-------------------------|
| 1. | Flow Rate | NLT 7300 L/Hr |
| 2. | Conductivity | <1.3 μ S/cm |
| 3. | Temperature (WFI) | NLT 80 Deg. C |
| 4. | Temperature (PW) | Ambient |
| 5. | UV Intensity | NLT 30 w/m ² |

NOTE: The Purified Water dump valves will be operated if the limits of conductivity & TOC are out of limit (mentioned above).

6.3 VELOCITY:

6.3.1 Velocity (M/Sec.) for 2" Dia Loop Line Pipe = Flow (LPH) X K /1000

Where, **K** = Conversion Factor (0.167)

6.3.2 Velocity (M/Sec.) for 1" Dia Loop Line Pipe = Flow (LPH) X K /1000

Where, **K** = Conversion Factor (0.802)

6.4 STAND BY PUMP CHANGE OVER:

Stop the distribution pump operation.

6.4.1 Check the electrical terminal connection and Switch ON manually for few seconds for the stand by pump rotation direction.

6.4.2 Switch OFF the pump and open the pump outer cover.

6.4.3 Wash properly the impeller blade and outer cover with Purified water.

6.4.4 Then again rinse impeller and outer cover with Purified water.

6.4.5 Fix the outer cover in a position.

6.4.6 Collect water in storage tank for pump final rinsing.

6.4.7 Close the tank outlet main valve and position the pipe accessories to connect the stand by Pump.

6.4.8 Tighten the all TC clamps for the suction line and keep open the discharge line end.

6.4.9 Connect the hose pipe in discharge end mouth and open the tank out let valve.

6.4.10 Keep the hose pipe end in main drain of the water system area.

6.4.11 Switch ON the pump and allow water to drain out for 5 Min.



Title: Operation of Purified Water Circulation Loop

| | | | |
|------------------------|--|-----------------------|--------|
| SOP No.: | | Revision No.: | 00 |
| Effective Date: | | Supersedes No. | Nil |
| Review Date: | | Page No. | 3 of 7 |

6.4.12 Switch Off the pump and connect the pump discharge line to main loop line system.

6.4.13 Record the details in stand by pump change over record as per **Annexure-III**.

6.4.14 Frequency: Before sanitizing the loop system / whenever required / if any pump breakdown / when ever plant shutdown.

(**Note:** During all the above mentioned frequency cases Sanitization process should be carried out.)

6.5 Record the Operation Details in “**Purified Water Circulation Loop Operation Record**” as shown in **Annexure-I**.

6.6 Record the Purified Water and Water for Injection (SCADA Printout) Details in **Annexure-II, Titled “Record for Purified Water”**.

7.0 ABBREVIATIONS:

| | |
|-----|------------------------------|
| SOP | Standard Operating Procedure |
| QA | Quality Assurance |
| No. | Number |
| PW | Purified Water |
| WFI | Water For Injection |
| UVD | Ultra Violet Detector |
| kg | Kilogram |
| cm | Centimeter |
| TOC | Total Organic Carbon |

8.0 ANNEXURES:

| ANNEXURE No. | TITLE OF ANNEXURE | FORMAT No. |
|---------------------|--|-------------------|
| Annexure-I | Purified Water Circulation Loop Operation Record | |
| Annexure-II | Record for Purified Water (SCADA Printout) | |
| Annexure-III | Stand by Pump Change Over Record | |



Title: Operation of Purified Water Circulation Loop

| | | | |
|------------------------|--|-----------------------|--------|
| SOP No.: | | Revision No.: | 00 |
| Effective Date: | | Supersedes No. | Nil |
| Review Date: | | Page No. | 4 of 7 |

9.0 DISTRIBUTION:

- Controlled Copy No. 01 Head Engineering
- Master Copy Quality Assurance

10.0 REFERENCES:

Not Applicable

11.0 REVISION HISTORY:

| Revision No. | Change Control No. | Details of Changes | Reason of Changes | Effective Date | Done By |
|--------------|--------------------|--------------------|-------------------|----------------|---------|
| | | | | | |

