

PHARMA DEVILS

Title: Passivation of Water System

SOP No.:	Revision No.:	00
Effective Date:	Supersedes No.	Nil
Review Date:	Page No.	1 of 4

1.0 OBJECTIVE:

To lay down a procedure for passivation of water system.

2.0 SCOPE:

This SOP is applicable for passivation of water system.

3.0 RESPONSIBILITY:

Officer/Executive – Engineering

4.0 ACCOUNTABILITY:

Head – Engineering

5.0 ABBREVIATIONS:

EREngineeringHNO3Nitric AcidLtd.LimitedPvt.PrivateQCQuality ControlQtyQuantitySOPStandard Operating Procedure

6.0 **PROCEDURE**:

- **6.1** Fill the Purified Water Storage tank up to the half level.
- **6.2** Ensure that Suction and Discharge Valve of pump are open and all user points & sampling valves are closed.
- 6.3 Start the recirculation pump and circulate water for 30min.
- 6.4 Drain the Water from the tank and Piping via user points & Sampling Valves.
- 6.5 Prepare Cleaning Solution by adding 2% of HNO₃ solution in Purified Water Storage tank.
- **6.6** Start the Circulation pump and circulate the Cleaning Chemical for 60min.
- 6.7 Solution Shall be Prepared in Purified Water, and Quantity of the Nitric acid shall be computed as:

Qty. of $HNO_3 = 100*2/\%$ assay of the Solution

- 6.8 Flush the Cleaning Chemical from the entire loop and tank via Sampling and User Points.
- **6.9** Prepare Passivation chemical by diluting Sodium Hydroxide (NaOH) of 100% purity with purified Water to make 2% of the solution tank.



Title: Passivation of Water System

SOP No.:	Revision No.:	00
Effective Date:	Supersedes No.	Nil
Review Date:	Page No.	2 of 4

- 6.10 Start the Recirculation pump & circulate the Passivation Chemical for 60min.
- **6.11** Drain the Passivation Chemical from the equipment's and Circulation loop via user Points and Sampling Points.
- 6.12 Transfer Fresh Purified Water to the Purified Water Storage Tank.
- 6.13 After this, flush the Purified Water till the pH is neutral (5.0 to 7.0) or to the in feed Water Quality.
- **6.14** Frequency: in case of any Modification, Discontinuation and New Installation.

7.0 ANNEXURES:

ANNEXURES No.	TITLE OF ANNEXURE	FORMAT No.
Annexure-I	Intimation Form	
	Log for Passivation of Purified Water Storage and Distribution System	

ENCLOSURES: SOP Training Record

8.0 **DISTRIBUTION:**

- Controlled Copy No. 01
- Controlled Copy No. 02
- Master Copy

Quality Assurance Engineering Quality Assurance

9.0 **REFERENCES**:

As per Manual.

10.0 REVISION HISTORY:

CHANGE HISTORY LOG

Revision	Change Control	Details of Changes	Reason for	Effective	Updated
No.	No.		Change	Date	By



PHARMA DEVILS ENGINEERING DEPARTMENT

Title: Passivation of Water System

SOP No.:	Revision No.:	00
Effective Date:	Supersedes No.	Nil
Review Date:	Page No.	3 of 4

ANNEXURE – I INTIMATION FORM

FROM,

Engineering Department

SUBJECT: Passivation of Purified Water Storage and Distribution System

Date of Passivation :

S.No.	Department	Name	Remark	Signature & Date			
1.	Quality Assurance						
2.	Production						
3.	Quality Control						





Title: Passivation of Water System

SOP No.:	Rev	vision No.:	00
Effective Date:	Su	persedes No.	Nil
Review Date:	Pag	ge No.	4 of 4

ANNEXURE – II LOG FOR PASSIVATION OF PURIFIED WATER STORAGE AND DISTRIBUTION SYSTEM

AREA / BLOCK:

HNO ₃ HNO ₃ puri		Qty. of purified	ied (Hrs.)		pH of purified water	pH of recirculated water (pH of recirculated water	Remark	Sign &	Review by Sign & Date	
		taken	water taken	From	То		should be ± 0.2 of pH of purified water)		Date	
L	1	<u> </u>	<u> </u>	1			1	I	I	<u> </u>