



STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Procedure for Operation and Cleaning of Fogger Machine	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

1.0 OBJECTIVE:

To lay down the procedure for Operation & Cleaning of Fogger Machine.

2.0 SCOPE:

This SOP is applicable for procedure for operation and cleaning of Fogger Machine used in the Microbiology Laboratory.

3.0 RESPONSIBILITY:

Microbiologist- Quality Control
Head – Quality Control

4.0 PROCEDURE:

4.1 Operation:

- 4.1.1 Connect the fogger machine to electric pin to main supply.
- 4.1.2 Keep the fogger machine 2 feet above the ground level.
- 4.1.3 Open the lid of tank of fogger machine
- 4.1.4 Add 5 ltrs of Gramicid 3% v/v solution and close the lid.
- 4.1.5 Adjust the blower at 45° angle.
- 4.1.6 Put on the MAIN switch
- 4.1.7 Switch on the blower, Fogger starts emitting fine fog.
- 4.1.8 After Completion of fumigation, switch off the fogger machine.

4.2 Precautions:

- 4.2.1 Make sure the area to be fumigated cleaned thoroughly.
- 4.2.2 Use proper safety apparel such as gowns, hand gloves and safety goggles during fumigation

4.3 Cleaning Procedure :

- 4.3.1 Fogger machine (Outer Surface) should be clean from lint free cloth Before use.

5.0 ANNEXURE (S):

Nil



STANDARD OPERATING PROCEDURE

Department: Microbiology	SOP No.:
Title: Procedure for Operation and Cleaning of Fogger Machine	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

6.0 REFERENCE (S):

SOP: Preparation, Approval, Distribution control, revision and Destruction of Standard operating Procedure (SOP).

7.0 ABBREVIATION (S) / DEFINITION (S) :

SOP : Standard Operating Procedure

REVISION CARD

S.No.	REVISION No.	REVISION DATE	DETAILS OF REVISION	REASON (S) FOR REVISION	REFERENCE CHANGE CONTROL No.
1.	00	---	---	New SOP	---