

PHARMA DEVILS MICROBIOLOGY DEPARTMENT

DESIGN QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR FLOW UNIT

Pre - Execution Approval

	Name	Designation	Signature	Date
Prepared By				
Reviewed By				
Reviewed By				
Approved By				



MICROBIOLOGY DEPARTMENT

DESIGN QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR FLOW UNIT

1.0 Objective:

The purpose of Design Qualification of the Horizontal Laminar Flow Unit is to ensure that all the critical aspects including safety, effective cleaning, maintenance, utility requirement, working space, have been considered while designing the equipment and they are properly documented.

2.0 Scope:

Scope is limited to the following

Equipment / System Name	Vertical Laminar Flow Unit	
ID Number		
Location	Cooling Zone	

3.0 Equipment / System Description:

KlenzFloTM Vertical Laminar flow clean air workstation provides a contained environment to protect the product. It provide clean zone and particle free conditions by unidirectional down ward airflow from the HEPA filters situated above the working area. This shall produce clean zone with air quality, which is equivalent or better than class 5 ISO 144644 - 1 / specifications.

The Vertical Laminar flow clean air workstation is required to provide local class 5 (ISO 144644 - 1) grade environment for aseptic processing. The equipment should be designed for continuous run and shall not experience any major breakdown because of workmanship. The equipment is intended to be in operation in all the three shift 7 days a week.

The Vertical Laminar flow clean air workstation shall be used for maintaining class 5 (ISO 14644 - 1) environment during unloading of sterilized equipments, accessories, materials, etc.

4.0 Checklist for Design Verification:

Design of Horizontal Laminar Flow unit shall be verified for the compliance with the critical parameters mentioned in the URS/ executed tender.



MICROBIOLOGY DEPARTMENT

DESIGN QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR FLOW UNIT

Parameters	User requirement / Technical Specification	Vendor Design / Specification	Remarks	
Construction				
Material of construction	SS 304	SS 304		
Curtain	Not Mentioned	Antistatic PVC curtain		
HEPA filter Screen	SS 304	SS 304		
HEPA filters	Not Mentioned	EU - 13, 610x915x75 mm (Antimicrobial)		
Pre filters	Not Mentioned	EU – 6, 765x298x45 mm (Biocidal)		
	Not Mentioned	EU – 6, 728x298x45 mm (Biocidal)		
	Not Mentioned	EU – 6, 840x298x45 mm (Biocidal)		
DOP Port	Should be provided	100 % DOP test port		
Electrical construction	 a) 300 lux along with one number of 5/15 amp, each switch socket (combined) b) Independent controls for lighting, blowers & sockets. 	 a) Fluorescent light: 4'L ; 40-Watts b) Independent controls for lighting, blowers & Sockets. 		



MICROBIOLOGY DEPARTMENT

DESIGN QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR FLOW UNIT				
Parameters	User requirement / Technical Specification	Vendor Design / Specification	Remarks	
	Acceptable Tolerand	ce for Control Parameter		
Air velocity	Between 0.50 m/sec to 0.45 m/sec	0.45 ± 0.05 m/sec		
Noise (Decibels)	Less than 65	Not Mentioned		
Differential pressure	Should be between 7 to 12 mm WC	Not Mentioned		
Particle count	As per ISO 14644 - 1	Complies		
	Type of Control System			
Relay / Switchgear	Should be provided	Provided		
Indicators	 Power ON Blower ON Light switch 	Provided		
Differential pressure gauge	0 to 25 mm WC	0 to 25 mm WC		
Accuracy of gauge	± 3 %	± 3 %		
Electrical power supply	440 - 3 phase 50 Hz / 220 1 phase 50 Hz	230 V; AC; 50 Hz		



MICROBIOLOGY DEPARTMENT

DESIGN QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR FLOW UNIT

Parameters	User requirement / Technical Specification	Vendor Design / Specification	Remarks
Special features	Not Mentioned	Soft touch key pad switch	
	<u> </u>	Safety	
Alarm	Audiovisual alarm for motor blower trip	Audiovisual alarm for motor blower trip	
МСВ	Make Hager / Merlin gerlin / Legrand	Provided	
Motor	Not Mentioned	0.5 HP; 1440 RPM;	
Blower (Impeller)	Not Mentioned	210 x 85 mm	
Overall Dimensions (WxDxH)	72 x 36 inches (2743 x 914 mm)	1065 x 2746 x 610 mm	



MICROBIOLOGY DEPARTMENT

DESIGN QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR FLOW UNIT

5.0 Any Changes identified towards equipment design / lay out.

.....

6.0 Recommendations and Conclusions:

7.0 References:

User requirement specification.

Copy of purchase order.

Risk assessment.

Document submitted by the supplier / vendor

8.0 Documents Attached

Document Title	Annexure No.	
List of changes / deviation	Ι	
User requirement specification.	II	
Copy of purchase order	III	
Risk assessment	IV	
Document submitted by the supplier / vendor	V	

PHARMA DEVILS MICROBIOLOGY DEPARTMENT



9.0 Abbreviations:

- : User requirement specification : Stainless Steel URS
- SS
- mm : millimeter
- WC : water content
- % : Percentage
- : Meter per second m/sec





MICROBIOLOGY DEPARTMENT

DESIGN QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR FLOW UNIT

Annexure - I

List of changes / Deviations

S.No.	Description of Change / Deviations	Justification

Verified By:

Approved By:



PHARMA DEVILS MICROBIOLOGY DEPARTMENT

DESIGN QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR FLOW UNIT

Post execution approval:

	Name	Designation	Signature	Date
Compiled By				
Reviewed By				
Reviewed By				
Approved By				