

QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

S No.	ITEM DESCRIPTION	PAGE No.
1.0	PROTOCOL APPROVAL	2
2.0	OVERVIEW:	3
2.1	Objective	3
2.2	Purpose	3
2.3	Scope	3
2.4	Responsibility	3-4
2.5	Execution Team	5
3.0	ACCEPTANCE CRITERIA	6
4.0	REQUALIFICATION CRITERIA	6
5.0	OPERATION QUALIFICATION PROCEDURE	7
5.1	System Description	7
5.2	Instruction of filling the checklist	8
5.3	Test instrument details	9
5.4	Verification of Calibrated component	10
5.5	Verification of functional checks	11-12
5.6	Verification of safety Feature(s)	13
5.7	Verification of supporting utilities	14
5.8	Verification of standard Operating Procedure	15
5.9	Training Record of Personnel(s)	16
5.10	Deficiency And Corrective Action(s) (Reports)	17
5.11	List of Annexure(s)	18
6.0	OPERATIONAL QUALIFICATION FINAL REPORT	19
6.1	Summary	19
6.2	Conclusion	19
6.3	Final Report Approval	20



1.0 PROTOCOL APPROVAL:

Signing of this approval page of Protocol indicates agreement with the qualification approach described in this document. If modification to the qualification approach becomes necessary, an addendum shall be prepared and approved. The protocol cannot be used for execution unless approved by the following signatories.

This Operational Qualification protocol of Automatic airjet bottle air and vacuum cleaning machine has been reviewed and approved by the following signatories:

FUNCTION	NAME	DESIGNATION	DEPARTMENT	SIGNATURE	DATE
PREPARED BY			QUALITY		
<u>Б</u> Т			ASSURANCE		
			QUALITY		
REVIEWED			ASSURANCE		
BY			ENGINEERING		
			PRODUCTION		
			HEAD		
APPROVED			OPERATION		
BY			QUALITY		
			ASSURANCE		



2.0 OVERVIEW:

2.1 OBJECTIVE:

The objective of the operational qualification is to prove that each operation proceeds as per design specification and the tolerances prescribed there in the document are the same at utmost transparency.

The Qualification of Automatic airjet bottle air and vacuum cleaning machine performed in view of Dry Syrup of Production Cepha Oral manufacturing facility.

2.2 PURPOSE:

The purpose of this protocol is to establish the documentary evidence to ensure that the installed Automatic airjet bottle air and vacuum cleaning machine shall operate reproducibly and consistently within its full dynamic range of operation according to manufacturer's specification.

2.3 SCOPE:

The Scope of this protocol is limited to the operational Qualification of Automatic airjet bottle air and vacuum cleaning machine, installed in Dry Syrup of Production Cepha Oral manufacturing facility.

2.4 **RESPONSIBILITY:**

In accordance with protocol, following functions shall be responsible for the qualification of system.

Execution Team (Comprising members from Production, Engineering and Quality Assurance) and their responsibilities are following:

- Prepares the qualification protocol.
- Ensures that the protocol is in compliance with current policies and procedures on system Qualification.
- > Distributes the finalized protocol for review and approval signatures.
- Execution of Qualification protocol.
- Review of protocol, the completed qualification data package, and the final report.



- The operational checks, calibration, SOP verification, verification of safety features, verification of utility supply shall be carried out by engineering persons and production person.
- The production operator / supervisor shall carry out the cleaning and operation of machine.

Head – Production/ Engineering:

- Review of protocol, the completed qualification data package, and the final report.
- > Assist in the resolution of validation deficiencies.

Head – Operation and Quality Assurance:

Review and approval of protocol, the completed qualification data package, and the final report.



2.5 EXECUTION TEAM:

The satisfactory operation of the Automatic airjet bottle air and vacuum cleaning machine shall be verified by executing the qualification studies described in this protocol .The successfully executed protocol documents that the Automatic airjet bottle air and vacuum cleaning machine is operational and is satisfactorily working.

Execution team is responsible for the execution of operation of Automatic airjet bottle air and vacuum cleaning machine. All executors involved with this protocol shall sign within the prescribed format given below.

NAME	DESIGNATION	DEPARTMENT	SIGNATURE	DATE



OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

3.0 ACCEPTANCE CRITERIA:

- 3.1 The equipment shall be operational as per its specified operating instructions
- 3.2 All SOPs for the equipment to be verified and checked
- 3.3 Training shall be given to all the concerned personnel
- 3.4 All the functionality of equipment components to be checked
- 3.5 The RPM of motor should be in the range of $\pm 5\%$.

4.0 **REQUALIFICATION CRITERIA:**

The machine shall be requalified if

- There are any major changes in system components which affect the performance of the system
- After major breakdown maintenance is carried out.
- As per revalidation date and schedule

5.0 OPERATIONAL QUALIFICATION PROCEDURE:

5.1	SYSTEM DESCRIPTION	DN:	
1.	Equipment Name		Automatic airjet bottle air and vacuum cleaning
		•	machine
2.	Supplier/Manufacturer	:	Perth Engineers & Consultant
3.	Model	:	
4.	Serial no.	:	
5.	Capacity	:	80-100 bottles per minute
6.	Dimension	:	2172 mm (L) X 600 mm (W) X 1045 mm (H)
7.	Location	:	Dry Syrup



OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

5.1.1 MACHINE DESCRIPTION:

Automatic airjet bottle air and vacuum cleaning machine is used for cleaning of bottles for further operation to subsequent machine.

- The drive unit is fitted with flange mounted motor on gear box.
- This reduction gear box is having closed helical gear drive with oil bath.
- Separate AC frequency drive is fitted with machine for speed control of conveyor.
- This helical reduction is the major speed-reducing unit.
- The drive assembly with gear box & motor is mounted on structure.
- The machine is provided with various interlocking features e.g. Roller type switch for clutch, Air Pressure Switch etc. for safety purpose.

5.2 INSTRUCTION FOR FILLING THE CHECKLIST:

- 5.2.1 Write down the actual observation in observation column as per design specification
- 5.2.2 Observation functional parameter should be write actual function in specified column.
- 5.2.3 Give the detailed information in the summary and conclusion part of the Operational Qualification report.
- 5.2.4 Whichever column is blank or not used 'NA' shall be used.



5.3 TEST INSTRUMENT DETAILS

This test is intended to describe the equipments/instruments and its complete details to have a traceability to the national standard which is to be used for the verification of the operation of the Automatic airjet bottle air and vacuum cleaning machine.

S.No.	Name Of Instrument	Inst. ID. Number	Calibration done on	Calibration Due date	Certificate Number

Checked by Date:

Remark: -----





5.4 Verification of Calibrated component :

This test is intended to describe the equipments/instruments and its complete details to have a traceability to the national standard, which is to be used for the verification of the operation of the Automatic airjet bottle air and vacuum cleaning machine.

S.No.	Name of Instrument	Inst. ID. Number	Calibration done on	Calibration valid up to	Certificate number

Checked by Date:

Remark: -----

QUALITY ASSURANCE DEPARTMENT



OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

5.5 VERIFICATION OF FUNCTIONAL CHECKS:

Procedure	Acceptance Criteria	Observation	Checked By Sign/Date
Verification of Key Functionality of	of Control Panel:		
Main Power ON/OFF switch	To supply power to the unit with Green indication 'ON'		
Start the machine by pressing Green Push Button	The machine will get started		
Stop the machine by pressing Red Push Button	The machine will get stopped		
Check the speed variation by Potentiometer	Speed of machine will get varied		
Run Test:			
Run the machine in empty	There should not be any abnormal sound		
Power Failure & Restoration Test:			J
Start the machine as per	The machine should not		
manufacturer's instruction. Trip the	start until & unless it is		
main incoming power supply. Wait	started manually.		
for sometime & switch ON the			
machine.			
Remark:			



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

5.6 VERIFI	CATION OF SAFET	TY FEATURES :		
Safety Features	Method of	Acceptance Criteria	Observation	Verified By
Description		Ston ploto & unit will		Sign/Date
Roller type switch for clutch	Feed the cross bottle at input	Star plate & unit will get stopped		
Air Pressure Switch	Close the compressed air supply while the machine in operation	Machine will get stopped & LOW AIR indication shall glow		
Earthing	Run machine and check with multimeter/ clamp meter at the body cover of the machine.	No current should be sensed over whole body		

Remark: -----



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

5.7 VERIFICATION OF SUPPORTING UTILITIES:

S.No.	Utility	Method Of Verification	Observation	Checked By
1.	Electrical Power Supply: 3 phase, 415V +/- 5%, 50Hz supply with neutral and proper earthing	Physically with clamp meter		Sign/Date
2.	Output Voltage from Transformer: 220 V +/- 5%	Physically with clamp meter		
3.	Compressed Air: 6 bar	Physically		

Remark: -----

Reviewed by (Sign/Date)

5.8 VERIFICATION OF STANDARD OPERATING PROCEDURE (SOP)

The following Standard Operating Procedures were verified as important for effective performance of Automatic airjet bottle air and vacuum cleaning machine operation.

S.No.	SOP TITLE	SOP NUMBER	VERIFIED BY	DATE

Remark: -----



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

5.9 TRAINING RECORD OF PERSONNEL (S):

S.No.	Name of Personnel	Designation	Sign & Date	Trained By	Remark

Remark: -----



OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

5.10 DEFICIENCY AND CORRECTIVE ACTION (S) REPORT (S):

Following deficiency was verified and corrective actions taken in consultation with the Engineering Department.

Description of deficiency:

Corrective action(s) taken:

Deviation accepted by (Sign/Date)

Deviation Approved by (Sign/Date)



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

5.11 Annexure (S)

Annexure No.	Details Of Annexure

Remarks (if any):

Done By & Date:

Verified By & Date:



OPERATIONAL QUALIFICATION AUTOMATIC AIRJET BOTTLE AIR AND VACUUM CLEANING MACHINE

6.0 OPERATIONAL QUALIFICATION FINAL REPORT:

6.1 SUMMARY:

6.2 CONCLUSION:

Prepared By Sign/ Date Checked By Sign/ Date



6.3 FINAL REPORT APPROVAL

It has been verified that all tests required by this protocol are completed, reconciled and attached to this protocol or included in the qualification summary report. Verified that all amendments and discrepancies are documented, approved and attached to this protocol (If applicable). Signature in the block below indicates that all items in this qualification report of Automatic airjet bottle air and vacuum cleaning machine have been reviewed and found to be acceptable and that all variations or discrepancies (If applicable) have been satisfactorily resolved. After the successful operational qualification of the Automatic airjet bottle air and vacuum cleaning machine, the equipment can be taken for performance qualification.

FUNCTION	NAME	DESIGNATION	DEPARTMENT	SIGNATURE	DATE
REVIEWED BY			QUALITY ASSURANCE		
			ENGINEERING		
			PRODUCTION		
APPROVED BY			HEAD OPERATION		
			QUALITY ASSURANCE		