

QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

S.No.	ITEM DESCRIPTION	PAGE No.
1.0	PROTOCOL APPROVAL	2
2.0	OVERVIEW:	3
2.1	Objective	3
2.2	Purpose	3
2.3	Scope	3
2.4	Responsibility	3-4
2.5	Execution Team	5
3.0	ACCEPTANCE CRITERIA	6
4.0	REQUALIFICATION CRITERIA	6
5.0	OPERATION QUALIFICATION PROCEDURE	7
5.1	System Description	7
5.2	Instruction of filling the checklist	8
5.3	Test instrument details	9
5.4	Verification of Calibrated component	10
5.5	Verification of functional checks	11-12
5.6	Verification of safety Feature(s)	13
5.7	Verification of supporting utilities	14
5.8	Verification of standard Operating Procedure	15
5.9	Training Record of Personnel(s)	16
5.10	Deficiency And Corrective Action(s) (Reports)	17
5.11	List of Annexure(s)	18
6.0	OPERATIONAL QUALIFICATION FINAL REPORT	19
6.1	Summary	19
6.2	Conclusion	19
6.3	Final Report Approval	20



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

1.0 PROTOCOL APPROVAL:

Signing of this approval page of Protocol indicates agreement with the qualification approach described in this document. If modification to the qualification approach becomes necessary, an addendum shall be prepared and approved .The protocol cannot be used for execution unless approved by the following authorities.

FUNCTION	NAME	DESIGNATION	DEPARTMENT	SIGNATURE	DATE
PREPARED BY			QUALITY ASSURANCE		
REVIEWED			QUALITY ASSURANCE		
BY			ENGINEERING		
			PRODUCTION		
APPROVED			HEAD OPERATION		
ВҮ			QUALITY ASSURANCE		



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

2.0 **OVERVIEW:**

2.1 **OBJECTIVE:**

The objective of the operational qualification is to prove that each operation proceeds as per design specification and the tolerances prescribed there in the document are the same at utmost transparency.

The Qualification of Automatic Measuring Cup Placement performed in view of Dry Syrup of manufacturing facility.

2.2 PURPOSE:

The purpose of this protocol is to establish the documentary evidence to ensure that the installed Automatic Measuring Cup Placement shall operate reproducibly and consistently within its full dynamic range of operation according to manufacturer's specification.

2.3 SCOPE:

This protocol shall define the test procedures, documentation, references and acceptance criteria to establish that the Automatic Measuring Cup Placement operates and performs as intended in accordance with the design qualification.

The Scope of this protocol is limited to the operational Qualification of Automatic Measuring Cup Placement installed in Dry Syrup in Production Cepha Oral Block of manufacturing,

Once the operational qualification of Automatic Measuring Cup Placement has been completed successfully, the equipment shall be preceded for the performance qualification procedure.

2.4 RESPONSIBILITY:

In accordance with protocol, following functions shall be responsible for the qualification of system.

Execution Team (Comprising members from Production, Engineering and Quality Assurance) and their responsibilities are following:

- > Prepares the qualification protocol.
- ➤ Ensures that the protocol is in compliance with current policies and procedures on system Qualification.



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

- ➤ Distributes the finalized protocol for review and approval signatures.
- > Execution of Qualification protocol.
- ➤ Review of protocol, the completed qualification data package, and the final report.
- ➤ The operational checks, calibration, SOP verification, verification of safety features, verification of utility supply shall be carried out by engineering persons and production person.
- ➤ The production operator / supervisor shall carry out the cleaning and operation of machine.

Head – Production/ Engineering:

- ➤ Review of protocol, the completed qualification data package, and the final report.
- Assist in the resolution of validation deficiencies.

Head – Operation and Quality Assurance:

➤ Review and approval of protocol, the completed qualification data package, and the final report.



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

2.5 EXECUTION TEAM:

The satisfactory operation of the Automatic Measuring Cup Placement shall be verified by executing the qualification studies described in this protocol .The successfully executed protocol documents that the Automatic Measuring Cup Placement is operational and is satisfactorily working.

Execution team is responsible for the execution of operation of Automatic Measuring Cup Placement. All executors involved with this protocol shall sign within the prescribed format given below.

NAME	DESIGNATION	DEPARTMENT	SIGNATURE	DATE



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

3.0 ACCEPTANCE CRITERIA:

- 3.1 The equipment shall be operational as per its specified operating instructions
- 3.2 All SOPs for the equipment to be verified and checked.
- 3.3 Training shall be given to all the concerned personnel.
- 3.4 All the functionality of equipment components to be checked.

4.0 REQUALIFICATION CRITERIA:

The machine shall be requalified if

- There are any major changes in system components which affect the performance of the system
- After major breakdown maintenance is carried out.
- As per revalidation date and schedule



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.0 OPERATIONAL QUALIFICATION PROCEDURE:

5.1 SYSTEM DESCRIPTION:

1 Equipment Name . Automatic Measuring Cup Placement

2 Supplier/Manufacturer . Parth Engineers & Consultant

3 Model :

4 Serial no. :

5 Capacity . 80-100 BPM

6 Dimension . 1800 mm (L) X 745 mm (W) X 1525 mm (H)

7 Location . Dry Syrup

5.1.1 MACHINE DESCRIPTION:

Automatic Measuring Cup Placement Machine is used to put cup on filled bottles & transfer it for further operation to subsequent machine.

- The drive unit is fitted with flange mounted motor (0.75 kw/ 1.87 A/ 280-400 V/ 50 Hz/ 1400 RPM) on reduction gear box.
- This reduction gear box is having closed helical gear drive with oil bath.
- Separate AC frequency drive is fitted with machine for speed control of conveyor.
- This helical reduction is the major speed- reducing unit.
- The drive assembly with gearbox & motor is mounted on structure.



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.2 INSTRUCTION FOR FILLING THE CHECKLIST

- 5.2.1 Write down the actual observation in observation column as per design specification
- 5.2.2 Observation functional parameter should be write actual function in specified column.
- 5.2.3 Give the detailed information in the summary and conclusion part of the Operational Qualification report.
- 5.2.4 Whichever column is blank or not used 'NA' shall be used.

5.3 TEST INSTRUMENT DETAILS

This test is intended to describe the equipments / instruments and its complete details to have traceability to the national standard which is to be used for the verification of the operation of the Automatic Measuring Cup Placement.

S.No.	Name of Instrument	Inst. ID. Number	Calibration done on	Calibration Due date	Certificate Number

Checked by Date:	
Remark:	
Reviewed by (Sign/Date)	



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.4 Verification of Calibrated component:

Reviewed by (Sign/Date)

This test is intended to describe the equipments/instruments and its complete details to have a traceability to the national standard, which is to be used for the verification of the operation of the Automatic Measuring Cup Placement.

S.No.	Name of Instrument	Inst. ID. Number	Calibration done on	Calibration valid up to	Certificate number

Checked by Dat	te:		
Remark:		 	



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.5 VERIFICATION OF FUNCTIONAL CHECKS:

5.5.1 VERIFICATION OF KEY FUNCTIONALITY:

NAME OF EQUIPMENT COMPONENT	SPECIFIED FUNCTION	OBSERVATION	Checked By Sign/Date
Conveyor assembly	To carry the filled & sealed		
	bottle for placing of measuring cup.		
Hopper assembly	To carry the measuring cup		
	for feed up to the bottle		
	through chute.		
Drive assembly	Main gear motor to drive		
	the feeder conveyor.		
Output capacity	80-100 bottle / Min.		
Machine Speed	Minbottle/ Min.		
	Maxbottle/Min.		

Remark:	 -



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.5.2 Verification of Functionality of Control Panel:

Procedure	Acceptance Criteria	Observation	Checked By Sign/Date
Machine START/ STOP Green	To start & stop Conveyor		
Actuator on Control panel	Green indication On/ Off		
Check the speed variation by	Speed of machine will get		
Potentiometer	varied		
Run Test:			
Run the machine in empty	There should not be any		
	abnormal sound		
Power Failure & Restoration Test:			
Start the machine as per	The machine should not		
manufacturer's instruction. Trip the	start until & unless it is		
main incoming power supply. Wait	started manually.		
for sometime & switch ON the			
machine.			

Remark:	 	 	



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.6 VERIFICATION OF SAFETY FEATURES:

Safety	Method of	Acceptance Criteria	Observation	Verified
Features	Verification			$\mathbf{B}\mathbf{y}$
Description				Sign/Date
Cup sensor	In operation, block	Machine shall get		
	the cup in the chute	stopped		
Earthing	Run machine and	No current should be		
	check with	sensed over whole		
	multimeter / clamp	body		
	meter at the body			
	cover of the			
	machine.			

Remark:		 	
Reviewed	by (Sign/Date)	 	

5.7 VERIFICATION OF SUPPORTING UTILITIES:

S.No.	Utility	Method Of	Observation	Checked
		Verification		$\mathbf{B}\mathbf{y}$
				Sign/Date
	Electrical Power Supply:	Physically with		
1.	3 phase, 230- 400 V, 50Hz	clamp meter		
	supply with neutral and proper			
	earthing			



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.8 VERIFICATION OF STANDARD OPERATING PROCEDURE (SOP)

The following Standard Operating Procedures were verified as important for effective performance of Automatic Measuring Cup Placement operation.

S.No.	SOP TITLE	SOP NUMBER	VERIFIED BY	DATE
Remark	ζ:			



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.9 TRAINING RECORD OF PERSONNEL (S):

S.No.	Name of Personnel	Designation	Sign & Date	Trained By	Remark

Remark:	



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.10 DEFICIENCY AND CORRECTIVE ACTION (S) REPORT (S):

Following deficiency was verified and corrective actions taken in consultation with the Engineering Department.

Corrective action(s) taken:

Deviation accepted by (Sign/Date)

Deviation Approved by (Sign/Date)



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

5.11 Annexure (S)

Annexure No.	Details Of Annexure
_	
Remarks (if any):	
Done By & Date:	Verified By & Date:



QUALITY ASSURANCE DEPARTMENT

		CONTRACTOR DELIVERATION	
C	PERAT	TIONAL QUALIFICATION FOR AUTOMATIC MEASURING	CUP PLACEMENT
	6.0	OPERATIONAL QUALIFICATION FINAL REPORT:	
	6.1	SUMMARY:	
	6.2	CONCLUSION:	
	Prepar Sign/	red By Date	Checked By Sign/ Date



QUALITY ASSURANCE DEPARTMENT

OPERATIONAL QUALIFICATION FOR AUTOMATIC MEASURING CUP PLACEMENT

6.3 FINAL REPORT APPROVAL:

It has been verified that all tests required by this protocol are completed, reconciled and attached to this protocol or included in the qualification summary report. Verified that all amendments and discrepancies are documented, approved and attached to this protocol, If applicable. Signature in the block below indicates that all items in this qualification report of Automatic Measuring Cup Placement have been reviewed and found to be acceptable and that all variations or discrepancies (If applicable) have been satisfactorily resolved. After the successful operational qualification of the Automatic Measuring Cup Placement the equipment can be taken for performance qualification.

FUNCTION	NAME	DESIGNATION	DEPARTMENT	SIGNATURE	DATE
DEVIEWED			QUALITY ASSURANCE		
REVIEWED BY			ENGINEERING		
			PRODUCTION		
APPROVED			HEAD OPERATION		
BY			QUALITY ASSURANCE		