



**INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER**

**Pre - Execution Approval**

	<b>Name</b>	<b>Designation</b>	<b>Signature</b>	<b>Date</b>
Prepared By				
Reviewed By				
Reviewed By				
Reviewed By				
Approved By				



## INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER

### 1.0 Objective:

- To determine that the equipment operates according to specifications, and to record all relevant information's and data to demonstrate its functions as intended for.

### 2.0 Scope:

Scope is limited to the following

<b>Equipment / System Name</b>	<b>M AIR T AIR SAMPLER</b>
<b>ID Number</b>	.....
<b>Location</b>	<b>Incubator Room</b>

### 3.0 Checklist for Operational verification:

Operation of M Air T air sampler is verified for the compliance with the critical parameters mentioned in the Functional Specification. Sequentially prepare the list of operating checks specifying critical parameters, perform them as per operating procedure as mentioned in the vender operating manual & document the same in the attached check list (Annexure - I).

### 4.0 Any Changes/Deviations identified during operating checks:

Refer Annexure – II

### 5.0 Identification & preparation of Standard Operating Procedures:

Prepared the SOP for Operation, Cleaning and Calibration of M Air T air sampler using the actual feedback from the operation checks and vendor-operating manual. The preventive maintenance procedure and schedule shall be applicable as per Annual Maintenance Contract of QC Instrument.

<u>S. No.</u>	<u>SOP Title</u>	<u>SOP Number</u>
1.	Operation, Cleaning and Calibration Of M Air T Air Sampler	.....



## INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER

### 6.0 Training

The List of QC (Microbiology) person trained during the operation qualification of M Air T air sampler is listed as per Annexure - III.

### 7.0 Recommendations and Conclusions:

### 8.0 References:

Installation Qualification

Operating Manual submitted by the Vendor.

### 9.0 Annexure

Annexure - I : Checklist for Operational Verification.

Annexure - II : List of Changes / Deviation.

Annexure - III : Training Detail

### 10.0 Abbreviations:

SOP : Standard Operating Procedure

QC : Quality Control

### Post execution approval:

	Name	Designation	Signature	Date
Compiled By				
Reviewed By				
Reviewed By				
Reviewed By				
Approved By				



**INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER**

**Annexure - I**

**Checklist for Operational Verification**

S.No.	Operating Parameters	Function /Acceptance criteria	Observation	Remarks
<b>A.</b>	<b>Verification of Cassette Sterility Qualification</b>			
1.	Cassette Catalogue Code	Not applicable		
2.	Incubation Time (hours)	Not applicable		
3.	Lot number of the cassette	Not applicable		
4.	Piece of cassette tested	Not applicable		
5.	M Air T Cassettes	Sterile/nonsterile at $32.5 \pm 2.5^{\circ}\text{C}$		
6.	M Air T Cassettes	Sterile/nonsterile at $22.5 \pm 2.5^{\circ}\text{C}$		
<b>B.</b>	<b>Verification of Cassette Growth Promotion Qualification</b>			
1.	Incubation Temperature ( $^{\circ}\text{C}$ )	Not applicable		
2.	Incubation Time (Hours)	Not applicable		
3.	Cassette Catalogue Number	Not applicable		
4.	Cassette Lot Number	Not applicable		



## INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER

S.No.	Operating Parameters	Function /Acceptance criteria	Observation	Remarks
<b>B.</b>	<b>Verification of Cassette Growth Promotion Qualification</b>			
5.	Name of the microorganisms use for Growth promotion test	Not applicable		
6.	<hr/> Cfu observed in Petriplates (90 mm)	Not applicable		
	<hr/> Cfu observed in M Air T Cassettes	Not applicable		
	Recovery	$\geq 70 \%$		
7.	<hr/> Cfu observed in Petriplates (90 mm)	Not applicable		
	<hr/> Cfu observed in M Air T Cassettes	Not applicable		
	Recovery	$\geq 70 \%$		



**INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER**

S.No.	Operating Parameters	Function /Acceptance criteria	Observation	Remarks
<b>B.</b>	<b>Verification of Cassette Growth Promotion Qualification</b>			
8.	Cfu observed in Petriplates (90 mm)	Not applicable		
	Cfu observed in M Air T Cassettes	Not applicable		
	Recovery	$\geq 70 \%$		
9.	Cfu observed in Petriplates (90 mm)	Not applicable		
	Cfu observed in M Air T Cassettes	Not applicable		
	Recovery	$\geq 70 \%$		
10.	Cfu observed in Petriplates (90 mm)	Not applicable		
	Cfu observed in M Air T Cassettes	Not applicable		
	Recovery	$\geq 70 \%$		



**INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER**

S.No.	Operating Parameters	Function /Acceptance criteria	Observation	Remarks
<b>C.</b>	<b>Verification of Functional Testing</b>			
1.	Clean Room	Class 100.000 for USP (D for E.P.)		
2.	Clean Room Volume	> 10 m <sup>3</sup>		
3.	Air born contamination	20 - 100 cfu (m <sup>3</sup> )		
4.	Clean Room Temperature	4 <sup>0</sup> - 45 <sup>0</sup> C		
5.	Clean Room Pressure	Atmospheric		
6.	Clean Room Humidity	5 % 95 % non condensing		
7.	Cassette Catalogue Code	Not applicable		
8.	Cassette Lot Number	Not applicable		
9.	Incubation Temperature (°C)	Not applicable		
10.	Incubation Time (Hours)	Not applicable		



**INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER**

S.No.	Operating Parameters	Function /Acceptance criteria	Observation	Remarks
<b>C.</b>	<b>Verification of Functional Testing</b>			
11	Sampling Position	Not applicable		
12.	Cfu observed in M Air T under validation	Not applicable		
	Cfu observed in Qualified M Air T	Not applicable		
13.	Cfu observed in M Air T under validation	Not applicable		
	Cfu observed in Qualified M Air T	Not applicable		
14.	Cfu observed in M Air T under validation	Not applicable		
	Cfu observed in Qualified M Air T	Not applicable		





**INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER**

**Annexure - II**

**List of Changes / Deviations**

S.No.	Description of Change / Deviations	Justification based on impact analysis

**Verified By:**

**Approved By:**



**INSTALLATION QUALIFICATION PROTOCOL FOR AIR SAMPLER**

**Annexure - III**

**Training Detail**

S.No.	Name of Trainee	Name of Trainer

**Verified By:**

**Approved By:**