

MICROBIOLOGY DEPARTMENT

### INSTALLATION QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR AIR FLOW

### **Pre - Execution Approval**

	Name	Designation	Signature	Date
Prepared By				
Reviewed By				
Reviewed By				
Reviewed By				
Approved By				

## 1.0 Objective:

- To establish documented evidence which will provide a High degree of assurance and reliability about the performance of the Vertical Laminar flow unit.
- To determine that the Vertical Laminar flow unit performs as intended by repeatedly running the system on its intended schedules and recording all relevant information and data. Results must demonstrate that performance consistently meets predefined specifications under normal conditions and where appropriate for worst case situations.

### **2.0 Scope:**

### Scope is limited to the following

Equipment / System Name	Vertical Laminar Flow Unit
ID Number	•••••
Location	Cooling Zone



MICROBIOLOGY DEPARTMENT

### INSTALLATION QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR AIR FLOW

### 3.0 Performance Verification Checklist:

Verify the performance of all the critical process/functions as per their respective procedure. Refer Performance Verification Checklist as per Annexure - I.

## 4.0 Any Changes/Deviations identified during operating checks:

Refer Annexure - II

### 5.0 Recommendations and Conclusions:

#### **6.0 References:**

**Design Qualification** 

Installation Qualification

Operational Qualification

Standard Operating Procedure

#### 7.0 Annexure

Annexure - I : Performance Verification Checklist

Annexure - II : List of Changes / Deviation.

Annexure - III: Performance Qualification Report.

#### 8.0 Abbreviations:

SOP : Standard Operating Procedure



MICROBIOLOGY DEPARTMENT

## INSTALLATION QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR AIR FLOW

## Post execution approval:

	Name	Designation	Signature	Date
Compiled By				
Reviewed By				
Reviewed By				
Reviewed By				
Approved By				



MICROBIOLOGY DEPARTMENT

## INSTALLATION QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR AIR FLOW

## Annexure - I Checklist for Performance Verification

S.No.	Test		Acceptance Criteria	Actual Result	Remarks	Verification Done by
		I <sup>st</sup> Location				
1.	HEPA Filter Integrity	2 <sup>nd</sup> Location	0. 001 %			
1.	(% Penetration)	3 <sup>rd</sup> Location	penetration			
		4 <sup>th</sup> Location				
		I <sup>st</sup> Location				
2a.	Average Air Velocity (fpm) (At supply)	2 <sup>nd</sup> Location	90 ± 10 %			
Za.		3 <sup>rd</sup> Location				
		4 <sup>th</sup> Location				
		I <sup>st</sup> Location				
2b.	Average Air Velocity (fpm) (At exhaust)	2 <sup>nd</sup> Location	Not applicable			
		3 <sup>rd</sup> Location				
		4 <sup>th</sup> Location				



MICROBIOLOGY DEPARTMENT

## INSTALLATION QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR AIR FLOW

S.No.	Test		Acceptance Criteria	Actual Result	Remarks	Verification Done by
		I <sup>st</sup>				
Pressure Drop (mm of WC)	2 <sup>nd</sup>	10 to 20				
3a.	(Across Supply HEPA Filter)	3 <sup>rd</sup>	mm of WC			
		4 <sup>th</sup>				
		I <sup>st</sup>				
21.	Pressure Drop (mm of WC)	2 <sup>nd</sup>	Not applicable			
3b. (Within Wo area)	(Within Work	3 <sup>rd</sup>				
		4 <sup>th</sup>				
				1.		
		> 0.5 µm	3520	2.		
4a. p	Particle Count (Cum. # of			3.		
	particles / cu.m of air) (If applicable)	of air) applicable)		1.		
	(I <sup>st</sup> Location)		29	2.		
				3.		



MICROBIOLOGY DEPARTMENT

## INSTALLATION QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR AIR FLOW

S.No.	Test		Acceptance Criteria	Actual Result	Remarks	Verification Done by
				1.		
		> 0.5 µm	3520	2.		
	Particle Count (Cum. # of particles / cu.m			3.		
4b.	of air) (If applicable)			1.		
	(2 <sup>nd</sup> Location)	> 5.0 µm	29	2.		
				3.		
		# of	3520	1.		
				2.		
	Particle Count (Cum. # of particles / cu.m			3.		
4c.	of air) (If applicable)			1.		
	(3 <sup>rd</sup> Location)	(3 <sup>rd</sup> Location) > 5.0 μm	29	2.		
				3.		
	Particle Count	m. # of les / cu.m f air) plicable) > 0.5 µm	3520	1.		
4d.	(Cum. # of particles / cu.m of air)			2.		
	(If applicable) (4 <sup>th</sup> Location)			3.		



MICROBIOLOGY DEPARTMENT

## INSTALLATION QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR AIR FLOW

S.No.	Test		Acceptance Criteria	Actual Result	Remarks	Verification Done by
4d.	Particle Count (Cum. # of particles / cu.m of air) (If applicable) (4 <sup>th</sup> Location)	> 5.0 µm	29	1.       2.       3.		
5.	Passive Air Sampling (Settle Plate Exposure) (Grade A)		1 (cfu/4 hrs)	1 <sup>st</sup> Day  2 <sup>nd</sup> Day  3 <sup>rd</sup> Day		
6.	Active Air Sampling (Grade A)		1 (cfu/m³)	1 <sup>st</sup> Day  2 <sup>nd</sup> Day  3 <sup>rd</sup> Day		

Note:

 $\textbf{Remarks:} \ \ \text{Performance Verification of Dynamic Pass Box - Complies} \ / \ \ \text{Does not comply}.$ 



MICROBIOLOGY DEPARTMENT

## INSTALLATION QUALIFICATION PROTOCOL FOR VERTICAL LAMINAR AIR FLOW

## Annexure - II

## **List of Changes / Deviations**

S.No.	Description of Change / Deviations	Justification based on impact analysis

Verified	By:
----------	-----

Approved By: