



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

| | |
|--------------------------------|----------------------------------|
| EQUIPMENT ID. No. | |
| LOCATION | Solution Preparation Room |
| DATE OF QUALIFICATION | |
| SUPERSEDES PROTOCOL No. | NIL |



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

CONTENTS

| S.No. | TITLE | Page No. |
|--------------|---|-----------------|
| 1.0 | Pre-Approval | 3 |
| 2.0 | Objective | 4 |
| 3.0 | Scope | 4 |
| 4.0 | Responsibility | 5 |
| 5.0 | Equipment Details | 6 |
| 6.0 | Equipment Description | 6 |
| 7.0 | Pre-Qualification Requirements | 7 |
| 8.0 | Critical Variables to be Met | 8-11 |
| 9.0 | References | 12 |
| 10.0 | Documents to be Attached | 12 |
| 11.0 | Deviation from Pre-Defined Specification, If Any | 13 |
| 12.0 | Change Control, If Any | 13 |
| 13.0 | Review (Inclusive of follow up action, If Any) | 13 |
| 14.0 | Conclusion | 14 |
| 15.0 | Recommendation | 14 |
| 16.0 | Abbreviations | 15 |
| 17.0 | Protocol Post Approval | 16 |



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

1.0 PRE – APPROVAL:

INITIATED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|--|-------------|------------------|-------------|
| OFFICER/EXECUTIVE (QUALITY ASSURANCE) | | | |

REVIEWED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|-------------------------------|-------------|------------------|-------------|
| HEAD (PRODUCTION) | | | |
| HEAD (ENGINEERING) | | | |

APPROVED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|-------------------------------------|-------------|------------------|-------------|
| HEAD (QUALITY ASSURANCE) | | | |



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

2.0 OBJECTIVE:

- To provide documented evidence for the Installation Qualification of Stirrer.
- To confirm that the equipment and its components are installed as per the Specifications mentioned in the design qualification document and other requirements given by supplier

3.0 SCOPE:

- The Scope of Qualification Document is limited to the Installation Qualification of **Stirrer** (**Make:**) to be installed in **Solution Preparation Room, Coating area**.
- The equipment shall be operated under the dust free environment and conditions as per the cGMP requirements.



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

4.0 RESPONSIBILITY:

The Validation Group, comprising of a representative from each of the following departments, shall be responsible for the overall compliance of this Protocol cum Report:

| DEPARTMENTS | RESPONSIBILITIES |
|--------------------------|---|
| Quality Assurance | <ul style="list-style-type: none">• Initiation, Review, Approval, Authorization, and Compilation of the Installation Qualification Protocol cum Report.• Co-ordination with Production and Engineering to carryout Installation Qualification.• Monitoring of Installation Qualification Activity.• Post Approval of Installation Qualification Protocol cum Report after Execution. |
| Production | <ul style="list-style-type: none">• Review & Pre Approval of Installation Qualification Protocol cum Report.• To Co-ordinate and support for Execution of Qualification study as per Protocol.• Post Approval of Installation Qualification Protocol cum Report after Execution. |
| Engineering | <ul style="list-style-type: none">• Review & Pre Approval of Installation Qualification Protocol cum Report.• Co-ordination, Execution and technical support in Stirrer Installation Qualification Activity.• Responsible for Trouble Shooting (if occurs during execution).• Post Approval of Installation Qualification Protocol cum Report after Execution. |



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

5.0 EQUIPMENT DETAILS:

| | |
|---------------------------------|-----------------------------|
| Equipment Name | Stirrer |
| Equipment ID. | |
| Manufacturer's Name | Techmach Engineering Works. |
| Model .No | GMP |
| Supplier's Name | Techmach Engineering Works. |
| Location of Installation | Solution Preparation Room |

6.0 EQUIPMENT DESCRIPTION:

Stirrer is suitable for emulsifying, dispersing, mixing and comminuting of liquid to Liquid products. It is based on rotor- stator principle. It is available in plain as well as water jacketed model which are suitable for heat sensitive products.

Three way cock system for drainage & recirculation of liquids provided as standard. Extra discharge spout provided as a standard for viscous products.

Special design facilitates adjustment of the grinding gap by an exterior screw by means of handle even during operation.

Stirrer is an important step in pharmaceutical manufacturing process; this equipment is a self Contained & portable unit for the process of size reduction. It uses the principle of impact of air. The product is dropped axially from the hopper in a communication chamber where it comes in contact with blades rotating at high speed.

Operation:

Product is fed to the operating area of a rotor, having a speed of 2800 RPM by specially designed feed device. The product is processed by high shear, pressure & friction between two Phase, and also, which exerts their force on it by means of pressing & releasing action..



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

7.0 PRE – QUALIFICATION REQUIREMENTS:

7.1 Verification of Documents:

- Executed and approved design qualification document
- Piping and instrumentation diagram (P& ID)
- Technical specification of equipment
- Certificate of material of construction of components.

7.1.1 Procedure:

- Verify the above mentioned documents for availability, completeness and approval status
- If any deviation is observed the same has to be recorded giving reasons for deviation and approved. Deviation should be approved by Authorized person.
- Approved Drawings and supporting documents would form a part of the IQ Protocol cum report.

7.1.2 Acceptance Criteria:

- All the documents should be available, complete and approved by respective authorities.



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

8.0 CRITICAL VARIABLES TO BE MET:

8.1 Installation Qualification Checklist:

| S.No. | Installation Check | Observation | Observed by (Engineering) Sign/ Date |
|-------|--|-------------|--|
| 1. | Check for the Dimensional accuracy | | |
| 2. | Check for the receipt of the consignment in good condition | | |
| 3. | Check for any scratches on the machine body | | |
| 4. | Check for the electrical panel. All Electrical connections should be as per the Circuit Diagram. | | |
| 5. | Check the Rotor Assembly Free Movements | | |
| 6. | Check The Grease in the Bearing Housing | | |
| 7. | Check the Direction of Rotation | | |

**Checked By
(Production)**

Sign/Date:

Verified By

(Quality Assurance)

Sign/Date:

Inference:

.....

.....

.....

.....

.....

.....

Reviewed By

(Manager QA)

Sign/Date:



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

8.2 Technical specification:

| S.No. | Component | Location | Observation | Observed by Engineering (Sign/Date) |
|-------|--------------|--|-------------|--------------------------------------|
| 1. | Model.No | GMP | | |
| 2. | Capacity | Std | | |
| 3. | Main Motor | Make : Hindustan Motor speed : 2800 RPM (±10%) Supply : 415 ,3Phase, 50 Hz Type : Flange mounted, TEFC Frame : 90 L KW/HP : 2.25/3 | | |
| 4. | FLP Starter | Make : FCG HP : 3 Relay : 4 to 6 amp | | |
| 5. | Castor Wheel | Make : Swift Size : 65 x 25mm Model : SSPU6525M | | |

**Checked By
(Production)**

Sign/Date:

Verified By

(Quality Assurance)

Sign/Date:

Inference:

.....

.....

.....

.....

Reviewed By

(Manager QA)

Sign/Date:



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

8.3 MOC Verification List:

| S.No. | Component | MOC | Observation | Observed by Engineering (Sign/Date) |
|-------|------------------|-------|-------------|-------------------------------------|
| 1. | Rotor | SS316 | | |
| 2. | Cap On Rotor | SS316 | | |
| 3. | Center Bolt | SS316 | | |
| 4. | Upper Stator | SS316 | | |
| 5. | Lower Stator | SS316 | | |
| 6. | Baffle | SS316 | | |
| 7. | Body Cover | SS304 | | |
| 8. | Top Cover | SS304 | | |
| 9. | Motor Housing | CI | | |
| 10. | Base For Housing | CI | | |

Checked By
(Production)
Sign/Date:

Verified By
(Quality Assurance)
Sign/Date:

Inference:

.....
.....
.....
.....

Reviewed By
(Manager QA)
Sign/Date:



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

8.4 SAFETY:

| Critical variables | Acceptance criteria | Observation | Observed by Engineering (Sign/Date) |
|------------------------|--|-------------|-------------------------------------|
| Mechanical Guard | Mechanical guard for all rotating parts. | | |
| Joints | Welding of joints without any welding burrs. | | |
| Metal Parts | All the metal parts should be properly grounded without any sharp Edges. | | |
| Leveling and Balancing | Equipment should be properly balanced & leveled. | | |
| Electrical Wiring | Electrical wiring should be as per approved drawings. | | |
| Noise Level | Below 80 db. | | |
| Emergency Switch | Provided easy access position. | | |

**Checked By
(Production)**

Sign/Date:

Verified By

(Quality Assurance)

Sign/Date:

Inference:

.....

.....

.....

.....

.....

**Reviewed By
(Manager QA)**

Sign/Date:



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

9.0 REFERENCES:

The Principle Reference is the following:

- Validation Master Plan
- Schedule-M – “Good Manufacturing Practices and Requirements of Premises, Plant and Equipment for Pharmaceutical Products.”
- WHO Essential Drugs and Medicines Policy, QA of Pharmaceuticals, Vol-2 – Good Manufacturing Practices and Inspection.

10.0 DOCUMENTS TO BE ATTACHED:

- Technical details for Equipment Requirement with Engineering Drawings.
- Certificate of MOC.
- Operation and Maintenance Manual.



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

11.0 DEVIATION FROM PRE-DEFINED SPECIFICATION, IF ANY:

.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....

12.0 CHANGE CONTROL, IF ANY:

.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....

13.0 REVIEW (INCLUSIVE OF FOLLOW UP ACTION, IF ANY):

.....
.....
.....
.....
.....
.....
.....
.....
.....
.....



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

14.0 CONCLUSION:

.....

.....

.....

.....

.....

.....

.....

.....

15.0 RECOMMENDATION:

.....

.....

.....

.....

.....

.....

.....

.....



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

16.0 ABBREVIATIONS:

| | | |
|------|---|--------------------------------------|
| AC | : | Alternative Current |
| cGMP | : | current Good Manufacturing Practices |
| DQ | : | Design Qualification |
| Hz | : | Hertz |
| IQ | : | Installation Qualification |
| Ltd | : | Limited |
| mm | : | Millimetre |
| MOC | : | Material of Construction |
| No. | : | Number |
| QA | : | Quality Assurance |
| SOP | : | Standard Operating Procedure |
| SRR | : | Stirrer |
| Std | : | Standard |
| V | : | Volt |
| WHO | : | World Health Organization |



PHARMA DEVILS

**INSTALLATION QUALIFICATION
PROTOCOL CUM REPORT
FOR
STIRRER**

PROTOCOL No.:

17.0 POST- APPROVAL:

INITIATED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|--|------|-----------|------|
| OFFICER / EXECUTIVE (QUALITY ASSURANCE) | | | |

REVIEWED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|-------------------------------|------|-----------|------|
| HEAD (PRODUCTION) | | | |
| HEAD (ENGINEERING) | | | |

APPROVED BY:

| DESIGNATION | NAME | SIGNATURE | DATE |
|-------------------------------------|------|-----------|------|
| HEAD (QUALITY ASSURANCE) | | | |