



PHARMA DEVILS
QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

OPERATIONAL QUALIFICATION
OF
PICO CONDUCTIVITY METER
(LABINDIA)



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

TABLE OF CONTENTS

1.0 Pre-Approval

2.0 Overview

- 2.1 Purpose
- 2.2 Scope
- 2.3 Responsibility
- 2.4 Requalification
- 2.5 Instrument Identification

3.0 Operational Qualification Procedure

- 3.1 Training
- 3.2 Key Functionality and Safety Features
- 3.3 SOP verification
- 3.4 Deficiency (if any) and Corrective Action Report

4.0 Acceptance Criteria

5.0 Summary

- 5.1 Conclusion
- 5.2 Post –Approval

6.0 Appendix

- 6.1 Abbreviations and Definitions



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

1.0 Pre-Approval:

Signing of this Approval page of Operational Qualification Protocol No. indicates agreement with the Operational Qualification approach described in this document. Should Modifications to the Operational Qualification become necessary, an addendum will be prepared and approved.

Written By	Signature	Date
Officer/Executive- Quality Control		

Checked By	Signature	Date
Executive/Manager – Quality Control		
Executive/Manager – Quality Assurance		

Approved By	Signature	Date
Manager - Quality Assurance		
General Manager - Works		



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

2.0 Overview:

2.1 Purpose:

The purpose of this protocol is:

- To verify the operational attributes of Pico Conductivity meter critical to serve the intended purpose.
- To establish the suitability of the draft SOP prepared for the operation of System.
- To document the observations for future reference.

2.2 Scope:

This protocol covers the operational qualification of Pico Conductivity meter.

2.3 Responsibility:

The validation group comprising of representatives from each of the following departments shall be responsible for the overall compliance with this protocol :

- ◆ Quality Control Department
- ◆ Quality Assurance Department
- ◆ Engineering Department

The Quality Control and Engineering shall be responsible for checking the operations and recording data as per the procedures outlined in this protocol.

Quality Control shall collect all the test data and shall compile the results to make the reports of qualification studies.

The Reports shall be checked by Production and Quality Assurance.

The post approval of the qualification shall be done by the Quality Assurance and plant Head.



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

2.4 Requalification:

Operational Qualification to be repeated incase of

- ◆ Replacement of any major component -
- ◆ Major modification in the existing instrument -
- ◆ During monitoring if instrument. is found to be malfunctioning -
- ◆ Shifting of the instrument from one location to another -

2.5 Instrument Identification

The Instrument is identified as: Pico Conductivity meter

Serial No. :

In-house Instrument No. :

Name of the Supplier : LABINDIA INSTRUMENT PVT LTD.

Purchase Order No. :



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

3.0 Operational Qualification Procedure

- 1) A draft SOP shall be prepared on the basis of manufacturer guide / instrument manual for operation before the Qualification testing.
- 2) Prior to the Qualification test, the Personnel shall be trained by the Engineer from the Manufacturer / supplier on the operational features of the instrument. This training shall be recorded in Section 3.1.
- 3) The trained personnel shall carry out the Operational Qualification along with the Service Engineer, following the Procedures mentioned under Section 3.2.1 for Key Functionality and Safety Features. Record the observations of Qualification Test in Test Data Sheet of Section 3.2.1. Checkpoints designed for the purpose of OQ are also aimed at verification of these draft SOP's.
- 4) Operate the instrument as per the draft SOP. Record the change if any and confirm the SOP. Report the confirmation of SOP in the Section 3.3.
- 5) Report the deficiency from the specified function, if any in the section 3.4



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

3.1 Training

Date:

Title: Pico Conductivity meter

Name of the Trainer(s): _____

S.No.	Name of the Trainee	Employee Number	Signature
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			

Signature of Trainer(s): _____

Date : _____



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

3.2. Key Functionality & Safety Features:

A. Purpose:

The purpose of this procedure is to demonstrate that the control panel and other manual operations (if any) of Conductivity meter (**Equipment No.**) function as specified by the manufacturer.

B. Testing:

1. Check all the displays on the panel are identified.
2. Turn on the power from the electrical panel.
3. Set the control(s) on the panel.
4. Verify functionality of each component on the panel against its Specified functions.
5. Observe and record the responses in the Test Data Sheet, under section 3.2.1.



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

3.2.1 TEST DATA SHEET (Confirmation of Services Connection):

S.No.	Test Particulars	Specified Function	Observations	Checked By
1.	Power on Check			
	Check the LCD display Status after Power on	Initialization Beep		
	Press the RESET Key and check the System Resets and initialises			
	Press active keys on the front panel & check for the buzzer beeps			
2.	ELECTRODE INPUT CHECK			
	To check the Electrode Input with Fixed Resistors source and Temperature sensor			
	<ul style="list-style-type: none">• Select Conductivity Measure function, Press TEMPCO Key, Enter TEMPCO =0• Temp. ATC = Manual @ 25 C , Press MEM CLEAR – To clear the CAL Memory.• Connect Fixed Resistor to Banana connectors of Electrode amplifier, and RTD Temp. Sensor to 5-pin shell connector.• Select the Fixed Resistors to input the signal and record/note the displayed value.			



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

Fixed Value Resistors @ Input		Displayed Value	Deviation Limit	Remarks Passed/Fail
Ohms	Siemens	Siemens	+/- 0.5 %	
	mS			
	mS			
	µS			

Verified By:

Name: _____ Signature: _____ Date: _____

3. Conductivity Electrode Calibration/Cell Constant Check

To check the conductivity Electrode Calibration using buffer 1.413 mS/cm, 147.0 µS/cm and distilled water having conductivity lesser than 2.0 µS/cm.

- Connect Conductivity Electrode (K=1) to Banana connector of Electrode amplifier, and RTD Temp. sensor to 5- pin shell connector.
- Press TEMCO key. Enter Temp. coefficient = 1.92% for KCL, ATC = AUTO @ 25 C,
- Press CCAL function key. Enter the value of cell constant and conductivity standards.
- Rinse electrode several times in distilled water between measurement and soak dry on tissue paper. (Do not rub the Electrode Pt-plates.)
- Dip it in fresh conductivity buffer as per the CAL procedure. Record/note the displayed value. After the calibration is over, Printer prints Calibration report.

Calibration Buffer Used mS.	CAL. Cell Constant. (Limit 80-120%) (From Calibration report)	Displayed Conductivity Value after Calibration	Remarks Passed/Failed
STD1=	K (0.1)=		
STD2=	K (1.0)=		



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

D) TEMPERATURE SENSOR (PT100):

Bath SET Temperature C	Reading on Thermometer C	Display Reading C	Deviation +- 0.5 C	REMARKS Passed/Failed

D) MEASUREMENT TEST:

- Check for the Printer Connection & make sure the Printer is switched ON.
- Load the paper in the printer and leave the printer in online – READY condition.
- Take the conductivity Standard which has lower value than the Calibration Standard used for taking reading and Storing.

REQUIRMENT CONDITIONS & OBSERVATION:

- Select the conductivity measure mode, Rinse the electrode with the distilled water and dip it in the Conductivity Standard as mentioned above.
- Display reads the Conductivity and cel constant . K value with reading no. Temperature and * for Standard (non GLP) measurement.
- Wait till the reading stabilities.
- Press ENTER Key to Store the reading.
- Now Press the Key with GLP written on it – the * indication from the display disappears (as GLP measurement)
- Press Enter again to Store the reading, and display reads and waits for parameter entries – Sample Name, I.D. No. of the Std. Sample and the reading is stored.
- Press ENTER Key to Store the new readings. Take Printout of the stored readings.

E) Measurement test

- Select Probe Check mode : Use Known Conductivity standard to estimate the electrodes Cell Constaant K (80 – 120 % limit)
- Display reads the cell Const. K, with message “ Incorrect Buffer or probe Faulty “
- Select Calibration Reminder Alarm mode, Check the alarm activation by setting the new time and date.
- Check for DATE & TIME

ERROR INDICATION CHECK

- Check for invalid entry by entering out of the range



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

S.No	Test Particulars	Specified Function	Observations	Checked By	
2	Timer mode setting : There are two mode of operation	Timer Mode Manual Mode			
3	Press the <table border="1" data-bbox="386 562 509 659"><tr><td>Timer</td></tr></table> key to change from Timer mode to Manual mode or from Manual mode to Timer mode	Timer	Initially when power supply is switched 'ON' the instrument is in Timer mode. If the Timer LED is ON then timer mode is selected, if Timer LED is OFF and Timer display shows '----' then MANUAL mode is selected.		
Timer					

Verified By:

Name: _____ Signature: _____ Date: _____

S. No	Test Particulars	Specified Function	Observations	Checked By
-------	------------------	--------------------	--------------	------------



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

S. No	Test Particulars	Specified Function	Observations	Checked By
4	<p>Timer Setting Range Setting</p> <p>Hours and Minutes</p> <p>b. Press <input type="button" value="▼"/> key</p> <p>followed by</p> <p><input type="button" value="Timer 1"/> / <input type="button" value="Timer 2"/> key</p> <p>c. Press <input type="button" value="▲"/> key</p> <p>followed by</p> <p><input type="button" value="Timer 1"/> / <input type="button" value="Timer 2"/> key</p>	<p>If this range is selected then the display will show '----' before time</p> <p>Time range in minutes and seconds should be selected</p> <p>Time range in hours and minutes should be selected</p>		

Verified By:

Name: _____ Signature: _____ Date: _____






S.No	Test Particulars	Specified Function	Observations	Checked By
------	------------------	--------------------	--------------	------------



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

S.No	Test Particulars	Specified Function	Observations	Checked By
5	<p>Time Setting</p> <p>a) Press the  key followed by timer key</p> <p>b) Set the Time by pressing  or  key</p> <p>c) Press  key</p> <p>d) Press </p>	<p>The display shows previously set time with extreme right digit flashing</p> <p>Time will get increased or decreased with every key pressed</p> <p>The position of the number to be changed will be shifted</p> <p>The value set will get registered</p>		

Verified By:










Name: _____ Signature: _____ Date: _____



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

S.No	Test Particulars	Specified Function	Observations	Checked By
6	<p>Temperature setting</p> <p>a) Press PROBE key </p> <p>b) Press  SET key followed by  key</p> <p>c) Set the temperature by pressing  or  key</p> <p>d) Press  SET key</p> <p>e) Press </p>	<p>Bath probe will get set</p> <p>Temperature display shows previous set temperature and the right digit on the display starts flashy</p> <p>Temperature will get increased or decreased with every key pressed</p> <p>The position of the number to be changed will be shifted</p> <p>The value set will be registered</p>		
7	<p>Confirm the set value of temperature by pressing  SET key followed By  key</p>	<p>Display should show the set value</p>		

Verified By:





Name: _____ Signature: _____ Date: _____



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

S. No	Test Particulars	Specified Function	Observations	Checked By
8	After pressing the  key start setting within 10 sec.	The mode gets terminated if starting setting exceeds 10 secs.		
9	Check the temperature of the bath and water in the Beaker is 37.9° and $37.0 \pm 2^\circ\text{C}$ respectively. By shifting the probe to required location.	The temperature display will show the temperature respective to the location of the probe after pressing  key		
10	TEMP  Press	Temperature controller will start		
B	OPERATION			
1	Put one tablet each in individual tube of the basket assembly (total six) followed by rises and assemble the basket on the extended arms. Press START  STOP	The instrument will start operating		

Verified By:

Name: _____ Signature: _____ Date: _____



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

S.No	Test Particulars	Specified Function	Observations	Checked By
2	Press the <input type="text"/> TIMER 1 and / or <input type="text"/> TIMER 2 <input type="text"/>	The up and down movement of the basket will start and continue till the set time		
3	Press <input type="text"/> TIMER 1 or <input type="text"/> TIMER 2 key before completion of the set time	The operation will halt		
4	Check for complete passing of the granules through the mesh at the bottom of the basket by lifting the basket out of water and note the time displayed on the board Press <input type="text"/> START key <input type="text"/> STOP .	The apparatus will stop its operation		
C	BOTTOM ILLUMINATION CHECK Put ON the switch provided at right side of the bottom of the apparatus	The green light shall be observed in the bath		

Verified By:

Name: _____ Signature: _____ Date: _____



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

3.3 SOP Verification:

Draft SOP No. :

Title : Disintegration Tester (USP) Model ED – 2AL

Operate the instrument. as per the draft SOP and record the details given below:

Operated By :

Checked By :

The operating personnel understand and follow the SOP description (Yes/No) : YES

Changes required in draft SOP (If any) : _____

NO _____

SOP to be revised (Yes/No) : NO

If yes, Review No. _____

Remarks: SOP Confirmed / Not Confirmed

Verified By:

Name : _____ Signature : _____ Date : _____



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

3.3 Deficiency (if any) and Corrective Action Report:

If there is no deficiency, then write N. A

Description of deficiency and date observed:

Person, responsible for corrective action and date assigned:

Corrective actions taken and date conducted:

Conducted By : _____ Approved By : _____

Date: _____ Date: _____

Comments (if any):

Verified By:

Name: _____ Signature: _____ Date: _____



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

4.0 Acceptance Criteria:

Operational Qualification shall be considered acceptable when all the conditions specified in various data sheets under section 3.0 have been met.

Any deviation from the acceptance criteria of the specific check point shall be reported and decision should be taken for the rejection, replacement or rectification of the instrument / component.



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

5.0 Summary:

Checks	Observations Yes / No	Remarks (if any)
Whether the acceptance criteria of the protocol and specific checkpoints are met.		

5.1 Conclusion:

The Disintegration Tester (USP) Model – **is / is not** qualifying the Operational Qualification tests as per Protocol No. The Instrument **can / cannot** be tested for its Performance Qualification as per Protocol No.



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

5.2 Post-Approval:

Name	Signature	Date
Manager - Engineering		
General Manager – Works		
Manager – Quality Assurance		

6.0 Appendix:

6.1 Abbreviations and Definitions

OQ	- Operation Qualification
mm	- Millimeter
Min	- Minutes
V	- Volt
Hz	- Hertz
cm	- Centimeter
N.A.	- Not Applicable
S. No.	- Serial Number
Sr.	- Senior
mV	- milli Volt
°C	- Degree Centigrade
AC	- Alternate Current
DC	- Direct Current
g	- Gram
RH	- Relative Humidity
USP	- United States Pharmacopoeia
SOP	- Standard Operating Procedure



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

OPERATIONAL QUALIFICATION FOR PICO CONDUCTIVITY METER

Acceptance criteria	The product, instrument., and / or process specifications and limits, such as acceptable quality level and unacceptable quality level, that are necessary for making a decision to accept or reject.
Operational qualification	The documented verification that all aspects of a facility, utility, or equipment that can affect product quality operate as intended throughout all anticipated ranges.
Validation	Establishing documented evidence that a system does what it purports to do.
Revalidation	Repetition of the validation process or a specific portion of it