



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

1.0 OBJECTIVE:

To lay down procedure for operation, cleaning and calibration of UV Spectrophotometer .

2.0 SCOPE:

This SOP is applicable for operation, cleaning and calibration of UV Spectrophotometer (Make: Perkin elmer, Model : LAMBDA 25).

3.0 RESPONSIBILITY – Execution- Executive QC
Checking -Assistant Manager QC

4.0 ACCOUNTABILITY - Manager Quality Control

5.0 PROCEDURE:

5.1 OPERATING PROCEDURE:

- 5.1.1 Switch on the instrument. Switch on the computer. Enter the login name and password to run the windows.
- 5.1.2 Click on "Perkin elmer UV win lab".
- 5.1.3 Enter the login name and password.
- 5.1.4 Go to METHOD, select method. Then go to task,data collection,Instrument.
Enter wavelengths, scan speed, data interval, ordinate mode, Slit width,lamp change as per the requirement.
- 5.1.5 Select SAMPLE INFORMATION,enter sample ID description.
- 5.1.6 Go to processing.
- 5.1.7 Go to ACTION press START. Make it auto zero to blank put blank solution click "OK" then put the sample solution click "OK".Instrument will start automatically.
- 5.1.8 After completing the process, go to reporting Press PRINT to print the result.
- 5.1.9 In case of any problem / clarification contact LAB INCHARGE.

5.2 CALIBRATION PROCEDURE

(Calibration Frequency: After every 6 months)

5.2.1 CONTROL OF ABSORBANCE

5.2.1.1 PREPARE THE FOLLOWING REAGENTS:

- 5.2.1.1.1 0.005M Sulphuric Acid : Add carefully 7.0 ml of conc. sulphuric acid in 100 ml of cold water and make up the volume to 250 ml. Further dilute 10 ml of this solution to 1000 ml with DI water in a 1 lt.volumetric flask.



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

- 5.2.1.1.2 Potassium dichromate Solution : Dissolve 57 to 63 mg of Potassium dichromate AR previously dried at 130 deg.C to constant weight, in 100 ml of 0.005 M sulphuric acid and dissolve by keeping in an ultrasonic bath. Make up the volume to 1000 ml with 0.005 M sulphuric acid.
- 5.2.1.1.3 Clean the cell chamber of the instrument. Soak the cuvettes in 2% v/v solution of detergent in water and wash with water.
- 5.2.1.1.4 Switch on the power supply of the instrument
- 5.2.1.1.5 Switch on the instrument. Switch on the computer. Enter the login name and password to run the windows.
- 5.2.1.1.6 Click on Perkin elmer UV win lab.
- 5.2.1.1.7 Enter the login name and password.
- 5.2.1.1.8 Go to METHOD select method. Then go to task,data collection,instrument.
Enter START and END wavelengths, lambda max. first and then lambda min. scan speed, data interval, ordinate mode, Slit width, lamp change.

5.2.1.2 INSTRUMENT SETUP

5.2.1.2.1

METHOD TYPE	SCAN
START WL	400 nm
STOP WL	200 nm
ORDINATE MODE	A
SLIT WIDTH	1.00 nm
SCAN SPEED	120 nm/min
DATA INTERVAL	1 nm
NUMBER OF CYCLE	1
CYCLE TIME	1 SEC.
UV LAMP	ON
VISIBLE LAMP	ON
LAMP CHANGE	326 nm



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

5.2.1.2.2 Then go to SAMPLE INFORMATION, enter sample ID description.

5.2.1.2.3 Go to processing, add peak table.

5.2.1.3 PROCEDURE

5.2.1.3.1 Fill the quartz cells with 0.005 M sulphuric acid and start the baseline scan by go to ACTION then START. After the completion of the scanning the instrument will display 0.000 against the absorbance.

5.2.1.3.2 Fill the sample quartz cell with potassium dichromate solution after thoroughly rinsing it with potassium dichromate solution. Press 'OK'.

5.2.3.3 Determine absorption at 235, 257, 313 and 350 nm. The permitted tolerance for Ultraviolet range is +/- 1 nm.

5.2.1.4 CALCULATIONS

5.2.1.4.1 Calculate the A(1%, 1 cm) as given below:

Absorbance

A (1%, 1 cm) = -----

Weight per 100 ml.

5.2.1.4.2 Check the absorbance using solution of potassium dichromate at the wavelengths indicated in the below table.

5.2.1.4.3 The values should conform to the values given in the following table.

WAVELENGTH(nm)	SPECIFIC ABSORBANCE	MAXIMUM TOLERANCE
A (1%, 1cm)		
235	124.5	122.9 to 126.2
257	144.5	142.8 to 146.2
313	48.6	47.0 to 50.3
350	107.3	105.6 to 109.0

5.2.2 CONTROL OF WAVELENGTHS

5.2.2.1 PREPARATION OF REAGENTS

5.2.2.1.2 Holmium Perchlorate solution:

Prepare 4% w/v solution of Holmium oxide in 1.4 M Perchloric acid.



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

5.2.2.2 INSTRUMENT SETUP

5.2.2.2.1 Set the following values against the menu displayed as done above :

METHOD TYPE	SCAN
START WL	600 nm
STOP WL	200 nm
ORDINATE MODE	A
SLIT WIDTH	1.00 nm
SCAN SPEED	120 nm/min
DATA INTERVAL	1 nm
NUMBER OF CYCLE	1
CYCLE TIME	1 SEC.
UV LAMP	ON
VISIBLE LAMP	ON
LAMP CHANGE	326 nm

5.2.2.3 PROCEDURE

5.2.2.3.1 Fill the quartz cells with 1.4 M Perchloric acid and start the baseline scan by go to ACTION then START. After the completion of the scanning the instrument will display 0.000 against the absorbance.

5.2.2.3.2 Fill the sample quartz cell with Holmium perchlorate solution after thoroughly rinsing it with Holmium perchlorate solution. press "OK".

5.2.2.3.3 Determine maxima.

5.2.2.3.4 The maxima shall be at the following wavelengths with permitted tolerance of ± 1 nm for ultraviolet range and ± 3 nm for visible range.

WAVELENGTH	TOLERANCE
241.15 nm	± 1 nm
287.15 nm	± 1 nm



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

361.5 nm ± 1 nm

536.3 nm ± 3 nm

5.2.2.3.5 Go to reporting,select templet save method and task ,print the result.

5.2.3. LIMIT OF STRAY LIGHT

5.2.3.1 PREPARATION OF REAGENTS

5.2.3.1.1 Potassium Chloride solution:

Prepare 12 gm/l solution of Potassium chloride in water.

5.2.3.2 INSTRUMENT SETUP

METHOD TYPE	WAVELENGTH PROGRAM
WAVELENGTH	200 nm
ORDINATE MODE	A
SLIT WIDTH	1.00 nm
RESPONSE	0 S
NUMBER OF CYCLE	1
CYCLE TIME	1 SEC.
UV LAMP	ON
VISIBLE LAMP	ON
LAMP CHANGE	326 nm

5.2.3.3 PROCEDURE

5.2.3.3.1 The absorbance of above Potassium Chloride solution in water. Check the absorbance of Potassium chloride solution at 200 nm in 1 cm cell using water as blank.

5.2.4 RESOLUTION

5.2.4.1 Record the spectrum of a 0.2 ml of Toluene AR in 1000 ml of Hexane AR from 260 nm to 275 nm.

5.2.4.2 Calculate ratio of minimum peak height to maximum peak height using Hexane as blank.

5.2.4.3 The minimum ratio of the absorbance at the maximum at 269nm to that at the minimum at 266nm shall be more than 1.5

5.2.4.4 INSTRUMENT SETUP



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

METHOD TYPE	SCAN
START WL	275 nm
STOP WL	260 nm
ORDINATE MODE	A
SCAN SPEED	120 nm/min
DATA INTERVAL	1 nm
SLIT WIDTH	1.00 nm
NUMBER OF CYCLE	1
CYCLE TIME	1 Sec.
UV LAMP	ON
VISIBLE LAMP	ON
LAMP CHANGE	326 nm

5.2.5. RESOLUTION POWER

- 5.2.5.1 Record the second derivative spectrum of 0.2 gm solution of Toluene AR in 1000 ml of methanol AR using methanol as the compensation liquid.
- 5.2.5.2 The spectrum shows the negative extreme located between two large negative extreme at 261 nm and 268 nm respectively.
- 5.2.5.3 The ratio of absorbance difference at 263 and 265 nm to the absorbance difference at 263 and 261 nm shall not be less than 0.2.

5.2.5.4 INSTRUMENT SETUP

METHOD TYPE	SCAN
START WL	290 nm
STOP WL	220 nm
ORDINATE MODE	A
SCAN SPEED	30 nm/min
DATA INTERVAL	0.1 nm
SLIT WIDTH	0.5 nm
NUMBER OF CYCLE	1



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

CYCLE TIME 1 Sec.
UV LAMP ON
VISIBLE LAMP ON
LAMP CHANGE 326 nm

5.3 CLEANING PROCEDURE

Frequency: Daily or after each use.

- 5.3.1 Open the front part of the instrument .
- 5.3.2 Wipe out any material in the quartz cell holder assembly by means of tissue paper.
- 5.3.3 Clean all the sampling accessories with tissue paper after analysis and keep them in proper place.
- 5.3.4 Clean the outer surface of the instrument with dry cotton cloth.
- 5.3.5 Record the details of cleaning in log card.

6.0 SAFETY & PRECAUTIONS:

Not Applicable

7.0 REVISION HISTORY:

Revision No.	Reason for Revision	Superseded from & Date

8.0 DISTRIBUTION:

Copy No.	Issuance Record				Withdrawal Record		Destruction Record	
	Date	Dept. issued	Name / Signature of receiver	Issued By Name / Signature	By	Sign/ Date	By	Sign/ Date

9.0 REFERENCES:



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

Not Applicable

10.0 ABBREVIATIONS & ANNEXURES:

SOP : Standard Operating Procedure

No. : Number

QC : Quality Control

Annexure I : Calibration Record of UV Spectrophotometer



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

ANNEXURE-I CALIBRATION RECORD OF UV SPECTROPHOTOMETER

PERFORMANCE DATE	DATE OF LAST PERFORMANCE DONE	NEXT DUE FOR PERFORMANCE

INSTRUMENT DETAILS		
INSTRUMENT NAME	INSTRUMENT MAKE	INSTRUMENT IDENTIFICATION NO.

1. Control of Absorbance:

Calculated Specific Absorbance A (1%, 1cm) = $\frac{\text{Absorbance} \times 1000}{W \text{ in gms}}$

W= weight of K₂Cr₂O₇ in mg = _____

S.No.	Wavelength (nm)	Observed Absorbance	Calculated Specific Absorbance
1	235		
2	257		
3	313		
4	350		

S.No.	Wavelength (nm)	Calculated Specific Absorbance A (1%, 1cm)	Maximum Tolerance
1	235		122.9 – 126.2
2	257		142.8 – 146.2
3	313		47.0 – 50.3
4	350		105.6 – 109.0

Remarks: Satisfactory/Not Satisfactory

2. Control of Wavelength:

Wavelength (nm)	Measured maxima at (nm)	Tolerance
241.15		± 1 nm
287.15		± 1 nm
361.5		± 1 nm
536.3		± 3 nm

Remarks: Satisfactory / Not Satisfactory



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE

Department: Quality Control	SOP No.:
Title: Operation, cleaning and calibration of UV Spectrophotometer	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

3. Limit of stray light:

Absorbance of KCl (12 gm/lit) solution in water (at 200nm) _____

Remarks: Satisfactory/Not Satisfactory

4. Resolution:

Ratio of absorbance at 269nm =
Absorbance at 266nm

Wavelength (nm)	Measured Absorbance
269	
266	

Remarks: Satisfactory/Not Satisfactory

5. Resolution Power:

Wavelength scan (290 – 220 nm)
Absorbance difference at 263nm -265nm (A)

Absorbance difference at 263nm -261nm (B)

Wavelength (nm)	Measured Absorbance
261	
263	
265	

= _____

= _____

Ratio of (A/B) = _____

Remarks : Satisfactory/Not Satisfactory

Calibration done by :

Date :

Checked by :

Date :

CONCLUSION

INSTRUMENT WORKING SATISFACTORY

INSTRUMENT NOT WORKING SATISFACTORY

PERFORMED BY			CHECKED BY		
NAME	SIGN	DATE	NAME	SIGN	DATE