



PHARMA DEVILS
QUALITY CONTROL DEPARTMENT

**ANALYTICAL METHOD VALIDATION OF RESIDUAL SOLVENT OF LEVOCETIRIZINE &
MONTELUKAST TABLETS**

**ANALYTICAL METHOD VALIDATION
PROTOCOL CUM REPORT
FOR
RESIDUAL SOLVENT OF LEVOCETIRIZINE &
MONTELUKAST TABLETS**

Name of the Product	:	LEVOCETIRIZINE & MONTELUKAST TABLETS
Department	:	QUALITY CONTROL
Protocol No.	:	
Effective Date	:	



PHARMA DEVILS
QUALITY CONTROL DEPARTMENT

ANALYTICAL METHOD VALIDATION OF RESIDUAL SOLVENT OF LEVOCETIRIZINE & MONTELUKAST TABLETS

1. PRE APPROVAL:

Document is Prepared, Reviewed and Approved by the following persons for Pre Approvals.

Prepared By :

Department	Name	Designation	Signature/ Date
Quality Control			

Reviewed By:

Department	Name	Designation	Signature/ Date
Quality Control			
Quality Assurance			

Approved By:

Department	Name	Designation	Signature/ Date
Quality Control			
Quality Assurance			



PHARMA DEVILS
QUALITY CONTROL DEPARTMENT

**ANALYTICAL METHOD VALIDATION OF RESIDUAL SOLVENT OF LEVOCETIRIZINE &
MONTELUKAST TABLETS**

INDEX

LIST OF CONTENT	PAGE NO.
Approval	02
Objective	04
Scope	04
Responsibility	04
References	04
Methodology	04-06
Parameters (Specificity and Precision)	07



ANALYTICAL METHOD VALIDATION OF RESIDUAL SOLVENT OF LEVOCETIRIZINE & MONTELUKAST TABLETS

- 2. OBJECTIVE:** The objective of this protocol cum report is to establish a documented evidence for the Residual Solvent Method Validation of Levocetirizine & Montelukast Tablets.
- 3. SCOPE:** This Protocol cum report is applicable for Residual Solvent test of Levocetirizine & Montelukast Tablets for Specificity and Precision, domestic market in Quality Control Department.
- 4. RESPONSIBILITY:** Officer/Executive-QC is responsible for the preparation of the protocol. The protocol shall be reviewed / approved by the Officer/Executive/Manager –QC/QA, approved by Head QA/QC.
- 5. REFERENCES:**
 - SOP: Operation and calibration of GC.
 - SOP: Procedure Validation Activities.
 - SOP: Operation and calibration of analytical balance.
 - SOP: Procedure for handling of reference/working standard.
 - SOP: Handling procedure of Gas chromatography column

6. METHODOLOGY:

RESIDUAL SOLVENT (BY GC):

Equipment /Instrument Required:

- Gas Chromatography
- Analytical Balance
- Calculator

Glass wares required:

- Volumetric Flask
- Beaker

Chromatographic conditions:

Column : A 0.53 mm x 30 m, 3.0µm fused-silica analytical column coated with 6% Cyanopropyl-phenyl 94% dimethylpolysiloxane (DB-624)

Detector : Flame-ionization detector

Carrier gas : Nitrogen



PHARMA DEVILS

QUALITY CONTROL DEPARTMENT

ANALYTICAL METHOD VALIDATION OF RESIDUAL SOLVENT OF LEVOCETIRIZINE & MONTELUKAST TABLETS

Flow rate : A Flow of about 2ml/minute
Make-up gas : Nitrogen, if recommended
Injection port temperature: 210°C
Detector temperature : 260°C
Column temperature : Programmed according to the following steps.

It is maintained at 40° C for 10 minutes, and then increased to 230°C at a rate of 30°C per minute, and hold for 5 minutes.

Head space parameter:

Vial incubation temperature: 90°C
Vial incubation time : 20 minutes
Transfer line temperature: 110°C
Syringe temperature : 100°C
Carrier gas : Nitrogen
Pressurization time : 2.0 minutes
Injection volume : 0.1 ml
Withdrawal time : 0.5 minutes
Carrier pressure : 15 psi
Split ratio : 2:1

Standard preparation: Prepare a solution, in DMSO. Dissolve 60 mg of Dichloromethane and 500 mg IPA and dilute to 100 ml with DMSO. Further dilute 10 ml of this solution to 100 ml with DMSO. [Note: Prepare fresh standard for daily use]

Test preparation: Take 20 tablets and crush to fine powder. Weigh 1.0 g of sample in 20 ml Headspace vial, add 5 ml of DMSO.

Procedure: Take 5 ml of standard solution in 6 different 20 ml headspace vial, Separately inject equal volumes of Blank (DMSO), Standard solution (Six injections) and sample solution (Duplicate) into the system, record the chromatogram and measure the peak responses. The relative standard deviation of the Dichloromethane (MDC) and Isopropyl alcohol (IPA) peak response from replicate injections is not more than 15%. Order of elution is MDC and IPA.



ANALYTICAL METHOD VALIDATION OF RESIDUAL SOLVENT OF LEVOCETIRIZINE & MONTELUKAST TABLETS

System suitability: The resolution between each peak more than 1.5 and percentage relative standard deviation for six replicate of standard is not more than 15%.

Calculation:

Dichloromethane in µg per g:

$$\text{Content in ppm} = \frac{A_t \cdot W_s \cdot 10 \cdot 5 \cdot P}{A_s \cdot 100 \cdot 100 \cdot W_t \cdot 100} \times 1000000$$

Where,

- A_t : Average area of the peak due to Dichloromethane in test preparation.
- A_s : Average area of the peak due to Dichloromethane in standard preparation.
- W_s : Weight of Dichloromethane working standard taken, in mg.
- W_t : Weight of sample taken, in mg.
- P : Potency of Dichloromethane working standard.

Isopropyl Alcohol in µg per g:

$$\text{Content in ppm} = \frac{A_t \cdot W_s \cdot 10 \cdot 5 \cdot P}{A_s \cdot 100 \cdot 100 \cdot W_t \cdot 100} \times 1000000$$

Where,

- A_t : Average area of the peak due to Isopropyl Alcohol in test preparation.
- A_s : Average area of the peak due to Isopropyl Alcohol in standard preparation.
- W_s : Weight of Isopropyl Alcohol working standard taken, in mg.
- W_t : Weight of sample taken, in mg.
- P : Potency of Isopropyl Alcohol working standard.

6.1 SPECIFICITY:

The specificity of an analytical method is its ability to measure unequivocally the analyte in the presence of components that may be expected to be present in the test matrix.

Note: For detailed of procedure, refer methodology



PHARMA DEVILS
QUALITY CONTROL DEPARTMENT

ANALYTICAL METHOD VALIDATION OF RESIDUAL SOLVENT OF LEVOCETIRIZINE & MONTELUKAST TABLETS

S.No.	Validation parameters	Acceptance Criteria	Observation	
1	Specificity	There should be no interference observed from Blank, mix Standard and sample to the Dichloromethane, Isopropyl Alcohol	Name of Analyte	RT
			Dichloromethane	--
			Isopropyl Alcohol	5.098
			No Interference observed from Blank, Mixed Standard and Sample	

6.2 PRECISION:

2	Precision	System Precision	The %relative standard deviation (% RSD) of the peak areas for Dichloromethane, Isopropyl Alcohol must not be more than 15.0%.	Name of Analyte	
				Dichloromethane	
				Isopropyl Alcohol	
		Method Precision	The Result of Test must comply with the Specification.	Name of Analyte	Area of Test
				Dichloromethane	Not detected
				Isopropyl Alcohol	22942602