## PHARMA DEVILS



PRODUCTION DEPARTMENT

#### STANDARD OPERATING PROCEDURE

Department: Production	SOP No.:
Title: Inventory and Storage of Sieves and Screens	Effective Date:
Supersedes: Nil	<b>Review Date:</b>
Issue Date:	Page No.:

#### **1.0 OBJECTIVE:**

To lay down a procedure for Inventory and Storage of Sieves, Screens.

## 2.0 SCOPE:

This procedure is applicable to Inventory Management and Storage of Sieves, Screens at Production Department.

## **3.0 RESPONSIBILITY:**

Technical Associate - for Execution Officer/ Executive Production Department- for verification and implementation of SOP Head Production Department- shall ensure compliance of the SOP

### 4.0 **DEFINITION(S):**

NA

#### 5.0 **PROCEDURE:**

#### 5.1 Inventory

- 5.1.1 On intimation of receipt from stores, collect all the sieves / screens.
- 5.1.2 Check the integrity of individual sieves / screens through vendor certificate and record the observations in Annexure-I.
- 5.1.3 Remark 'OK' for pass sieves / screens, and 'Not OK' for sieves / screens which fails the integrity test. Return back the defective sieves / screens to stores.
- 5.1.4 Allot a code no. as per the respective departmental SOP, to the pass sieves / screens and transfer the respective storage area.
- 5.1.5 Record the observations in Annexure-I for sieves and Annexure-II for screens.
- 5.1.6 Clean the sieves / screens as per SOP.

#### 5.2 Storage

- 5.2.1 Cover the clean and dry sieve / screen in polybag and put cleaned label having previous product name, cleaned by and checked by signature and date as per SOP.
- 5.2.2 Keep the covered sieves on the backside gap of cleaned and dried rectangular trolley or Screens can be kept in the rack.

# PHARMA DEVILS



PRODUCTION DEPARTMENT

#### STANDARD OPERATING PROCEDURE

Department: Production	SOP No.:
Title: Inventory and Storage of Sieves and Screens	<b>Effective Date:</b>
Supersedes: Nil	<b>Review Date:</b>
Issue Date:	Page No.:

5.2.3 After completion of activity clean and place sieve / screen in polybag, label and store as per procedure mentioned in step no. 5.2.1 to 5.2.2.

#### 6.0 ABBREVIATION(S):

SOP: Standard operating procedure.

#### 7.0 **RERERENCE**(S):

SOP: Cleaning of new equipments / machineries / in-process product containers / accessories SOP: Status labeling

#### 8.0 ANNEXURE(S):

Annexure-I	:	Sieves Inventory record
Annexure-II	:	Screens Inventory record

#### 9.0 **DISTRIBUTION:**

- 9.1 Master Copy : Quality Assurance
- 9.2 **Controlled Copy (s) :** Production department, Quality Assurance
- 9.3 **Reference Copy (s)** : Production Department

PHARMA DEVILS									
PRODUCTION DEPARTMENT									
STANDARD OPERATING PROCEDURE									
Department: Production     SOP No.:									
			eves and Screens			Effective Da	nte:		
Superse	des: Nil					Review Date	e:		
Issue Da	ite:					Page No.:			
ANNEXURE I SIEVES INVENTORY RECORD									
Month/Year: Sieve No.:									
S.No.	Quantity	Received	Supplier	Date of Receipt	Sieve Code	No.		Checked By	
S.No.	Date		Sieve Code no.	Visual	Inspection	Scrappe	d On	Remarks	
				ОК	Not OK				
Done By: Checked By:									

PHARMA DEVILS PRODUCTION DEPARTMENT																
STANDARD OPERATING PROCEDURE																
Department: Production     SOP No.:																
Title: In	ventory and	Storage of Sieves an	nd Screens				]	Effective Dat	te:							
Supersedes: Nil Review Date:									:							
Issue Da	ate:						]	Page No.:								
ANNEXURE II SCREENS INVENTORY RECORD Aperture Size:																
S.No.	Quant	ity Received	Supplier		Date of Receipt Screen G			ode no.	Checked By							
S. No.	Date				pped On	Remarks										
				OK	Not OK											
Done B	v:							(	Done By:							