



VALIDATION PROTOCOL FOR COATING OF TABLETS

- 1. Objective:**
Coating of tablets using Neocota/Coating pan.
To validate the process of coating by ensuring that the tablets when coated, as per the Batch Manufacturing Record are uniformly coated and confirm to the specifications.
- 2. Scope:**
Applicable to all products undergoing coating (Aqueous, Film, Enteric and Seal).
- 3. Principle:**
Coating of tablets using Neocota / Coating pan.
- 4. Site of the Study:**
Hormone Coating Department.
- 5. Responsibility:**
Representatives from: Production
Quality Control
Engineering
Quality Assurance
(Individuals to name to be recorded in the report)
- 6. Description of the Equipment to be used:**
Equipment: Neocota/ Coating pan.
CODE No.: *Details to be recorded in the report*
CAPACITY: *Details to be recorded in the report.*
Date of equipment qualification to be recorded in the validation report.
- 7. BMR and SOP's to be followed:**
 - 7.1 **SOP for Operating Neocota:** SOP No. to be recorded in the report.
 - 7.2 **SOP for Operating Vernier caliper:** SOP number to be recorded in the report.
 - 7.3 **SOP for Operating Hardness tester:** SOP number to be recorded in the report.
 - 7.4 **SOP for Operating DT apparatus:** SOP number to be recorded in the report.
 - 7.5 **SOP for Operating Analytical balance:** SOP number to be recorded in the report.
 - 7.6 **Batch Manufacturing Record:** *Manufacturing Code No. to be recorded in the report.*
- 8. Controls:**
 - 8.1 Requirement:**
Validated Analytical Methods for estimation of active ingredient/s: *Details to be recorded in Report.*
 - 8.2 Calibration:**
Calibration of Equipment/Testing Apparatus. *Calibration details to be recorded in the report.*
 - 8.3 Training:**
Availability of training records of personnel involved in the validation exercise should be recorded in the validation report.




VALIDATION PROTOCOL FOR COATING OF TABLETS

8.4 Precautions:

- 8.4.1 Ensure instrument Qualification is been done before use.
- 8.4.2 Ensure proper instructions are followed as laid down in BMR and SOP.
- 8.4.3 Ensure Calibration status of testing apparatus before use.
- 8.4.4 Safety Aspects while operation of equipment and process must be ensured.

9. Validation Procedure:

- 9.1 Operate the equipment as per the SOP.
- 9.2 Record the following in the validation report.
 - 9.3.1 Product Name
 - 9.3.2 Batch Number
 - 9.3.3 Equipment
 - 9.3.4 Quantity loaded
 - 9.3.5 Batch size in Kgs.
 - 9.3.6 Volume occupied
 - 9.3.7 Number of lots
 - 9.3.8 Shape of pan
 - 9.3.9 Number of baffles
 - 9.3.10 Working capacity of pan
 - 9.3.11 Number of spray guns
 - 9.3.12 Ambient temperature of the coating cubicle/area
 - 9.3.13 Date/s of experiment
- 9.4 Take a composite sample of 50 tablets from 5 positions as shown below from the tablet bed during the process, after time intervals of 1 hour.



Check the followings:

- a. Average weight of 20 tabs. for weight rise.
- b. Individual Weight and Thickness variation on 20 tabs.
- c. Any visual defects.

9.5 Sample 100 tablets before and after coating, check the following:

- a. Average weight of 50 tablets for weight rise
 - b. Individual Weight variation on 50 tabs.
 - c. Individual Thickness variation on 50 tabs.
 - d. Hardness (on 5 tabs).
 - e. Disintegration time of 6 tablets.
 - f. Any visual defects.
- 9.6 Check the dissolution of tablets. (Refer QC reports).
 - 9.7 Record the process parameters.
 - 9.8 Check and record the spray rate at a regular interval of one hour till the end of any one Lot.



VALIDATION PROTOCOL FOR COATING OF TABLETS

- 9.10 Take a composite sample of 1000 tablets and do a rejection analysis to identify the types of defects.
- 9.11 Record the following process parameters at time intervals of every 1 hour (or) as per the intervals given in the BMR throughout the coating process for the Batch or lot.
- 9.11.1 Inlet, Exhaust and bed temperature
 - 9.11.2 Pan RPM
 - 9.11.3 Atomizing air pressure
 - 9.11.4 Visual inspection for over wetting or under spraying or over drying.
 - 9.11.5 Check the spray gun to tablet bed distance initially at the start of the batch/Lot.
 - 9.11.6 Check and record the spray rate at a regular interval of one hour (or) as per the intervals given in the Batch Manufacturing Record till the end of the batch or lot.
- 9.12 Attach all data records to the validation report.
- 9.13 Compare the data with Acceptance Criteria.
- Note:** Tablets collected for in-process tests to be destroyed after Completion of in-process tests as per SOP.

10. Acceptance criteria:

- 10.1 The Average weight, average weight rise and thickness of coated tablets should comply to the standard limits given in the BMR.
- 10.2 Hardness and disintegration time of coated tablet should comply to the standard limits given in the Batch Manufacturing Record.
- 10.3 Visual inspection after every interval should not show signs of over wetting or under spraying or over drying.
- 10.4 The process parameters throughout the coating should comply to the standard limits given in the Batch Manufacturing record.
- 10.5 Dissolution results should comply to the QC specifications.
- 10.6 The quantity of rejection after inspection should not be more than 1.0% out of the total quantity sorted.

11. Non Compliance:

11.1 Deviation:

Details of deviations (including justification of acceptance) done to successfully carry out the validation exercise and should be investigated. (Attach the details in the Validation report).

11.2 Out of Specification:

Any out of specification result observed should be recorded in the validation report and investigation should be done.

12. Type of validation:

Concurrent Validation/Revalidation.

13. Frequency of validation:

- 12.1 New Product: Three successful validation exercises.
- 12.2 Revalidation (Existing Product): One validation exercise within five years.
- 12.3 Revalidation (After major change): Three successful validation exercises



VALIDATION PROTOCOL FOR COATING OF TABLETS

14. Results/Observations:

Record the observations during the study and results obtained from Quality Control Department in the Validation Report.

15. Summary of findings:

Summarize the findings of the Validation Study to draw an inference.

16. Recommendation:

Record the recommendations based on the interpretation of the results of the Validation Report.

17. Team Approval:

The individuals who have performed the Validation Study, supervised the validation, completed the records, performed the testing of the product should approve the validation report.

18. Review:

The Validation Report should be reviewed by Unit Quality Assurance and Unit Head. The report should include any follow up action if required.

19. Approval by:

Validation Report should be finally approved by Unit Quality Assurance and Unit Head.

20. Attachment:

Annexure (if any) attached to the Validation Report should be recorded.

21. Abbreviations:

OOS : Out of specification
BMR : Batch Manufacturing Record
T.I Sheet : Technical Information Sheet
SOP : Standard Operating Procedure
A.R No. : Analytical Reference Number