



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

1. OBJECTIVE:

To demonstrate that the critical operation involved in container/bottle packing of _____ Tablets/ Capsules, B.No._____, B. Size_____ & Pack Size_____, are capable of consistently packing batch / batches which meet all quality parameters and the limits of acceptance criteria.

2. JUSTIFICATION FOR SELECTION OF ITEM/ EQUIPMENT/ PROCESS/ PRODUCT/ SYSTEM:

3. SCOPE:

Applicable to container/bottle packing line of XXXX tablets / capsules of XX / XXX count container/bottle.

4. SITE OF STUDY: _____

5. RESPONSIBILITY:

Representatives from,

Packaging : _____

Quality Assurance : _____

Engineering : _____

6. ** SOP AND BPR TO BE FOLLOWED:

Master Batch Packing Record of Manufacturing Code: _____ ,
Pack Code: _____ and Version No. _____

6. SOP for Bottle un-scrambler and cleaning machine, SOP No.: _____

1 SOP for Tablet/Capsule Counting Machine, SOP No.: _____

SOP for Tablet/Capsule count verification Machine, SOP No.: _____

SOP for Silica gel bag insertion machine, SOP No.: _____

SOP for Cotton insertion machine, SOP No.: _____

SOP for weighing balance, SOP No.: _____

SOP for Capping Machine, SOP No.: _____

SOP for Bottle un-scrambler and cleaning machine, SOP No.: _____

SOP for Tablet/Capsule Counting Machine, SOP No.: _____

SOP for Bottle un-scrambler and cleaning machine, SOP No.: _____

SOP for Tablet/Capsule Counting Machine, SOP No.: _____

SOP for Tablet/Capsule count verification Machine, SOP No.: _____

SOP for Silica gel bag insertion machine, SOP No.: _____



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

SOP for Cotton insertion machine, SOP No.: _____

SOP for weighing balance, SOP No.: _____

SOP for Capping Machine, SOP No.: _____

SOP for Torque Testing Apparatus, SOP No.: _____

SOP for Induction Sealing Machine, SOP No.: _____

SOP for Bottle un-scrambler and cleaning machine, SOP No.: _____

SOP for Tablet/Capsule Counting Machine, SOP No.: _____

SOP for Tablet/Capsule count verification Machine, SOP No.: _____

** Machines as applicable to individual products should be mentioned.

CONTROLS:

** Qualification and Calibration details of machine:

Equipment Name	Code No.	Qualification done on	Calibration done on	Calibration due on
Bottle un-scrambler and cleaning machine				
Tablet/Capsule Counting Machine				
Tablet/Capsule Count Verification Machine				
Silica gel bag insertion machine				
Cotton insertion machine				
Weighing balance				
Capping Machine				
Torque Testing Apparatus				



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Equipment Name		Code No.	Qualification done on	Calibration done on	Calibration due on
Induction Sealing Machine	% Power display				
	Conveyer speed display				
Leak test Apparatus	Timer				
	Vacuum Indicator				
Retorquer Machine					
Sticker labeling Machine					
Problue adhesive melter machine					
Leaflet (Pack insert) placement machine					
Cartonator					
Print check system					
Shrink wrapping machine					
Checkweigher					
Shipper sealing machine					
Shipper strapping machine					
Shipper weighing balance					

** Machines as applicable to individual products should be mentioned.

Identify each sampled container/bottle with interval no., date and time



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

****Packing materials/ components :**

Packaging Configuration	Item Code	Item Name	Description
xx- Count Container/bottles		Container/bottle	
		Cap	
		Silica gel bag	
		Cotton / Rayon Sani coil	
		Sticker label	
		Pack insert	
		Patient information leaflet / Medication pad	
		Inner Carton	
		Outer carton	
		Shipper	
		Strap	
Xxx -Count Container/bottles		Container/bottle	
		Cap	
		Silica gel bag	
		Cotton / Rayon Sani coil	
		Sticker label	
		Pack insert	
		Patient information leaflet / Medication pad	
		Inner Carton	
		Outer carton	
		Shipper	
		Strap	

** Packing materials/ components as applicable to individual product should be mentioned.

Training:

Name	Training reports availability	Status

Precautions: Safety aspects while operation of equipment and process should be ensured.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

VALIDATION PROCEDURE:

Note : Draw samples at frequencies given for individual test. Ensure sample is taken at the end of batch run also. In case of small batch size, ensure that minimum 3 intervals at start, middle and end at optimum batch run are considered for sample collection.

Speed setting of container/ bottle packing line will change as per the count, container/ bottle size and shape and tablet / capsule size and shape.

Same marked container/ bottle and cap (if applicable for that test) should be used throughout the batch for the following tests :

- No foil sensor challenge test
- Label absence challenge
- Functioning of Checkweigher

Upon completion of the above tests, the container/ bottle and cap (if applicable for that test) should be destroyed and recorded in the packing material reconciliation and destruction.

Batch packing should be carried out as per BPR :

Validation should be done as per validation protocol no. _____

Version no.: _____

Date(s) of Validation		
Batch No.		
Batch size (Tablets /Capsules)		
Quantity to be packed (containers/bottles)		
Fill count (Tablets /Capsules)		
Container/bottle	Size	
	Code No.	
Cap	Size	
	Code No.	

Packing process flow diagram : As per reference protocol no. : _____

Version no.:



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Silica gel bag **::

Sampling Time	Containers/bottles										Sampled and Checked by	
	1		2		3		4		5			
	Quality	Qty.	Quality	Qty.	Quality	Qty.	Quality	Qty.	Quality	Qty.		

Enter V if acceptable (intact bag) and x if Not acceptable

Weighing of cotton / rayon sani coil **: :

Weighing balance code no. _____ .

Sampling Time	Containers/bottles					Sampled and Checked by
	1 (in g)	2 (in g)	3(in g)	4 (in g)	5 (in g)	

Closing and opening torque (for capped container/bottle) **: :

Sampling Time	Closing torque					Opening torque					Sampled and Checked by	
	1 (in-lb)	2 (in-lb)	3 (in-lb)	4 (in-lb)	5 (in-lb)	1 (in-lb)	2 (in-lb)	3 (in-lb)	4 (in-lb)	5 (in-lb)		

No foil sensor challenge:

Sampling Time	No foil sensor challenge	Sampled and Checked by

For no foil sensor challenge: X Rejected



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Induction sealing :

Starting time _____, **Completion time :** _____

	Sealer Height (Between the induction sealer head and cap of container/bottle)	Induction sealing quality** (5 Containers/ bottles)					Leak Test (containers/bottles)	Checked by
		1	2	3	4	5		
Low conveyer speed:								
Low Dower percentage:								
Low conveyer speed:								
High power percentage:								

For induction sealing: Enter √ if acceptable and x if Not acceptable

- Note :**
1. Induction sealing height should be kept constant.
 2. After completion of the leak test, container/bottle and cap of each parameter should be destroyed and recorded in the packing material reconciliation and destruction.

Closing and opening torque (For re-torquing):**

Sampling Time	Closing Torque.					Opening Torque					Sampled and Checked by	

Challenge and quality inspection of labels by vision system :

Challenge time	Pharmacode/code No. /issue date verification (label)	Overprinting verification (label)	Challenged by

Enter X if machine stopped / rejected

Note : After completion of test, label or container/bottle with label should be destroyed and recorded in the packing material reconciliation and destruction.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Quality of label **::

Sampling Time	Container Bottles					Sampled and Checked by
	1	2	3	4	5	

Enter √ if acceptable and X if not acceptable.

Label absence challenge :

Challenge time	Containers/Bottles					Challenged by
	1	2	3	4	5	

Inspection of dispensing of hot melt glue on cap**:

Sampling Time	Containers/Bottles					Checked by
	1	2	3	4	5	

Position of pack insert (leaflet)**:

Sampling time	Container/bottle					Sampled and Checked by
	1	2	3	4	5	

Enter √ if acceptable and X if Not acceptable

Pharmacode verification of pack insert (leaflet):

Challenge Time	Pharmacode verification (Pack insert)	Challenged by

Enter X if machine stopped / rejected

Note: After completion of pharmacode verification test, pack insert should be destroyed and recorded in packing material reconciliation and destruction.

Functioning of Checkweigher

Challenge Time	Functioning of checkweigher	Challenged by



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE
VALIDATION (TABLETS/CAPSULES)**

Enter √ if acceptable and X if Not acceptable

Challenge test for carton and pack insert (leaflet) in cartonator

In process Time	In process challenge test		Result	Done by
	Challenge test for carton sensor (carton missing)			
	Challenge test for Leaflet sensor (Leaflet missing)			
	Challenge test for pharmacode of Pack insert (Leaflet) Back side and front side			
	Challenge test for pharmacode of carton			
	Visual inspection	Embossing/ overprinting quality on carton		
		Carton fill value		
		Carton quality		
		Pack insert (Leaflet) folding quality		

Enter X if machine stopped/rejected.

Enter √ if acceptable and x if Not acceptable.

Note: After completion of pharmacode verification test for carton and pack insert (leaflet), the same should be destroyed and recorded in packing material reconciliation and destruction.

Inspection of defects (AQL inspection, Ref. SOP: _____):
For Container/bottle filling operation (Sample size (n): 50 Containers/Bottles):
Sampling Time:



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Test	Critical AQL: 0.10 Accept: 0 Reject: 1	Checked By
Fill count out of range		
Improper Silica gel bag quality		
Improper Cotton / Rayon Sani coil quality		
Incorrect product		
Hazardous foreign material (For example, glass, metal)		
Incorrect component (container/bottle and cap)		
Extraneous foreign material		

Test	Critical AQL: 0.65 Accept: 1 Reject: 2	Checked By
Damaged container/bottles and caps		
Closing and opening torque out of range		
Improper capping		
Missing ribs closers/cap		
Missing aluminium layer in cap		

Test	Major AQL: 1.0 Accept: 1 Reject: 2	Checked By
Excess powder		
Improper Tablet / Capsule quality		

Test	Minor AQL: 2.5 Accept: 3 Reject: 4	Checked By
Scratches on containers/bottles and caps		

Test	Critical AQL: 0.65 Accept: 1 Reject: 2	Checked By
Incomplete seal		
Burnt Induction sealing		

For Labeling operation (Sample size (n): 50 container/bottles)

Sampling Time:



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE
VALIDATION (TABLETS/CAPSULES)**

Test	Critical AQL: 0.10 Accept: 0 Reject: 1	Checked By
Incorrect, Missing or Damaged Label		
Incorrect, Missing or Illegible Lot no.		
Incorrect, Missing or Illegible Expiration date		
Incorrect NDC no.		
Missing / Wrong Pack inserts and Patient information leaflet		
Other product labels mix-up		

Test	Major AQL: 2.5 Accept: 3 Reject: 4	Checked By
Smudging of ink and ink lifting		

Test	Minor AQL: 4.0 Accept: 5 Reject: 6	Checked By
Non-Conforming Print		
Presence of stains		
Minor Scratches on Label		

**** Note:** Upon completion of the physical testing, container/bottles and tablet/ capsule used for non destructive testing should be returned to the line officer for verification and addition to packing line.

HIGH SPEED SETTING :

Set up parameters for high speed

Equipment Setting time from _____ to _____

Note : Write NA wherever not applicable

Conduct following test on high speed setting :

**** Note :** Upon completion of the physical testing, container/bottles and tablet / capsule used for non destructive testing should be returned to the line officer for verification and addition to packing line.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Enter √ if acceptable (intact bag) and x if Not acceptable

Weighing of cotton / rayon sani coil:

Weighing balance code no. _____ .

Sampling Time	Container/bottle					Sampled and Checked by
	1 (in g)	2 (in g)	3 (in g)	4 (in g)	5 (in g)	

Closing and opening torque (for capped container/bottle)**:

Sampling Time	Closing torque					Opening torque					Sampled and Checked by	
	1 (in-lb)	2 (in-lb)	3 (in-lb)	4 (in-lb)	5 (in-lb)	1 (in-lb)	2 (in-lb)	3 (in-lb)	4 (in-lb)	5 (in-lb)		

No foil sensor challenge:

Sampling Time	No foil sensor challenge	Sampled and Checked by

For no foil sensor challenge: X Rejected

Induction sealing :

Starting time: _____ , Completion time : _____

	Sealer Height (Between the induction sealer head and cap of container/bottle)	Induction sealing quality** (5 Containers/ bottles)					Leak Test (containers/ bottles)	Checked by
		1	2	3	4	5		
High conveyer speed:								
High power percentage:								
High conveyer speed:								
Low power percentage:								

Note : 1. Induction sealing height should be kept constant.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

2. After completion of the leak test, container/bottle and cap of each parameter should be destroyed and recorded in the packing material reconciliation and destruction.

Closing and opening torque (For re-torquing)**:

Sampling Time	Closing torque					Opening torque					Sampled and Checked by
	1 (in-lb)	2 (in-lb)	3 (in-lb)	4 (in-lb)	5 (in-lb)	1 (in-lb)	2 (in-lb)	3 (in-lb)	4 (in-lb)	5 (in-lb)	

Challenge and quality inspection of labels by vision system :

Challenge time	Pharmacode/ code No. /issue date verification (label)	Overprinting verification (label)	Challenged by

Enter X if machine stopped / rejected

Note: After completion of test, label/Container/bottle with label should be destroyed and recorded in the packing material reconciliation and destruction.

Quality of label **::

Sampling time	Containers/bottles					Sampled and Checked by
	1	2	3	4	5	

Enter √ if acceptable and X if Not acceptable

Label absence challenge:

Challenge time	Containers/bottles					Challenged by
	1	2	3	4	5	

Enter X if machine stopped / rejected

Inspection of dispensing of hot melt glue on cap**:

Sampling time	Containers/bottles					Checked by
	1	2	3	4	5	



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

--	--	--	--	--	--	--

Enter √ if acceptable and x if Not acceptable

Position of pack insert (leaflet)** :

Sampling time	Containers/bottles					Sampled and Checked by
	1	2	3	4	5	

Enter √ if acceptable and x if Not acceptable

Pharmacode verification of pack insert (leaflet):

Challenge Time	Pharmacode verification (Pack insert)	Challenged by

Enter X if machine stopped / rejected

Note: After completion of pharmacode verification test, pack insert should be destroyed and recorded in packing material reconciliation and destruction.

Functioning of checkweigher:

Challenge Time	Functioning of checkweigher	Challenged by

Enter √ if acceptable and x if Not acceptable

Challenge test for carton and pack insert (leaflet) in cartonator

In-process time	In-process challenge test	Result	Done by	
	Challenge test for carton sensor (carton missing)			
	Challenge test for Leaflet sensor (Leaflet missing)			
	Challenge test for pharmacode of carton			
	Challenge test for pharmacode of pack insert (leaflet) Back side and front side			
	Visual inspection	Embossing/ overprinting quality		
		Carton fill value		
		Carton quality		
Pack insert (leaflet) folding quality				

Enter X if machine stopped/rejected. Enter √ if acceptable and x if Not acceptable.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Note: After completion of the test, carton and pack insert (leaflet) should be destroyed and recorded in packing material reconciliation and destruction.

Inspection for defects (AQL inspection):

For container/bottle filling operation (Sample size (n): 50 container/bottles, Ref.

SOP: _____)

Sampling Time:

Test	Critical AQL: 0.10 Accept: 0 Reject: 1	Checked By
Fill count out of range		
Improper Silica gel bag quality		
Improper Cotton / Rayon Sani coil quality		
Incorrect product		
Hazardous foreign material (For example, glass, metal)		
Incorrect component (container/bottle and cap)		
Extraneous foreign material		

Test	Critical AQL: 0.65 Accept: 1 Reject: 2	Checked By
Damaged containers/bottles and caps		
Closing and opening torque out of range		
Improper capping		
Missing ribs, closers/cap		
Missing aluminium layer in cap		

Test	Major AQL: 1.0 Accept: 1 Reject: 2	Checked By
Excess powder		
Improper Tablet / Capsule quality		

Test	Minor AQL: 2.5 Accept: 3 Reject: 4	Checked By
Scratches on containers/bottles and caps		

Sampling Time:

Test	Critical AQL: 0.65 Accept: 1 Reject: 2	Checked By
Incomplete seal		
Burnt Induction sealing		



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

For Labeling operation (Sample size (n): 50 container/bottles)
(Ref. SOP: _____)

Sampling Time:

Test	Critical AQL: 0.10 Accept: 0 Reject: 1	Checked By
Incorrect, Missing or Damaged Label		
Incorrect, Missing or Illegible Lot no.		
Incorrect, Missing or Illegible Expiration date		
Incorrect NDC no.		
Missing / Wrong Pack inserts and Patient information leaflet		
Other product labels mix-up		

Test	Major AQL: 2.5 Accept: 3 Reject: 4	Checked By
Smudging of ink and ink lifting		

Test	Minor AQL: 4.0 Accept: 5 Reject: 6	Checked By
Non-Conforming Print		
Presence of stains		
Minor Scratches on Label		

** Note : Upon completion of the physical testing, container/bottles and tablet/capsule used for nondestructive testing should be returned to the line officer for verification and addition to _____ packing line.

OPTIMUM SPEED SETTING:

Note: In case of small batch size, ensure that minimum 3 intervals at start, middle and end at Optimum batch run are considered for AQL sample collection.

Set up parameters for optimum speed



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Equipment Setting time from _____ to _____

Date	Equipment Name	Speed (Container/bottles per Minute)	Done by	Checked by
	Tablets/capsule counting machine	Optimum speed		
	Silica gel bag insertion machine			
	Cotton/ Rayon Sani coil insertion machine			
	Capping machine			
	Induction Sealing machine			
	Retorquer machine			
	Sticker Labeling machine			
	Problue adhesive melter machine			
	Leaflet placement machine			
	Checkweigher			
	Cartonator (if applicable)			
	Print check system			

Note : Write NA wherever not applicable

Conduct following test on optimum speed setting :

**** Note :** Upon completion of the physical testing, container/bottles and tablet/ capsule used for nondestructive testing should be returned to the line officer for verification and addition to packing line.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Fill Count**:

Sr. No.	Date	Sampling Time	Container/bottle										Sampled and Checked by	
			Side A1	Side B1	Side A2	Side B2	Side A3	Side B3	Side A4	Side B4	Side A5	Side B5		
1														
2														
3														
4														
5														
6														
7														
8														
9														
10														
11														
12														
13														
14														
15														
16														
17														
18														
19														
20														
21														
22														
23														
24														
25														
26														
27														
28														
29														
30														
31														
32														
33														
34														
35														

Total No. of sample (Includes slow and High speed) :



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Silica gel bag **

Sr. No.	Date	Sampling Time	Container/bottle										Sampled and Checked by
			1		2		3		4		5		
			Quality	Qty	Quality	Qty.	Quality	Qty.	Quality	Qty.	Quality	Qty.	
1													
2													
3													
4													
5													
6													
7													
8													
9													
10													
11													
12													
13													
14													
15													
16													
17													
18													
19													
20													
21													
22													
23													
24													
25													
26													
27													
28													
29													
30													
31													
32													
33													
34													
35													

Total No. of Sample (Includes Slow and High Speed):

Enter √ if acceptable (intact bag) and X if not acceptable.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Weighing of cotton / rayon sani coil **: _____

Weighing balance code no. _____.

Sr. No.	Date	Sampling Time	Container/bottle					Sampled and Checked by
			1 (in g)	2 (in g)	3 (in g)	4 (in g)	5 (in g)	
1								
2								
3								
4								
5								
6								
f								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								
21								
22								
23								
24								
25								
26								
27								
28								
29								
30								
31								
32								
33								
34								
35								

Total no. of sample (Includes slow and High speed): _____



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Closing and opening torque (for capped container/bottle) :**

Sr. No.	Date	Sampling Time	Closing torque					Opening torque					Sampled and Checked by
			1 [in-Ib]	2 [in-Ib]	3 [in-Ib]	4 [in-Ib]	5 [in-Ib]	1 [in-Ib]	2 [in-Ib]	3 [in-Ib]	4 [in-Ib]	5 [in-Ib]	
1													
2													
3													
4													
5													
6													
7													
8													
9													
10													
11													
12													
13													
14													
15													
16													
17													
18													
19													
20													
21													
22													
23													
24													
25													
26													
27													
28													
29													
30													
31													
32													
33													
34													
35													

Total no. of sample (Includes slow and High speed):



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

No foil sensor challenge :

Sr. No.	Date	Sampling Time	No foil sensor challenge	Sampled and Checked by
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				
32				
33				
34				
35				

For no foil sensor challenge: X Rejected



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Induction sealing :

** Induction sealing machine parameter and **induction sealing quality:

Sr. No.	Date	Time	Sealer Height (Between the induction sealer head and cap of the container)	Optimum % Power	Optimum Speed of conveyor belt (mpm)	Induction sealing quality					Checked by
						1	2	3	4	5	
1											
2											
3											
4											
5											
6											
7											
8											
9											
10											
11											
12											
13											
14											
15											
16											
17											
18											
19											
20											
21											
22											
23											
24											
25											
26											
27											
28											
29											
30											
31											
32											
33											
34											
35											

Total no. of sample (Includes slow and High Speed)



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Leak test:

Sr. No.	Date	Sampling Time	Leak test on containers /bottles	Sampled and Checked by
1 / Initial				
2 / Middle				
3/End				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				
32				
33				
34				
35				

Total no. of samples (Includes slow and High speed):

Note: Upon completion of the leak test, container/bottle and cap should be destroyed and recorded in the packing material reconciliation and destruction.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Closing and opening torque (For re-torquing)**:

Sr. No.	Date	Sam-pling Time	Closing torque					Opening torque					Sampled and Checked by
			1 (in-lb)	2 (in-lb)	3 (in-lb)	4 (in-lb)	5 (in-lb)	1 (in-lb)	2 (in-lb)	3 (in-lb)	4 (in-lb)	5 (in-lb)	
1													
2													
3													
4													
5													
6													
7													
8													
9													
10													
11													
12													
13													
14													
15													
16													
17													
18													
19													
20													
21													
22													
23													
24													
25													
26													
27													
28													
29													
30													
31													
32													
33													
34													
35													

Total no. of Sample (Includes Slow and High Speed)



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Challenge and quality inspection of labels by vision system:

Sr. No.	Date	Challenge time	Pharmacode/code No. /issue date verification (label)	Overprinting verification (label)	Challenged by
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					
28					
29					
30					
31					
32					
33					
34					
35					
Total no. of Sample (Includes Slow and High Speed)					

Enter X if Machine Stopped/ Rejected

Note : After completion of test, label or Container/bottle with label should be destroyed and)\ recorded in the packing material reconciliation and destruction.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Quality of label **::

Sr. No.	Date	Sampling time	Container/bottle					Sampled and Checked by
			1	2	3	4	5	
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								
21								
22								
23								
24								
25								
26								
27								
28								
29								
30								
31								
32								
33								
34								
35								

Total no. of sample (includes Slow and High Speed)

Enter √ if acceptable and X if not acceptable.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Label absence challenge:

Sr. No.	Date	Challenge time	Container/bottle					Challenged by
			1	2	3	4	5	
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								
21								
22								
23								
24								
25								
26								
27								
28								
29								
30								
31								
32								
33								
34								
35								

Total No. of Samples:

Enter X if machine stopped / rejected



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE
VALIDATION (TABLETS/CAPSULES)**

Inspection of dispensing of hot melt glue on cap :**

Sr. No.	Date	Sampling Time	Container/bottle					Checked by
			1	2	3	4	5	
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								
21								
22								
23								
24								
25								
26								
27								
28								
29								
30								
31								
32								
33								
34								
35								

Total no. of samples (includes Slow and high Speed)

Enter √ if acceptable and x if not acceptable



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Position of pack insert (leaflet)** :

Sr. No.	Date	Sampling time	Container/bottle					Sampled and Checked by
			1	2	3	4	5	
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								
21								
22								
23								
24								
25								
26								
27								
28								
29								
30								
31								
32								
33								
34								
35								

Total no. of samples (includes Slow and high Speed)

Enter √ if acceptable and x if not acceptable



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Pharmacode verification of pack insert (leaflet):

Sr. No.	Date	Challenge Time	Pharmacode/ verification (Pack insert)	Challenged by
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				
32				
33				
34				
35				

Total no. of samples (includes Slow and high Speed)

Enter X if machine stopped / rejected

Note: After completion of Pharmacode verification test pack insert should be destroyed and recorded in packing material reconciliation and Destruction.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Functioning of checkweigher:

Sr. No.	Date	Challenge time	Functioning of checkweigher	Challenge by
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				
32				
33				
34				
35				

Total no. of samples (includes Slow and high Speed)

Enter √ if acceptable and x if not acceptable



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Sr. No.	Date	Challenge time	Challenge Test For				Visual inspection			Done by
			Carton Sensor (Carton Missing)	Leaflet sensor (Leaflet Missing)	Pharmacode of Carton	Pharmacode of Leaflet (Pack Insert) Back side and Front Side	Embossing /Overprinting Quality	Carton Fill Value	Carton Quality	
31										
32										
33										
34										
35										

Total no. of samples (includes Slow and high Speed)

For challenge: Enter X if machine stopped / rejected

For visual inspection: Enter √ if acceptable and x if Not acceptable.

Note: After completion of test, carton and pack insert (leaflet) should be destroyed and recorded in packing material reconciliation and destruction.

Sampling and inspection of defects (AQL inspection, Ref. SOP: _____):

Note : In case of small batch size, ensure that minimum 3 intervals at start, middle and end of Optimum batch run are considered for AQL sample collection.

Container/bottle filling :

As per reference protocol No. _____ , Version no. _____

Sampling Interval	Date	Interval Time	Sample Quantity	Sampled by
Optimum Speed 1.				
2.				
3.				
4.				
5.				
6.				
7.				
8.				
9.				
10.				
Total no. of Samples:				

****Note:** Upon completion of the physical testing, container/bottles used for nondestructive testing should be returned to the line officer for verification and addition to packing line.



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

Time											
Interval No.		1	2	3	4	5	6	7	8	9	10
Number of container/bottles inspected											
Sample size (n): Major AQL: 2.5 Accept: Reject: (Ref. SOP: _____)	Smudging of ink and ink lifting										
Sample size (n): Minor AQL: 4.0 Accept: Reject: (Ref. SOP: _____)	Non-Conforming Print										
	Presence of stains										
	Minor Scratches on Label										
Inspected By/ Date:											

Tests	Date		
	Time		
		Initial	Middle
Shipper fill value			
Shrink quality			
Carton quality			
Shipper labeling quality			
Transcription accuracy			
Strapping quality			
Checked by			



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

**VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE
VALIDATION (TABLETS/CAPSULES)**

REQUIRED TESTING AND ACCEPTANCE/ EVALUATION CRITERIA:

As per reference protocol No. _____, Version no. _____

DETAILS OF DEVIATIONS/ NON CONFORMANCE:

TYPE OF VALIDATION:

Concurrent Validation / Revalidation.

FREQUENCY:

Three consecutive batches for each container/bottle configuration.

Three consecutive batches in case of change in any machine from the line.

Revalidation: Every three years on one batch for each configuration.

RESULT AND OBSERVATIONS:

RISK MANAGEMENT STUDY :



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE VALIDATION (TABLETS/CAPSULES)

SUMMARY OF FINDINGS OF VALIDATION ACTIVITY:

Fill count	:	<input type="text"/>
Silica gel bag quality and quantity	:	<input type="text"/>
Cotton/ Rayon sani coil weight	:	<input type="text"/>
Opening torque (after capping)	:	<input type="text"/>
Closing torque (for capping)	:	<input type="text"/>
No foil sensor challenge	:	<input type="text"/>
Induction Sealing machine parameters	:	<input type="text"/>
Induction sealing quality	:	<input type="text"/>
Leak test	:	<input type="text"/>
Opening torque (after retorquing)	:	<input type="text"/>
Closing torque (for retorquing)	:	<input type="text"/>
Pharmacode/code No. /issue date verification (label)	:	<input type="text"/>
Overprinting verification (label)	:	<input type="text"/>
Quality of labeling	:	<input type="text"/>
Label absence challenge	:	<input type="text"/>
Dispensing of hot melt glue on cap	:	<input type="text"/>
Position of pack insert (leaflet) on cap	:	<input type="text"/>
Pharmacode verification of pack insert (leaflet)	:	<input type="text"/>
Functioning of checkweigher	:	<input type="text"/>
Carton sensor in cartonator challenge	:	<input type="text"/>



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

**VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE
VALIDATION (TABLETS/CAPSULES)**

Leaflet sensor in cartonator challenge	:	<input type="text"/>
Pharmacode verification of carton in cartonator	:	<input type="text"/>
Pharmacode verification of pack insert (leaflet)	:	<input type="text"/>
Back side and front side in cartonator	:	<input type="text"/>
Overprinting/Embossing quality of carton in cartonator	:	<input type="text"/>
AQL for container/bottle filling operation	:	<input type="text"/>
AQL for labeling operation	:	<input type="text"/>

Within Acceptance Criteria

Not within Acceptance Criteria

NA

Not Applicable

RECOMMENDATIONS:

TEAM APPROVAL

User

Date:

Quality Assurance

Date:

Engineering

Date:

REVIEW:



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

**VALIDATION REPORT FOR GUIDELINE FOR CONTAINER/BOTTLE PACKING LINE
VALIDATION (TABLETS/CAPSULES)**

APPROVED BY:

Unit Quality Assurance Head

Date:

NOTED BY:

Unit Head

Date:

ANNEXURES:



ABBREVIATIONS

mg	:	milligram
SOP	:	Standard Operating Procedure
BPR	:	Batch Packing Record
No.	:	Number
%	:	Percent
AQL	:	Acceptable Quality Level
Ref.	:	Reference
NDC	:	National Drug Code
NA	:	Not applicable
Qty.	:	Quantity
g	:	gram
In	:	inch
lb	:	pound
Sr.	:	Serial
mpm	:	Meter per minute