



**PHARMA DEVILS**

**PROTOCOL FOR CLEANING VALIDATION -  
PRODUCTION BULK  
(BLENDING)**

**PROTOCOL No.:**

**PROTOCOL**

**FOR**

**CLEANING VALIDATION**

**PRODUCTION BULK**  
**(BLENDING)**



**PHARMA DEVILS**

**PROTOCOL FOR CLEANING VALIDATION -  
PRODUCTION BULK  
(BLENDING)**

**PROTOCOL No.:**

**TABLE OF CONTENTS**

<b>S.No.</b>	<b>Details</b>	<b>Page No.</b>
<b>1.</b>	<b>OVERVIEW:</b>	<b>3</b>
1.1	Introduction:	3
1.2	Objective:	3
1.3	Scope:	3
1.4	Responsibility:	3
<b>2.</b>	<b>EXECUTION TEAM:</b>	<b>4</b>
<b>3.</b>	<b>TRAINING RECORD:</b>	<b>4</b>
<b>4.</b>	<b>REQUIREMENT FOR CLEANING VALIDATION:</b>	<b>4</b>
<b>5.</b>	<b>SYSTEM/EQUIPMENT DESCRIPTION:</b>	<b>5</b>
<b>6.</b>	<b>CLEANING VALIDATION METHODOLOGY OR PROCEDURE:</b>	<b>6</b>
<b>7.</b>	<b>ACCEPTANCE CRITERIA:</b>	<b>8</b>
<b>8.</b>	<b>REVALIDATION CRITERIA:</b>	<b>9</b>
<b>9.</b>	<b>OBSERVED NON-CONFORMANCES (IF ANY):</b>	<b>9</b>
<b>10.</b>	<b>VALIDATION REPORT:</b>	<b>9</b>
<b>11.</b>	<b>ABBREVIATIONS:</b>	<b>10</b>
<b>12.</b>	<b>LIST OF ANNEXURE:</b>	<b>10</b>
<b>13.</b>	<b>REFERENCE DOCUMENTS (IF ANY):</b>	<b>10</b>
<b>14.</b>	<b>REVISION HISTORY:</b>	<b>11</b>



**PHARMA DEVILS**

# **PROTOCOL FOR CLEANING VALIDATION - PRODUCTION BULK (BLENDING)**

**PROTOCOL No.:**

## **1. OVERVIEW:**

### **1.1 Introduction:**

The purpose of this protocol for assurance of cleaning methodology to detect previous product residue from equipment/equipments parts. To provide the necessary documentary evidence that the cleaning procedure can reduce the residual contamination of previous product below the established limit so that such residue do not affect the quality and safety of the subsequent product to be manufactured in the same equipment.

### **1.2 Objective:**

The objective of the cleaning validation protocol is to assure that the cleaning procedures of equipments, in-process containers in Production Bulk (Blending) facility can reproducibly remove residue of the product to levels below the established acceptance limits after manufacturing.

### **1.3 Scope:**

The scope of this validation activity is limited to validate the cleaning process of equipments in blending area of production Bulk department of Auronext Pharma Pvt. Ltd., Bhiwadi (Rajasthan).

**Reason for revision:** Periodic cleaning verification is included in revalidation criteria.

### **1.4 Responsibility:**

To conduct the cleaning validation study, a team shall be formed. The team shall contain the members from the Quality Control, Engineering, Production and Quality Assurance Departments. The Validation team is described through the following responsibility:

Quality Assurance	To prepare the protocol & report To provide the training To execute & supervise the study
Quality Control, Engineering, Production and Quality Assurance	To review the protocol To conduct the study To collect and analyze samples
<b>Reviewer-1 :</b> Production In-Charge/Designee	To review validation protocol. To review validation report.
<b>Reviewer-2:</b> CQA- Validation	To review validation protocol. To review validation report.
<b>Approver-1:</b> Head Production /Designee	To approve validation protocol. To approve validation report.
<b>Approver-2:</b> Head-QA/Designee	To approve validation protocol. To approve validation report.

## **2. EXECUTION TEAM:**

Following personnel shall be responsible for the execution of validation study:

Production : To conduct the validation study as per protocol.



- Quality Assurance : To monitoring the activity and collect the sample for chemical analysis as per protocol
- Quality Control : To conduct the microbiology monitoring, analysis of samples and reporting of results.
- Engineering : To provide utility and maintenance support.

**3. TRAINING RECORD:**

**3.1 Purpose:**

The purpose of the training is to familiarize the trainees with the purpose and procedure of cleaning validation activity.

**3.2 Scope:**

This training is applicable to the protocol for cleaning validation.

**3.3 Topics:**

The following topics shall be covered during training: Identifying the responsibility of involved person.

4.3.1 Purpose & procedure of cleaning validation.

4.3.2 Documentation practices to be followed.

4.3.3 General precautions / guidelines to be followed during validation.

**4. REQUIREMENT FOR CLEANING VALIDATION:**

**4.1 Documental Requirements:**

S. No.	Document No.	Title of SOP
01	PB048	Operation and cleaning of Blending and Filling machine
02	AN-QC-CL-CAL-0016	Procedure for calibration of high performance liquid chromatography.
03	AN-QC-CL-OPI-0033	Procedure for operation of high performance liquid chromatography (Model-Shimadzu).
04	AN-QC-CL-OPI-0046	Procedure for operation of high performance liquid chromatography (Model-Waters).
05	QC/GTP015	Bacterial endotoxin test procedure.
06	QC/GTP018	Particulate matter testing.
07	QC/GTP017	Bio burden testing by filtration method.

**4.2 Apparatus / Instrument Requirements:**

Polypropylene shaft (TEXWIPE),10×10 cm<sup>2</sup> stainless steel coupon, Test tube (Glass bottle), WFI, Rinse Sample bottle.

**5. SYSTEM/EQUIPMENT DESCRIPTION:**

**5.1 API Solubility:**

5.1.1 Meropenem: Sparingly soluble in water.

**5.2 Solubility Matrix:**

DESCRIPTION	SOLUBILITY (1g in listed ml)	mg/ml
-------------	------------------------------	-------



Very soluble (VS)	Less than 1	>1000
Freely soluble (FS)	From 1 to 10	100 to 1000
Soluble (S)	From 10 to 30	33.33 to 100
Sparingly soluble (SPS)	From 30 to 100	10 to 33.33
Slightly soluble (SLS)	From 100 to 1000	1 to 10
Very slightly soluble (VSLS)	From 1000 to 10 000	0.1 to 1
Practically insoluble (PI)	More than 10 000	<0.1

**5.3 Equipment Description:**

The list of equipments and their part is given below with their surface area:

S. No.	Name of Equipment	Calculated Surface Area (dm <sup>2</sup> )
1.	Blending Bin	171.51
2.	Distance piece	3.14
3.	Screw feeder	7.85
Total Surface Area		182.51
		18251 cm <sup>2</sup>

*☞ If different capacity of equipments used in product, higher contact surface area shall be considered for calculation.*

**6. CLEANING VALIDATION METHODOLOGY OR PROCEDURE:**

**6.1 General Recording Instructions:**

- 6.1.1 Read the contents of the document thoroughly before proceeding for Execution of the activity (in case of doubts / contradictions / contact the approvers of the document for clarifications).
- 6.1.2 Recording of all the observations and data shall be done as per **SOPAN-QA-GEN-0013 “Good Documentation Practices”**.

**6.2 Cleaning Methodology:**

- 6.2.1 The cleaning validation study shall be performed for three consecutive batches.
- 6.2.2 Cleaning procedure shall be evaluated for first three batches of all the new products, which is introduced in the premises.
- 6.2.3 Clean the equipments as per respective SOP (**Operation and cleaning of Blending and Filling machine**).
- 6.2.4 Collect the swab and rinse sample as per sampling procedure.
- 6.2.5 The swab sampling location is defined in **Annexure – 01**(Sampling locations of equipment).



**6.3 Product matrix:**

**Table-1**

S. No.	Name of product	API	Therapeutic Daily Dose Maximum (TDD)(mg)	Minimum Batch Size (kg)	Solubility (in water)
1	Meropenem for Injection	Meropenem	6000	24	Sparingly soluble

**6.3.1 Maximum Allowable Residue (MAR) calculation:**

**6.3.1.1 10 ppm Criteria:**General limit for maximum allowed concentration (mg/mg or ppm) of “previous” Substance in the next batch. = **0.00001mg** or 10 ppm.

☞ *Cleaning validation study the MACO (Maximum allowable carry over) limit based on 10ppm criteria and Dose criteria is not considered due to single product manufacturing facility.*

**6.4 Sampling Procedure and Test Methodology:**

**6.4.1 Physical verification:**

6.4.1.1 After completion of cleaning Production and QA personnel shall physically verify the equipment for cleanliness. It should be visually clean then only further activity shall be performed.

**6.4.2 Rinse Sample:**

6.4.2.1 Rinse the individual equipment as per respective SOP.

6.4.2.2 Send the entire samples along with blank sample to QC department for analysis of following test:

- Visual particles
- Sub visual particles
- Residue of previous product
- BET
- Bio burden

6.4.2.3 Depyrogenated container shall be used for BET sample collection.

6.4.2.4 Use the pre validated analytical method for analysis.

☞ *If samples are not analyzed immediately then samples are stores at 2 to 8°C and Samples shall be analyzed within 24 hours for chemical analysis.*

☞ *If samples are not analyzed immediately then samples are stores at 2 to 8°C and Samples shall be analyzed within 12 hours for microbial analysis.*

☞ *For bio burden sample use autoclaved bottle.*

**6.4.3 Swab Sample:**

6.4.3.1 Sampling locations (points) are selected based on the worst case criteria:

- Hard to clean area
- Inaccessible area

The sampling locations are defined in **Annexure – 01**.

6.4.3.2 Take the clean swab made up of Polyurethane foam having Polypropylene shaft (TEXWIPE) and put it in the test tube containing diluents and squeeze the swab along with the sides of test tube to remove the excess diluents from it.



6.4.3.3 Take out the wet swab from the test tube without touching the tip of swab.

6.4.3.4 Place the one side tip of swab at the identified location and apply it on the 10 x 10 cm<sup>2</sup> areas of the equipment/equipments parts (locations as defined in **Annexure – 01**). Swab samples are taken as shown in the **Figure – 01**.

☞ *In case of non-regular shape (e.g. cylindrical piping) simulate to the possible accessible area. Use appropriate swab holding devices to get the proper access of sampling points as described in the swab location.*

6.4.3.5 Put the swab stick into the test tube without touching the tip.

6.4.3.6 Label the test tube with the information of Swab sample, equipment part name, Sampling location, product name, batch no., sampled by and date.


6.4.3.7 Send the test tube to the QC department for analysis of following test:

- Residue of previous product
- BET
- Bio burden

6.4.3.8 Use the pre validated analytical method for analysis.

☞ *If samples are not analyzed immediately then samples are stores at 2 to 8°C and Samples shall be analyzed within 24 hours for chemical analysis.*

☞ *If samples are not analyzed immediately then samples are stores at 2 to 8°C and Samples shall be analyzed within 12 hours for microbial analysis.*



**Figure – 01**

## 6.5 Recovery, LOQ and LOD Study:

Recovery study, LOQ and LOD study is performed as per the analytical cleaning method validation protocol.

6.5.1 Analytical Method Validation protocol for Meropenem for Injection by HPLC.

### 6.5.2 Testing Plan:

6.5.2.1 **Rinse sample:** Visually inspect to particulate matter and analyze the rinse sample by suitable HPLC method for residue of previous product in the sample. BET shall be analyzed as per GTP QC/GPT015 and bio burden shall be analyzed as per GTP QC/GPT017. Use the pre validated analytical method for analysis.

6.5.2.2 **Swab sample:** Analyze the swab sample by suitable HPLC method for residue of previous product in the sample. BET shall be analyzed as per GTP QC/GPT015 and bio burden shall be analyzed as per GTP QC/GPT017. Use the pre validated analytical method for analysis.

☞ *Particulate matters to be check firstly, if it is complying then proceed for next analysis.*

☞ *For sampling and dispensing different accessories to be used for each different product.*



**7. ACCEPTANCE CRITERIA:**

---

- 7.1 Should be visually clean.
- 7.2 Should be free from visible particle and sub visible particle  $\geq 10\mu$ : NMT 25/ml  
 $\geq 25\mu$ : NMT 3/ml
- 7.3 Residue of previous product in next product should be not more than 10 ppm (0.001%).
- 7.4 The acceptance limit for microbial load is as follows:
  - Bacterial Endotoxin Test shall be  $< 0.125$  EU/ml.
  - Bio burden shall be less than  $< 10$  CFU / 100ml.

**8. REVALIDATION CRITERIA:**

---

**Process to be revalidated in any of the following case:**

- 8.1 Equipment change with different design specification.
- 8.2 Major modifications / changes to the equipments.
- 8.3 Change in the cleaning procedure.
- 8.4 Change in the manufacturing process.
- 8.5 Addition of new molecule.
- 8.6 MACO value decreased.
- 8.7 **Periodic verification:**
  - 8.7.1 Periodic verification shall be done once in a year ( $\pm 30$  days) with product, precisely the worst case product to cover all the equipments for those equipment on which worst case product is not taken, consider the worst case product from the group of products sharing that equipment.
    - ☞ *Periodic verification may consider the next worst product, in case the actual worst case has not been planned during the review period.*
    - ☞ *Periodic verification shall be carried out as per approved protocol.*

**9. OBSERVED NON-CONFORMANCE (IF ANY):**

---

All the non-conformance shall be addressed as per “Reporting and Monitoring of Process Non-Conformance in the Automated Quality Management System Software (SOP: CQA-CP-GEN-014)”.

**10. VALIDATION REPORT:**

---

On completion of the cleaning validation, evaluation of the test results shall be carried out for achievement of the acceptance criteria. A summary report shall be prepared clearly stating the outcome of the cleaning validation against the predetermined acceptance criteria.

**11. ABBREVIATIONS:**





**PHARMA DEVILS**

**PROTOCOL FOR CLEANING VALIDATION -  
PRODUCTION BULK  
(BLENDING)**

**PROTOCOL No.:**

- 11.1 MACO : Maximum Allowable Carry Over
- 11.2 MAR : Maximum Allowable Residue
- 11.3 HPLC : High Performance Liquid Chromatography
- 11.4 NMT : Not More Than
- 11.5 BET : Bacterial Endotoxin Test
- 11.6 TDD : Therapeutic Daily Dose
- 11.7 MBS : Minimum Batch Size
- 11.8 SF : Safety Factor
- 11.9 cm<sup>2</sup> : Centimeter Square
- 11.10 CFU : Colony Forming Unit
- 11.11 EU : Endotoxin Unit
- 11.12 PPM : Parts Per Million
- 11.13 LOQ : Limit of Quantitation
- 11.14 LOD : Limit of Detection
- 11.15 WFI : Water for Injection
- 11.16 API : Active Pharmaceutical Ingredient

**12. LIST OF ANNEXURE:**

Annexure No.	Annexure Title
01	Sampling location of equipment
02	Visual inspection for cleaning validation / verification
03	Equipment / equipment parts cleaning record

**13. REFERENCE DOCUMENT (IF ANY):**

- 13.1 Management of Validation/Qualification documents in DMS:CQA-CP-GEN-062
- 13.2 Quality Policy: CQA-CP-GEN-040
- 13.3 Cleaning Validation Programme: AN-QA-GEN-0030
- 13.4 Preparation of Validation and Qualification Protocol and its Control: AN-QA-GEN-0008

**14. REVISION HISTORY:**

Superseded Protocol		S. No.	Step No.	Changes made
Protocol No. / Version No.	Effective Date			



**PHARMA DEVILS**

**PROTOCOL FOR CLEANING VALIDATION -  
PRODUCTION BULK  
(BLENDING)**

**PROTOCOL No.:**

<b>Superseded Protocol</b>		<b>S. No.</b>	<b>Step No.</b>	<b>Changes made</b>
<b>Protocol No. / Version No.</b>	<b>Effective Date</b>			
CLV/QA/003/ Version No. 02	17 <sup>th</sup> Nov. 2014	Changes included as per change request: CRF/QA17/046		
		1	NA	Editorial changes as per SOP - AN-QA-GEN-0030
		2	8.7	Periodic verification is included



**PHARMA DEVILS**

**PROTOCOL FOR CLEANING VALIDATION -  
PRODUCTION BULK  
(BLENDING)**

**PROTOCOL No.:**

**Annexure – 01: Sampling location of equipments**

<b>Equipment Name:</b>	Blending Machine	<b>ID No.:</b>	PB/BLM-01
<b>Location Photo</b>		<b>Location Identification &amp; rational for selection</b>	
		<p>L-01 (Bottom surface of bin): Hard to clean area and physical verification is not feasible.</p> <p>L-02 (Bottom surface of bin): Hard to clean area and physical verification is not feasible.</p> <p>L-03 (Outlet of bin): Hard to clean area and physical verification is not feasible.</p> <p>L-04 (Distance piece): Hard to clean area and physical verification is not feasible.</p> <p>L-05 (Screw feeder): Hard to clean, maximum chance to deposition of product on mesh pores.</p>	

<b>Done by (Signature &amp; Date)</b>	<b>Checked by (Signature &amp; Date)</b>

**Annexure – 02: Visual inspection for cleaning validation / verification**

S. No.	Product Name/ Batch No.	Equipment Name	Equipment ID No.	Part Inspected	Visually clean (Yes/No)	Checked by Sign/Date



