



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

Cleaning Validation of Lyophilizer

**CLEANING VALIDATION
PROTOCOL
OF
LYOPHILIZER**



Cleaning Validation of Lyophilizer

TABLE OF CONTENTS

1.0 APPROVAL 3

2.0 OVERVIEW:..... 4

 2.1 INTRODUCTION:..... 4

 2.2 OBJECTIVE: 4

3.0 SCOPE AND PURPOSE 4

 3.1 SCOPE..... 4

 3.2 PURPOSE..... 4

5.0 CLEANING VALIDATION APPROACH..... 4

6.0 CLEANING VALIDATION SAMPLE TEST AND LIMIT 5

7.0 DEVIATION: 5

8.0 VALIDATION ACCEPTANCE CRITERIA 5

9.0 REVALIDATION CRITERIA 5

10.0 SUMMARY AND CONCLUSION:..... 5

11.0 REFERENCE DOCUMENT (IF ANY): 5

12.0 LIST OF ANNEXURE: 6



Cleaning Validation of Lyophilizer

1.0 APPROVAL

The cleaning validation protocol is prepared by the Validation team for the project “Integrated Sterile Bulk and Formulations Facility”, under the authority QA Head. Hence this document before being effective shall be approved by the QA Head.

PREPARED BY		
NAME/ FUNCTIONAL AREA	DESIGNATION	SIGNATURE /DATE
Validation & QA		

CHECKED BY		
NAME/ FUNCTIONAL AREA	DESIGNATION	SIGNATURE /DATE
Production		
Quality Control		
Quality Assurance		

APPROVED BY		
NAME/ FUNCTIONAL AREA	DESIGNATION	SIGNATURE /DATE
Quality Assurance – Head		



Cleaning Validation of Lyophilizer

2.0 OVERVIEW:

2.1 Introduction:

The purpose of this protocol is to establish that the Cleaning procedure for Lyophilizer.

2.2 Objective:

The objective of the cleaning validation protocol is to establish & assure with documented evidence that the cleaning procedure for Lyophilizer, remove residue of the riboflavin to levels below the established acceptance limits and verify the visually for removal of Riboflavin.

3.0 SCOPE AND PURPOSE:

3.1 Scope:

This CLV is prepared to describe the validation requirement and plan for:

- All cleaning procedures used for cleaning of lyophilizer, which are used in processing of Sterile API and lyophilized vial.

3.2 Purpose:

Cleaning Validation Protocol (CLV) serves to specify and to co-ordinate the validation activity required for all cleaning processes which are to be followed for lyophilizer to be used for processing of sterile API product and lyophilized vial.

CVP will serve to ensure that cleaning procedures used for lyophilizer at Integrated Sterile Bulk and Formulation Facility (which may adversely affect product quality if found to not efficiently clean the equipment as per the requirement) are validated as per cGMP requirements.

The PVP is subordinate to the Validation Master Plan (VMP) of Integrated Sterile Bulk and Formulation Facility of Rajasthan Antibiotics Limited, Bhiwadi.

4.0 CLEANING VALIDATION APPROACH:

This cleaning validation protocol is prepared based on project cleaning validation plan or originated based on the guidance & instruction of VMP. Lyophilizer has to be cleaned before the start of a production process. The cleaning is done to remove the traces of the product.

4.1 Execution Method:

Spraying the riboflavin (0.02% w/v) solution on shelf and inner surface of lyophilizer, after that kept or drying then start the cleaning procedure as per SOP and collect the sample and send to Quality control.



Cleaning Validation of Lyophilizer

S.No.	Name of Sample	Test Required	Limit
1.	Rinse Sample	pH	5.0 to 7.0
		Previous Product Content	NMT 10 PPM
		Conductivity	NMT 1.3 $\mu\text{s/cm}$
		Particulate matter	
		Bioburden	10 cfu/ 100 ml
		BET	0.25 Eu/ml
2.	Check the visually for removal of riboflavin		Not seen

5.0 CLEANING VALIDATION SAMPLE TEST AND LIMIT:

6.0 DEVIATION: There is any deviation observed during the execution of cleaning validation of lyophilizer record as per SOP "Handling of Deviation" and attach with validation report.

7.0 VALIDATION ACCEPTANCE CRITERIA: All the parameters are meeting with the predefined specifications.

8.0 REVALIDATION CRITERIA:

Change in Manufacturing Equipments and Introduction of New Products.

9.0 SUMMARY AND CONCLUSION:

Summary and Conclusion concluded on the basis of three cleaning validation batches sample analysis results.

10.0 REFERENCE DOCUMENT (IF ANY):

S.No.	Document Name
1.	Cleaning Validation Plan
2.	Validation Master Plan (VMP)



Cleaning Validation of Lyophilizer

11.0 LIST OF ANNEXURE:

Annexure No.	Annexure Title
01	Training Record
02	Cleaning Validation – Sampling Plan
03	List of Equipment