



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

PERFORMANCE QUALIFICATION

FOR

COMPUTER SYSTEM

OF

QC (Viscometer)

System Name	QC (Viscometer)	
System ID		
Location	Instrument Lab	
Effective Date		

Document No.:

Page 1 of 15



QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

TABLE OF CONTENT

1.	PRE APPROVALS	3
2.	SIGNATURE OF EXECUTOR	4
3.	REVISION HISTORY	4
4.	OBJECTIVE	5
5.	SCOPE	5
6.	SYSTEM DESCRIPTION	
7.	ROLE AND RESPONSIBILITY	6
8.	REFERENCES	7
9.	DOCUMENTATION PROCEDURE	8
10.	QUALIFICATION COMPLETION AND APPROVAL	8
11.	ACCEPTANCE CRITERIA	8
12.	PERFORMANCE VERIFICATION TEST	9
1	2.1 Verification of Control Loops Test	9
13.	PERFORMANCE QUALIFICATION TEST STATUS	11
14.	DISCREPANCIES HANDLING DURING COMPUTER QUALIFICATION	11
	DISCREPANCY AND CORRECTIVE ACTION FORM	
16.	ABBREVIATION	13
17.	ATTACHMENT SUMMARY	14
18.	PERFORMANCE QUALIFICATION SUMMARY & CONCLUSION	14
19.	POST APPROVALS	15



QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

1. PRE APPROVALS:

The signature listed below indicates the pre-approval of this performance qualification. This approval is joint responsibility of listed functional areas.

	DOCUMENT DEVELOPMENT	SIGN / DATE
Name	:	
Designation	:	

DOCUMENT REVIEW AND APPROVAL (M/S)			
Sign / Date	:		
Name	:		
Designation	:		
	Engineering		
Sign / Date	:		
Name	:		
Designation	:		
	IT		
Sign / Date	:		
Name	:		
Designation	:		
	Quality Control		

	DOCUMENT AUTHORIZATION (M/S)			
Sign / Date	:			
Name	:			
Designation	:			
	Quality Assurance			



QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

2. SIGNATURE OF EXECUTOR:

All the executer involved in this document have to sign within prescribed format given below.

M/s

Name	Designation	Signature	Initial	Date

M/s

Name	Designation	Signature	Initial	Date

3. REVISION HISTORY:

Date	Supersedes	Reason for Revision



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

4. OBJECTIVE:

5. SCOPE:

This document is applicable to validation of Hardware and Software of computer system installed and performs at M/s. Quality control department. This system requirement specification shall define the documentation, references and acceptance criteria to establish that the validation of Hardware and Software of Computer system after modification is installed, Operated and performs in accordance with the guidelines laid down by the manufacturer of the system.

6. SYSTEM DESCRIPTION:

Computer system of QC_VISCO defines the controlling of analytical instrument connected to the system. It controlled the process of instrument connected to it by the software installed in the CS. The operator interface is carried out by CS screen. The CS is used to feed required parameters and set points in the system during operation. The system is connected to data server for printing and data backup. The system is secured by IT through password only within the system.





QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

7. ROLE AND RESPONSIBILITY:

The validation team comprising of representative from each of the following departments should be responsible for overall compliance with this validation plan.

Department	Responsibilities
Validation Agency	To collect the necessary data for performance qualification activities.
()	> To prepare and execute the performance qualification in coordination with
	engineering, validation and quality assurance team.
	Comply with regulatory / Guidelines / Standards / validation plan requirements
	throughout the validation life cycle.
	> To submit performance qualification for approval.
Engineering	> To provide the necessary data for performance qualification activities.
(M /s)	> To review the performance qualification.
IT	> To provide the necessary data for performance qualification activities.
(M /s)	> To review the performance qualification.
Quality Control (M/s)	> To provide the necessary data for performance qualification activities.
()	> To review the performance qualification.
Quality Assurance (M/s)	> To approve and authorized the performance qualification.





QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

8. REFERENCES:

The publication listed below form part of this reference documents. Each publication shall have latest revision in effect on the date of this document is approved for execution.

GAMP 5	Good Automated Manufacturing Practices, Version 5, Guideline Document for Automated Systems from International Society of Pharmaceutical Engineering
21 Code of Federal	Current Good Manufacturing Practice in Manufacturing,
Regulations (CFR), Part 210	Processing, Packing, or Holding off Drugs; General
21 Code of Federal Regulations (CFR), Part 211	Current Good Manufacturing Practice for finished Pharmaceuticals
21 Code of Federal	21 Code of Federal Regulations (CFR), Part 11
Regulations (CFR), Part	Electronic Records, Electronic Signatures, Final Rule Electronic
11	Submissions; Establishment of Public Docket, Notice
ICH Q9	International Conference of Harmonization (ICH) quality risk assessment Q9
EU GMP	Laying down the principles and guidelines of GMP in respect of medicinal products for human use.
WHO	Appendix 5, validation of computerized systems.



PHARMA DEVILS QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

9. DOCUMENTATION PROCEDURE:

- Qualification activities will be performed as defined in the approved document.
- All documentation will be completed during the execution of the qualification.
- Recording of information will be made in permanent ink.
- Fill out complete information in the verification table provided.
- Do not keep any space blank. Mark blank space with a single line throughout the appropriate space with mentioning NA (Not Applicable) and put initial and date.
- Correct the mistakes by drawing a single line through the incorrect data, recording the correct information and then initial sign and date the change.

10. QUALIFICATION COMPLETION AND APPROVAL:

- Verify that all tests required by qualification are completed and attached.
- Verify that all amendments and discrepancies are documented, approved and attached.
- If all items in the qualification for the Computer System Based system (QC_VISCO) have been reviewed and found to be acceptable, sign the corresponding block in the qualification completion and approval form.

11. ACCEPTANCE CRITERIA:

- Installation completion as per manufacturer's recommendations & cGMP requirements.
- Installation of major components as per the design specifications.
- The supply of all necessary documentation from manufacturer.
- The system is operating as intended and is under state of control.
- Performance features meet system requirements and system specifications.



QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

12. PERFORMANCE VERIFICATION TEST:

12.1 Verification of Control Loops Test:

Objective	:	To Verify the performance of Process.
Tools Required	:	Not Applicable
Procedure	:	1. Start the equipment in normally.
		2. Login with user level id.
		3. Set require method/sequences for the test.
		4. Start process and observe the set process method.
		5. If printing facility available, attached the printout of whole integrated
		control loop test.
Acceptance	:	Computer system should able to control the set process method within the
Criteria		specified limit.

Verification Table:

Description	Specified	Actual result	Discrepancy? (Y/N)
Check all pre- requirement condition for system start and set the required method.	All condition should be healthly and method should be set in within range.		
Login with user level id.	Login shall be Successfully.		
Set the method/sequences.	Machine should start and control process of method/sequences.		
Start the method/sequences.	Observe the set process.		
Take the print report of control loop test	Printed and set method/sequences should be same and attached the		

Document No.:



PHARMA DEVILS QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

Description	Spe	cified		Actual result	Discrepancy? (Y/N)
	printout in attachment.				
Remarks:					
Meet the acceptanc	e Criteria [] Yes [] No	Reference Att	achment No. []
Checked by :				Date:	
Verified by :				Date:	





QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

13. PERFORMANCE QUALIFICATION TEST STATUS:

The performance qualification test status is as per below mentioned table.

Test Description	Status (Pass / Fail)	Discrepancy? (Y/N)
Verification of Control Loops Test		

14. DISCREPANCIES HANDLING DURING COMPUTER QUALIFICATION:

- In case of discrepancy observed during qualification, document in the defined column in each table and document the details of the observation in the discrepancy log sheet.
- Inform to engineering IT quality control and quality assurance about discrepancy.
- Investigate the discrepancy and ensure the possible impact.
- If discrepancy does not have potential to impact on operation as well as performance of the system, close the discrepancy with proper justification.
- The engineering IT QC and QA will decide whether discrepancy is acceptable or not.
- If discrepancy is acceptable, provide conclusion and recommendation if any into respective column.



QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

15. DISCREPANCY AND CORRECTIVE ACTION FORM:

Protocol Reference	
Discrepancy Number	

DISCREPANCY:

Describe the Discrepancy	
Reported by	Date

CORRECTIVE ACTION:

Describe corrective action taken (Attach additional sheets if necessary)		
Reported by	Date	

DISPOSITION ACTION:

Acceptable?	Yes	No		
Discussion				
Approved by			Date	

OMPLETION:

Completed by	Date
--------------	------

Document No.:



QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

16. ABBREVIATION:

Abbreviations	Description	
GMP	Good Manufacturing Practices	
SRS	System Requirement and Specification	
IQ	Installation Qualification	
OQ	Operation Qualification	
PQ	Performance Qualification	
QA	Quality Assurance	
SOP	Standard Operating Procedure	
NA	Not Applicable	
ICH	International Conference of Harmonization	
VAC	Alternate Current Voltage	
VDC	Direct Current Voltage	
WHO	World Health Organization	



QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

17. ATTACHMENT SUMMARY:

Attachment No.	Description

18. PERFORMANCE QUALIFICATION SUMMARY & CONCLUSION:

Compiled by: _____

Date:

Document No.:

Page 14 of 15



QUALITY ASSURANCE DEPARTMENT

PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

19. POST APPROVALS:

The signature listed below indicates the post approval of this performance qualification. This approval is joint responsibility of listed functional areas.

	DOCUMENT DEVELOPMENT	SIGN / DATE
Name	:	
Designation	:	

	DOCUMENT REV	TEW AND APPROVAL (M/S)
Sign / Date	•	
Name	:	
Designation	:	
	Engineering	
Sign / Date	:	
Name	:	
Designation	:	
	IT	
Sign / Date	:	
Name	:	
Designation	:	
	Quality control	

	DOCUMENT AUTHORIZATION (M/S)
Sign / Date	:
Name	:
Designation	:
	Quality Assurance

Document No.:

Page 15 of 15