



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

**PERFORMANCE QUALIFICATION
FOR
COMPUTER SYSTEM
OF
GAS CHROMATOGRAPHY (PC)**

System Name	GAS CHROMATOGRAPHY (PC)
System ID	
Location	Instrument Lab
Effective Date	



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

TABLE OF CONTENT

1. PRE APPROVALS	3
2. SIGNATURE OF EXECUTOR.....	4
3. REVISION HISTORY.....	4
4. OBJECTIVE	5
5. SCOPE.....	5
6. SYSTEM DESCRIPTION	5
7. ROLE AND RESPONSIBILITY.....	6
8. REFERENCES.....	7
9. DOCUMENTATION PROCEDURE.....	8
10. QUALIFICATION COMPLETION AND APPROVAL	8
11. ACCEPTANCE CRITERIA.....	8
12. PERFORMANCE VERIFICATION TEST.....	9
12.1 Verification of Control Loops Test.....	9
13. PERFORMANCE QUALIFICATION TEST STATUS.....	11
14. DISCREPANCIES HANDLING DURING COMPUTER QUALIFICATION	11
15. DISCREPANCY AND CORRECTIVE ACTION FORM	12
16. ABBREVIATION	13
17. ATTACHMENT SUMMARY	14
18. PERFORMANCE QUALIFICATION SUMMARY & CONCLUSION.....	14
19. POST APPROVALS	15



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

1. PRE-APPROVALS:

The signature listed below indicates the pre-approval of this performance qualification. This approval is joint responsibility of listed functional areas.

DOCUMENT DEVELOPMENT		SIGN / DATE
Name :	_____	
Designation :	_____	

DOCUMENT REVIEW AND APPROVAL (M/S.....)	
Sign / Date :	_____
Name :	_____
Designation :	_____
	Engineering
Sign / Date :	_____
Name :	_____
Designation :	_____
	IT
Sign / Date :	_____
Name :	_____
Designation :	_____
	Quality Control

DOCUMENT AUTHORIZATION (M/S.....)	
Sign / Date :	_____
Name :	_____
Designation :	_____
	Quality Assurance



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

2. SIGNATURE OF EXECUTOR:

All the executor involved in this document have to sign within prescribed format given below.

M/s

Name	Designation	Signature	Initial	Date

M/s

Name	Designation	Signature	Initial	Date

3. REVISION HISTORY:

Date	Supersedes	Reason for Revision



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

4. OBJECTIVE:

The objective of the system requirement specification defines the requirements of hardware software, functions and documentation for the Computer System (GC_PC) installed, and performs at M/s. This document will be cross-referenced with the relevant system specification to ensure that the system specifications have been met according to M/s. requirements.

5. SCOPE:

This document is applicable to validation of Hardware and Software of computer system installed and performs at M/s. Quality control department. This system requirement specification shall define the documentation, references and acceptance criteria to establish that the validation of Hardware and Software of Computer system after modification is installed, Operated and performs in accordance with the guidelines laid down by the manufacturer of the system.

6. SYSTEM DESCRIPTION:

Computer system of GC_PC defines the controlling of analytical instrument connected to the system. It controlled the process of instrument connected to it by the software installed in the CS. The operator interface is carried out by CS screen. The CS is used to feed required parameters and set points in the system during operation. The system is connected to data server for printing and data backup. The system is secured by IT through password only within the system.



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

7. ROLE AND RESPONSIBILITY:

The validation team comprising of representative from each of the following departments should be responsible for overall compliance with this validation plan.

Department	Responsibilities
Validation Agency (.....)	<ul style="list-style-type: none">➤ To collect the necessary data for performance qualification activities.➤ To prepare and execute the performance qualification in coordination with engineering, validation and quality assurance team.➤ Comply with regulatory / Guidelines / Standards / validation plan requirements throughout the validation life cycle.➤ To submit performance qualification for approval.
Engineering (M/s.)	<ul style="list-style-type: none">➤ To provide the necessary data for performance qualification activities.➤ To review the performance qualification.
IT (M/s.)	<ul style="list-style-type: none">➤ To provide the necessary data for performance qualification activities.➤ To review the performance qualification.
Quality Control (M/s.)	<ul style="list-style-type: none">➤ To provide the necessary data for performance qualification activities.➤ To review the performance qualification.
Quality Assurance (M/s.)	<ul style="list-style-type: none">➤ To approve and authorized the performance qualification.



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

8. REFERENCES:

The publication listed below form part of this reference documents. Each publication shall have latest revision in effect on the date of this document is approved for execution.

GAMP 5	Good Automated Manufacturing Practices, Version 5, Guideline Document for Automated Systems from International Society of Pharmaceutical Engineering
21 Code of Federal Regulations (CFR), Part 210	Current Good Manufacturing Practice in Manufacturing, Processing, Packing, or Holding off Drugs; General
21 Code of Federal Regulations (CFR), Part 211	Current Good Manufacturing Practice for finished Pharmaceuticals
21 Code of Federal Regulations (CFR), Part 11	21 Code of Federal Regulations (CFR), Part 11 Electronic Records, Electronic Signatures, Final Rule Electronic Submissions; Establishment of Public Docket, Notice
ICH Q9	International Conference of Harmonization (ICH) quality risk assessment Q9
EU GMP	Laying down the principles and guidelines of GMP in respect of medicinal products for human use.
WHO	Appendix 5, validation of computerized systems.



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

9. DOCUMENTATION PROCEDURE:

- Qualification activities will be performed as defined in the approved document.
- All documentation will be completed during the execution of the qualification.
- Recording of information will be made in permanent ink.
- Fill out complete information in the verification table provided.
- Do not keep any space blank. Mark blank space with a single line throughout the appropriate space with mentioning NA (Not Applicable) and put initial and date.
- Correct the mistakes by drawing a single line through the incorrect data, recording the correct information and then initial sign and date the change.

10. QUALIFICATION COMPLETION AND APPROVAL:

- Verify that all tests required by qualification are completed and attached.
- Verify that all amendments and discrepancies are documented, approved and attached.
- If all items in the qualification for the Computer System Based system (GC_PC) have been reviewed and found to be acceptable, sign the corresponding block in the qualification completion and approval form.

11. ACCEPTANCE CRITERIA:

- Installation completion as per manufacturer's recommendations & cGMP requirements.
- Installation of major components as per the design specifications.
- The supply of all necessary documentation from manufacturer.
- The system is operating as intended and is under state of control.
- Performance features meet system requirements and system specifications.



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

12. PERFORMANCE VERIFICATION TEST:

12.1 Verification of Control Loops Test

Objective : To Verify the performance of Process.

Tools Required : Not Applicable

Procedure : 1. Start the equipment in normally.
2. Login with user level id.
3. Set require method/sequences for the test.
4. Start process and observe the set process method.
5. If printing facility available, attached the printout of whole integrated control loop test.

Acceptance Criteria : Computer system should able to control the set process method within the specified limit.

Verification Table:

Description	Specified	Actual result	Discrepancy? (Y/N)
Check all pre- requirement condition for system start and set the required method.	All condition should be healthy and method should be set in within range.		
Login with user level id.	Login shall be Successfully.		
Set the method/sequences.	Machine should start and control process of method/sequences.		
Start the method/sequences.	Observe the set process.		
Take the print report of control loop test	Printed and set method/sequences should be same and attached the printout		



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

Description	Specified	Actual result	Discrepancy? (Y/N)
	in attachment.		

Remarks:

Meet the acceptance Criteria [] Yes [] No Reference Attachment No. []

Checked by : _____

Date: _____

Verified by : _____

Date: _____



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

13. PERFORMANCE QUALIFICATION TEST STATUS:

The performance qualification test status is as per below mentioned table.

Test Description	Status (Pass / Fail)	Discrepancy? (Y/N)
Verification of Control Loops Test		

14. DISCREPANCIES HANDLING DURING COMPUTER QUALIFICATION:

- In case of discrepancy observed during qualification, document in the defined column in each table and document the details of the observation in the discrepancy log sheet.
- Inform to engineering IT quality control and quality assurance about discrepancy.
- Investigate the discrepancy and ensure the possible impact.
- If discrepancy does not have potential to impact on operation as well as performance of the system, close the discrepancy with proper justification.
- The engineering IT QC and QA will decide whether discrepancy is acceptable or not.
- If discrepancy is acceptable, provide conclusion and recommendation if any into respective column.



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

15. DISCREPANCY AND CORRECTIVE ACTION FORM:

Protocol Reference	
Discrepancy Number	

DISCREPANCY:

Describe the Discrepancy	
Reported by	Date

CORRECTIVE ACTION:

Describe corrective action taken (Attach additional sheets if necessary)	
Reported by	Date

DISPOSITION ACTION :

Acceptable?	Yes	No
Discussion		
Approved by	Date	

COMPLETION:

Completed by	Date
--------------	------



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

16. ABBREVIATION:

Abbreviations	Description
GMP	Good Manufacturing Practices
SRS	System Requirement and Specification
IQ	Installation Qualification
OQ	Operation Qualification
PQ	Performance Qualification
QA	Quality Assurance
SOP	Standard Operating Procedure
NA	Not Applicable
ICH	International Conference of Harmonization
VAC	Alternate Current Voltage
VDC	Direct Current Voltage
WHO	World Health Organization



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

17. ATTACHMENT SUMMARY:

Attachment No.	Description

18. PERFORMANCE QUALIFICATION SUMMARY & CONCLUSION:

Compiled by: _____

Date: _____



PERFORMANCE QUALIFICATION FOR COMPUTER SYSTEM

System Name: GAS CHROMATOGRAPHY (PC)

System ID:

19. POST APPROVALS:

The signature listed below indicates the post approval of this performance qualification. This approval is joint responsibility of listed functional areas.

DOCUMENT DEVELOPMENT	SIGN / DATE
Name : _____ Designation : _____	

DOCUMENT REVIEW AND APPROVAL (M/S.....)
Sign / Date : _____ Name : _____ Designation : _____ Engineering
Sign / Date : _____ Name : _____ Designation : _____ IT
Sign / Date : _____ Name : _____ Designation : _____ Quality control

DOCUMENT AUTHORIZATION (M/S.....)
Sign / Date : _____ Name : _____ Designation : _____ Quality Assurance