



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

**INSTALLATION QUALIFICATION
FOR
COMPUTER SYSTEM OF
HPLC**

| | |
|-----------------------|-----------------------|
| System Name | HPLC |
| System ID | |
| Location | Instrument Lab |
| Effective Date | |



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

TABLE OF CONTENTS

| | |
|--|-----------|
| 1. PRE APPROVALS | 3 |
| 2. SIGNATURE OF EXECUTOR..... | 4 |
| 3. REVISION HISTORY | 4 |
| 4. OBJECTIVE | 5 |
| 5. SCOPE..... | 5 |
| 6. SYSTEM DESCRIPTION | 5 |
| 7. ROLE AND RESPONSIBILITY..... | 6 |
| 8. REFERENCES..... | 7 |
| 9. DOCUMENTATION PROCEDURE..... | 7 |
| 10. QUALIFICATION COMPLETION AND APPROVAL | 8 |
| 11. ACCEPTANCE CRITERIA..... | 8 |
| 12. INSTALLATION VERIFICATION TEST | 9 |
| 12.1 Identification of System Details | 9 |
| 12.2 Verification of Master Documents | 10 |
| 12.3 Verification of Capacity Requirement | 11 |
| 12.4 Verification of Hardware Components | 12 |
| 12.5 Verification of Software Components | 14 |
| 12.6 Verification of Physical and Logical Security Control..... | 15 |
| 12.7 Verification of Test Instruments Calibration and Traceability | 17 |
| 12.8 Verification of Power Utility..... | 18 |
| 12.9 Verification of Environment Condition | 19 |
| 12.10 Verification of Communication Link Between Server To Computer System | 20 |
| 12.11 Verification of General System Installation | 21 |
| 12.12 Verification of Standard Operating Procedure | 22 |
| 13. INSTALLATION QUALIFICATION TEST STATUS | 23 |
| 14. DISCREPANCIES HANDLING DURING COMPUTER SYSTEM QUALIFICATION | 24 |
| 15. DISCREPANCY AND CORRECTIVE ACTION FORM | 25 |
| 16. ABBREVIATION | 26 |
| 17. ATTACHMENT SUMMARY | 27 |
| 18. INSTALLATION QUALIFICATION SUMMARY & CONCLUSION | 27 |
| 19. POST APPROVALS..... | 28 |



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

1. PRE-APPROVALS:

The signature listed below indicates the preapproval of this installation qualification. This approval is joint responsibility of listed functional areas.

| DOCUMENT DEVELOPMENT | SIGN / DATE |
|---|-------------|
| Name : _____ Designation : _____ | |

DOCUMENT REVIEW AND APPROVAL (M/S)

| |
|---|
| Sign / Date : _____ Name : _____ Designation : _____ Engineering |
| Sign / Date : _____ Name : _____ Designation : _____ IT |
| Sign / Date : _____ Name : _____ Designation : _____ Quality Control |

DOCUMENT APPROVAL (M/S)

| |
|---|
| Sign / Date : _____ Name : _____ Designation : _____ Quality Assurance |
|---|



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

2. SIGNATURE OF EXECUTOR:

All the executor involved in this document have to sign within prescribed format given below.

M/s

| Name | Designation | Signature | Initial | Date |
|------|-------------|-----------|---------|------|
| | | | | |
| | | | | |

M/s

| Name | Designation | Signature | Initial | Date |
|------|-------------|-----------|---------|------|
| | | | | |
| | | | | |

3. REVISION HISTORY:

| Date | Supersedes | Reason for Revision |
|------|------------|---------------------|
| | | |
| | | |



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

4. OBJECTIVE:

The objective of the system requirement specification defines the requirements of hardware software, functions and documentation for the Computer System (HPLC) installed at M/s. This document will be cross-referenced with the relevant system specification to ensure that the system specifications have been met according to M/s. requirements.

5. SCOPE:

This document is applicable to validation for Hardware and Software system of HPLC installed after modification at M/s..... This installation qualification shall define the documentation, references and acceptance criteria for validation of Hardware and Software system of HPLC is installed in accordance with the guidelines laid down by the manufacturer of the system.

6. SYSTEM DESCRIPTION:

Computer system of HPLC defines the controlling of analytical instrument connected to the system. It controlled the process of instrument connected to it by the software installed in the CS. The operator interface is carried out by CS screen. The CS is used to feed required parameters and set points in the system during operation. The system is connected to data server for printing and data backup. The system is secured by IT through password only within the system.



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

7. ROLE AND RESPONSIBILITY:

The validation team comprising of representative from each of the following departments should be responsible for overall compliance with this validation plan.

| Department | Responsibilities |
|---|---|
| Validation Agency (.....) | <ul style="list-style-type: none">➤ To collect the necessary data for installation qualification activities.➤ To prepare and execute the installation qualification in coordination with engineering, validation and quality assurance team.➤ Comply with regulatory / Guidelines / Standards / validation plan requirements throughout the validation life cycle.➤ To submit installation qualification for approval. |
| Engineering (M/s.) | <ul style="list-style-type: none">➤ To provide the necessary data for installation qualification activities.➤ To review the installation qualification. |
| IT (M/s.) | <ul style="list-style-type: none">➤ To provide the necessary data for installation qualification activities.➤ To review the installation qualification. |
| Quality Control (M/s.) | <ul style="list-style-type: none">➤ To provide the necessary data for installation qualification activities.➤ To review the installation qualification. |
| Quality Assurance (M/s.) | <ul style="list-style-type: none">➤ To approve and authorized the installation qualification. |



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

8. REFERENCES:

The publication listed below form part of this reference documents. Each publication shall have latest revision in effect on the date of this document is approved for execution.

| | |
|---|--|
| GAMP 5 | Good Automated Manufacturing Practices, Version 5, Guideline Document for Automated Systems from International Society of Pharmaceutical Engineering |
| 21 Code of Federal Regulations (CFR), Part 210 | Current Good Manufacturing Practice in Manufacturing, Processing, Packing, or Holding off Drugs; General |
| 21 Code of Federal Regulations (CFR), Part 211 | Current Good Manufacturing Practice for finished Pharmaceuticals |
| 21 Code of Federal Regulations (CFR), Part 11 | 21 Code of Federal Regulations (CFR), Part 11 Electronic Records, Electronic Signatures, Final Rule Electronic Submissions; Establishment of Public Docket, Notice |
| ICH Q9 | International Conference of Harmonization (ICH) quality risk assessment Q9 |
| EU GMP | Laying down the principles and guidelines of GMP in respect of medicinal products for human use. |
| SRS | System Requirement Specification |
| WHO | Appendix 5, validation of computerized systems. |

9. DOCUMENTATION PROCEDURE:

- Qualification activities will be performed as defined in the approved document.
- All documentation will be completed during the execution of the qualification.
- Recording of information will be made in permanent ink.
- Fill out complete information in the verification table provided.
- Do not keep any space blank. Mark blank space with a single line throughout the appropriate space with mentioning NA (Not Applicable) and put initial and date.
- Correct the mistakes by drawing a single line through the incorrect data, recording the correct information and then initial sign and date the change.



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

10. QUALIFICATION COMPLETION AND APPROVAL:

- Verify that all tests required by qualification are completed and attached.
- Verify that all amendments and discrepancies are documented, approved and attached.
- If all items in the qualification for the validation to Hardware and Software system of HPLC been reviewed and found to be acceptable, sign the corresponding block in the qualification completion and approval form.

11. ACCEPTANCE CRITERIA:

- Installation completion as per manufacturer's recommendations & cGMP requirements.
- Installation of major components as per the design specifications.
- The supply of all necessary documentation from manufacturer.



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12. INSTALLATION VERIFICATION TEST:

12.1 Identification of System Details:

Objective : This test sheet is intended to verification of equipment details.

Tools Required : Not Applicable

Procedure : 1. Record Equipment Name
2. Record Identification No.
3. Record Equipment Location

Acceptance : Data recorded from the equipment shall match with the data specified in

Criteria verification table.

Verification Table:

| Equipment Details | Specified As | As observed | Discrepancy? (Y/N) |
|--------------------|----------------|-------------|-----------------------|
| Equipment Name | HPLC | | |
| Identification No. | | | |
| Location | Instrument Lab | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.2 Verification of Master Documents:

Objective : To verify the availability of related master documents.

Tools Required : Not Applicable

Procedure : 1. Verify Documents Name.
2. Verify Documents Reference.
3. Verify Documents Availability.

Acceptance : Documents should be available.

Criteria

Verification Table:

| Documents Name | Documents Reference | Availability (Yes/No) | Verified (Yes/No) | Discrepancy? (Y/N) |
|--------------------|------------------------|-----------------------|-------------------|--------------------|
| SRS | | | | |
| Operational Manual | Refer attachment No. 1 | | | |
| BOM | Refer attachment No. 2 | | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date: _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.3 Verification of Capacity Requirement:

Objective : To verify the processing capacity of Computer System

Tools required : Not Applicable

Procedure : Physical verification of Capacity Requirement as per SRS

1. Not more than 50% of the installed hard disk capacity in PC components should be consumed by installed software.
2. Historical data storage capacity should allow for online retrieval of at forever of any historical data.

Acceptance criteria : 1. Capacity Requirement of the control system shall match with SRS.

Verification Table:

| S.No. | Item Name | Expected | Actual (Yes/No) | Discrepancy ? (Y/N) |
|-------|--------------------------|---|-----------------|---------------------|
| 1. | Local Electronic Storage | Not more than 50 % should be consumed by installed software | | |
| 2. | Historical data storage | Not more than 50 % should be consumed by historical data | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.4 Verification of Hardware Components:

Objective : To verify the installed hardware components as per the SRS.

Tools Required : Not Applicable

Procedure : 1. Verify Hardware Name.
2. Verify Hardware Make/ Assemble By
3. Verify Hardware Model No./Specification

Acceptance : Installed hardware component should match with SRS.

Criteria

Verification Table:

| Name | Make/ Assemble By | Model No./Specification | Qty. | Installation (Yes/No) | Discrepancy? (Y/N) |
|-----------------------|---------------------------------|----------------------------|------|--------------------------|-----------------------|
| Monitor | Acer | V196HQL | 01 | | |
| CPU | Acer | Veriton-IC6404 | 01 | | |
| Keyboard | Acer | NA | 01 | | |
| Mouse | Acer | NA | 01 | | |
| UPS | Emerson Network System | 20 kVA | 01 | | |
| RAM | Acer | 4 GB | 01 | | |
| Processor | Intel | I3 3.60 GHz | 01 | | |
| Printer | Canon | 3300 | 01 | | |
| Analytical Instrument | HPLC Agilent Technologies | 1260 Infinity II | 01 | | |



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

Verification Table: Communication Port:

| S.No. | Port Type | Qty. | Installation (Yes/No) | Discrepancy? (Y/N) |
|-------|-----------|------|-----------------------|--------------------|
| 1. | USB | 2 | | |
| 2. | Ethernet | 2 | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.5 Verification of Software Components:

Objective : To verify the installed software components as per the SRS.

Tools Required : Not Applicable

Procedure : 1. Verify Software Name.
2. Verify Software Version
3. Verify operating system.
4. Verify software backup availability.

Acceptance Criteria : Installed software component should match with SRS.

Verification Table:

A. For Software Components:

| S.No. | Software Name | Version | Installation (Yes/No) | Discrepancy? (Y/N) |
|-------|----------------|--------------|-----------------------|--------------------|
| 1. | Adobe Reader | 11.0.00 | | |
| 2. | SAM CoDeC Pack | 5.72 | | |
| 3. | Windows | 7 64 bit SP1 | | |

B. Operating system details:

| S.No. | Window | Product key/ Licence key | Discrepancy? (Y/N) |
|-------|-----------|--------------------------|--------------------|
| 1. | Windows 7 | | |

C. Software Backup Availability:

| S.No. | Available (Yes/ No) | Discrepancy? (Y/N) |
|-------|---------------------|--------------------|
| | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.6 Verification of Physical and Logical Security Control:

Objective : Verify the physical and logical security of Computer System.

Tools Required : Not Applicable

Procedure : 1. Verify physical Security.
2. Verify logical security of Application Window and Software
3. Verify User for Application access

Acceptance : Physical security should be maintained. Logical

Criteria

- Window Login password Should be available.
- Application Software should have multiple numbers of user's role with user name.

A. Verification Table for Physical Security:

| System | Security | Availability (Yes/No) | Discrepancy? (Y/ N) |
|-----------------|--------------------------------|-----------------------|---------------------|
| Computer system | PC should be physically secure | | |

B. Verification Table for Logical Security Window:

| Specified user | Logical security available (Yes/No) | Discrepancy? (Y/ N) |
|------------------|-------------------------------------|---------------------|
| Logical security | | |



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

C. Verification Table for Logical Security Application Software:

| Specified user | Logical security available (Yes/No) | Discrepancy? (Y/ N) |
|----------------|-------------------------------------|---------------------|
| Analyst | | |
| Reviewer | | |
| Admin | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.7 Verification of Test Instruments Calibration and Traceability:

Objective : To verify the test instruments traceability.

Tools Required : Not Applicable

Procedure : 1. Verify Certificate No.
2. Verify Traceability.
3. Verify Calibration Date.
4. Verify Calibration Due Date.

Acceptance : Test instruments should be calibrated at the execution.

Criteria

Verification Table:

| Certificate No. | Traceability (Attachement No.) | Calibration Done On | Calibration Due On | Verified (Yes/No) | Discrepancy? (Y/N) |
|-----------------|-----------------------------------|------------------------|-----------------------|----------------------|-----------------------|
| | Refer attachment No. 3 | | | | |
| | Refer attachment No. 4 | | | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.8 Verification of Power Utility:

Objective : To verify the installed power utility is as per specifications.

Tools Required : Digital Multimeter

Procedure : 1. Switch ON the Power Supply.
2. Put the Multimeter in AC/DC range.
3. Record the supply voltage.

Acceptance : Measured voltage shall match with the specified voltage.

Criteria

Verification Table:

Supply Voltage Measurement:

| Name | Specified Voltage | Measured Voltage | Discrepancy? (Y/N) |
|------------------------------|-------------------|------------------|--------------------|
| Computer System Power Supply | 230 VAC | | |
| Analytical instrument | 220-240 VAC | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.9 Verification of Environment Condition:

Objective : To verify the environment conditions.

Tools Required : Digital Thermo Hygrometer

Procedure : 1. Switch ON the thermo hygrometer.
2. Record maximum temperature.
3. Record maximum relative humidity.

Acceptance : Test instruments should be calibrated at the execution.

Criteria

Verification Table:

| Name | Temperature | Relative Humidity | Measured Results | Discrepancy? (Y/N) |
|---|-------------|-------------------|------------------|--------------------|
| Computer System Environmental Condition | NMT 25 °C | NA | | |
| Analytical instrument | NMT 25 °C | NA | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.10 Verification of Communication Link Between Server To Computer System:

Objective : To Verify the communication link with server to client

Tools Required : Not Applicable

Procedure : 1. Verify and record the communication link between PC to file server
2. Verify and record the communication link between PC to Printer.

Acceptance : Communication link ping with PC to server should be executed and report

Criteria should be proper.

Verification Table:

| S.No. | Source | Destination | Ping Executed (Yes/No) | Discrepancy? (Y/N) |
|-------|-----------------|-----------------------|------------------------|--------------------|
| 1. | Computer system | File Server | | |
| 2. | Computer system | Printer | | |
| 3. | Computer system | Analytical Instrument | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No Refer Attachment No. []

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.11 Verification of General System Installation:

- Objective : To verify the general system installation.
- Tools Required : Not Applicable
- Procedure : 1. Check all the test given in verification table.
2. Record the result in verification table.
- Acceptance : All the test result should match with expected result.
- Criteria

Verification Table:

| Description | Expected Result (Yes/No) | Discrepancy? (Y/N) |
|--|-----------------------------|-----------------------|
| Major components should be protected from shock. | | |
| No visible physical damage should be available. | | |
| Sufficient space should be available for maintenance. | | |
| System identification nameplate should be available. | | |
| System should be installed with all necessary instruments. | | |
| Earthing should be connected properly. | | |
| Power and signal cable should be separate. | | |
| Unterminated and broken wire should not be open. | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

12.12 Verification of Standard Operating Procedure:

Objective : To verify the availability of related standard operating procedure.

Tools Required : Not Applicable

Procedure : 1. Verify SOP Name.
2. Verify SOP No.
3. Verify SOP Availability.

Acceptance : Documents should be available.

Criteria

Verification Table:

| SOP Name | SOP No. | Availability (Yes/No) | Discrepancy? (Y/N) |
|---|---------|-----------------------|--------------------|
| Standard Operating Procedure of Backup / Restoration of Analytical Instrument Data | | | |
| Standard Operating Desktop Policy for Computer Operated Analytical Instrument | | | |
| Standard Operating Procedure on Password Policy for Software in Laboratory and GMP System | | | |

Remarks:

Meet the acceptance Criteria [] Yes [] No

Checked by : _____

Date : _____

Verified by : _____

Date : _____



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

13. INSTALLATION QUALIFICATION TEST STATUS:

The installation qualification test status is as per below mentioned table.

| Test Description | Status (Pass / Fail) | Discrepancy? (Y/N) |
|--|---------------------------------|-------------------------------|
| Identification of System Details | | |
| Verification of Master Documents | | |
| Verification of Capacity Requirement | | |
| Verification of Hardware Components | | |
| Verification of Software Components | | |
| Verification of Physical and Logical Security Control | | |
| Verification of Test Instruments Calibration and Traceability | | |
| Verification of Power Utility | | |
| Verification of Environmental Condition | | |
| Verification Of Communication Link Between Server To Computer System | | |
| Verification of General System Installation | | |
| Verification of Standard Operating Procedures | | |



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

14. DISCREPANCIES HANDLING DURING COMPUTER SYSTEM QUALIFICATION:

- In case of discrepancy observed during qualification, document in the defined column in each table and document the details of the observation in the discrepancy log sheet.
- Inform to User, engineering and quality assurance about discrepancy.
- Investigate the discrepancy and ensure the possible impact.
- If discrepancy does not have potential to impact on operation as well as performance of the system, close the discrepancy with proper justification.
- The User, engineering and QA will decide whether discrepancy is acceptable or not.
- If discrepancy is acceptable, provide conclusion and recommendation if any into respective column.



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

15. DISCREPANCY AND CORRECTIVE ACTION FORM:

| | |
|--------------------|--|
| Protocol Reference | |
| Discrepancy Number | |

DISCREPANCY:

| | |
|--------------------------|------|
| Describe the Discrepancy | |
| | |
| | |
| | |
| Reported by | Date |

CORRECTIVE ACTION:

| | |
|--|------|
| Describe corrective action taken (Attach additional sheets if necessary) | |
| | |
| | |
| | |
| Reported by | Date |

DISPOSITION ACTION:

| | | |
|-------------|------|----|
| Acceptable? | Yes | No |
| Discussion | | |
| | | |
| | | |
| | | |
| Approved by | Date | |

COMPLETION:

| | |
|--------------|------|
| Completed by | Date |
|--------------|------|



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

16. ABBREVIATION:

| Abbreviations | Description |
|----------------------|---|
| GMP | Good Manufacturing Practices |
| IQ | Installation Qualification |
| OQ | Operation Qualification |
| QA | Quality Assurance |
| SOP | Standard Operating Procedure |
| NA | Not Applicable |
| ICH | International Conference of Harmonization |
| mA | Mili Ampere |
| VAC | Alternate Current Voltage |
| VDC | Direct Current Voltage |
| RH | Relative Humidity |
| CS | Computer System |
| NMT | Not More Than |
| BOM | Bill of material |
| WHO | World Health Organization |



INSTALLATION QUALIFICATION FOR COMPUTER SYSTEM

System Name: HPLC

System ID:

19. POST APPROVALS:

The signature listed below indicates the post approval of this installation qualification. This approval is joint responsibility of listed functional areas.

| DOCUMENT DEVELOPMENT | SIGN / DATE |
|---|-------------|
| Name : _____ Designation : _____ | |

DOCUMENT REVIEW AND APPROVAL (M/S)

Sign / Date : _____
Name : _____
Designation : _____
Engineering

Sign / Date : _____
Name : _____
Designation : _____
IT

Sign / Date : _____
Name : _____
Designation : _____
Quality Control

DOCUMENT APPROVAL (M/S)

Sign / Date : _____
Name : _____
Designation : _____
Quality Assurance