



## STANDARD OPERATING PROCEDURE

**Department:** Formulation, Research and Development

**Title:** Operation and Cleaning of Octagonal Blender

**Supersedes:** Nil

**Issue Date:**

- 1.0 OBJECTIVE:**  
To lay down a procedure for operation and cleaning of Octagonal Blender.
- 2.0 SCOPE:**  
This SOP is applicable to Formulations Research and Development.
- 3.0 RESPONSIBILITY:**
- 3.1 Executive – FR&D for operation and cleaning.
  - 3.2 Group Leader – FR&D will ensure the maintenance.
  - 3.3 Head – FR&D for implementation and compliance.
- 4.0 DEFINITION(s):**  
Not Applicable
- 5.0 PROCEDURE:**
- 5.1 OPERATION:**
- 5.1.1 Check for the cleanliness of the machine
  - 5.1.2 Remove the lid of the blender for loading the material.
  - 5.1.3 Load the material and replace the lid and tighten the screws.
  - 5.1.4 Attach the blender to the motor shaft.
  - 5.1.5 Close the pipe guard and ensure that the limit switch is pressed. The pipe guard should be in horizontal position.
  - 5.1.6 Press the main switch 'ON' on the front panel.
  - 5.1.7 Press 'ESC' key on the program keyboard module.
  - 5.1.8 Press , or + key to increase or decrease the parameter values respectively.
  - 5.1.9 Set the parameter value to 58, the corresponding frequency is shown as 70 on the screen. (This frequency corresponds to approximately 25 rpm)
  - 5.1.10 Press the red button, the values stop flashing.
  - 5.1.11 Press 'SEL' key, the red dot adjacent to 58 flashes.
  - 5.1.12 Then press , or + key to increase or decrease the frequency.
  - 5.1.13 Press ← ( ENTER ) key, the dot stops flashing.
  - 5.1.14 Then finally press 'ESC' key.
  - 5.1.15 Turn the A/M SWITCH on the front panel to auto mode. The panel view touch sensitive screen displays 'SCREEN 1' (as per FT2RD003/A01).
  - 5.1.16 Press the 'SELECT RECIPE' icon, the panel view displays 'SCREEN 2' (as per FT2RD003/A02) with options to enter the product code.
  - 5.1.17 Press the arrow keys , + , the pointer < moves to the desired product code.
  - 5.1.18 Press enter key then press 'PRE' icon.
  - 5.1.19 Press the space provided to enter the product name, it flashes then press enter.



# PHARMA DEVILS

Formulation, Research and Development Department

## STANDARD OPERATING PROCEDURE

**Department:** Formulation, Research and Development

**Title:** Operation and Cleaning of Octagonal Blender

**Supersedes:** Nil

**Issue Date:**

- 5.1.20 The screen displays the alphabet table.
- 5.1.21 Select the letters by using = key. Press 'SEL' icon after selecting every letter and then finally press enter. <
- 5.1.22 Similarly enter the operator's name.
- 5.1.23 Press 'NEXT' icon to view 'SCREEN 3'(as per FT2RD003/A03).
- 5.1.24 Press the space provided for the 'SET TIME' then press enter.
- 5.1.25 Enter the value in 'seconds' and press enter.
- 5.1.26 Press 'PRC START' icon to start the blending process.
- 5.1.27 Press 'PAUSE' icon to pause in between for sometime, then press 'RELEASE PAUSE' to continue the process from the stopped time.
- 5.1.28 Press 'PRC STOP' to stop the process completely.
- 5.1.29 Press 'PREV' icon to return to the previous screen.
- 5.1.30 'SCREEN 4' (as per FT2RD003/A04) displays 'PROCESS END' after completion of blending for the set time.
- 5.1.31 Finally press 'ACK' icon to acknowledge the completion of process.
- 5.1.32 Press the red emergency knob on the front panel to stop the process entirely.
- 5.1.33 Press the 'INCH' button on the panel to turn the blender to the desired position while removing the material.

## 5.2 CLEANING INSTRUCTIONS

- 5.2.1 Switch off the machine and clean the motor guard with a wet cloth followed by dry cloth.
- 5.2.2 Separate the blender from the motor shaft and clean with potable water followed by DM water.
- 5.2.3 Use suitable solvents for dissolving the drug in case of water insoluble materials and then rinse the blender thoroughly with DM water.
- 5.2.4 Wipe the blender with a clean dry cloth.

## 6.0 REFERENCE(s)

Not Applicable

## 7.0 ABBREVIATION(s)

Abbreviation	Full Description
SOP	Standard Operating Procedure
FR&D	Formulation, Research and Development
DM	Demineralised

## 8.0 FLOWCHART(s)

Not Applicable



# PHARMA DEVILS

Formulation, Research and Development Department

## STANDARD OPERATING PROCEDURE

**Department:** Formulation, Research and Development

**Title:** Operation and Cleaning of Octagonal Blender

**Supersedes:** Nil

**Issue Date:**

### 9.0 ANNEXURE(s)

Annexure No.	Details/Title of Annexure	Format No. (Current version)
FT2RD003/A01	Screen 1	-
FT2RD003/A02	Screen 2	-
FT2RD003/A03	Screen 3	-
FT2RD003/A04	Screen 4	-
FT2RD003/A05	Equipment log sheet	FT2RD003/F01







# PHARMA DEVILS

Formulation, Research and Development Department

## STANDARD OPERATING PROCEDURE

**Department:** Formulation, Research and Development


**Title:** Operation and Cleaning of Octagonal Blender

**Supersedes:** Nil

**Issue Date:**

### Annexure 03

#### Screen 03





# PHARMA DEVILS

Formulation, Research and Development Department

## STANDARD OPERATING PROCEDURE

**Department:** Formulation, Research and Development

**Title:** Operation and Cleaning of Octagonal Blender

**Supersedes:** Nil

**Issue Date:**

Annexure 04

Screen 04

PROCESS END

ACK

