


PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 1 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE





PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 2 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

CONTENTS

1.0.0	PRE APPROVAL SIGNATURES:	3
2.0.0	GENERAL:	4
2.1.0	PURPOSE:	4
2.2.0		
2.3.0	BACKGROUND:	4
2.4.0	REVISION HISTORY:	4
2.5.0		
2.6.0	VALIDATION TEAM:	5
2.7.0		
2.8.0	TRAINING RECORD:	7
2.9.0		8
2.10.		9
2.11.	0 PLC SYSTEM SCHEMATIC DIAGRAM	10
3.0.0	OPERATIONAL QUALIFICATION TEST	
0.010		
3.1.0	VERIFICATION OF MASTER TEST INSTRUMENTS	
3.2.0		
3.3.0		
3.4.0		
3.5.0	VERIFICATION OF PLC INPUTS/ OUTPUTS	
3.6.0	VERIFICATION OF SECURITY LEVELS	
3.7.0	VERIFICATION OF HMI SCREENS	
3.8.0	VERIFICATION OF RANGE OF SET PARAMETERS	
3.9.0		
3.10.	0 VERIFICATION OF COMMUNICATION FAILS RECOVERY CONDITION	
3.11.	0 VERIFICATION OF ALARMS AND INTERLOCKS	
3.12.	0 VERIFICATION AND RECORD AN INTEGRATED CONTROL LOOP TEST	45
4.0.0	DEVIATION REPORT AND DISCREPANCY REPORT	47
5.0.0	SUMMARY REPORT	48
6.0.0	TERMINOLOGIES	49
7.0.0	LIST OF ATTACHMENTS	50
8.0.0	LIST OF ABBREVIATIONS	51
9.0.0	POST APPROVAL SIGNATURES	52



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 3 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

1.0.0 PRE APPROVAL SIGNATURES:

The signatures below indicate pre approval of this operational qualification document and it is ready for execution. Any changes or modifications to the intent or the acceptance criteria of this operational qualification document, following approval, requires the generation of an amendment which must be approval prior to execution.

OPERATIONAL QUALIFICATION PRE APPROVAL

Function	Name	Department	Designation	Signature/Date
			4 1	
Prepared by	An	Engineering		
Reviewed by		Engineering		
Reviewed by		Production		
Reviewed by		Quality Assurance		
	Final Approval: Fin	al approval has been	n given by the followin	ng
Function	Name	De	esignation	Signature/Date
Approved by			ad Quality ssurance	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 4 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

2.0.0 <u>GENERAL:</u>

2.1.0 <u>PURPOSE</u>:

The purpose of this operational qualification document is to verify and document that the PLC system of "Automatic Vertical Round Bottle Sticker Labelling Machine" has been operated and fulfill its intended use when placed in its intended environment.

The purpose of the Operational Qualification is to provide documented evidence to demonstrate that the PLC system is operated and performed as per the manufacturer specifications.

2.2.0 <u>SCOPE</u>:

This Operational Qualification will be performed on "Automatic Vertical Round Bottle Sticker Labelling Machine" which is located in "Dry Syrup".

This operational qualification document describes the PLC system hardware and software, equipment details, test procedures, documentation, references and acceptance criteria used to establish that "Wash Area" has been operated in accordance with the master documentations.

2.3.0 BACKGROUND:

The **"Automatic Vertical Round Bottle Sticker Labelling Machine"** is a new system purchase specifically for use at

2.4.0 **REVISION HISTORY:**

Version No.	Effective Date	Reason for Change
00		New Document
00		New Document



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 5 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

2.5.0 <u>REFERENCES:</u>

The test and execution procedure within the scope of the qualification document are consistence with the following reference.

Guideline	Details	
GAMP-5	Good Automated Manufacturing Practices	
21 CFR Part 210	Code of Federal Regulations, Current Good Manufacturing Practices in Manufacturing Processing, Packing.	
21 CFR Part 211	Code of Federal Regulations, Current Good Manufacturing Practices for finished Pharmaceuticals.	
EU GMP Annex-11	European Union Good Manufacturing Practices Annexure-11	

2.6.0 VALIDATION TEAM:

Validation team is responsible for the execution of operational qualification of PLC system. Validation team comprises.

Name	Department Designa	ation Sign & Date
	Engineering	
	Production	
	QA	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 6 of 52


OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

2.7.0 <u>RESPONSIBILITY:</u>

- Collect all manuals, electrical wiring diagram and documentary or any other data necessary for the preparation, execution of operational qualification document from M/S.
- > Preparation and execution of Operational Qualification document.
- > Initiate Qualification study in coordination with Production, Quality Assurance and Engineering.
- Provide training to the persons, who present during execution, of this study.

En	gineering	Pr	oduction	Qu	uality Assurance
	Co-ordinate during execution of Qualification activities.	A	Co-ordinate during execution of Qualification activities.	A	Co-ordinate during execution of Qualification activities.
>	To provide utilities for Qualification activity.	~	Provide personnel for facilitating the execution of Qualification activity.	A	To review and approve the Qualification document.
	To review the operational qualification document.	AA	Check that test requirements To Review the operational qualification document.		


No Name of Trainee Designation of Signature of OK/ To be of	ÂU	PHARMA QUALITY ASSU OPERATIONAL QU TOMATIC VERTICAL		YE DATE: 7 of 52		
Sr. Name of Trainee Designation of Trainee Signature of Trainee Evaluation OK/ To be retrained Signature of evaluation Image: Im	Follow executi	ing persons have been trained on of this qualification docu	ment.	llification document	and will execute	e/help in
St. Name of Trainee Designation of Trainee Signature of Trainee OK/ To be retrained of evaluato Image: Imag	Venue	of Training:		Date of Training:		
Name: Name:		Name of Trainee			OK/ To be	Signature of evaluator
Name: Name:		3		3		
Name: Name:		5				
Name: Name:			<u>nus</u>	minh		
Name: Name:			<u> </u>			
Name: Name:				-		
Name: Name:						
Name: Name:						
Name: Name:						
Designation: Designation:	Nam	e:				
Signature: Signature:		-		-		



2.9.0 PLC DESCRIPTION:

The dictionary defines automation as "the creation and application of technology to monitor and control the production and delivery of products and services."

"Programming Logic Controller" (PLC) is an industrial computer control system that continuously monitors the state of input devices and makes decisions based upon a custom program, to control the state of devices connected as outputs.



PLC consists of input modules or points, a Central Processing Unit (CPU), and output modules or points. An input accepts a variety of digital or analog signals from various field devices (sensors) and converts them into a logic signal that can be used by the CPU. The CPU makes decisions and executes control instructions based on program instructions in memory.

Output modules convert control instructions from the CPU into a digital or analog signal that can be used to control various field devices (actuators). A programming device is used to input the desired instructions. These instructions determine what the PLC will do for a specific input. An operator interface device allows process information to be displayed and new control parameters to be entered.

The PLC is used many inputs or modules to sense and measure physical quantities of equipment, such as motion, temperature, level, current, voltage, position, and pressure etc. Depending on the status of inputs which sensed by inputs or modules, processor controls various output module to energize or drive the field devices such as valves, motor starters and contactors etc that apply power circuit voltages to the control devices.



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 9 of 52


OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

Digital or discrete input/output has only two states, one is ON and another is OFF. Input and output have light emitting diode (LED) to indicate the state of each input/output. Analog input/output allow to monitor and controlling analog voltage and control.

"Human Machine Interface" (HMI) is platform which is assist the operator to supervise and control the equipment. Operator has displayed information from the HMI and gives the command to PLC then PLC will execute the command.

2.10.0 SYSTEM DESCRIPTION

The PLC Based control system should have capability to control Automatic Vertical Round Bottle Sticker Labelling Machine. Purpose of the PLC system is to monitor, operate and control the machine. PLC System gets Digital and analog signal from Field devices like Emergency Stop, No Bottle, Label Gap Sensor, And Other Sensors. The data is processed in CPU of PLC and give Output as per the logics of Servo Pass, Servo Lock and etc.





PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:


PAGE No.: 11 of 52


OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.0.0 OPERATIONAL QUALIFICATION TEST

<u>Sr. No.</u>	<u>Test Details</u>
1	VERIFICATION OF MASTER TEST INSTRUMENTS
2	VERIFICATION OF CALIBRATION CERTIFICATES OF FIELD INSTRUMENTS
3	VERIFICATION OF MAINTENANCE QUALIFICATION DOCUMENTS
4	VERIFICATION OF LED INDICATION OF PLC SYSTEM
5	VERIFICATION OF PLC INPUTS/ OUTPUTS
6	VERIFICATION OF HMI SCREENS
7	VERIFICATION THE RANGE OF SET PARAMETERS
8	VERIFICATION AND TESTING OF POWER LOSS RECOVERY CONDITION
9	VERIFICATION OF COMMUNICATION FAILS RECOVERY CONDITION
10	VERIFICATION OF ALARMS AND INTERLOCKS
11	VERIFICATION OF INTEGRATED CONTROL LOOP

	PROTOCOL No.: QUALITY ASSURANCE DEPARTMENT REVISION No: 00 EFFECTIVE DATE: PAGE No.: 12 of 52 OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE							
3.1.0 <u>VE</u>	RIFICATIO	OF MASTER TEST INSTRUMENTS						
թւ	ırpose	: This test is verified to master instrument which is us	sed for testing.					
Sc	cope	: Recording of master test instrument name, ID no date.	and calibration date and due					
Pı	rocedure	 Record the details of Master Instrument as be Instrument Name - Instrument Tag Plate Instrument ID No - Instrument Tag Plate Calibration date and Due date – Calibration 						
Di	iscrepancy	: If any discrepancy is encountered which prevent originally intended, document the discrepancy report						
	cceptance riteria	: Recorded data from manufacturer tag plates and m calibrated.	naster instruments should be					








OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE


3.2.1 DATA TABLE OF CALIBRATION CERTIFICATES OF SENSORS

Sr. No.	Sensor Name	Sensor ID No	Calibration Date	Calibration Due Date	Sign. & date
				<u></u>	
	0]
			72		
			3		
			R	 	
		~ ~ ~]

	PHARMA SCH	ALADS	PROTOCOL No.:
	QUALITY ASSURANCE DE		REVISION No: 00
- V	CONTLA REPORTANCE D		EFFECTIVE DATE:
			PAGE No.: 16 of 52
OPERA	TIONAL QUALIFI	CATION DOCUMENT O	
		FOR	
AUTOMATIC	VERTICAL ROUN	ID BOTTLE STICKER LA	ABELLING MACHINE
Comments/ Remarks:			
)
		bach	
Function	Name	Department	Sign. & Date
Tested by		J Contraction	
Tested by		Engineering	
Verified by		Engineering	
Reviewed by		QA	
Kevieweu by		QA	

		PHARMA SCHOLARS	PROTOCOL No.:
		QUALITY ASSURANCE DEPARTMENT	REVISION No: 00
			EFFECTIVE DATE:
			PAGE No.: 17 of 52
	OPE	RATIONAL QUALIFICATION DOCUMENT	
		FOR	
	AUTOMAT	FIC VERTICAL ROUND BOTTLE STICKER	LABELLING MACHINE
3.3.0	VERIFICATIO	ON OF MAINTENANCE QUALIFICATION DOCI	UMENTS
	Purpose	: Verify and review the documents of maintenance	e.
	Scope	: Verify and record the maintenance qualification	document.
	Procedure	 Check the availability of documents. Check the every entry for equipment maint 	enance is listed in log-book.
	Discrepancy	: If any Discrepancy is encountered which prevoriginally intended, document the Discrepancy F	
	Acceptance Criteria	• All documents should available and maintain.	







QUALITY ASSURANCE DEPARTMENT

PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 20 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE


3.4.1 DATA TABLE OF PLC SYSTEM LED INDICATION IN OFF CONDITION

LED	Expected state of LED	Actual state of LED	Meets acceptance criteria	Sign.& date
	PLC PROCES	SOR (Mitsubishi FX5	U-80MR-ES)	
Power	OFF		Yes () No ()	
Run	OFF		Yes () No ()	
Error	OFF		Yes () No ()	
Comm 1	OFF		Yes () No ()	

3.4.2 DATA TABLE OF PLC SYSTEM LED INDICATION IN ON CONDITION

LED	Expected state of LED	Actual state of LED	Meets acceptance criteria	Sign.& date
	PLC PROCES	SOR (Mitsubishi FX5)	U-80MR-ES)	
Power	ON	. .	Yes () No ()	
Run	ON		Yes () No ()	
Error	OFF		Yes () No ()	
Comm 1	Blink		Yes () No ()	

	PHARMA SCH		REVISION No: 00
V	QUALITY ASSURANCE D	EPARTMENT	EFFECTIVE DATE:
4 F			
			PAGE No.: 21 of 52
OPERA	TIONAL QUALIFI	CATION DOCUMENT OF FOR	PLC SYSTEM
AUTOMATIC	VERTICAL ROUN	ID BOTTLE STICKER LAI	BELLING MACHINE
nments/ Remarks:			
linents/ Kemarks.			
	C	9	
Function	Name	Department	Sign. & Date
	Name	Department	Sign. & Date
Tested by	Name	Department Engineering	Sign. & Date
	Name	Department	Sign. & Date
Tested by	Name	Department Engineering	Sign. & Date
Tested by Verified by	Name	Department Engineering Engineering	Sign. & Date
Tested by Verified by	Name	Department Engineering Engineering	Sign. & Date
Tested by Verified by	Name	Department Engineering Engineering	Sign. & Date
Tested by Verified by	Name	Department Engineering Engineering	Sign. & Date
Tested by Verified by	Name	Department Engineering Engineering	Sign. & Date
Tested by Verified by	Name	Department Engineering Engineering	Sign. & Date
Tested by Verified by	Name	Department Engineering Engineering	Sign. & Date





PROTOCOL No.:


REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 23 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

- **Discrepancy** : If any discrepancy is encountered which prevent completion of the report as originally intended, document the discrepancy report.
- Acceptance: > All inputs and outputs must be verified to meet wiring diagram of PLC systemCriteriaand function as per design document.
 - All analog inputs/outputs range must be verified to meet calibration certificates of individual device.
 - Verification of all analog inputs/ outputs range should be within 2% of set value.





QUALITY ASSURANCE DEPARTMENT

PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 24 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.5.1 DATA TABLE OF PLC DIGITAL INPUTS

		PLC Digital Input ((DVP 12SA211T)		
PLC Address	Description	Expected state of LED	Actual state of LED	Meets acceptance criteria	Sign. & date
VO	No Bottle	ON		Yes () No ()	
X0		OFF		Yes () No ()	
V 1	Label Gap Sensor	ON		Yes () No ()	
X1		OFF		Yes () No ()	
X2	Spare	NA		Yes () No ()	
N2	SOAP Pump O/L	ON		Yes () No ()	
X3		OFF		Yes () No ()	
X4	Spare	NA		Yes () No ()	
X5	Spare	NA		Yes () No ()	
X6	Spare	NA		Yes () No ()	



QUALITY ASSURANCE DEPARTMENT

PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 25 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.5.2 DATA TABLE OF PLC DIGITAL OUTPUTS

PLC Digital Outputs (DVP 12SA211T)								
PLC Address	DescriptionExpected state of LEDActual state of LED		Meets acceptance criteria	Sign. & date				
Y0	Comes Dulos	ON		Yes () No ()				
	Servo Pulse	OFF		Yes () No ()				
X 71	Correct Locale	ON		Yes () No ()				
Y1	Servo Lock	OFF		Yes () No ()				
Y2	Spare	NA	COOL S	Yes () No ()				
V2	Delay 1	ON		Yes () No ()				
Y3	Relay 1	OFF		Yes () No ()				

Comments/ Remarks:

Function	Name	Department	Sign. & Date
Tested by		Engineering	
Verified by		Engineering	

	PHARMA SCHOLARS GUALITY ASSURANCE DEPARTMENT ERATIONAL QUALIFICATION DOCUMENT O FOR FIC VERTICAL ROUND BOTTLE STICKER L	
	QA QA	
3.6.0 <u>VERIFICATIO</u>	ON OF SECURITY LEVELS	
Purpose	: Verify and testing of different security levels to access.	prevent the unauthorized user
Scope	: Verify and record the user access to PLC.	2
Procedure	 Switch ON the PLC. Challenge all the parameters with entering diff Record the warning messages from PLC (if an 	
Discrepancy	: If any Discrepancy is encountered which preven originally intended, document the Discrepancy Re	
Acceptance Criteria	System must prevent the unauthorized users. Accer particular level as per configured.	ess Rights should be limited to

\frown	DUNDAN COUNT	A DC	PROTOCOL No.:
	PHARMA SCHOL		
	QUALITY ASSURANCE DEPART	14/EN]	REVISION No: 00 EFFECTIVE DATE:
T			EFFECTIVE DATE:
			PAGE No.: 27 of 52
OPERA 7	FIONAL QUALIFICA	TION DOCUMENT OF P	LC SYSTEM
		FOR	
AUTOMATIC	VERTICAL ROUND B	OTTLE STICKER LABE	LLING MACHINE
Comments/ Remarks:			
Function	Name	Department	Sign. & Date
Tested by	- Mar	Engineering	
		Engineering	
Verified by	V		
Verified by Reviewed by		QA	
		ÇA ÇA	
		QA	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 28 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.7.0 **VERIFICATION OF HMI SCREENS** Purpose : To verify each function of HMI screens, function keys. Scope : Verify and record all HMI screens and function keys. Procedure : \succ Take program backup of HMI and verify with the actual function screens. > If backup not possible then take screens from manual and verify with the actual function screens. > Check all programmable function keys for their actual response in each screen. Record the results. > Check Function Key Command are properly programmed by operating output devices. > Verify the display on HMI with actual machine conditions and record the actual results. > Open a screen; check values displayed on HMI for each display object. Record observations. Discrepancy : If any discrepancy is encountered which prevent completion of the report as originally intended, document the discrepancy report. Acceptance All programmable keys and displays should perform as per define function : Criteria



QUALITY ASSURANCE DEPARTMENT

PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 29 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.7.1 DATATABLE OF HMI SCREENS

Sr. No.	Screen Name	Available Yes/No	Meets acceptance criteria	Sign. & date
01	Main Menu Screen		Yes () No ()	
02	Manual Screen		Yes () No ()	
03	Setting Parameter Screen		Yes () No ()	
04	Delay Setting Screen		Yes () No ()	
05	Auto Screen		Yes () No ()	



QUALITY ASSURANCE DEPARTMENT

PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 30 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.7.2 DATATABLE OF HMI DISPLAY

Parameter	Span	Actual Result	Meet Acceptance Criteria	Sign. & Date
Auto Mode Screen				
Machine Speed (BPM)	1 to 200		Yes() No()	
Counter	Numeric		Yes() No()	
Setting Parameter Screen				
Servo Speed Factor (MPM)	0.01 to 9.99		Yes() No()	
Delay Setting Screen	C		9	
Delay Time (100ms) 0 to 200			Yes () No ()	
On Time (100ms) 0 to 200			Yes() No()	
Label Start (sec) 0.00 to 9.9			Yes() No()	
Manual Screen	Aus	J VILL	22	
Machine Speed (BPM)	1 to 200		Yes () No ()	
Counter	Numeric	<u> </u>	Yes () No ()	
nments/ Remarks:				
nments/ Remarks:	Name	Department	Sig	gn. & Date
	Name	Department Engineering	-	gn. & Date
Function	Name			gn. & Date



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 31 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.8.0 **VERIFICATION OF RANGE OF SET PARAMETERS** Purpose : To check and verify the range of set parameters in HMI screens. Scope : Verify and record the minimum and maximum values of set parameter. : \succ Enter minimum value for the given span and record it. Procedure Enter maximum value for the given span and record it. Enter value above and below the acceptable span and observe the response. > If value cannot be entered **"above"** maximum and **"below"** minimum then record the message as "Value cannot be entered". Discrepancy : If any discrepancy is encountered which prevent completion of the report as originally intended, document the discrepancy report. Discrepancy : If any discrepancy is encountered which prevent completion of the report as originally intended, document the discrepancy report. Acceptance : All programmable keys and displays should perform as per define function Criteria



PROTOCOL No.:

••••

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 32 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.8.1 DATATABLE OF HMI SCREENS

Parameter	Span	Min span	Value is set in system	Max span	Value is set in system	Below Range setting	Value is not set in system	Upper Range setting	Value is not set in system	Meet Acceptanc e Criteria	Sign. & Date
Auto Mode Scre	een				A 2	S man	3				
Machine Speed (BPM)	1 to 200		Yes () No ()	6	Yes () No ()		Yes () No ()		Yes () No ()	Yes () No ()	
Setting Paramet	ter Screen										
Servo Speed Factor (MPM)	0.01 to 9.99		Yes () No ()		Yes () No ()	v	Yes () No ()		Yes () No ()	Yes () No ()	
Delay Setting So	creen										
Delay Time (100ms)	0 to 2000		Yes () No ()		Yes () No ()		Yes () No ()		Yes () No ()	Yes () No ()	
On Time (100ms)	0 to 2000		Yes () No ()	C	Yes () No ()		Yes () No ()		Yes () No ()	Yes () No ()	
Label Start (sec)	0.00 to 9.99		Yes () No ()		Yes () No ()		Yes () No ()		Yes () No ()	Yes () No ()	



PROTOCOL No.:


•••••

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 33 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE





PROTOCOL No.:

••••

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 34 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

Parameter	Span	Min span	Value is set in system	Max span	Value is set in system	Below Range setting	Value is not set in system	Upper Range setting	Value is not set in system	Meet Acceptanc e Criteria	Sign. & Date
Manual Screen							77				
Machine Speed (BPM)	1 to 200		Yes () No ()	E	Yes () No ()		Yes () No ()		Yes () No ()	Yes () No ()	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 35 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

	e		
Function	Name	Department	Sign. & Date
Tested by		Engineering	
Verified by	And	Engineering	
Reviewed by		QA	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 36 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.9.0 VERIFICATION OF POWER LOSS RECOVERY CONDITION

Purpose	: Verify and testing of power loss recovery condition.
Scope	: Record value of process set-parameters after the power fails.
Procedure	 Switch "ON" the PLC system. Enter value in set parameter field in HMI screen and record it. Start the machine. Cut power supply during machine is in "ON" condition. Restore power supply after 240 seconds. Restart the machine and again record the value which is entered before. Ensure that machine should not start unless and until start command is given.
Discrepancy	: If any Discrepancy is encountered which prevent completion of the report as originally intended, document the Discrepancy Report.
Acceptance Criteria	 After Power restore the value of process set-parameter should not changed. The process time should start from previous value when the power loss. The machine should not start without command by user when power failure condition occurs. Machine works normally after Power loss recovery condition.



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 37 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.9.1 DATA TABLE OF POWER LOSS RECOVERY CONDITION

Parameter Description	Value Before power loss	Value after power recovery	Expected Result	Meets acceptance criteria	Sign. & date
Auto Mode Scree	n				
Machine Speed (BPM)			Value remains same	Yes () No ()	
Setting Paramete	r Screen				
Servo Speed Factor (MPM)	9		Value remains same	Yes () No ()	
Delay Setting Scr	een 🧹	and I			
Delay Time (100ms)	(Value remains same	Yes () No ()	
On Time (100ms)			Value remains same	Yes () No ()	
Label Start (sec)			Value remains same	Yes () No ()	
Manual Screen					
Machine Speed (BPM)			Value remains same	Yes () No ()	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 38 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

Comments/ Remarks:

Function	Name	Department	Sign. & Date
Tested by		Engineering	
Verified by		Engineering	
Reviewed by	A COLOR	QA	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 39 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.10.0 VERIFICATION OF COMMUNICATION FAILS RECOVERY CONDITION

Purpose	: Verify and testing of communication fails recovery.
Scope	: Record the value of process set-parameter after communication fails.
Procedure	 Switch "ON" the PLC system. Enter value in set parameter field in HMI screen and record it. Start the machine in auto mode. While the machine is in running condition, unplug the communication link cable between PLC and HMI. Take a screen shot or photograph of communication link failure message. Try to change recipe in HMI screen. Records the result or error message, if no message appears then write "recipe cannot be changed". Reconnect the communication link cable after 240 seconds. Record the values again.
Discrepancy	: If any Discrepancy is encountered which prevent completion of the report as originally intended, document the Discrepancy Report.
Acceptance Criteria	 : In communication link failure condition, message must appear in HMI screen. > When communication link fails the recipe should not be changed or saved. > After communication link reconnect, the set-parameters should not be changed. > Machine works normally after communication failure condition.



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 40 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.10.1 DATA TABLE OF COMMUNICATION FAILURE CONDITION

Procedure	Expected result	Actual Result	Meets acceptance criteria	Sign. & date
Start the system in auto mode and disconnect the PLC – HMI communication cable	"Err" Message should display in HMI screen.		Yes () No ()	
Reconnect the PLC – HMI communication cable	Message should remove from HMI screen.		Yes () No ()	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 41 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

Parameter Description	Value Before Communicati on unplug	Value after Communicat ion link reconnect	Expected Result	Meets acceptance criteria	Sign. & date
Auto Mode Scree	en				J
Machine Speed (BPM)			Value remains same	Yes () No ()	
Setting Paramete	er Screen				
Servo Speed Factor (MPM)			Value remains same	Yes () No ()	
Delay Setting Scr	reen				
Delay Time (100ms)	8		Value remains same	Yes () No ()	
On Time (100ms)	Ś		Value remains same	Yes () No ()	
Label Start (sec)		his	Value remains same	Yes () No ()	
Manual Screen					
Machine Speed (BPM)			Value remains same	Yes () No ()	
uments/ Remarks:					
iments/ Remarks:	: Nam	e	Department	S	ign. & Date
		e	Department Engineering	S	ign. & Date
		e		S	ign. & Date



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:


PAGE No.: 42 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.11.0 VERIFICATION OF ALARMS AND INTERLOCKS

Purpose	: Verify that all alarms and interlocks of the PLC system and its functions.
Scope	: Check and record the alarms and interlocks.
Procedure	 Start the equipment in normally. Find cause for particular alarm and record it. Record the alarm message, which indicate on HMI. Check and observe the interlock of equipment. If printing facility available, attached the printout of generated alarms otherwise take photograph of each.
Discrepancy	: If any Discrepancy is encountered which prevent completion of the report as originally intended, document the Discrepancy Report.
Acceptance Criteria	: All alarms and interlocks of machine should work properly.



Cause	Expected Result	Actual Result	Correction	Meets acceptance criteria	Sign. & date
Emergency Switch ()perated				
If emergency push button is pressed	"Emergency Switch Operated" alarm message should display in HMI screen and Auto Cycle Deactivated And The Total System Should Go In Hold Mode.		Release the emergency push button	Yes () No ()	



PROTOCOL No.:

REVISION No: 00


EFFECTIVE DATE:

PAGE No.: 44 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

Comments/ Remarks:



Function	Name	Department	Sign. & Date
Tested by		Engineering	
Verified by		Engineering	
Reviewed by		QA	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 45 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.12.0 VERIFICATION AND RECORD AN INTEGRATED CONTROL LOOP TEST

Scope	: Check and record of an integrated control loop test.
Procedure	 Start the equipment in normally. Login with higher level id and password. Set require recipe for the test for process start Start process and observe the set process parameters. Record the reading of set process parameters until the completion of process If printing facility available, attached the printout of whole integrated control loop test.
Discrepancy	: If any Discrepancy is encountered which prevent completion of the report as originally intended, document the Discrepancy Report.
Acceptance Criteria	PLC system should able to control the set process parameter within the specified limit



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 46 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM

FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

3.12.1 DATA TABLE OF CONTROL LOOP TEST

Set P	arameters	<u>Set Value</u>	
Machine Spe	ed (RPM)		
Sr. No	Specified	Actual Observation	Meet Acceptance Criteria
1	Process Should start as per set parameter		Yes () No ()
nents/ Rema	rks:		

Function	Name	Department	Sign. & Date
Tested by		Engineering	
Verified by		Engineering	
Reviewed by		QA	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 47 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

4.0.0 DEVIATION REPORT AND DISCREPANCY REPORT

•	ficiency and its classification*			Category	
Sr.No.	Defici	Deficiency			
	prrective action, Responsible person				
Sr. No.	Recommended corrective action	n	Responsibility	Assigned date	
			0		
Engineerir (Sign and d			As <mark>surance</mark> and date)		
-	s taken (For Category C deficiency)				
Sr.	Corrective action taken	-	Sign	Date	
		a Maria			
Closure remarks	: Allowed / Not allowed to proceed furt	ner			
Reviewed and ap	pproved by Engineering:				
	proved by Quality Assurance:				

Follow-up Compliance (For category C deficiency):

Recomr	nended corrective actions taken (Action taken v	within stipulated period)	
Sr.	Corrective action taken	Sign	Date
Closure	e remarks:		
Review	ed and approved by Engineering:		

*Category A: Equipment/instrument/system accepted with deficiency

Category B: Conditional acceptance of equipment, deficiency to be corrected within stipulated period

Category C: Deficiency to be rectified before proceeding further



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 48 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

5.0.0 SUMMARY REPORT

	C	2	
Function	Name	Department	Sign. & Date
Tested by		Engineering	
Verified by		Engineering	
Reviewed by		QA	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 49 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

6.0.0 <u>TERMINOLOGIES</u>

✤ Access security :

For protection that ensures system access only to authorized persons on their assigned access level.

***** Automated system :

A system that automatically, without human intervention, controls or monitors a specific set of sequential activities; such as a plant process, laboratory function, or data processing operation.

♦ Installation Qualification (IQ) :

Document evidence that verify the equipment and its sub-system has been installed properly as per the specification.

✤ Operational Qualification (OQ):

Document evidence that the equipment related system or subsystem has been operated properly as per specification.

✤ Validation :

Documented evidence which provides a high degree of assurance that a specific process will consistently produce a product meeting its predetermined specifications and quality attributes.

* PLC:

PLC is an industrial computer control system that continuously monitors the state of input devices and makes decisions based upon a custom program, to control the state of devices connected as outputs.

♦ CPU:

The **central processing unit** (**CPU**) is the portion of a computer system that carries out the instructions of a computer program, to perform the basic arithmetical, logical, and input/output operations of the system.

SOP :Standard Operating Procedure



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 50 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

7.0.0 LIST OF ATTACHMENTS

<u>Sr.</u> <u>No.</u>	<u>Reference</u>	Description Of Attachment Master Test Instrument Calibration Certificate Field Instrument Calibration Certificate	
1	Attachment-1		
2	Attachment-2		
3	Attachment-3	HMI Screen	
4	Attachment-4	Communication Fail	



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 51 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

8.0.0 LIST OF ABBREVIATIONS

<u>Acronym</u>		Description
CPU	\rightarrow	Central Processing Unit
CGMP	\rightarrow	Current Good Manufacturing Practices
GAMP	\rightarrow	Good Automated Manufacturing Practices
GMP	\rightarrow	Good Manufacturing Practices
ID	\rightarrow	Identification Number
ΙΟ	\rightarrow	Input Output
RQ	\rightarrow	Re- Qualification
PLC	\rightarrow	Programmable Logic Controller
SOP	\rightarrow	Standard Operating Procedure
UPS	\rightarrow	Uninterruptible Power Supply
VAC	\rightarrow	Volts Alternating Current
VDC	\rightarrow	Volts Direct Current
VMP	\rightarrow	Validation Master Plan
HMI	\rightarrow	Human Machine Interface



PROTOCOL No.:

REVISION No: 00

EFFECTIVE DATE:

PAGE No.: 52 of 52

OPERATIONAL QUALIFICATION DOCUMENT OF PLC SYSTEM FOR

AUTOMATIC VERTICAL ROUND BOTTLE STICKER LABELLING MACHINE

9.0.0 POST APPROVAL SIGNATURES

The signatures below indicate post approval of this Operational Qualification document and it is executed properly. All variations or discrepancies have been satisfactorily resolved.

Executed by	9	Engineering	b ?	
Reviewed by	5	Engineering	1	
Reviewed by		Production	2	
Reviewed by		Quality Assurance		
Function	Name	De	esignation	Signature/Date
Function	Name			Signature/Date
Approved by			ad Quality ssurance	