



# PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Department:** Microbiology

**SOP No.:**

**Title:** Antimicrobial Effectiveness Testing

**Effective Date:**

**Supersedes:** Nil

**Review Date:**

**Issue Date:**

**Page No.:**

### 1.0 PURPOSE

To lay down the procedure to perform Antimicrobial effectiveness testing to be performed using standard culture.

### 2.0 SCOPE

It is applicable to the microbiology lab for Antimicrobial effectiveness testing.

### 3.0 RESPONSIBILITY

Microbiology personnel

### 4.0 PROCEDURE

#### 4.1 Precaution taken during antimicrobial effectiveness testing.

4.1.1 Glassware to be used shall be sterilized.

4.1.2 Media to be used shall be pre incubated

4.1.3 Standard culture to be used shall not be more than 4 passage

#### 4.2 Test organisms used.

Use cultures of the following microorganisms *Candida albicans* ATCC No. 10231, *Aspergillus niger* ATCC No. 16404, *Escherichia coli* ATCC No. 8739, *Pseudomonas aeruginosa* ATCC No. 9027, and *Staphylococcus aureus* ATCC No. 6538. The viable microorganisms used in the test must not be more than five passages removed from the original ATCC culture or any other equivalent cultures.

#### 4.3 Preparation of inoculums

4.3.1 Prepare the inoculums as per the SOP No. SOP/QCM/016.

4.3.2 To harvest the bacterial and *Candida albicans* cultures, use sterile peptone saline, wash the surface growth, collecting it in a suitable glassware, and adding sufficient sterile peptone saline to obtain a microbial count of about  $1 \times 10^8$  colony-forming units (CFU) per ml. To harvest the cells of



# PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Department:** Microbiology

**SOP No.:**

**Title:** Antimicrobial Effectiveness Testing

**Effective Date:**

**Supersedes:** Nil

**Review Date:**

**Issue Date:**

**Page No.:**

*Aspergillus niger*, use sterile peptone saline containing 0.05% of polysorbate 80, and add sufficient sterile peptone saline to obtain a count of about  $1 \times 10^8$  CFU per ml.

4.3.3 Determine the number of CFU per ml in each suspension, using the conditions of media and microbial recovery incubation times listed in *Table below* to confirm the initial CFU per ml. This value serves to calibrate the size of inoculum used in the test. The bacterial and yeast suspensions are to be used within 24 hours of harvest, but the fungal preparation may be stored under refrigeration for up to 7 days.

4.3.4 Culture Conditions for Inoculum Preparation as below table

| Organism                | Suitable Medium | Incubation Temperature        | Incubation Time | Microbial Recovery Incubation Time |
|-------------------------|-----------------|-------------------------------|-----------------|------------------------------------|
| Escherichia coli        | SCM /SCA        | $32.5 \pm 2.5^\circ \text{C}$ | 18 to 24 hours  | 3 to 5 days                        |
| Pseudomonas aeruginosa  | SCM /SCA        | $32.5 \pm 2.5^\circ \text{C}$ | 18 to 24 hours  | 3 to 5 days                        |
| Staphylococcus aureus   | SCM /SCA        | $32.5 \pm 2.5^\circ \text{C}$ | 18 to 24 hours  | 3 to 5 days                        |
| <i>Candida albicans</i> | SCM /SDA        | $22.5 \pm 2.5^\circ \text{C}$ | 44 to 52 hours  | 3 to 5 days                        |
| Aspergillus niger       | SCM /SDA        | $22.5 \pm 2.5^\circ \text{C}$ | 6 to 10 days    | 3 to 7 days                        |

### 4.4 Test Procedure

4.4.1 The test can be conducted either in five original containers if sufficient volume of product is available in each container or in five sterile capped bacteriological containers of suitable size into which a sufficient volume (Not less than) of product has been transferred.

4.4.2 Inoculate each container with one of the prepared and standardized inoculum, and mix. The volume of the suspension inoculum used is between 0.5% and 1.0% of the volume of the product.

4.4.3 The concentration of test microorganisms that is added to the product are such that the final concentration of the test preparation after inoculation is between  $1 \times 10^5$  and  $1 \times 10^6$  CFU per ml of the product.



# PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Department:** Microbiology

**SOP No.:**

**Title:** Antimicrobial Effectiveness Testing

**Effective Date:**

**Supersedes:** Nil

**Review Date:**

**Issue Date:**

**Page No.:**

- 4.4.4 The initial concentration of viable microorganisms in each test preparation is estimated based on the concentration of microorganisms in each of the standardized inoculum as determined by the plate-count method.
- 4.4.5 Incubate the inoculated containers at  $22.5 \pm 2.5^{\circ}\text{C}$ . Sample each container at the appropriate intervals specified in point number 4.5. Record any changes observed in appearance at these intervals.
- 4.4.6 Determine by the plate-count procedure the number of CFU present in each test preparation for the applicable intervals
- 4.4.7 Incorporate an inactivator (neutralizer) of the specific antimicrobial in the plate count or in the appropriate dilution prepared for plating (if required).
- 4.4.8 Using the calculated concentrations of CFU per ml present at the start of the test, calculate the change in  $\log_{10}$  values of the concentration of CFU per ml for each microorganism at the applicable test intervals, and express the changes in terms of log reductions.

### 4.5 Frequency and Acceptance criteria

As per EP:

| Category of products                   | Case | Bacterial log reduction |        |        |          |           |           | Fungal log reduction |           |           |
|--|------|-------------------------|--------|--------|----------|-----------|-----------|----------------------|-----------|-----------|
|  |      | 6 hrs                   | 24 hrs | 2 days | 7 th day | 14 th day | 28 th day | 7 th day             | 14 th day | 28 th day |
| Parenteral and ophthalmic preparations | A    | 2                       | 3      | -      | -        | -         | NR        | 2                    | -         | NI        |
|  | B    | -                       | 1      | -      | 3        | -         | NI        | -                    | 1         | NI        |
| Topical preparations                   | A    | -                       | -      | 2      | 3        | -         | NI        | -                    | 2         | NI        |
|  | B    | -                       | -      | -      | -        | 3         | NI        | -                    | 1         | NI        |
| Oral preparations                      | A    | -                       | -      | -      | -        | 3         | NI        | -                    | 1         | NI        |



# PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Department:** Microbiology

**SOP No.:**

**Title:** Antimicrobial Effectiveness Testing

**Effective Date:**

**Supersedes:** Nil

**Review Date:**

**Issue Date:**

**Page No.:**

### As per USP:

| Category  | Bacteria  | Yeast and Molds  |
|---|---|--|
| Injections, other parenterals including emulsions, otic products, sterile nasal products, and ophthalmic products made with aqueous bases or vehicles | Not less than 1.0 log reduction from the initial calculated count at 7 days, not less than 3.0 log reduction from the initial count at 14 days, and no increase from the 14 days' count at 28 days. | No increase from the initial calculated count at 7, 14, and 28 days. |
| Topically used products made with aqueous bases or vehicles, nonsterile nasal products, and emulsions, including those applied to mucous membranes.   | Not less than 2.0 log reduction from the initial count at 14 days, and no increase from the 14 days' count at 28 days.  | No increase from the initial calculated count at 14 and 28 days.     |
| Oral products other than antacids, made with aqueous bases or vehicles.   | Not less than 1.0 log reduction from the initial count at 14 days, and no increase from the 14 days' count at 28 days.  | No increase from the initial calculated count at 14 and 28 days.     |
| Antacids made with an aqueous base  | No increase from the initial calculated count at 14 and 28 days.  |  |

### As per Harmonized criteria:



# PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Department:** Microbiology

**SOP No.:**

**Title:** Antimicrobial Effectiveness Testing

**Effective Date:**

**Supersedes:** Nil

**Review Date:**

**Issue Date:**

**Page No.:**

| Category of products                   | Case | Bacterial log reduction |        |        |         |          |          | Fungal log reduction |          |          |
|--|------|-------------------------|--------|--------|---------|----------|----------|----------------------|----------|----------|
|  |      | 6 hrs                   | 24 hrs | 2 days | 7th day | 14th day | 28th day | 7th day              | 14th day | 28th day |
| Parenteral and ophthalmic preparations | A    | 2                       | 3      | -      | NI      | NI       | NR       | 2                    | NI       | NI       |
|  | B    | -                       | 1      | -      | 3       | NI       | NI       | NI                   | 1        | NI       |
| Topical preparations                   | A    | -                       | -      | 2      | 3       | NI       | NI       | -                    | 2        | NI       |
|  | B    | -                       | -      | -      | -       | 3        | NI       | -                    | 1        | NI       |
| Oral preparations                      | A    | -                       | -      | -      | -       | 3        | NI       | -                    | 1        | NI       |
| Antacids                               | -    | -                       | -      | -      | -       | NI       | NI       | -                    | NI       | NI       |

NR: No recovery      NI: No increase

A: Recommended efficacy expected

B: In justified cases where the A criteria cannot be attained, the B criteria must be satisfied.

### 5.0 ABBREVIATIONS AND DEFINITIONS

|      |                                     |
|------|-------------------------------------|
| SOP  | Standard Operating Procedure        |
| QCM  | Quality Control Microbiology        |
| QAD  | Quality assurance Department        |
| Rev. | Revision                            |
| No.  | Number                              |
| ATCC | American Type of Culture Collection |
| CFU  | Colony Forming Unit                 |
| SCM  | Soyabean Casein Digest Medium       |
| SCA  | Soyabean Casein Digest Agar         |
| SDA  | Sabouraud Dextrose Agar             |
| °C   | Degree Centigrade                   |
| %    | Percentage                          |



# PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

## STANDARD OPERATING PROCEDURE

**Department:** Microbiology

**SOP No.:**

**Title:** Antimicrobial Effectiveness Testing

**Effective Date:**

**Supersedes:** Nil

**Review Date:**

**Issue Date:**

**Page No.:**

### 8.0 REVISION LOG

| Revision Number | Effective Date | Reason for Revision |
|-----------------|----------------|---------------------|
|                 |                |                     |
|                 |                |                     |
|                 |                |                     |
|                 |                |                     |
|                 |                |                     |
|                 |                |                     |

### 6.0 REFERENCE DOCUMENTS

SOP: Procedure for Storage and Preparation of microbiological Culture media”.

### 7.0 ANNEXURE/ATTACHMENTS

None