



**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Microbiology	<b>SOP No.:</b>
<b>Title:</b> Maintenance of Cultures	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

1.0 **OBJECTIVE**

To lay down a procedure for maintaining of Bacterial cultures.

2.0 **RESPONSIBILITY**

Microbiologist / Executive.

3.0 **ACCOUNTABILITY**

Head - Quality Control

3.0 **PROCEDURE**

- 4.1 Bacterial cultures are very sensitive and if not sub-cultured they can change their morphological & biochemical characteristics. Sub-culturing or periodic transfer is a very delicate technique which has to be carried out by avoiding chances of contamination.
- 4.3 Prepare the media as mentioned in SOP.
- 4.2 Prepare slants in clean 18 mm diameter rimless test tubes and pre-incubate for 48 hrs at 30-35°C to check any contamination.
- 4.3 Boil the media till get completely dissolve. Cool to about 42 - 45°C and adjust the required pH. Fill each test tube with about 9 to 10ml of the medium & plug the tubes with non absorbent cotton & sterilize by autoclaving.
- 4.4 After autoclaving, prepare the solutions by slanting the test tubes to desirable angle.
- 4.5 After the medium get solidified, keep the prepared slants into the incubators (37°C for 48 hrs.) for pre incubation.
- 4.8 Prior to sub culturing incubation is necessary because to check for the contamination of slants occurs in process.
- 4.9 After incubation, clean the exterior surface of test tube with IPA 70% solution.
- 4.10 Mark each test tube with the name of organism & date of subculture.
- 4.11 At every passage from mother culture check the Gram character and morphology of the culture. In any case the passage from mother culture shall not exceed 5.
- 4.12 Subculture should be done in the Laminar flow clean air station & in between the two gas burners to avoid the chances of cross contamination. Take a loop full of culture from previous stock cultures or mother culture which ever is applicable, and inoculate in freshly prepared sterile slant



**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Microbiology	<b>SOP No.:</b>
<b>Title:</b> Maintenance of Cultures	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

- 4.13 Inoculate each organism in two slants (prepare two slants of each organism). After inoculating, incubate all the tubes at required temperature (i.e. 30-35°C for bacteria and 20-25°C for fungus).
- 4.13 After incubation check the cultures visibly for any contamination.
- 4.15 Mark the Newly prepared lot of cultures tubes as **WORKING CULTURES** and one as **STOCK CULTURE**.
- 4.16 After observing the growth, keep the newly prepared cultured tubes in clean double plain polythene bags and preserved the tubes in the refrigerator at 2 – 8 °C.
- 4.17 Prepare new slants from the **STOCK CULTURE** of previous month (at each month, as a new passage) and should proceed from more than four passage.
- 4.18 After completion of two passages from a single **STOCK CULTURE**, take the passage from **MOTHER CULTURE** slant, proceed for the next passage, by preparing **STOCK CULTURE** and **WORKING CULTURE**.
- 4.19 Enter the date of sub-culturing and passage number in the format as mentioned in Annexure I for sub-culturing of micro-organism.
- 4.20 After every one year procure the new culture from any national recognized institution with certificate for authenticity of the cultures.
- 4.21 The working cultures shall be used for preparing suspensions for positive control; for Growth Promotion test etc.
- 4.22 The **STOCK CULTURE** shall be used only for sub culturing.
- 4.23 After one month when the new cultures are ready for use destroy the old cultures.
- 4.0 **DISPOSAL OF CULTURES :**
- 5.1 Aseptically pipette 10ml of disinfectant solution in the culture tubes.
- 5.2 Sterilise at 121 °C for 20 mins.
- 5.3 After sterilization collect the media in the double polybag, pour sufficient disinfectant solution and tie it properly. Dispose this bag in Incinerator.
- 6.0 **MEDIA USED FOR PREPARING SLANTS:**
1. For Bacteria : Soyabean Casein digest Agar.
  2. For Fungi : Sabaroud Dextrose Agar.



**PHARMA DEVILS**  
MICROBIOLOGY DEPARTMENT

**STANDARD OPERATING PROCEDURE**

<b>Department:</b> Microbiology	<b>SOP No.:</b>
<b>Title:</b> Maintenance of Cultures	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

**7.0 ABBREVIATIONS:** NIL

**8.0 REFERENCE**  
USP and Recommendations by culture supplier.

**9.0 ANNEXURE :** Annexure I

