



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

Title	:	Comparative Study of Media for Microbiological Limit Test of Water			
Protocol No.	:		Page No.	:	1 of 11
Report No.	:				

VALIDATION PROTOCOL
COMPARATIVE STUDY OF MEDIA
FOR MICROBIOLOGICAL LIMIT TEST
OF WATER

Location	:	FORMULATION UNIT
Protocol No.	:	
Supersedes	:	NONE
Issue Date	:	
No. of Pages	:	7 Pages and 4 Annexures



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

Title : **Comparative Study of Media for Microbiological Limit Test of Water**

Protocol No. : **Page No.** : **2 of 11**

Report No. :

CONTENTS

Point No.	Section Title	Page No.
	Protocol Approval Sheet	2
1.0	Objective	3
2.0	Scope	3
3.0	Responsibility	3
4.0	Validation Team Members	3
5.0	Abbreviations	3
6.0	Safety Considerations	4
7.0	Pre-requisites for validation	4
8.0	Acceptance Criteria	5
9.0	Procedure	5
10.0	Revalidation	6
11.0	Deviations and Investigations	6
12.0	Validation Report	6
13.0	List of Annexures/formats attached	6
14.0	References	7



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

Title : **Comparative Study of Media for Microbiological Limit Test of Water**

Protocol No. : **Page No.** : **3 of 11**

Report No. :

PROTOCOL APPROVAL SHEET

Prepared by

Functional Area	Name	Designation	Signature	Date
Microbiology				

Checked by

Functional Area	Name	Designation	Signature	Date
In charge, Microbiology				
Quality Assurance				

Approved by

Functional Area	Name	Designation	Signature	Date
Head, Quality Control				
Head, Quality Assurance				



Title	:	Comparative Study of Media for Microbiological Limit Test of Water			
Protocol No.	:		Page No.	:	4 of 11
Report No.	:				

1.0 Objective

The objective of this study is to determine the suitable media for microbiological limit test of water, with respect to use three different types of basic media to ensure that the capability to support the microbial recovery from water samples, when operated as per the pre approved analysis method.

2.0 Scope

The protocol shall be applicable at formulation Plant of The validation study of the comparative study of media for microbiological limit test of water shall be carried out for one time to check the efficiency of media and analysis method is suitable for recovery of microbial load.

3.0 Responsibility

Microbiology : Preparation and review of validation protocol; conducting the experiment as per the approved protocol and compilation of data; preparation of validation report.

Quality Assurance : Review of validation protocol, results and report.

Head, QC : To check the protocol with respect to its intended purpose and to make evaluations on compiled data from the test; Final approval of protocol and report.

Head, QA : Review of protocol for the correctness and adequacy of the text and the experiment, regulatory compliance; final approval of protocol. Review of report for the correctness of the data and compliance of the protocol and approval of the report.

4.0 Validation Team members

Validation team shall comprise of the representatives from following functions:

- Microbiology
- Quality Assurance



Title	:	Comparative Study of Media for Microbiological Limit Test of Water		
--------------	---	---	--	--

Protocol No.	:		Page No.	:	5 of 11
---------------------	---	--	-----------------	---	----------------

Report No.	:	
-------------------	---	--

5.0 Abbreviations

SOP	: Standard Operating Procedure
SCS	: Standardized Cell Suspension
LAF	: Laminar Air Flow
CFU	: Colony Forming Unit
TSA	: Tryptone Soy Agar
PCA	: Plate Count Agar

6.0 Safety Considerations

- 6.1 All the procedures should be carried out aseptically.
- 6.2 All the materials that are to be used should be sterile.
- 6.3 All the materials used for handling the microorganisms should be decontaminated.

7.0 Pre-requisites for Validation

- 7.1 The following pre-requisites are to be used for purified water testing with respect to the sample size for Microbiological evaluations.
 - 7.1.1 Standardized cell suspension of *Escherichia coli*, *Salmonella*, *Pseudomonas aeruginosa*, *Staphylococcus aureus* and *Bacillus subtilis* (SCS)
 - 7.1.2 Tryptone Soy Agar
 - 7.1.3 R2A Agar
 - 7.1.4 Plate Count Agar
 - 7.1.5 Glass Bottles (100 ml & 500 ml capacity)
 - 7.1.6 Measuring cylinders (10 & 100 ml capacity)
 - 7.1.7 Normal saline
 - 7.1.8 70 % Iso Propyl Alcohol
 - 7.1.9 Petri plates (90 mm)
 - 7.1.10 Colony counter
 - 7.1.11 Calibrated Incubators
 - 7.1.12 Calibrated Compound Microscope
 - 7.1.13 Validated Autoclave
 - 7.1.14 Validated LAF unit



Title	:	Comparative Study of Media for Microbiological Limit Test of Water		
--------------	---	---	--	--

Protocol No.	:		Page No.	:	6 of 11
---------------------	---	--	-----------------	---	----------------

Report No.	:	
-------------------	---	--

7.1.15 Membrane Filtration Apparatus

7.1.16 0.45 μ membrane filters

7.1.17 Forceps

7.2 Documents

SOPs of the relevant instruments and testing procedures.

8.0 Acceptance Criteria

8.1 The estimated number of cells in positive control shall be less than 100 for the Microbial Limit Test.

8.2 The recovery shall be not less than 70% for the CFU obtained from positive control.

8.3 All the three media should pass 'Growth Promotion Test for Microbiological Culture Media'.

8.4 The negative control should not show any growth.

9.0 Procedure

9.1 Collect the samples as per routine sampling plan and follow the procedure as per SOP on 'sampling procedure for water samples'.

9.2 Transfer the samples to the microbiology laboratory for analysis in a closed condition.

9.3 Analyze the samples in duplicate mode followed by Membrane filtration technique with a sample quantity of 10.0 ml, in triplicate as follows;

9.3.1 Filter the sample (Aseptically collected in sterile measuring cylinder from sample bottle), through sterilized 0.45 μ membrane filter using sterilized filtration assembly with vacuum pump.

9.3.2 Rinse the membrane with 100 ml of sterile saline solution (Sterile normal saline shall be prepared by dissolving sodium chloride in purified water in the ratio of 9.0 g to 1000 ml respectively).

9.3.3 Aseptically remove the membrane with the help of a sterile forceps and put on a sterile pre incubated Tryptone Soy Agar plate.

9.4 Analyze the same samples followed by point no. 9.3 to 9.3.2, aseptically remove the membrane with the help of a sterile forceps and put on a sterile pre incubated R2A Agar plate & follow the same procedure for PCA.

9.5 Incubate these Petri plates at $32.5 \pm 2.5^{\circ}\text{C}$ for 72 hours and observe for every 24 hours.



Title	:	Comparative Study of Media for Microbiological Limit Test of Water			
Protocol No.	:		Page No.	:	7 of 11
Report No.	:				

- 9.6 Pipette out 1.0 ml of each SCS (10 – 100 cfu) and add into 10 ml of sterile normal saline & filter this solution through sterilized 0.45 μ membrane filter using sterilized filtration assembly with vacuum pump.
- 9.7 Rinse the membrane with 100 ml of sterile saline solution.
- 9.8 Aseptically remove the membrane with the help of a sterile forceps and put on a sterile pre incubated Tryptone Soy Agar plate for Positive Control on TSA & incubate along with samples & observe for every 24 hours.
- 9.9 Follow the same procedure from point no. 9.6 to 9.8 for R2A Agar and PCA.
- 9.10 Keep one no. of each TSA, R2A Agar and PCA plates as negative control and incubate along with samples & observe for every 24 hours.
- 9.11 Record the details in given Data Sheet as per Annexure – 3.
- 9.12 Follow the same procedure for consecutive 30 days and monitor the trend.
- 9.13 After completion of this study, report shall be prepared on the basis of results obtained at 72 hours incubation.
- 9.14 Based on the test results the optimum media usage shall be confirmed and shall be recommended for the routine practice.

10.0 Re-validation

- 10.1 In the view of following listed, the test shall be re-conducted if there is:
 - 10.1.1 Any major change in the specification of purified water.
 - 10.1.2 Any major change in the testing procedure of purified water.

11.0 Deviations and Investigations

- 11.1 Any deviation to this protocol and thereupon investigation shall be recorded as per SOP.

12.0 Validation Report

- 12.1 Based on the outcome from this validation study, a report shall be prepared by Microbiology and Quality Assurance. The validation report shall be reviewed and then approved by all functional heads of all the concerned departments. Validation Report shall include following:
 - 12.1.1 Cover page of the Report (as per Annexure-1).
 - 12.1.2 Validation Report Approval Sheet (as per Annexure-2).



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

Title : **Comparative Study of Media for Microbiological Limit Test of Water**

Protocol No. : **Page No.** : **8 of 11**

Report No. :

12.1.3 Data Sheets (as per Annexure-3).

12.1.4 Validation Report summary (as per Annexure-4).

13.0 List of Annexures/Formats attached

1.	Cover page of the Report	Annexure-1
2.	Validation Report Approval Sheet	Annexure-2
3.	Data Sheet	Annexure-3
4.	Validation Report Summary & Conclusion	Annexure-4

14.0 References

14.1 In – House Method.



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

Title	:	Comparative Study of Media for Microbiological Limit Test of Water				
Protocol No.	:		Page No.	:		
Report No.	:					

DATA SHEET

Sampling Point Name:

Sampled on:

Analyzed on:

Method of Analysis: Membrane Filtration Technique

Autoclave I.D. No.:

Media Lot No.:

Sampling Point No.:

Sampled By:

Analyzed By:

Analyzed Sample Quantity: 10 ml

LAF I.D. No.:

Incubator I.D. No.:

Date	Incubation Time	Observation (Sample)															Observed By	Remarks
		R2A Agar					Tryptone Soy Agar					Plate Count Agar						
		Plate-1	Plate-2	Avg.	+ ve	- ve	Plate-1	Plate-2	Avg.	+ ve	- ve	Plate-1	Plate-2	Avg.	+ ve	- ve		
	After 24 hours																	
	After 48 hours																	
	After 72 hours																	

Checked By:

Date:



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

Title : **Comparative Study of Media for Microbiological Limit Test of Water**

Protocol No. : **Page No.** :

Report No. :

Summary & Conclusion:

The comparative study of media for Microbial Analysis of Purified Water was carried out for consecutive 30 days as per the routine sampling schedule, using R2A agar, Plate Count Agar and Soyabean Casein Digest Agar by Membrane Filtration Technique.

During the study it has been observed that the microbial counts are within the minimum and maximum range as given below.

Name of Medium	No. of Counts Observed (cfu)	
	Maximum	Minimum
Soyabean Casein Digest Agar (TSA)	25	54
R2A Agar	22	49
Plate Count Agar	22	55

The results observed in all three mediums are almost equal. Hence it is concluded that Soyabean Casein digest Medium is the suitable medium for analysis of water so Soyabean Casein Digest Agar shall be continued as earlier it being used.

Prepared By Date	Checked By Date
(Microbiologist)	(Quality Assurance)



PHARMA DEVILS

MICROBIOLOGY DEPARTMENT

Title	:	Comparative Study of Media for Microbiological Limit Test of Water		
Protocol No.	:		Page No.	:
Report No.	:			

REPORT APPROVAL SHEET

Prepared by

Functional Area	Name	Designation	Signature	Date
Microbiology				

Checked by

Functional Area	Name	Designation	Signature	Date
Microbiology				
Quality Assurance				

Approved by

Functional Area	Name	Designation	Signature	Date
Head, Quality Control				
Head, Quality Assurance				