PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE			
Department: Production	SOP No.:		
Title: Cleaning and Filling of Isopropyl Alcohol Bottles	Effective Date:		
Supersedes: Nil	Review Date:		
Issue Date:	Page No.:		

1.0 OBJECTIVE:

To lay down a procedure for Cleaning and Filling of Isopropyl Alcohol Bottles.

2.0 SCOPE:

This SOP is applicable for Cleaning and Filling of Isopropyl Alcohol Bottles used in to the Production Area at

3.0 RESPONSIBILITY:

Officer / Executive – Production

4.0 ACCOUNTABILITY:

Head – Production

5.0 **DEFINITIONS:**

Not Applicable

6.0 PROCEDURE:

- **6.1** Isopropyl Alcohol Bottles shall be cleaned with Purified Water, if required using scrubber or brush for cleaning the Isopropyl Alcohol Bottles.
- **6.2** Finally rinse the Isopropyl Alcohol Bottles with Water for Injection.
- **6.3** The Isopropyl Alcohol Bottles shall be allowed to air dry.
- **6.4** Finally mop the dried Isopropyl Alcohol Bottles with 70% IPA from outside.
- 6.5 After cleaning, area supervisor shall check the Isopropyl Alcohol Bottles for satisfactory cleaning.
- **6.6** After satisfactory cleaning fill the Isopropyl Alcohol in IPA Bottles.
- **6.7** Record the cleaning details in "**IPA Bottles Cleaning and Filling Record**" as per format shown in **Annexure-I.**
- **6.8** Frequency: Daily



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCEDURE		
Department: Production	SOP No.:	
Title: Cleaning and Filling of Isopropyl Alcohol Bottles	Effective Date:	
Supersedes: Nil	Review Date:	
Issue Date:	Page No.:	

7.0 ABBREVIATIONS:

WFI Water for Injection

Ltd. Limited
No. Number

QA Quality Assurance.

IPA Isopropyl Alcohol

8.0 ANNEXURES:

ANNEXURE No.	TITLE OF ANNEXURE	FORMAT No.
Annexure-I	IPA Bottles Cleaning and Filling Record	

9.0 **DISTRIBUTION**:

• Controlled Copy No.01 Production

• Controlled Copy No.02 Quality Assurance

• Master Copy Quality Assurance

10.0 REFERENCES:

Not Applicable

11.0 REVISION HISTORY:

	Revision No.	Change Control	Details of	Reason of	Effective Date	Done By
l		No.	Changes	Changes		
	00	Not Applicable	Not Applicable	New SOP		



PHARMA DEVILS

PRODUCTION DEPARTMENT

STANDARD OPERATING PROCE	DURE
Department: Production	SOP No.:
Title: Cleaning and Filling of Isopropyl Alcohol Bottles	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:
ANNEXURE-I	
IPA BOTTLE FILLING AND CLEA	ANING RECORD

Department: Frequency: Daily

Date	Area	Bottle Cleaning Date	Cleaned By Sign & Date	IPA Filled By Sign & Date	Checked By Sign & Date